

PHILIPS

Allegro

Enhancing PET to a new level of performance

Johann Fernando, Ph.D.
Director, PET Science

ANZSNM Physics Workshop 2002




Outline

- Allegro technology
 - Clinical goals of PET
 - 3D Imaging to meet clinical goals
 - GSO ideal for 3D imaging
 - Detector design to optimize performance
- Cs-137 attenuation correction benefits
- 3D RAMLA benefits


**If you cannot see the lesion
you cannot localize it**

Clinical Goal

***Image the smallest possible lesion with a low contrast
in a background of radioactivity***





Large number of counts in the image – high sensitivity



3D imaging maximizes Sensitivity

Allegro is the only PET System designed for 3D imaging from the beginning

- 3D imaging maximizes the sensitivity



3D Imaging provides 5 times Sensitivity vs. 2D

- 2D Imaging is limited by the Septa for FDG Whole-Body scans

	2D	3D
Sensitivity	1x	5x
Count rate	1x	5x
Throughput	1 hour	~ 30 mins.
Scatter	20%	40%

Source: Moszynski et al. IEEE NS-45: 472-477, 1998

3D Imaging Requirements

Feature

- Fast decay time
- Good energy resolution
- High stopping power
- Good timing resolution

Benefit

- Permits high singles rate
- Limits randoms and scatter
- Gives high coinc. fraction
- Reduces randoms

These 3D imaging requirements are fulfilled by the



GSO - Fast, High Density Scintillator


GSO is the Material of Choice for 3D PET

	<u>BGO</u>	<u>NaI</u>	<u>LSO</u>	<u>GSO</u>
Decay Time	300	230	40	60
Deadtime (nsec)	700	240	100	120
Energy Res. (%)	15	8	12	8
Stopping Power	11.6	30.7	12.3	16.0
Light Output	15%	100%	75%	25%

- GSO combines high stopping power of BGO with excellent energy resolution of NaI
- Fast decay time of GSO is necessary to keep up with the increased photon flux in 3D and reduce detector deadtime

Source: Moszynski et al. IEEE NS-45: 472-477, 1998

GSO Improves Systems Reliability due to High Temperature Stability



- Temperature stability of GSO is much better than LSO or BGO
- LSO and BGO based scanners require special water cooling systems
- Makes GSO the material of choice for combined PET/CT tomographs

Moszynski et al. IEEE NS-45: 472-477, 1998

System down time for detector calibrations is minimized

LSO systems have compromised Transmission capabilities

- Coincidence transmission imaging is a poor second choice
 - Ge-68 source
 - Source needs to be replaced about once/year
 - Maximum data rate limited to maximum true rate, i.e ~100-400 Kcps
- Long acquisition times for Transmission scans due to coincidence mode acquisitions
- Poor image quality of the transmission images due to radioactive background from LSO
 - Radioactivity increases the background counts in Transmission images
 - particularly high attenuation regions

Source: Radioisotope decay data; SNM proceedings (2001); internal analysis

GSO Advantages versus LSO/BGO

- Excellent energy resolution
- Short decay time
- Good stopping power
- Very good temperature stability
- Not radioactive

Optimum for
3D Imaging
and PET/CT

The design of the system is as important as the selection of the scintillator

Detector Designs


- Different crystal/PMT schemes are used in designing PET scanners
- Trade-off between count rate capability, system energy resolution and complexity
 - Continuous Positioning Detector Design
 - Block Design

Continuous Positioning Detector Technology optimizes use of GSO for 3D imaging

Combining single crystals with continuous light guide

- Much better system **energy resolution**
- High **count rate capability**
- Small crystal size for **excellent spatial resolution**
- Better **3D imaging capabilities**


PIXELAR™
Continuous pixelated PET detectors



Sources: Adam, L.E., Watson, C.C., Nuklearmedizin 1999,38:61-65; Cherry S.R. et.al. IEEE NS Vol. 42 (4),1064 – 1068, 1995; Philips measurements


Allegro PIXELAR Detector vs Competitor Block Detectors

Continuous Positioning PIXELAR Detector - Allegro



- New technology
- Designed for 3D imaging
- Variation of light collection <20%
- Excellent energy resolution: <18%

Block Detector All competitor systems




- > 20 year old technology
- Designed for 2D imaging
- Variation of light collection ~300%
- Poor energy resolution: 25-35%

Crystal & detector design determines system energy resolution

Siemens Accel

LSO

10-13 %^{1,2)}



Block

~300 %⁴⁾

Crystal Type

Crystal Energy Resolution

+

Detector Design


=

Light collection variation

Allegro

GSO

8 %⁵⁾



Continuous pixelated PET detectors

PIXELAR

< 20 %⁵⁾

>25 %³⁾

System Energy Resolution


< 18 %⁵⁾*

1) A. Saoudi et al., IEEE NS 46(6), 1925-1928,1999 , 2)Moszynski et al. IEEE NS-45(3), 472-477, 1998
3) Adam, L.E., Watson, C.C., Nuklearmedizin 1999,38:61-65; 4) Cherry S.R. et.al. IEEE NS Vol. 42 (4),1064 – 1068, 1995


⁵⁾ Philips measurements

Good Energy Resolution is Essential to Reduce Randoms

Singles rates are very high in 3-D mode



Singles and coincidence energy spectra for 20 cm Ø, 70 cm long phantom




Relative Randoms Fraction

Randoms Fraction decreases exponentially with increased LLD

An LLD setting of 410 keV on the Allegro enables excellent image quality

Good Energy Resolution is Essential to Reduce Scatter




Scatter fraction as a function of LLD for 20 cm Ø, 20 cm long phantom with line source

Scatter Fraction decreases linearly with increased LLD

An LLD setting of 410 keV on the Allegro enables excellent image quality

Allegro is Designed for 3D Imaging to Provide High Image Quality




Allegro design allows use of Cs-137 for transmission for **better throughput, image quality and cost**

Allegro (Cs-137 source)	LSO & BGO Block Detectors (Ge-68 coincidence source)
<p>High count rates in excess of 3Mcps</p> <p>Fast Acquisition times (45 sec)</p> <p>High quality transmission images</p> <p>No need to replace source (30 yr half life)</p>	<p>10X lower count rates than Cs-137</p> <p>Long acquisition times</p> <p>Need for image segmentation</p> <p>Need to replace source annually (271 day half life -> \$10K/year)</p>

ALLEGRO™

Example of Cs-137 transmission image

- 43 sec. acquisition per step
- Transmission count rates > 3 Mcps



Ability to see bone structure and the Pelvis provide the anatomical landmarks required to perform image fusion accurately

Images courtesy of J. Karp, Ph.D.
University of Pennsylvania

Why 3D RAMLA?

Acquire in 3D -> Reconstruct in 3D for best image quality



3D RAMLA Benefits

- Accurate reconstruction for 3D imaging
- Improves image quality
- Reduces artifacts
- Minimizes noise

Fully 3D iterative recon improves image quality

FORE + OSEM

3D RAMLA




University of Pennsylvania PET Center


- No bladder artifact
- Sharper edges and definition of the Liver
- Uniform distribution
- Much better definition at the edge of the FOV

Fully 3D iterative recon improves image quality

FORE + OSEM



3D RAMLA



University of Pennsylvania PET Center


- No bladder artifact
- No artifacts
- Better definition of structure

Fully 3D iterative recon improves spatial resolution

University of Pennsylvania PET Center

2D

FORE + OSEM




3D

3D RAMLA

Much better definition of structure

Allegro High Performance PET System



Allegro Key Benefits

- **High quality PET images** through 3D mode data acquisition that optimizes sensitivity
- Excellent energy resolution of the PIXELAR detector technology provides optimum contrast resolution to **detect small low contrast lesions**
- High Density, fast scintillator **GSO for high patient throughput**
- **Accurate full 3D image reconstruction** to visualize low contrast lesions provided by the unique 3D RAMLA software
- **Seamless integration with Radiation Oncology** through **Syntegra** image fusion and multi modality connectivity