



Acid Labile Subunit ELISA

For the quantitative determination of acid labile subunit (ALS) in human serum and plasma (EDTA, heparin, citrate).

For Research use Only. Not For Use In Diagnostic Procedures.

Catalog Number: 22-ALSHU-E01

Size: 96 Wells

Version: 10.03.2016 Version 15- ALPCO September 6, 2018

1. Intended Use

This enzyme immunoassay kit is suited for measuring Acid Labile Subunit (ALS) in human serum or EDTA/heparin/citrate plasma for research use only. **For research use only. Not for use in diagnostic procedures.**

2. Introduction

The Insulin-like Growth Factors (IGF) – I and II are bound to specific binding proteins in circulation (IGFBP). Until today seven different proteins have been identified IGFBP-1 to 7 [1, 2]. IGF bioavailability, transport and storage is regulated or facilitated by these binding proteins which are expressed differentially according physiological and developmental requirements. The most abundant IGFBP in circulation is IGFBP-3. Together with IGFBP-5 it is able to form the so called ternary complex with IGF and the acid labile subunit (ALS) [3-5]. In the circulation nearly all IGF is bound in this ternary complex and thus not able to cross the endothelial barrier. Only very small amounts of IGF or IGFBP-3 exist outside this complex [6, 7]. The acid-labile subunit is an important part of the IGF-storage mechanism in circulation. Research demonstrates that in ALS deficiency or in ALS knock-out mice the concentration of IGF and IGFBP-3 in the circulation is significantly decreased resulting in impaired growth [10].

The acid-labile subunit, is a synthesized as propeptide of 605 amino acids. The signal peptide, necessary for ALS secretion (AA 1-27) is cleaved off during the transport process (Swiss-Prot P35858 Version 82). The mature protein consists of 578 amino acids and contains about 20 leucin rich sequence repeats. Beside the leucin-rich repeats several potential N-linked glycosylation sites have been described. Miller BS et al. were able to demonstrate that incomplete glycosylation of IGFs, ALS and IGFBP-3 results in a decreased serum concentration of these proteins. Mutations in or the complete knock out of the ALS gene result in IGF/IGFBP-3 deficiency and therewithin disturbances of growth [9,10].

The first ALS immunoassay was described by Baxter RC in 1990 [6]. Research with this in-house radioimmunoassay showed that ALS is present in high concentrations in serum of healthy subjects, but not detectable in other body fluids like amniotic fluid, cerebrospinal fluid or seminal plasma, in spite of the fact that these body fluids contain high levels of IGFBP-3.

3. Principle of The Assay

The ALS ELISA is a so-called Sandwich-Assay. It utilizes two specific and high affinity antibodies for this protein. These antibodies were created by immunization of rabbits with specific peptides as previously described by Khosravi and Stadler [16, 17].

The ALS in the sample binds to the immobilized first antibody on the microtiter plate. In the following step, the biotinylated second anti-ALS-Antibody binds in turn to the immobilized ALS. Next the Streptavidin-POD-Conjugate binds to the biotinylated antibody and in the closing substrate reaction the turn of the color is catalyzed, quantitatively depending on the ALS-level of the samples.

Initially the assay was calibrated against an internal serum and measurement results were expressed as mU/mL. After successful production of eukaryotic recombinant ALS the calibration was transferred to mass units (see calibration/traceability). Additionally, recombinant material was used to quantify the ALS content of the calibrators in mass units. Analysis revealed that 1mU ALS is equivalent to 5 ng ALS and all previous assay data describing the assay performance was accordingly transferred to ng/mL.

4. Warnings and Precautions

For Research Use Only. Not for Use in Diagnostic Procedures.

Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood. ALPCO is not liable for any loss or harm caused by noncompliance with the instructions provided. A Safety Data Sheet is available on ALPCO's website and by request.

Do not use damaged or microbially contaminated or spilled material.

Caution: This kit contains material of human and/or animal origin. Therefore, all components and samples should be treated as potentially infectious.

Appropriate precautions and good laboratory practices must be used in the storage, handling, and disposal of the kit reagents. The disposal of the kit components must be made according to the applicable regulations.

Human Serum

The following components contain human serum: **Control Serum KS1, KS2 and Standards A-F**. The sources of human sera for the control sera provided in this kit was tested and found non-reactive for Hepatitis-B surface antigen (HBsAg), Hepatitis C virus (HCV), and Human Immunodeficiency Virus 1 and 2 (HIV) antibodies. No known test methods can offer total assurance of the absence of infectious agents. Therefore, all components and samples should be treated as potentially infectious.

Reagents AK, VP, WP

Contain as preservative a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (<0.015%).

H317	May cause an allergic skin reaction.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.
P272	Contaminated work clothing should not be allowed out of the workplace.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.
P333+P313	If skin irritation or rash occurs: Get medical advice/ attention.
P302+P352	IF ON SKIN: Wash with plenty of soap and water.
P501	Dispose of contents/container in accordance with local/regional/national/ international regulations.

Substrate Solution (S)

The TMB-Substrate (S) contains 3,3',5,5' tetramethylbenzidine (<0.05%).

H315	Causes skin irritation.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
P261	Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Stop Solution (SL)

The Stop solution contains 0.2 M acid sulfuric acid (H₂SO₄).

H290	May be corrosive to metals.
H314	Causes severe skin burns and eye damage.
P280	Wear protective gloves/ protective clothing/ eye protection/ face protection.

P301+P330+	IF SWALLOWED: rinse mouth.
P331	Do NOT induce vomiting.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P309+P310	IF exposed or feel unwell: Immediately call a POISON CENTER or doctor/physician.

4.1. General first aid procedures:

Skin contact: Wash affected area, rinse immediately with plenty of water for at least 15 minutes. Remove contaminated clothes and shoes.

Eye contact: In case of contact with eyes, rinse immediately with plenty of water for at least 15 minutes; while rinsing spread the eyelids.

Ingestion: If swallowed and if the affected person is conscious, rinse out the mouth with plenty of water; seek medical advice immediately.

5. Sample

5.1. Sample Type

Serum and Plasma

30 IU/mL Sodium Heparin, 3.8 g/L Sodium Citrate or 0.0068 mol/L EDTA did not interfere with ALS measurement.

5.2. Sample Collection

The blood sample for serum preparation should be obtained according to standardized venipuncture procedures. Hemolytic reactions have to be avoided.

5.3. Required sample volume:

- 10 µL

5.4. Sample Stability

In firmly closable sample vials:

- Storage at 20-25°C: 3 days
- Storage at -20° C: min. 2 years
- Freeze-thaw cycles: max. 5

The storage of samples over a period of 2 years at -20°C, showed no influence on the reading. Freezing and thawing of samples should be minimized, 5 freeze-thaw cycles showed no effect on the measured ALS concentration.

5.5. Interference

Hemoglobin, triglyceride and bilirubin in the sample do not interfere to a concentration of 1 µg/mL, 100 mg/mL and 200 µg/mL, respectively. However, the use of hemolytic, lipemic or icteric samples should be validated by the user.

5.6. Sample Dilution


- Dilution: 1:150 with Sample Buffer PP
- Pipette 1490 μ L Sample Buffer PP (red colored) in PE-/PP-Tubes (application of a multi-channel pipette is recommended in larger series), add 10 μ L Serum or Plasma (dilution 1:150) and mix each tube immediately. After mixing use 50 μ L of this solution within 1 hour per determination in the assay (pipetting control = red coloring of the solution in the wells).
- Sample stability after dilution of the sample: maximum 1 hour at 20-25°C.
- In most determinations (e.g. Serum or Plasma samples and no extreme values expected), the dilution of 1:150 with Sample Buffer PP is suitable, respectively the assay covers the range from 0.53 ng/mL - 30 μ g/mL where 1:50 is the minimal tested sample dilution.
- If required, the dilution with Sample Buffer PP could be performed at lower or higher dilutions.

6. Materials

6.1. Materials Provided

The reagents listed below are sufficient for 96 wells including the standard curve.

MTP	Microtiter plate , ready for use, coated with rabbit-anti-ALS-antibody. Wells are separately breakable.	(8x12) wells
PP	Sample Buffer PP , ready for use, red colored, please use for the reconstitution of Standards A-F and Controls KS1/KS2 and for the dilution of Samples and Controls KS1/KS2 .	1 x 125 mL
A-F	Standards , lyophilized (native human ALS), concentrations are given on vial labels and on quality certificate in ng/mL.	6 x 1 mL
KS1	Control Serum 1 , lyophilized (human serum), concentration is given on quality certificate in ng/mL.	1 x 250 μL
KS2	Control Serum 2 , lyophilized (human serum), concentration is given on quality certificate in ng/mL.	1 x 250 μL
VP	Dilution Buffer , ready for use, please use for the dilution of Antibody Conjugate AK .	1 x 7 mL
AK	Antibody Conjugate, 50-fold concentrate , contains the biotinylated anti-human rabbit ALS Antibody. Dilute before use 1:50 in Dilution Buffer VP and use 50 μL for each well in the assay. Attention: Please dilute Antibody Conjugate AK freshly according to daily requirements.	1 x 140 μL
EK	Enzyme Conjugate EK , ready-to-use, contains HRP (Horseradish-Peroxidase)-labeled Streptavidin.	1 x 12 mL
WP	Washing Buffer , 20-fold concentrated solution, dilute 1:20 in distilled or deionized water.	1 x 50 mL
S	Substrate , ready for use, horseradish-peroxidase-(HRP) substrate, stabilized tetramethylbenzidine.	1 x 12 mL
SL	Stop Solution , ready for use, 0.2 M sulfuric acid.	1 x 12 mL
-	Sealing Tape , for covering the microtiter plate .	2 x

	Instructions for use	1 x
--	Quality Certificate	1 x

6.2. Materials Required, but Not Provided

- Distilled or deionized water for dilution of the Washing Buffer **WP (distilled water)**, 950 mL.
- Precision pipettes and multichannel pipettes with disposable plastic tips
- Polyethylene PE/Polypropylene PP tubes for dilution of samples
- Vortex-mixer
- Microtiter plate shaker (350 rpm)
- Microtiter plate washer (recommended)
- Microplate reader ("ELISA-Reader") with filter for 450 and ≥ 590 nm
- Timer

7. Technical Notes

Storage Conditions

Store the kit at 2-8°C after receipt until its expiry date. The lyophilized reagents should be stored at -20 °C after reconstitution. Avoid repeated thawing and freezing.

Storage Life

The shelf life of the components **after initial opening** is warranted for **4 weeks**. Store the unused strips and microtiter wells **airtight** together with the desiccant at 2-8°C in the clip-lock bag, use in the frame provided. The **reconstituted components** standards **A-F** and Control Sera **KS1 and KS2** must be stored at -20°C (max. 4 weeks). For further use, thaw quickly but gently (avoid temperature increase above room temperature and avoid excessive vortexing). Up to 3 of the freeze-thaw cycles did not influence the assay. The 1:20 diluted Washing Buffer **WP** is 4 weeks stable at 2-8°C

Preparation of reagents

Bring all reagents to room temperature (20 – 25°C) before use. Possible precipitates in the buffers have to be resolved before usage by mixing and/or warming. Reagents with different lot numbers cannot be mixed.

Reconstitution

The Standards **A – F** and Control Sera **KS1** and **KS2** are reconstituted with the Sample Buffer **PP**. It is recommended to keep reconstituted reagents at room temperature for 15 minutes and then to mix them thoroughly but gently (no foam should result) with a vortex mixer.

Dilution

After reconstitution **dilute** the Controls **KS1** and **KS2** with the Sample Buffer **PP** in the same ratio (**1:150**) as the sample.

The required volume of **Antibody Conjugate AK** is prepared by **1:50** dilution of the provided 50-fold concentrate with **Dilution Buffer VP**. Please dilute Antibody Conjugate freshly according to daily requirements.

The required volume of Washing Buffer **WP** is prepared by 1:20 dilution of the provided 20-fold concentrate with distilled or deionized water.

Assay Procedure

When performing the assay Standards **A-F**, Controls **KS1**, **KS2** and the samples should be pipetted as fast as possible (e.g. <15 minutes). To avoid distortions due to differences in incubation times, Enzyme Conjugate **EK** as well as the succeeding Substrate Solution **S** should be added to the plate in the same order and in the same time interval as the samples. Stopping Solution **SL** should be added to the plate in the same order as Substrate Solution **S**.

All determinations (Standards **A-F**, Controls **KS1**, **KS2** and samples) should be assayed in duplicate. For optimal results, accurate pipetting and adherence to the protocol are recommended.

Incubation

Incubation at room temperature means: Incubation at 20 - 25°C. The Substrate Solution **S**, stabilized H₂O₂-tetramethylbenzidine, is photosensitive – store and incubate in the dark.

Shaking

The incubation steps should be performed at mean rotation frequency of a particularly suitable microtiter plate shaker. We are recommending 350 rpm. Due to certain technical differences deviations may occur, in case the rotation frequency must be adjusted. Insufficient shaking may lead to inadequate mixing of the solutions and thereby to low optical densities, high variations and/or false values, excessive shaking may result in high optical densities and/or false values.

Washing

Proper washing is of basic **importance** for a secure, reliable and precise performance of the test. Incomplete washing is common and will adversely affect the test outcome. Possible consequences may be uncontrolled unspecific variations of measured optical densities, potentially leading to false results of the examined samples. Effects like high background values or high variations may indicate washing problems.

All washing must be performed with the provided Washing Buffer **WP** diluted to usage concentration. Washing volume per washing cycle and well must be 300 µL at least.

The danger of handling with potentially infectious material must be taken into account.

When using an **automatic microtiter** plate washer, the respective instructions for use must be carefully followed. Device adjustments, e.g. for plate geometry and the provided washing parameters, must be performed. Dispensing and aspirating manifold must not scratch the inside well surface. Provisions must be made that the remaining fluid volume of every aspiration step is minimized. Following the last aspiration step of each washing cycle, this could be controlled, and possible remaining fluid could then be removed, by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

Manual washing is an adequate alternative option. Washing Buffer may be dispensed via a multistepper device, a multichannel pipette, or a squirt bottle. The fluid may be removed by dynamic swinging out the microtiter plate over a basin. If aspirating devices are used, care has to be taken that

the inside well surface is not scratched. Subsequent to every single washing step, the remaining fluid should be removed by inverting the plate and repeatedly tapping it dry on non-fuzzy absorbent tissue.

8. Assay Procedure

Preparation of reagents		Reconstitution:	Dilution
AK	Antibody-conjugate	-	Dilute before use with 1:50 Dilution Buffer VP
A-F	Standards	in 1 mL Sample Buffer PP	-
KS1 and KS2	Control Sera	In 250 µL Sample Buffer PP	1:150 with Sample Buffer PP
WP	Washing Buffer	-	1:20 with distilled water
Sample and Control Sera KS1 &KS2 Dilution: 1:150 in Sample Buffer PP (red colored; e.g. 10 µL in 1490 µL PP). Mix directly and use within max. 60 min . Use 50 µL per determination (pipetting control= red coloration)			
Before assay procedure bring all reagents to room temperature 20-25°C .			
Assay Procedure in Double Determination:			
Pipette	Reagents	Well Position	
50 µL	1:50 diluted Antibody Conjugate	Pipette in <u>all</u> required number of wells	
50 µL	Standard A (0 ng/mL)	A1/A2	
50 µL	Standard B (7.5 ng/mL)	B1/B2	
50 µL	Standard C (31.25 ng/mL)	C1/C2	
50 µL	Standard D (62.5 ng/mL)	D1/D2	
50 µL	Standard E (125 ng/mL)	E1/E2	
50 µL	Standard F (200 ng/mL)	F1/F2	
50 µL	Control Serum KS1 (1:150 diluted)	G1/G2	
50 µL	Control Serum KS2 (1:150 diluted)	H1/H2	
50 µL	Sample (1:150 diluted)	in the rest of the wells according the requirements	
Cover the wells with the sealing tape.			
Sample Incubation: 2 hours at 20-25°C, 350 rpm			
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well	In each well	
100 µL	Enzyme Conjugate EK	In each well	
Cover the wells with the sealing tape.			
Incubation: 30 Minutes at 20-25°C, 350 rpm			
5 x 300 µL	Aspirate the contents of the wells and wash 5 x with 300 µL each Washing Buffer WP/ well	In each well	

100 µL	Substrate Solution S	In each well
Incubation: 30 Minutes in the Dark at 20-25°C		
100 µL	Stop Solution SL	In each well
	Measure the absorbance within 30 min at 450 nm with ≥ 590 nm as reference wavelength.	

9. Quality Control

Good laboratory practice requires that controls are included in each assay. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance. The test results are only valid if the test has been performed following the instructions. Moreover, the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state or local standards/laws. All standards and kit controls must be found within the acceptable ranges as stated on the QC Certificate. If the criteria are not met, the run is not valid and should be repeated. Each laboratory should use known samples as further controls.

9.1. Quality criteria

For the evaluation of the assay it is required that the absorbance values of the Standard A should be below 0.25, and the absorbance of Standard F should be above 1.00.

Samples, which yield higher absorbance values than Standard F, should be re-tested with a higher dilution.

10. Evaluation of Results

10.1 Establishing the Standard Curve

Standard	A	B	C	D	E	F
ng/mL	0	7.5	31.25	62.5	125	200
mU/mL	0	1.5	6.25	12.5	25	40

- 1) Calculate the **mean absorbance** value for the Standard A from the duplicated determination (well A1/A2).
- 2) Subtract the mean absorbance of the Standard A from the mean absorbance of all other samples and standards.
- 3) Plot the standard concentrations on the x-axis versus the mean value of the absorbance of the standards on the y-axis.
- 4) Recommendation: Calculation of the standard curve should be done by using a computer program, because the curve is in general (without respective transformation) not ideally described by linear regression. **A higher-grade polynomial, or four parametric logistic (4-PL) curve fit or non-linear regression** are usually suitable for the evaluation (as might be spline or point-to-point alignment in individual cases).

- 5) The ALS concentration in ng/mL of the samples and controls **KS1** and **KS2** can be calculated by **multiplication** with the respective **dilution factor**.

10.2 Example of a typical standard curve

The following data is for demonstration only and cannot be used in place of data generation at the time of assay.

	A	B	C	D	E	F
ng/mL	0	7.5	31.25	62.5	125	200
OD (450-620 nm)	0.049	0.282	1.248	1.929	2.54	2.934




Figure 1: Example standard curve

The example shown in Figure 1 **cannot** be used for calculation of your test results. A standard curve must be conducted for each assay run conducted.

10.3 Example calculation of ALS concentrations

Sample dilution: 1:150

Measured extinction of your sample	1.5
Measured extinction of the Standard A	0.049

Your measurement program will calculate the ALS concentration of the diluted sample automatically by using the difference of extinction values of sample and Standard A for the calculation. You only have to determine the most suitable curve fit.

In this example, the following equation is solved by the program to calculate the ALS concentration in the sample:

$$1.451 = (-0.184 - 3.09) / (1 + (x / 51.3)^{0.939}) + 3.09$$

$$51.201 = x$$

If the dilution factor (1:150) is taken into account the ALS concentration of the undiluted sample is

$$51.201 \times 150 = 7680 \text{ ng/mL}$$

10.4 Limitations of Procedure

The ALS ELISA is based on antibodies. Generally, this technique could be sensitive to heterophilic antibodies or rheumatic factors in the sample. Their influence is reduced by assay design but cannot be excluded completely.

11. Analytical Performance Characteristics

11.1 Calibration – Traceability

No international standard or reference preparation of ALS is available. Initially, the ALS ELISA was calibrated against a human serum standard. In a second step, the test system was recalibrated with eukaryotic, recombinant ALS. The recombinant ALS was measured in three different kit lots. A comparison of the measured results is shown in Figure 2. The analysis revealed a factor of 5 to convert mUnits into mass units (ng/mL).

According to the function $y = 4.997x$, the factor of 5 is used in the conversion of mUnits (mU/mL) into mass units (ng/mL):

1 mUnit ALS \approx 5 ng recombinant ALS




Figure 2: Assay Calibration, relation of mUnits and mass units

A previously conducted comparative analysis of serum samples demonstrates that the ALS ELISA provides results that are comparable to an in-house assay used by an academic group (see Figure 3).




Figure 3: Comparative analysis of a competitive immunoassay with the ALS ELISA (serum samples: n=25)

11.2 Analytical Sensitivity

The analytical sensitivity was assessed by measuring the blank and calculating the theoretical concentration of the blank + 2SD. In three measurements a range of 0.15 – 1.15 ng/mL with a mean sensitivity of 0.53 ng/mL was detected.

11.3 Precision Data

Intra-Assay Variance

Two samples have been measured 22 times in the same assay. The results are shown in Table 2. The measured coefficient of variation (CV) is 6.7% on average.

Table 2 Intra-Assay Variation

	Number of determinations	Mean value [ng/mL]	Standard deviation [ng/mL]	VC [%]
Sample 1	22	4556	298	6.55
Sample 2	22	6694	458	6.84

Inter-Assay Variance

Serum samples were measured in independent assays. On average the coefficient of variation was 8.96% (SD 6.11). Results are shown in Table 3.


Table 3 Inter-Assay Variation

	Number of single determinations	Mean value [ng/mL]	Standard deviation [ng/mL]	VC [%]
Sample 1	39	4980	485	10
Sample 2	45	5530	525	9
Sample 4	12	3225	230	7

11.4 Linearity

Linearity was tested by dilution of three native serum samples with high ALS content. The optical density of each dilution was measured and the results are shown in Figure 4. Serial dilution of three samples within a range of 1:50 – 1:500 revealed a good linearity measured by linear regression analysis ($R^2 > 0.95$).

Figure 4: Linearity, measured signal intensity [OD450] of differentially diluted samples. The recommended dilution is 1:150.



11.5 Interference

Interference of bilirubin and triglycerides was tested by adding different amounts of these substances to human serum containing ALS. For comparison the same amount of buffer without any substance was also added to the serum. Table 4 demonstrates that none of the tested substances exert any influence on the measurement of ALS in human serum.

Table 4 Interference

Triglyceride [mg/mL]	0	12,5	25	50	100
ALS [ng/mL]	5809	5403	5780	5383	5813
Bilirubin [µg/mL]	0	25	50	100	200
ALS [ng/mL]	5809	5283	5431	5771	5439
Hemoglobin [µg/mL]	0	0.125	0.25	0.5	1
ALS [ng/mL]	5809	5667	5315	6015	6100

11.6 Species Cross-Reactivity

Several commercially available animal sera have been diluted 1:10 and the diluted samples were used in this assay. Very low values were detected in serum samples of chicken, cattle, dog, rat, donkey, mouse, goat, sheep, guinea pig, fetal calf serum. On average the signal intensity was about 0.1 (corresponding Standard A value: 0.04).

12. References

1. Ballard. J.. et al.. *On the nomenclature of the IGF binding proteins*. Acta Endocrinol (Copenh). 1989. 121(5): p. 751-2.
2. Wilson. E.M.. Y. Oh. and R.G. Rosenfeld. *Generation and characterization of an IGFBP-7 antibody: identification of 31kD IGFBP-7 in human biological fluids and Hs578T human breast cancer conditioned media*. J Clin Endocrinol Metab. 1997. 82(4): p. 1301-3.
3. Baxter. R.C.. *Characterization of the acid-labile subunit of the growth hormone-dependent insulin-like growth factor binding protein complex*. J Clin Endocrinol Metab. 1988. 67(2): p. 265-72.
4. Baxter. R.C. and J.L. Martin. *Structure of the Mr 140.000 growth hormone-dependent insulin-like growth factor binding protein complex: determination by reconstitution and affinity-labeling*. Proc Natl Acad Sci U S A. 1989. 86(18): p. 6898-902.
5. Holman. S.R. and R.C. Baxter. *Insulin-like growth factor binding protein-3: factors affecting binary and ternary complex formation*. Growth Regul. 1996. 6(1): p. 42-7.
6. Baxter. R.C.. *Radioimmunoassay for insulin-like growth factor (IGF) II: interference by pure IGF-binding proteins*. J Immunoassay. 1990. 11(4): p. 445-58.
7. Blum. W.F. and M.B. Ranke. *Use of insulin-like growth factor-binding protein 3 for the evaluation of growth disorders*. Horm Res. 1990. 33 Suppl 4: p. 31-7.
8. Miller. B.S.. et al.. *The insulin-like growth factor system in children with congenital disorders of glycosylation*. Clin Endocrinol (Oxf). 2009.
9. Fofanova-Gambetti. O.V.. et al.. *Three novel IGFALS gene mutations resulting in total ALS and severe circulating IGF-I/IGFBP-3 deficiency in children of different ethnic origins*. Horm Res. 2009. 71(2): p. 100-10.
10. Hwa. V.. et al.. *Total absence of functional acid labile subunit. resulting in severe insulin-like growth factor deficiency and moderate growth failure*. J Clin Endocrinol Metab. 2006. 91(5): p. 1826-31.
11. Heath. K.E.. et al.. *Primary acid-labile subunit deficiency due to recessive IGFALS mutations results in postnatal growth deficit associated with low circulating insulin growth factor (IGF)-I, IGF binding protein-3 levels and hyperinsulinemia*. J Clin Endocrinol Metab. 2008. 93(5): p. 1616-24.
12. Domene. H.M.. et al.. *Acid-labile subunit deficiency: phenotypic similarities and differences between human and mouse*. J Endocrinol Invest. 2005. 28(5 Suppl): p. 43-6.
13. Fischer. F.. et al.. *Associations of insulin-like growth factors. insulin-like growth factor binding proteins and acid-labile subunit with coronary heart disease*. Clin Endocrinol (Oxf). 2004. 61(5): p. 595-602.
14. Morrison. K.M.. et al.. *Sample pre-treatment determines the clinical usefulness of acid-labile subunit immunoassays in the diagnosis of growth hormone deficiency and acromegaly*. Eur J Endocrinol. 2007. 156(3): p. 331-9.
15. Tzanela. M.. et al.. *Growth hormone binding protein and acid labile subunit levels in the assessment of acromegaly treatment*. Hormones (Athens). 2005. 4(3): p. 148-54.
16. Stadler. S.. et al.. *Monoclonal anti-acid-labile subunit oligopeptide antibodies and their use in a two-site immunoassay for ALS measurement in humans*. J Immunol Methods. 2001. 252(1-2): p. 73-82.
17. Khosravi. M.J.. et al.. *Acid-labile subunit of human insulin-like growth factor-binding protein complex: measurement, molecular, and clinical evaluation*. J Clin Endocrinol Metab. 1997. 82(12): p. 3944-51.