

Custom Shoes: Indications and Construction

What are custom shoes?

- Custom shoes are often the best option for patients whose condition cannot be effectively treated with readily available footwear options, modifications or customization.
- When the need for such footwear arises and the therapeutic goals have been defined, molds and measurements of the feet are taken in their entirety. These molds are turned into the actual lasts over which the shoes are constructed.



Custom Shoes

- Often, people who think they need custom shoes really just need properly fitting shoes.
- Pros
 - Shoe should (in theory) fit the foot perfectly
 - Will accommodate any degree of deformity
 - Should be quite comfortable for the patient
 - Design and feature possibilities are endless
- Cons
 - Cost
 - Turnaround time of 1-3 months
 - Appearance
- Consider structural and functional demands of foot and patient when configuring shoe and choosing materials.
 - E.g. Don't order a featherweight construction for a 300-pound patient.



The Key to Success

- The most important thing when it comes custom shoes and your patient is *expectation management*.
- If your patient has a severely deformed foot and the shoe is made from a cast of that severely deformed foot and is made to fit that severely deformed foot perfectly...guess what the shoe is going to look exactly like.
- Don't simply show them a catalog or the shoemaker's website. Have a sample or pictures of actual patients' custom shoes to show them.



EXPECTATION



REALITY

Objectives

- Accommodate deformity
- Relieve pain
- Provide support, control and stability
- Improve mobility
- Restore safe, energy-efficient gait
- Prevent further deformity or worsening of condition

Indications

- In most cases, custom shoes are used when there are no factory-made footwear options that will work for a particular patient. Most of the time, they are necessitated by severe deformity.
- Indications:
 - Extreme size differences in feet such as in a post-polio patient
 - Congenital foot deformities and conditions such as clubfoot, cerebral palsy or spina bifida
 - Charcot arthropathy
 - Rheumatoid or psoriatic arthritis
 - Charcot-Marie-Tooth disease
 - Lymphoedema
 - Severe forefoot deformities
 - Partial foot amputation
 - Severe adult acquired flatfoot deformity
 - Trauma such as crush injuries, burns or gun shot wounds
 - Accommodation of bulky AFO



Features and Benefits

- Custom shoes offer many material choices that simply aren't available (anymore) in factory-made orthopedic/therapeutic footwear.
 - Deerskin
 - Plastazote® lining
- Custom shoes offer features that are incredibly difficult or impossible to incorporate into, or add onto, a factory-made shoe.
 - Extra, extra, extra depth
 - Internal heel lifts to accommodate severe equinus
 - Accommodation for posterior heel deformities
- Custom shoes allow the pedorthist to adequately care for patients with:
 - Chronic and severe fluctuating edema
 - Asymmetric forefoot shapes and widths
 - Severe rocker bottom feet



Alternatives

- Custom shoe on affected foot only
- Mismates (split sizes)
- Relasting to accommodate deformities



Foot Capture for Custom Shoes

Casts

- Casts must be semi-weight-bearing and closed-toe. Patient should be seated with hip, knee and ankle at 90° (if possible). If patient is not able to achieve 90° at the ankle – or that is not the most functional position for the patient – then the mold should be taken with the foot and ankle in the alignment of best function.
- Cast should be taken on a soft, conformable casting pad that is at least 1" thick – unless shoemaker requests/requires use of a casting board with a heel rise and toe spring.
- Preferable casting techniques are:
 - Plaster bandage
 - Bi-valve
 - Circumferential wrap
 - Circumferential wrap using fiberglass casting tape
- Plaster wraps must be wrapped loosely, and the toes and metatarsals must not be restricted. Bring the cast to at least 1" above the height of the shoe being made. Explain to patient what you are trying to achieve so that they don't change or move during the casting, and they maintain the 90° angle. Check the final cast and see that:
 - The cast stands on its own.
 - The lower leg is vertical (unless there is a rigid deformity).
 - Any areas of concern are clearly marked.



Molding Techniques

- Plaster Bi- Valve




Casts

- STS® socks are not recommended. It is difficult to get a tight, accurate mold around the Achilles tendon. Also, they have a tendency to shrink during transit, especially in the back of a hot shipping trailer.
- When using fiber-glass casting material, be aware of the following issues:
 - The most common problem with fiber-glass casts is that the cast does not show the end of the longest toe. Proper wrapping around the end of the longest toe will eliminate this kind of problem. Most practitioners wrap the toes first, then the rest of the foot.
 - The other common problem with fiber-glass casts is that the wrapping is too tight, which makes the cast narrower and longer than the actual foot.
- When shipping plaster casts, be sure wait until they are completely dry before packing them. Wet casts may collapse or change shape under pressure of packing materials or packing process.



Foot Tracing

- Always provide weight-bearing foot tracings so that the custom last maker will know how much the foot splays and elongates under weight-bearing as compared to the semi-weight-bearing cast.
- Weight-bearing tracing is done while patient is in standing position.
- In the case that the patient has difficulties standing, a semi-weight bearing tracing will still help the last maker determine if the toes were curled or compressed by the casting material.
- Hold the pen or pencil in perpendicular position to the paper. Remember to indicate whether the tracing is weight-bearing or semi-weight-bearing.



Foam Impressions

- Sending foam impressions are optional. However, they can be extremely helpful, especially if the patient needs additional offloading in a specific area. It is easier to indicate this more precisely in the foam.
- They are highly recommended as they are a very helpful reference to the cast. Foot impressions provide better definition of the plantar surface of the foot than the plaster cast. They show *much* better definition than a fiberglass cast.


Photographs

- Foot pictures are also optional. But in case of any ulcerations/callusing, amputations, or severe foot deformities, pictures can provide better understanding of the foot condition, shape and shoe requirements.

Measurements


Circumferential measurements can be quite helpful to the last maker when it comes to accurately reproducing the foot.

1. Circumference at the ball of the foot
2. Circumference at the instep
3. Circumference around foot at the back of the heel
4. Circumference at the ankle



Order Forms

- Perhaps the most important factor in getting the shoes you think you're ordering is the order form.
- Use the order form specific to the manufacturer.
- Fill the form out completely and accurately.
- Contact the lab if you aren't sure about something.

Apis		APIS CUSTOM SHOES ORDER FORM	
Apis Footwear, 2239 Tyler Ave, S. El Monte, CA 91733 1-888-937-2747			
Account No.: _____		Business Name: _____	
P. O. No.: _____		Ship Address: _____	
Contact: _____		_____	
Phone: _____		_____	
Email: _____		City: _____ State: _____ Zip: _____	
Patient Name: _____		Country: _____	
Gender: ☐ Male ☐ Female Weight: _____ Height: _____		Ship Via: ☐ Ground ☐ Next Day Air ☐ 2 nd Day Air	
CUSTOM SHOE SPECIFICATIONS		DIAGNOSES AND FOOT MEASUREMENTS	
Shoe Style: _____ Color: _____ Qty: _____		Diagnosis: ☐ Diabetes ☐ Neuropathy	
Shoe Size: L _____ R _____ ☐ Match Shoe Length		Ulcer: ☐ L ☐ R Metatarsalgia: ☐ L ☐ R	
Boot Height Including Sole: _____ Toe Elongation: _____		Charcot: ☐ L ☐ R Amputated Toes: ☐ L ☐ R	
Extra Toe Box Height: L _____ R _____		Heel Spur: ☐ L ☐ R Hammer Toes: ☐ L ☐ R	
Closure: ☐ Laces ☐ Velcros ☐ Semi-Surgical ☐ Surgical			
Lining: ☐ Leather ☐ Mesh ☐ Plastazote ☐ Soft Cloth		Circumferences:	
ORTHOTIC SPECIFICATIONS		Foot Measurements: Ball: L _____ R _____	
A5513: _____ pairs or L _____ R _____		Foot Length: L _____ R _____	
TOE FILLERS: _____ pairs or L _____ R _____		Ball Width: L _____ R _____	
Tri-lam: ☐ Plastazote + PPT + EVA		Instep: L _____ R _____	
Top Cover: ☐ Plastazote ☐ Spenco ☐ Leather ☐ Vinyl		Heel: L _____ R _____	
☐ Other: _____		Toe Height: L _____ R _____	
Base: ☐ EVA ☐ Thermocork ☐ Cork		SOLE MODIFICATIONS	
☐ Other: _____		Forefoot Rocker ☐ L ☐ R Mild Rocker ☐ L ☐ R	
Base Density: ☐ Soft ☐ Medium ☐ Hard		Heel-to-toe Rocker ☐ L ☐ R Sever Rocker ☐ L ☐ R	
Heel Cup: ☐ Flat ☐ Medium ☐ Deep		Double Rocker ☐ L ☐ R Rocker Bar ☐ L ☐ R	
Met Pad: ☐ L ☐ R		Negative Rocker ☐ L ☐ R LOP Rocker ☐ L ☐ R	
Met Bar: ☐ L ☐ R		$\frac{3}{4}$ Length Shank ☐ L ☐ R Bevel Heel ☐ L ☐ R	
Arch Pad: ☐ L ☐ R		Full Length Shank ☐ L ☐ R SACH Heel ☐ L ☐ R	
Heel Pad: ☐ L ☐ R		Medial Buttress ☐ L ☐ R Snug Heel ☐ L ☐ R	
Lateral Flange: ☐ L ☐ R		Lateral Buttress ☐ L ☐ R Detach Sole ☐ L ☐ R	
Medial Flange: ☐ L ☐ R		Medial Wedge L _____ R _____ Medial Flare L _____ R _____	
Offload: ☐ L ☐ R		Lateral Wedge L _____ R _____ Lateral Flare L _____ R _____	
Mark areas for offloading		Heel-to-toe Lift L _____ R _____ Heel Lift L _____ R _____	
		Build Lift: ☐ On outsole ☐ Inside shoe	
Select dits for toe-fillers:		Soling: ☐ Light-weight ☐ Heavy Duty ☐ Rigid	
Left: ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ TMA		MISC. OPTIONS	
Right: ☐ 1 ☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ TMA		☐ Extra tongue padding ☐ Extra collar padding	
ADDITIONAL INSTRUCTIONS _____		☐ Composite Toe Box ☐ Reinforce Heel Counter	
_____		_____	
_____		_____	

*In compliance with federal guidelines, for any custom shoe re-orders, please submit new foot impressions, casts, or scans.
*Custom products are not refundable but we will remake at no additional charge.

Rev. 2023.11

Account Name _____ Phone (____) _____
 Address _____ Date _____
 City _____ State _____ Zip _____
 Practitioner Name _____ e-mail _____

Patient Name _____ M ☐ F ☐
 Patient Weight _____ lbs Diagnosis _____
 Principal Reason for the Shoe _____
 Is Patient currently ☐ Ambulatory OR ☐ Non-Ambulatory?

PLEASE INCLUDE WEIGHT BEARING TRACINGS OF PATIENT'S FEET

SHOE STYLE

- ☐ LOW TOP ☐ HIGH TOP
☐ CHUKKA ☐ OTHER _____

We recommend Chukka or High Top for those patients with: midfoot collapse, charcot deformity, severe edema, transect amputations, severe pes plano valgus, and any type of brace.

UPPER SPECIFICATIONS

- OPENING ☐ Regular OR ☐ Surgical
 CLOSURE ☐ Lace OR ☐ Velcro
☐ PADDED COLLARS ☐ PADDED TONGUES
☐ PLASTAZOTE LINING
☐ OTHER _____

COLOR

- ☐ BLACK ☐ DARK BROWN ☐ OTHER _____

CAST MODIFICATION

- ☐ EXTRA HIGH TOE BOX
☐ SNUG HEEL FIT
☐ EXTRA TOE ELONGATION
☐ DEPRESS AS MARKED ON CAST/DIAGRAM
☐ DUPLICATE AND RETURN CAST
☐ OTHER _____

SHOE WEIGHT

- ☐ LIGHT ☐ REGULAR ☐ HEAVY DUTY

INSOLE

- ☐ 1/4" PINK AND WHITE PLASTAZOTE
☐ CORK INSOLE (for AFO)
☐ TOE FILLER _____ Left _____ Right
☐ 1 EXTRA pair CUSTOM INSOLES
☐ 2 EXTRA pairs CUSTOM INSOLES
☐ OTHER _____

LIFTS

- ☐ INTERNAL OR ☐ EXTERNAL

LEFT RIGHT

HEEL _____ in. _____ in.
 BALL _____ in. _____ in.
 TOE _____ in. _____ in.

OUTSOLE

- ☐ REGULAR ☐ ROCKER
☐ MID-SOLE FOR CALIPER _____ Left _____ Right
☐ LEAVE SOLE OFF FOR ADJUSTMENT
☐ SOLE STIFFENER _____ Left _____ Right
☐ FLARES ☐ Medial OR ☐ Lateral
 _____ Left _____ Right
☐ WEDGES ☐ Medial OR ☐ Lateral
 _____ Left _____ Right
☐ OTHER _____


IF YOU HAVE ANY QUESTIONS - GIVE US A CALL!

OVER

FOOT TYPE

- ARE THE PATIENT'S FEET; ☐ RIGID OR ☐ FLEXIBLE?
 WHEN WEIGHT BEARING ARE TOES; ☐ UP OR ☐ DOWN?
 WHEN NON-WEIGHT BEARING ARE TOES; ☐ UP OR ☐ DOWN?
 BRING LOWER LEG TO 90° BY; ☐ CORRECTING CAST TO NEUTRAL
 OR ☐ BUILDING EXTERNAL LIFT UNDER HEEL
 OR ☐ LEAVE AS IS / DO NOT CORRECT
 ARE PATIENT'S FEET SUBJECT TO EDEMA; ☐ YES OR ☐ NO?

SPECIAL INSTRUCTIONS/NOTES



PLEASE INCLUDE WEIGHT BEARING TRACINGS OF PATIENT'S FEET

VERY IMPORTANT

Please provide a weight-bearing tracing
of your patient
on the inside pages of this order form.

TRU-MOLD SHOES

Orthopedic Shoes - Custom Designed Since 1955

42 Breckenridge St. Buffalo, NY 14213
(716) 881-4484 | Fax (716) 881-0406
info@trumold.com

SHIPPING:

☐ Two Day ☐ Ground ☐ Overnight

Please Fill this form completely

COMPANY NAME

Date _____

P.O.# _____ Phone _____

ORDERED BY _____

Fax # _____ E-mail _____

SHIP TO

Address _____

City _____ State _____ Zip _____

BILL TO

Address _____

City _____ State _____ Zip _____

Patient Name _____

Address _____

City _____ State _____ Zip _____

Sex _____ Weight _____ Age _____

Occupation _____

Had Tru-Mold Shoes before (Approx. Date) _____ Invoice CAD/CAM # _____

Repeat CAD/CAM order as before, no changes _____

Repeat CAD/CAM order but with the following changes shown on order form _____

Rx. DIAGNOSIS

Diabetic _____ Post Polio _____ Amputee (where): _____

Charcot _____

Casted over AFO: _____ LF _____ RF _____ Arthritic _____

Casted Prosthetic Foot _____ LF _____ RF _____

Other _____

Style _____ Color _____

Regular Opening _____ Wide Vent _____

Semi-Surgical _____ Scalloped Heel _____

Surgical _____

Other _____

CAST MODIFICATION


Reg. Elongation (5/8") _____ 1/8" Extra (5/8" + 1/8") _____

High Toe Box _____ Extra High Toe Box _____

Other _____

INDICATE ANY LESIONS

(I.E. AMPUTATIONS OR ULCERATIONS)

**PLASTAZOTE® MOLDED INSERT**

☐ One Pair ☐ Two Pairs ☐ Three Pairs

1/4" Pink + 1/4" White _____ 1/4" Poron Bottom _____

1/2" White (firm) _____ Thermocork _____

1/4" White (firm) _____ Spenco® Cover _____

1/2" Pink (soft) _____ Leather Cover _____

1/4" Pink (soft) _____ Met Pads _____

1/8" Poron Middle _____ Other _____

1/8" Poron Bottom _____

Recess insert and base of shoe ☐

LINING

Regular (full leather) _____ Flannel _____

Spenco® Heel _____ Vamp _____

Plastazote® Heel _____ Vamp _____

Tan Athletic _____ Black Athletic _____

Other _____ Fleece _____

CONSTRUCTION

Regular _____ Reinforce Counters _____

Heavy Duty _____

SHORT LIMB RAISE

L _____ R _____

Heel _____

Ball _____

Toe _____

Raise inside shoe cover with leather _____

Raise outside shoe with soling _____

Other _____

BASE MODIFICATIONS

L _____ R _____

Lateral Flare _____

Medial Flare _____

Lateral Wedge _____

Medial Wedge _____

Wide Base _____

Reinforced for Brace _____

Other _____

EXTERNAL ADDITIONS

Amputation filler - match shape and length _____ L _____ R _____

No amputation filler _____

Steel shank 3/4 length _____

Steel shank full length _____

Attach Tru-Mold's caliper plate _____

Attach customer's caliper plate/stirrup _____

Thermo Plastic toe cap _____

Other _____

SOLING

Standard (15 Iron) _____ Mini Rib _____ Silvano _____

Heavy Duty Work _____ Infinity _____ Lug _____

Rockersole Forefoot L _____ R _____

Rockersole Heel to toe L _____ R _____

Leave outside wedge and final outsole off L _____ R _____

Leave final outsole off L _____ R _____

Other _____

UPPERS

Lace _____ Tongue Loops _____

Velcro Strap _____ Speed-lace (top 3 eyelets) _____

Velcro D-Ring _____ Ortho-lace (combination lace/velcro) _____

Perforate Arch _____ Medial T-strap L _____ R _____

Perforate Top _____ Lateral T-strap L _____ R _____

Pull Loops _____ Attach T-strap yes _____ no _____

Other _____

Special "Foot Notes"

"MOLDED SHOES ARE ONLY AS GOOD AS THE CASTS RECEIVED"

Revised 4/28/22

TECHNICAL ASSISTANCE 1-800- "THE MOLD" (next page please)



ORDER / SPECIFICATION FORM

36 Mease Street / Buffalo, NY 14213 / 1-800-435-0965 / Fax 716-881-4349 / www.jerry miller shoes.com

VERY IMPORTANT

Please provide a weight-bearing tracing of your patient on the inside pages of this order form.

Date: _____ PO#(if needed) _____

Shipping: ☐ Two Day ☐ Ground ☐ Overnight

Ship to:

Name _____

Address _____

City _____ State _____ Zip _____

E Mail _____

Phone _____

Fax _____

Contact Person Regarding Order: _____

Bill to:

Name _____

Address _____

City _____ State _____ Zip _____

Patient Name _____

Sex _____ Weight _____ Age _____

Occupation _____

Had Jerry Miller Shoes before - Approx. Date _____

Invoice/CAD/CAM # _____

Repeat CAD/CAM order as before, no changes _____

Repeat CAD/CAM order but with the following changes shown on order form. _____

Diagnosis

☐ Diabetic ☐ Club Foot ☐ Charcot ☐ Arthritic

☐ Post Polio

☐ Amputee - which areas? _____

Special _____

Last Modification:

☐ Regular Elongation (5/8")

☐ 1/8" Extra Elongation (5/8"+1/8")

☐ High Toe Box

☐ Extra High Toe Box

☐ Special _____

Removable Plastazote® Molded Insert

☐ None ☐ 1/2" White (firm)

☐ 1/4" Pink + 1/4" White ☐ 1/2" Pink (soft)

☐ 1/4" White (firm) ☐ Thermo-cork

☐ 1/4" Pink (soft) ☐ Spenco Top Cover

☐ 1/8" Poron Middle ☐ Leather Top Cover

☐ 1/8" Poron Bottom

☐ 1/4" Poron Bottom

Special _____

☐ One Pair ☐ Two Pairs ☐ Three Pairs

Type of Casting

☐ Casted over AFO ☐ LF ☐ RF

☐ Casted over Prosthetic Foot ☐ LF ☐ RF

Style

Color _____

Please note that the colors may vary from swatches as printed

in catalog or shown on website

Type of Opening

☐ Regular Opening

☐ Opening To Ball of Foot (Semi-Surgical)

☐ Opening To End of Toes (Surgical)

Type of Closure

☐ Lace

☐ Speed Lace (Top three cycles unless specified)

☐ Velcro Strap

☐ Velcro with D-Ring

☐ Special _____

Lining

☐ Regular (glove leather) ☐ Spenco

☐ Fleece ☐ Black Sports

☐ Plastazote ☐ Tan Sports

☐ Flannel

☐ Check here if lining is for HEEL ONLY

Finishing

☐ Perforate Top ☐ Met Pads

☐ Perforate Arch ☐ Stitch Tongue to Medial Side

☐ Pull Loops ☐ Stitch Tongue to Lateral Side

☐ Tongue Loops ☐ Special _____

Soling

☐ 15 Iron (at no additional cost)

☐ Leather Heel (outside) ☐ Infinity

☐ Leather Sole ☐ Topy Cover

☐ Mini Rib ☐ Work Cover

☐ Soft Spike ☐ Lug

☐ Soft Vibram ☐ Outside Heel

☐ Neoline Heel ☐ Silvano

Special _____

Rocker Sole

Left Foot: ☐ Heel to Toe ☐ Forefoot ☐ Aggressive

Right Foot: ☐ Heel to Toe ☐ Forefoot ☐ Aggressive

Construction

☐ Regular

☐ Heavy Duty

☐ Reinforce Counters

Short Limb Raise

Heel ☐ L _____ ☐ R _____

Ball ☐ _____ ☐ _____

Toe ☐ _____ ☐ _____

Raise inside shoe ☐ _____ ☐ _____

covered with leather ☐ _____ ☐ _____

Raise outside shoe ☐ _____ ☐ _____

on soling ☐ _____ ☐ _____

Special: _____

Base Modifications

Lateral Flare ☐ L ☐ R

Medial Flare ☐ L ☐ R

Wide Base ☐ L ☐ R

Lateral Wedge ☐ L ☐ R

Medial Wedge ☐ L ☐ R

External Additions

Amputation Filler — ☐ L ☐ R

match length and shape ☐ L ☐ R

No Amputation Filler ☐ L ☐ R

Steel Shank 3/4 length ☐ L ☐ R

Steel Shank full length ☐ L ☐ R

Reinforce for Brace ☐ L ☐ R

Leave Outside Wedge and ☐ L ☐ R

Final Outsole off ☐ L ☐ R

Leave Final Outsole off ☐ L ☐ R

Attach Jerry Miller's 7/8" Caliper ☐ L ☐ R

Attach Customer's 7/8" Caliper/Stump ☐ L ☐ R

T-Strap Lateral ☐ L ☐ R ☐ Unattach?

T-Strap Medial ☐ L ☐ R ☐ Unattach?

Thermo Plastic Toe Cap ☐ L ☐ R

Steel Toe ☐ L ☐ R

If Standard Steel Toes cannot accommodate certain

deformities, shapes and sizes, may we substitute a

Thermo Plastic Toe Cap? ☐ Yes ☐ No

Special Instructions (please print):

Indicate any Lesions (ie. Amputations or Ulcerations)

☐ Recess on Insert Only

☐ Recess on Insert and Base of Shoe

LEFT FOOT



Left Lateral (Outside) Left Medial (Inside) Dorsum Plantar

RIGHT FOOT



Plantar Dorsum Right Medial (Inside) Right Lateral (Outside)

Revised 7/2022

3825 Investment Ln, Suite # 10
Riviera Beach, FL 33404
Phone (561) 840-6792
Fax (561) 840-6799
www.rightwaycms.com



CUSTOM SHOE WORKORDER

To serve properly,
please fill form
out completely and legibly

COMPANY NAME

Phone Fax

Purchase Order #

Date Account #

PRACTITIONER NAME

PATIENT NAME

Male Female Age Weight

Occupation

Had Rightwaycms before Date

CAST MODIFICATION

Correct Cast to 90°

Standard Toe Elongation (¼)

Extra Toe Elongation (¼)

High Toe Box Extra High Toe Box

Duplicate Cast* Met-Pad* L R

Base Depression* L R

SHOE STYLE Color

Heavy Duty Leather*

OPENING

Regular Semi-Surgical Surgical

Velcro® 1 2 Laces

Velcro D-ring Normal Direction (Standard) L R

Velcro D-ring Reverse Direction L R

Velcro Flat Reverse Direction L R

Velcro Flat Normal Direction L R

T-Strap* L R

Hook* Speed Lace*

SHIP TO:

Address

City State Zip

DIAGNOSIS

Diabetic Toes Overlapped L R

Post Polio Hammered L R

Amputee Arch No Deformity L R

Charcot Flexible Deformity L R

Arthritic Rigid Deformity L R

Other

REMOVABLE CUSTOM MOLDED INSERTS

Standard: ¼" Pink + ¼" White Plastazote®

¼" Pink + 1/8" Poron® + ¼" White Plastazote®

Cover* Leather* Spenco®* Poron®*

Extra Inserts 1 pair* 2 pairs** L R

Build-Up Inserts* L R

Other:

SPECIAL LINING

B-foam Cushion (standard)

Full Leather*

Plastazote®*

Vamp*

Heel*

*Collar and Tongue come padded Standard Unless

Otherwise Noted

Collar (no padding) Tongue (no padding)

Other:

CONSTRUCTION

Regular (light)

Elastic Counters (Standard) Double Counter*

Other:

SOLING

L W (12 I) Standard (18 I)

(24 I)* Heavy Duty Ribbed Soling*

Other:

SHORT LEG BUILD-UP

Heel* L R Inch

Ball* L R Inch

Toe* L R

Outsole* Insole*

Other:

SOLE MODIFICATIONS

Heel Flare* Lateral Medial L R

Sole Flare* Lateral Medial L R

Heel Wedge* Lateral Medial L R

Sole Wedge* Lateral Medial L R

Wide Base* Rocker sole Rocker Heel

Other:

PLEASE MARK AREAS OF SPECIAL ATTENTION ON THE PICTURES BELOW AND ON THE CAST

Right

Left



SPECIAL INSTRUCTIONS

Include A Weight-bearing Tracing

(HOLD A PENCIL VERTICALLY WHEN TRACING)

EXTERNAL ADDITIONS

Toe Filler* L R

To Match (length) L R

Short L R

Steel Shank* L R

7/8" Caliper Plate* L R

Metatarsal Bar* L R

Reinforce for Brace L R

Other:

* THERE IS AN ADDITIONAL CHARGE FOR
THIS MODIFICATION.

PLEASE SEE PRICE LIST.

THE BASIC INFORMATION REQUIRED
OR YOU MAY INCUR EXTRA COST AND DELAYS

SEND LITERATURE ORDER FORMS

NOTICE: Molded shoes are only as good as the cast received

Construction



Match factory shoe on contralateral side

- Most custom shoemakers will match a factory shoe for you. Keep in mind, however that it will not look exactly like a matched pair as the one was made from a mold of a foot with severe deformities.



Shoe Modifications and Footbeds

- Any shoe modification you can think of can be incorporated into a custom shoe.
- Likewise, any sort of foot orthosis, posting, offloading or correction may also be built in to the footbed of a custom shoe.
- In both cases, the correction, support or offloading can be utilized to a much higher degree and better effect in a custom shoe without the inherent restraints of a factory-made shoe.

Tips and Reminders

- Make sure you use the manufacturer's own order forms
- Make sure you have the most up-to-date price list
- Make sure you understand what features are included in the base price of the style you selected, and which are add-ons. You need an accurate estimate to be able to determine what to charge for the shoes.
- Always be sure to send weight-bearing tracings with the casts. Remember, measurements and photos, while not required, are helpful.
- Foam impressions are a good idea if you need special offloading or the plantar surface of the foot is abnormal or unusual.
- If the shoe doesn't fit as expected, it's best to send it back to be adjusted rather than adjust it yourself as this may void the warranty.

Resources

Hersco Ortho Labs

<https://hersco.com/>

Apis

<https://apisfootwear.com/>

Jerry Miller ID Shoes

<https://jerrymillershoes.com/>

Davis Shoe Therapeutics

<https://www.davisshoes.com/>

Tru Mold Shoes

<https://trumold.com/>

Crary

<https://www.craryshoes.com/>

Eneslow Shoes & Orthotics

<https://eneslow.com/>

Rightway CMS

<https://rightwaycms.com/>

Medicare and Custom Shoes

- Under the Therapeutic Shoes for Persons with Diabetes benefit:
 - A5501 - For diabetics only, fitting (including follow-up), custom preparation and supply of shoe molded from cast(s) of patient's foot (custom molded shoe), per shoe
 - A custom-molded shoe is one that:
 1. Is constructed over a positive model of the beneficiary's foot; and
 2. Is made from leather or other suitable material of equal quality; and
 3. Has removable inserts that can be altered or replaced as the beneficiary's condition warrants; and
 4. Has some form of shoe closure.
 - A5501 includes two additional pairs of inserts.
 - New casts must be taken every year.
- As an integral component of a brace:
 - L3230 - Orthopedic footwear, custom shoe, depth inlay, each

Case Studies


Case Study #1



Case Study #1



Case Study #2



Case Study #2



Case Study #2



Case Study #3



Case Study #3



Case Study #4



Case Study #4



Case Study #5



Case Study #5



Case Study #6



Case Study #7



The End