Osseous Technologies of America
10 Centre Dr.
Orchard Park, NY 14127
716-667-9211

Medical Device and QA Services
76 Stockport Road
Timperley, Cheshire, UK WA157SN

REF BC 100

THE BONE COLLECTOR

The bone collector technique was developed to prevent the waste of valuable bone that can be used to correct small defects around dental implants and periodontal compromised teeth. The procedure is based on scientific data which shows that the bone collected is vital and has the potential to promote bone formation.

Product characteristics:

- Easy and fast way to harvest autogenous bone
- Disposable
- Supplied with one filters
- Shut off valve

If you have a Hivac suction, an autoclavable adapter should be used.

Note: Adapter catalog number: Hi-Vac Adp

The Bone Collector Technique:

To obtain the best results with the disposable bone collector use two suction hoses. One for blood and saliva and the other attached to the bone collector that must work only suctioning at the donor site. This will prevent contamination of the collected bone with other fluids. After drilling and collecting the bone, allow copious wash with saline solution of the obtained tissue. Keep this material wet in saline until it will be placed on the recipient site.

Do not reuse

WARNINGS:


Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.

Re-use of this device may result in serious cross infection or be harmful to the patient. Re-sterilization may adversely affect material properties and device performance.

STERILE R Sterilized by gamma radiation

Contents sterile unless package is opened or damaged. Examine package to ensure it has not been opened and has not been damaged. If the sterile packaging is opened or damaged, do not use and contact Osseous Technologies of America immediately. After use, this product may be a potential

biohazard. Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

	Do not reuse
	Consult IFU
	Sterilization Using Gamma
	Use By
	Lot Number
	Catalogue Number
	Manufacturer
	European Authorized Representative

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