

Application Instructions




Bio Skin® supports should be applied to dry skin only. Oils, lotions and/or perspiration activate "SkinLok™" making correctly sized supports seem small and difficult to apply. This response to moisture should not cause confusion with respect to proper sizing. Always dry the body segment thoroughly prior to application.

1. Take the support with both hands and stretch the top and bottom of the support two to three times. This helps break-in the material.
2. Pull support on over elbow, positioning the elbow in the center of the support.


LATERAL EPICONDYLITIS-Tennis Elbow

Position the tear drop pad 2" from lateral epicondyle head, toward the wrist.



MEDIAL EPICONDYLITIS-Golfer's Elbow

Position the tear drop pad 2" from medial epicondyle head, toward the wrist.



Easy Care Instructions

Wash Bio Skin® regularly to prevent body salt and oil build-up that naturally occurs after the product has been worn a few times. Attach the hook to the loop to prevent lint build-up. This will give your Bio Skin® a longer life span.



Hand Wash

Hang Dry

Do Not Bleach

Do Not Dry Clean

Cropper Medical, Inc.
custserv@bioskin.com
800.541.2455
www.bioskin.com

