

# SLEEP SAT

## Patient Guide

This Oximetry test has been ordered by your clinician to monitor your oxygen levels and heart rate while you sleep.

**(!) There is no alarm on the SleepSat device, and it is not intended for real time monitoring.**

---

### For Best Results


- Remove nail polish.
- Avoid laying on the hand with the sensor, putting it above your head, or tucking it underneath a pillow.
- Turn the oximeter on **ONLY** when you are **IN BED** and ready to fall asleep.
- If the finger sensor irritates your finger, switch to a different finger.
- Charge device fully between nights.

# When You Are Ready For Bed:

1

---


Plug finger sensor into SleepSat wrist unit.



2

---

Put wrist band on your non-dominant arm. Adjust to fit.



### 3

---

Insert one of your fingers into the sensor.



(!) Make sure your finger is inserted **COMPLETELY!** Your fingertip should touch the end of the sensor.

### 4

---

Attach the sensor cable loosely to the back of the hand.



Press in on the power button firmly until the display lights up. Lie still for 2 minutes before falling asleep



**(!) IF THE FINGER SENSOR COMES OFF OR IS REMOVED DURING THE NIGHT, OR, IF THE WRIST UNIT CONTINUALLY FLASHES "SEN OFF":**

Put the sensor back on, making sure your finger is inserted completely. Press the power button firmly until the display lights up.

6

---

### IN THE MORNING:

Take the oximeter off.  
The oximeter will turn OFF  
automatically after 2 minutes.

7

---

### IF TESTING MULTIPLE NIGHTS:

Remove probe and plug  
SleepSat into USB Charger.

Charge the  
SleepSat  
for 4-6  
hours, then  
repeat  
steps 1-6.



## AT THE SCHEDULED DROP-OFF TIME:

As instructed by your healthcare provider, return all equipment including:

- ✓ SleepSat Wrist Unit
- ✓ Finger Sensor
- ✓ This Instruction Card
- ✓ Wall Charger & Cable (If provided)