## Digital Forehead and Ear Thermometer

Please read this guidebook  
completely before operating this unit.

Includes: 1 Thermometer, 1 User's Manual, 1 Cleaning Guide, 2 AAA Batteries

- Thank you for purchasing the Digital Forehead and Ear Thermometer.
- Remove safety tab from battery. Open battery case and gently pull out the plastic safety tab from underneath the battery. Thermometer will not activate until this has been performed.
- For safe and proper use of this product, be sure to read and fully understand Safety Precautions contained in this instruction manual.
- Keep this Manual with you for a convenient reference.

## ***Safety Precautions***

---

- Warning marks and symbols are contained for your safe and proper use of this product and prevention of any injury to you and others.
- See table below for description of warning marks and symbols:

### **WARNING**

Means a possibility of personal injury in case of improper use.

---

### **NOTICE**

Means a possibility of personal injury or property damage in case of improper use.

**\*Property damage covers any damage to house, family property, domestic animal and pet.**

## ***Marks or Symbols***

---

### **PROHIBITION**

Means Forbidden with detailed items expressed in words or figures within or beside the mark.

### **MUST OBSERVE**

Means Obligatory with detailed items expressed in words or figures within or beside the mark.

### **IMPLICATION OF SYMBOL**

Means Applied part of type BF.

## ***Intended for Use***

---

The device is intended for the intermittent measurement and monitoring of human body temperature from ear canal or forehead by consumers in the home.

## Warning

---

It is dangerous for patients to perform a self-evaluation and self-treatment based on the measuring results. Be sure to follow doctors' instruction.

**\*A self-evaluation may cause deterioration of diseases.**

Don't touch or blow infrared sensor.

**\*A dirty infrared sensor may cause inaccuracy.**

Clean a dirty infrared sensor with a soft, wet cloth in a gentle manner, or with a cotton swab slightly moistened with alcohol.

**\*Do not clean with a toilet tissue or paper towel because scratching of the infrared sensor could occur, causing inaccuracy.**

Keep the thermometer out of children's reach

**\*A self-measurement forced by children may cause ear damage.**

**For accidental swallowing of battery, please consult the doctor at once.**

Measurement shall be suspended until ear is warmed up after any outdoor activities

**\*Using a water-filled pillow or bag, or any immediate measurement after an outdoor activity in winter may cause a lower result.**

For temperature difference between storage area and measurement site, condition the thermometer for about 30 minutes in room temperature (measurement site).

**\*Failing to do so may cause inaccuracy.**

If the user's body temperature exceeds 100.4°F/38°C and / or feel uncomfortable, the user must contact medical practitioner.

Stop using the product in occurrence of any pain.

**\*It is likely to damage the ear canal.**

It is not recommended to use the product by users suffering from ear disease including otitis external and tympanitis.

**\*These conditions may worsen.**

Don't use the product in a wet ear resulting from swimming or bathing.

**\*It is possible to damage the ear or infrared sensor.**


Don't throw batteries into fire.

**\*Battery is likely to explode in fire.**

## Recommendations

- Don't use the product for other purposes.
- This device is NOT waterproof. Clean it with a soft, wet cloth or cotton swab.
- Do NOT clean with any chemical solvent or expose to direct sunlight or high temperatures.
- Don't drop or impose any vibration or impact on the product.
- Don't use a mobile phone nearby when the thermometer is working.
- Please act according to the native law in handling and disposing of the batteries.
- Take out the batteries if you are not going to use the unit for a long period of time.

## Product Structure


## Body Temperature Chart

A person's normal temperature will vary and trends to decrease with age. The following table shows normal temperature ranges by age:

Age	Temperature (°F)	Temperature (°C)
0 - 2 yrs	97.5 - 100.4 °F	36.4 - 38.0 °C
3 - 10 yrs	97.0 - 100.0 °F	36.1 - 37.8 °C
11 - 65 yrs	96.6 - 99.7 °F	35.9 - 37.6 °C
> 65 yrs	96.4 - 99.5 °F	35.8 - 37.5 °C

The range of normal temperature varies from person to person and can be influenced by many factors such as time of day, measure from different sites, level of activity, medications, emotion, etc. We recommend that you practice with the thermometer when you or family members are healthy, so you are familiar with the readings when they are ill.

## LCD Display

1. Day sign in time mode
2. A.M. and P.M. sign in time mode
3. Memory sign means the memory mode
4. Indicator of the series of records in memory mode
5. Year in time mode or the series of records in memory mode
6. Year sign
7. Month and day in time mode or memory mode
8. Month sign
9. Night sign in time mode
10. Sleep sign
11. Low battery indicator
12. Temperature in measuring mode or time in time mode
13. State of measurement sign; it displays while measuring temperature
14. Sign of temperature measured in ear
15. Sign of temperature measured on forehead
16. Measuring temperature sign
17. SET sign means the setting mode
18. Speaker sign
19. Day sign

## ***Temperature Taking Hints***

---

1. The right ear reading may differ from the reading taken at the left ear. Therefore, always take the temperature in the same ear.
2. The ear must be free from obstructions or excess earwax build-up to take an accurate reading
3. External factors may influence ear temperatures, including when an individual has:
  - Been lying on one ear or the other
  - Had his/her ears covered
  - Been exposed to very hot or very cold temperatures
  - Been recently swimming or bathingIn these cases, remove the individual from the situation and wait 20 minutes prior to taking a temperature
4. Use the untreated ear if prescription eardrops or other ear medications have been placed in the ear canal

## ***Preparations***

---

1. Replace any low battery to ensure full power supply.
2. For proper measurement, condition the thermometer for 30 minutes in measurement site.
3. An unexpected fluctuation of temperature measurement is likely to occur when reading is taken in front of an air conditioner or heater.
4. Be sure to keep ear or forehead clean for an accurate reading.
5. To achieve accurate measurements, it is very important to check that the infrared sensor is clean before using.
6. Measurement shall be done at a rest time to ensure an accurate result. Temperature is likely to increase immediately after exercise or bathing.
7. Fahrenheit/Celsius switch: Press **[MEM]** and Forehead measurement button at the same time to change between Fahrenheit and Celsius.

## How To Take Measurement


### 1. To measure temperature in the **EAR**:

Remove the dust cover from the thermometer and press any button to turn the unit on. Ensure the probe tip is clean before measuring. Place the probe tip into the ear canal. Press the Ear Button to start the measurement. You'll hear 1 beep and continue to hold. After the completion of measurement you'll hear two short beeps. Now you can read the result on the display or hear it through a speaker, if the TALK function is on. You can measure again after "000" disappears.

### 2. To measure temperature on the **FOREHEAD**:

Remove the dust cover from the thermometer and press any button to turn the unit on. Ensure the probe tip is clean before measuring. **Hold the thermometer so the probe tip is in the center of the forehead, midway between the eyebrows and hairline, but NOT touching the skin. Keep the probe tip approximately 0.25 in (1 cm) from the forehead.** Press the Forehead Button to start the measurement. You'll hear 1 beep and continue to hold. After the completion of measurement you'll hear two short beeps. Now you can read the result on the display or hear it through a speaker, if the TALK function is on.

- When body temperature is over 99.6°F (37.6°C), the Forehead measurement button will turn red color indicating a fever.
- The result will be stored into memory automatically. If 30 sets of memory are stored already, the first set of data will be removed, and the current measurement will be stored as the last record.
- Measure 3 times maximum in a consecutive way. It is possible that a proper body temperature measure is obstructed despite of a normal recovery. For continuous measurement up to 4 times or above, at least wait for 10 minutes before new measurement.
- Displaying the result of measurement:
  - Displaying the true value of the result for measurements between 89.6°F – 109.3°F (32.0°C and 42.9°C);
  - "HI" displayed if the results are higher than 109.3°F/42.9°C;
  - "LO" displayed if the results are lower than 89.6°F/32.0°C.

- After the measurement is completed, after 30 seconds, the thermometer automatically returns to the time mode.
- The display will turn off automatically after 1 minute if there are no buttons pressed. (When the sleep mode is activated).

## Setup Modes

---

Clock/Room Temperature mode:

- Display time in a sequence of year, month, date, hour, minute with a flickering ":" stands for second. The LCD displays current time and room temperature in alternate 5 seconds.
- Press **[MEM]** button to report the room temperature and current time.

Clock mode  
(Clock)


Clock mode  
(Room Temperature)

Clock mode  
(Battery low)

- Remove the battery cover and press and hold **[TALK]** button for 2 seconds to toggle the language which you want or non talking function.
- When the TALK function is on, a speaker symbol is on. If it is turned off, the speaker symbol will not appear on the display. The default language in the TALK function is English.

Clock  
(Non-talking)

Clock  
(Non-sleep mode)

<b>EN</b>	= English
<b>SP</b>	= Spanish
<b>FA</b>	= French
<b>GE</b>	= German
<b>I7</b>	= Italian
<b>AS</b>	= Russian
<b>--</b>	= Voice Off


- 5 Press **[MEM]** button to retrieve the recorded temperature or press and hold **[MEM]** button for 2 seconds to enter into set up mode.

Setup mode:

Press hold **[MEM]** button for 2 seconds to enter into setup mode.

1. Press **[MEM]** button to increase data or press hold **[MEM]** button to increase the data rapidly. Unit returns to clock mode automatically 30 seconds after last key is pressed.
2. Remove the battery cover and press **[TALK]** to shift between a setting of HOUR, MINUTE, YEAR, MONTH, DAY and SLEEP enable/ disable.

For Sleep mode:

3. The default setting is on. System enters into sleep mode after last key is pressed for 1 minute. **Notice:** If sleep mode is set to off, the thermometer will stay on and ready to use until batteries are drained. For saving power, we recommend setting sleep mode on.
4. To clear memory and all settings, remove batteries for 1 minute and re-install to begin set up.

## Memory mode

1. Press **[MEM]** button in clock mode to enter into memory mode;
2. “-” will be displayed if no data is stored;
3. Maximum storage is 30 sets of data with date and memory number shown;
4. Unit returns to clock mode automatically 5 seconds after last key pressed.

Note: Press any key to stop the speech during it is announcing.

Memo. no date

Memo. date


## Recommendations

1. For new users, it is recommended to measure the same ear 3 times for the first time. The maximum reading shall be recorded in case of different results since any thermometer suffering a continuous operation is deemed to see a slight error of ( $\pm 0.5^{\circ}\text{F}$ ,  $\pm 0.3^{\circ}\text{C}$ ).
2. Wait about 30 seconds before measuring the same patient again to avoid excessive cooling of the skin.
3. For users in sound health, both ears are almost the same in terms of body temperature. Please select the ear that shows a stable and relatively higher temperature.
4. Please clean the ear before measurement.

## Battery Replacement

---

Remove the battery cover. Use only new alkaline batteries (2x AAA). Put the battery cover back tightly.

## Maintenance, Storage and Calibration

---

### Maintenance:

1. Remove all stains from the device with a soft dry cloth.
2. Clean the end of the probe as follows: wipe its surface gently with a swab of cotton wool or soft cloth moistened with alcohol. The thermometer can be used after the alcohol dries.
3. Do not wash the thermometer with water, abrasive detergent or detergent containing benzene.

### Storage:

1. After each use clean the thermometer.
2. Replace the dust cover and always store with the cover on.
3. Do not leave the thermometer exposed to direct sun, high temperatures, humidity, fire, flames, shocks or hits.
4. If you do not use the thermometer for a long time, take out the batteries.

### Calibration:

1. The thermometer is factory-calibrated. If the thermometer is used according to the user manual, periodic calibration is not required. Do not try to modify or reassemble the thermometer.

## Troubleshooting

---

<b>Problem</b>	<b>Check lists</b>	<b>Action</b>
No response/ Automatic reset when pull out insulator	Battery used up?	Change new battery
	Battery in wrong polarity?	Take out battery, reinsert battery correctly
	Poor battery contact	
Battery symbol on LCD	Low Battery	Change new battery
"Lo" on screen (result below 89.6°F or 32.0°C)	Thermometer correctly pointed at the eardrum / forehead?	Follow user's manual to point at the eardrum / forehead.
"Hi" on screen (result over/exceed 109.3°F or 42.9°C)	Please check the operation method	Follow user's manual for proper measurement.

## Specifications

---

- Product Name: Digital Forehead and Ear Thermometer
- Power supply: DC3V (2 x AAA battery)
- Power consumption: < 50mA (for voice prompt)
- Temperature sensing Part: IR sensor
- Measuring range: 89.6°F – 109.3°F (32.0°C – 42.9°C)  
Room Temperature: 41.0°F – 139.8°F (5.0°C – 59.9°C)
- Measuring accuracy:  $\pm 0.4^{\circ}\text{F}$  (or  $\pm 0.2^{\circ}\text{C}$ )  
from 95.9°F – 107.6°F (or 35.5°C – 42.0°C)  $\pm 0.5^{\circ}\text{F}$  ( $\pm 0.3^{\circ}\text{C}$ )  
from 89.6°F – 95.8°F (or 32.0°C – 35.4°C)  
from 107.7°F – 109.3°F (or 42.1°C – 42.9°C)  
Room Temperature:  $\pm 2^{\circ}\text{F}$  ( $\pm 1^{\circ}\text{C}$ )
- Outer dimension: 120(L) x 40(W) x 65(D) mm
- Weight: 60g (without batteries)
- Working environment: 50.0°F – 104.0°F (10°C – 40°C) For Ear  
59.0°F – 104.0°F (15°C – 40°C) For Forehead
- Relative humidity: 20% RH – 85% RH
- Atmosphere Pressure: 700hPa-1060hPa
- Transportaion and storage condition:  
Temperature: 14°F - 131°F (-10°C to 55°C)  
Relative Humidity: 20%-90%  
Atmosphere Pressure: 700hPa-1060hPa

## Normalized Symbols

---

Refer to instruction manual

Applied part of type BF

Electrical devices are recyclable material and should not be disposed of with household waste after their useful life! Help us to protect the environment and save resources and take this device to the appropriate collection points. Please contact the organization which is responsible for waste disposal in your area if you have any questions.

## Guidance and Manufacturer's Declaration

Guidance and manufacturer's declaration – electromagnetic emissions		
<b>The device is intended for use in the electromagnetic environment specified below. The customer or the user assures that it is used in such an environment.</b>		
Emissions test Compliance Electromagnetic environment - guidance		
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR11	Class B	The device is suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage <b>fluctuations / flicker</b> emissions IEC 61000-3-3	Not applicable	
RF emissions CISPR 14-1	Complies	The device is not suitable for interconnection with other equipment.
RF emissions CISPR 15	Complies	The device is not suitable for interconnection with other equipment.

**Guidance and manufacturer's declaration – electromagnetic immunity**

**The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.**

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential mode ±2kV common mode	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 Cycles 70% UT (30% dip in UT) for 25 Cycles <5% UT (>95% dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

**Guidance and manufacturer's declaration – electromagnetic immunity**

**The device is intended for use in the electromagnetic environment specified below. The customer or the user should assure that it is used in such an environment.**

			Portable and mobile RF Communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance.
Conducted R F IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	$d = \left[ \frac{3.5}{V1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[ \frac{3.5}{E1} \right] \sqrt{P}$ 80 MHz to 800 MHz $d = \left[ \frac{7}{E1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz
			<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the Transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range. b Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

NOTE 1: At 80 MHz ends 800 MHz. the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**a)** Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To

assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**b)** Over the frequency range 150 kHz to 80 MHz, field strengths should be less than  $[V_i]$  V/m.

### Recommended separation distances between portable and mobile RF communications equipment and the device

**The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the as recommended below, according to the maximum output power of the communications equipment.**

Rated maximum transmitter output power of W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{\sqrt{P}} \right] \sqrt{P}$	150 kHz to 80 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	150 kHz to 80 MHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0,01	0.12	0.12	0.23
0,1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

## ***Warranty***

---

One-year warranty is available from purchasing date, excluding user-caused failures listed below:

- (1) Failure resulted in unauthorized disassembly and modification
- (2) Failure resulted in unexpected drop during application or transportation
- (3) Failure resulted in operation away from proper instruction in user's manual

## ***Declaration of Conformity***

---

IEC60601-1, IEC60601-1-2, ASTM E1965-2009 revised, ISO10993-5, ISO10993-10, ISO10993-1, ISO 14971, EN 12470-5.

**NOTES:**

1) Under the environment with electrostatic discharge, the unit may malfunction and may require user to reset the unit.

2) Dispose with care, please consult the retailer for details.

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©2017 HealthSmart® International