



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

Terrats Medical SL  
% Linda Schulz  
Regulatory Affairs  
PaxMed International, LLC  
12264 El Camino Real, Suite 400  
San Diego, California 92130

August 8, 2017

Re: K170588  
Trade/Device Name: DESS Dental Smart Solutions  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: Class II  
Product Code: NHA  
Dated: June 6, 2017  
Received: July 7, 2017

Dear Linda Schulz:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely,

  
**Mary S. Runner -S**

for  
Lori Wiggins, MPT, CLT  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)

K170588

Device Name

DESS Dental Smart Solutions

Indications for Use (Describe)

DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.

All digitally designed custom abutments for use with TiBase or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.

Compatible Implant Systems

Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)
3i Certain <sup>®</sup>	3.25, 4.0, 5.0	3.4, 4.1, 5.0
3i OSSEOTITE <sup>®</sup>	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0
OsseoSpeed <sup>™</sup>	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5
NobelActive <sup>®</sup>	3.5, 4.3, 5.0	NP, RP
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP
Straumann <sup>®</sup> Bone Level	3.3, 4.1, 4.8	NC, RC
Straumann <sup>®</sup> Tissue Level	3.3, 4.1, 4.8	RN, WN
Tapered Screw-Vent <sup>®</sup>	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

**510(k) Summary**  
**Terrats Medical SL**  
**DESS Dental Smart Solutions**  
K170588  
August 8, 2017

**ADMINISTRATIVE INFORMATION**

Manufacturer Name	Terrats Medical SL Avenida La Ferrería (Pol Ind La Ferrería) 62 Montcada i Reixac, 08110 Spain Telephone +34 93 564 60 06 Fax +34 93 564 73 17
Official Contact	Roger Terrats, COO
Representative/Consultant	Linda Schulz, BSDH, RDH Kevin Thomas, PhD PaxMed International, LLC 12264 El Camino Real, Suite 400 San Diego, CA 92130 Telephone: +1-858-792-1235 Fax: +1-858-792-1236 Email: LSchulz@paxmed.com KThomas@paxmed.com

**DEVICE NAME AND CLASSIFICATION**

Trade/Proprietary Name	DESS Dental Smart Solutions
Common Name	Dental implant abutment
Classification Name	Endosseous dental implant abutment
Classification Regulations	21 CFR 872.3630, Class II
Product Code	NHA
Classification Panel	Dental Products Panel
Reviewing Branch	Dental Devices Branch

**PREDICATE DEVICE INFORMATION**

Primary Predicate  
K120414, OsseoSpeed™ Plus, Astra Tech AB

Reference Predicates  
K072878, Modification to: Locator Implant Anchor, Zest Anchors, Inc.  
K092341, Low Profile Abutment, Biomet 3i  
K150203, Medentika CAD/CAM Abutments, Medentika GmbH  
K150367, Neodent Implant System, JJGC Indústria e Comércio de Materiais Dentários SA

Compatible Implant System Predicates

K063341	3i OSSEOTITE Certain® Dental Implants	Implant Innovations, Inc.
K063286	OSSEOTITE® Dental Implants	Implant Innovations, Inc.
K101732	OsseoSpeed™	Astra Tech AB
K073075	FRIADENT Implant Systems	DENTSPLY International, Inc.
K142260	NobelActive®	Nobel Biocare AB
K073142	NobelReplace Hexagonal Implants	Nobel Biocare AB
K050705	TiUnite® Implants	Nobel Biocare AB
K050406	NOBELSPEEDY™ Implants	Nobel Biocare USA LLC
K022562	Various Brånemark System Implants – Immediate Function Indication	Nobel Biocare AB
K140878	Straumann® Bone Level Tapered Implants	Straumann USA, LLC
K062129	P.004 Implants	Institut Straumann AG
K130222	Straumann® Dental Implant System SLActive and Roxolid Product Families	Straumann USA, LLC
K112160	Tapered Screw-Vent® X Implant	Zimmer Dental, Incorporated

INDICATIONS FOR USE

DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.

All digitally designed custom abutments for use with TiBase or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.

Compatible Implant Systems

Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5
NobelActive®	3.5, 4.3, 5.0	NP, RP
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7

DEVICE DESCRIPTION

DESS Dental Smart Solutions is a dental implant abutment system that includes seven abutment design types (Healing, Temporary, Straight, TiBase, Pre-milled Blank, DESS LOC, Multi-Unit), and ten abutment connections compatible with eleven implant systems. Platform diameters range from 3.3 mm to 5.7 mm. Corresponding implant body diameters range from 3.25 mm to 6.0 mm. The following table outlines the body and platform diameters by abutment design and corresponding implant line.



Abutments are offered in a variety of connection types to enable compatibility with a large number of currently marketed implants. Straight, Temporary and TiBase abutments have a SelectGrip® surface. DESS LOC Abutments have a ZrN coating. Selected DESS screws include DLC coating. DESS Dental Smart Solutions abutments are straight abutments. All abutments are provided non-sterile.

**PERFORMANCE DATA**

Non-clinical testing data submitted, referenced, or relied upon to demonstrate substantial equivalence included: sterilization validation to an SAL of 10<sup>-6</sup> according to ISO 17665-1 and ISO 17665-2 to ensure sterilization of the final finished device; biological evaluation according to ISO 10993-1 and cytotoxicity testing according to ISO 10993-5 for all surfaces to demonstrate that all devices are non-cytotoxic; SEM evaluation and measurement of the ZrN and SelectGrip surfaces to demonstrate suitability of the surface; and engineering and dimensional analysis of OEM implant bodies, OEM abutments, and OEM abutment fixation screws to confirm compatibility.

No clinical data were included in this submission.

**EQUIVALENCE TO MARKETED DEVICE**

The subject device is substantially equivalent in indications and design principles to the predicate devices shown above. Below are summary tables comparing the Indications for Use and the technological characteristics of the subject device and the predicate devices.

Comparison of Indications for Use Statements

Indications for Use Statement																																					
<b>Subject Device</b>																																					
DESS Dental Smart Solutions	DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations. All digitally designed custom abutments for use with TiBase or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.																																				
Terrats Medical SL	<p style="text-align: center;">Compatible Implant Systems</p> <table border="1"> <thead> <tr> <th>Implant System Compatibility</th> <th>Implant Diameter (mm)</th> <th>Platform Diameter (mm)</th> </tr> </thead> <tbody> <tr> <td>3i Certain®</td> <td>3.25, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>3i OSSEOTITE®</td> <td>3.25, 3.75, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>OsseoSpeed™</td> <td>3.5, 4.0, 5.0</td> <td>3.5/4.0, 4.5/5.0</td> </tr> <tr> <td>FRIADENT XiVE</td> <td>3.4, 3.8, 4.5</td> <td>3.4, 3.8, 4.5</td> </tr> <tr> <td>NobelActive®</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>NobelReplace Conical</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>Nobel Replace Trilobe</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Brånemark</td> <td>3.5, 3.75/4.0, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Straumann® Bone Level</td> <td>3.3, 4.1, 4.8</td> <td>NC, RC</td> </tr> <tr> <td>Straumann® Tissue Level</td> <td>3.3, 4.1, 4.8</td> <td>RN, WN</td> </tr> <tr> <td>Tapered Screw-Vent®</td> <td>3.7, 4.1, 4.7, 6.0</td> <td>3.5, 4.5, 5.7</td> </tr> </tbody> </table>	Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)	3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0	3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0	OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0	FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5	NobelActive®	3.5, 4.3, 5.0	NP, RP	NobelReplace Conical	3.5, 4.3, 5.0	NP, RP	Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP	Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP	Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC	Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN	Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7
Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)																																			
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0																																			
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0																																			
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0																																			
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5																																			
NobelActive®	3.5, 4.3, 5.0	NP, RP																																			
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP																																			
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP																																			
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP																																			
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC																																			
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN																																			
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7																																			

<b>Indications for Use Statement</b>																																													
<b>Primary Predicate Device</b>																																													
K120414 OsseoSpeed™ Plus Astra Tech AB	<p>Implants: The Astra Tech Dental Implants are intended for both one- and two-stage surgical procedures in the following situations and with the following clinical protocols:</p> <ul style="list-style-type: none"> <li>• replacing single and multiple missing teeth in the mandible and maxilla,</li> <li>• immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge,</li> <li>• especially indicated for use in soft bone applications where implants with other implant surface treatments may be less effective,</li> <li>• immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm, or in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.</li> </ul> <p>The intended use for OsseoSpeed™ Plus 3.0S is limited to replacement of maxillary lateral incisors and mandibular incisors.</p> <p>Abutments: Astra Tech Implant System Plus abutments are intended to be used in conjunction with Astra Tech Implant System Plus in fully edentulous or partially edentulous maxillary and/or mandibular arches to provide support for crowns, bridges or overdentures.</p> <p>Atlantis Abutments: The Atlantis™ Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prosthesis can be cemented, screw retained or friction fit to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p> <p>The Atlantis™ Crown Abutment in Zirconia is intended for use with an endosseous implant to function as a substructure that also serves as the final restoration, in partially or completely edentulous patients. The prosthesis is screw retained. The abutment screw is intended to secure the crown abutment to the endosseous implant.</p>																																												
<b>Reference Predicate Devices</b>																																													
K072878 Modification to: Locator Implant Anchor Zest Anchors, Inc.	Not available																																												
K092341 Low Profile Abutment Biomet 3i, Inc.	Biomet 3i Low Profile Abutments are intended for use as accessories to endosseous dental implants to support a prosthetic device in a partially or completely edentulous patient. A dental abutment is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis is screw retained to the abutment.																																												
K150203 Medentika CAD/CAM Abutments Medentika GmbH	<p>Medentika Preface CAD/CAM Abutments are intended for use with dental implants as a support for single or multiple tooth prostheses in the maxilla or mandible of a partially or fully edentulous patient.</p> <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Implant System Compatibility</th> <th style="text-align: left;">Series</th> <th style="text-align: left;">Implant Diameter (mm)</th> <th style="text-align: left;">Platform Diameter (mm)</th> </tr> </thead> <tbody> <tr> <td>Nobel Biocare Replace™ Select</td> <td>E</td> <td>3.5, 4.3, 5.0, 6.0</td> <td>3.5, 4.3, 5.0, 6.0</td> </tr> <tr> <td>Nobel Biocare NobelActive™</td> <td>F</td> <td>3.0, 3.5, 4.3, 5.0</td> <td>3.0, 3.5, 3.9 (4.3), 3.9 (5.0)</td> </tr> <tr> <td>Biomet 3i Osseotite® Certain®</td> <td>H</td> <td>3.25, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>Biomet 3i Osseotite®</td> <td>I</td> <td>3.25, 3.75, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>Nobel Biocare Brånemark</td> <td>K</td> <td>3.3, 3.75, 4.0, 5.0</td> <td>3.5, 4.1, 4.1, 5.1</td> </tr> <tr> <td>Straumann Bone Level</td> <td>L</td> <td>3.3, 4.1, 4.8</td> <td>3.3, 4.1, 4.8</td> </tr> <tr> <td>Straumann Standard</td> <td>N</td> <td>3.3, 4.1, 4.8</td> <td>3.5 (NNC), 4.8, 6.5</td> </tr> <tr> <td>Zimmer Tapered Screw-vent®</td> <td>R</td> <td>3.3, 3.7, 4.1, 4.7, 6.0</td> <td>3.5, 4.5, 5.7</td> </tr> <tr> <td>Astra Tech OsseoSpeed™</td> <td>S</td> <td>3.0, 3.5, 4.0, 4.5, 5.0</td> <td>3.0, 3.5, 4.0, 4.5, 5.0</td> </tr> <tr> <td>Dentsply Friadent® Frialit/XiVE®</td> <td>T</td> <td>3.4, 3.8, 4.5, 5.5</td> <td>3.4, 3.8, 4.5, 5.5</td> </tr> </tbody> </table> <p>Medentika PreFace is intended for use with the Straumann® CARES® System. All digitally designed abutments for use with Medentika CAD/CAM Abutments are intended to be manufactured at a Straumann® CARES® validated milling center.</p>	Implant System Compatibility	Series	Implant Diameter (mm)	Platform Diameter (mm)	Nobel Biocare Replace™ Select	E	3.5, 4.3, 5.0, 6.0	3.5, 4.3, 5.0, 6.0	Nobel Biocare NobelActive™	F	3.0, 3.5, 4.3, 5.0	3.0, 3.5, 3.9 (4.3), 3.9 (5.0)	Biomet 3i Osseotite® Certain®	H	3.25, 4.0, 5.0	3.4, 4.1, 5.0	Biomet 3i Osseotite®	I	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0	Nobel Biocare Brånemark	K	3.3, 3.75, 4.0, 5.0	3.5, 4.1, 4.1, 5.1	Straumann Bone Level	L	3.3, 4.1, 4.8	3.3, 4.1, 4.8	Straumann Standard	N	3.3, 4.1, 4.8	3.5 (NNC), 4.8, 6.5	Zimmer Tapered Screw-vent®	R	3.3, 3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7	Astra Tech OsseoSpeed™	S	3.0, 3.5, 4.0, 4.5, 5.0	3.0, 3.5, 4.0, 4.5, 5.0	Dentsply Friadent® Frialit/XiVE®	T	3.4, 3.8, 4.5, 5.5	3.4, 3.8, 4.5, 5.5
Implant System Compatibility	Series	Implant Diameter (mm)	Platform Diameter (mm)																																										
Nobel Biocare Replace™ Select	E	3.5, 4.3, 5.0, 6.0	3.5, 4.3, 5.0, 6.0																																										
Nobel Biocare NobelActive™	F	3.0, 3.5, 4.3, 5.0	3.0, 3.5, 3.9 (4.3), 3.9 (5.0)																																										
Biomet 3i Osseotite® Certain®	H	3.25, 4.0, 5.0	3.4, 4.1, 5.0																																										
Biomet 3i Osseotite®	I	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0																																										
Nobel Biocare Brånemark	K	3.3, 3.75, 4.0, 5.0	3.5, 4.1, 4.1, 5.1																																										
Straumann Bone Level	L	3.3, 4.1, 4.8	3.3, 4.1, 4.8																																										
Straumann Standard	N	3.3, 4.1, 4.8	3.5 (NNC), 4.8, 6.5																																										
Zimmer Tapered Screw-vent®	R	3.3, 3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7																																										
Astra Tech OsseoSpeed™	S	3.0, 3.5, 4.0, 4.5, 5.0	3.0, 3.5, 4.0, 4.5, 5.0																																										
Dentsply Friadent® Frialit/XiVE®	T	3.4, 3.8, 4.5, 5.5	3.4, 3.8, 4.5, 5.5																																										
K150367 Neodent Implant System JJGC Indústria e Comércio de Materiais Dentários SA	<p>Titanium Base Abutment is a titanium base placed onto Neodent dental implants to provide support for customized prosthetic restorations. It is used with a coping and crown, or crown alone, and is indicated for cement-retained single or multi-unit restorations, or screw-retained single restorations.</p> <p>PreFace Abutment is a titanium abutment to be used in fabricating a full custom abutment and placed onto Neodent dental implants to provide support for customized prosthetic restorations. PreFace Abutments are indicated for screw-retained single restorations or cement-retained single or multi-unit restorations.</p>																																												



Comparison of Technological Characteristics

	Subject Device	Primary Predicate	Reference Predicates			
	DESS Dental Smart Solutions  Terrats Medical SL	K120414  OsseoSpeed™ Plus  Astra Tech AB	K072878  Modification to: Locator Implant Anchor  Zest Anchors, Inc.	K092341  Low Profile Abutment  Biomet 3i, Inc.	K150203  Medentika CAD/CAM Abutments  Medentika GmbH	K150367  Neodent Implant System  JJGC Indústria e Comércio de Materiais Dentários SA
<b>Design</b>						
Prosthesis Attachment	Cement-retained Screw-retained	Cement-retained Screw-retained	Overdenture attachment	Screw-retained	Cement-retained	Cement-retained Screw-retained
Restoration	Single-unit Multi-unit	Single-unit Multi-unit	Overdenture	Single-unit Multi-unit	Single-unit Multi-unit	Single-unit Multi-unit
Abutment Platform Diameter (mm)	3.4 – 5.7	3.0 – 5.4	2.5 – 6.5	3.4 – 5.0	3.0 – 7.0	3.0 – 6.0
Abutment Angle	Straight	Straight to 30°	Straight	Straight to 30°	Straight to 30°	Straight to 30°
Abutment/Implant Interface	Internal, External	Internal	Internal, External	Internal, External	Internal, External	Internal
<b>Material</b>						
Abutment	Ti-6Al-4V	Ti-6Al-4V Zirconia, Gold, PEEK	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V
Screw	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V	Ti-6Al-4V

DESS Dental Smart Solutions abutments are substantially equivalent in design, function, material, size, and Indications for Use to OsseoSpeed Plus (K120414) and Medentika CAD/CAM Abutments (K150203). All are intended for use with endosseous dental implants in the maxilla and mandible to provide prosthetic support. Digital files for DESS Dental Smart Solutions abutments and for Medentika CAD/CAM Abutments are to be sent to a validated milling center for manufacture. Differences in the type of restoration named or specific milling center stated in the Indications for Use statement do not affect the intended use.

Subject device abutment designs and function are substantially equivalent to design and function of abutments included in K120414, K072878, K092341, K150203 and K150367. Implant/abutment interface compatibility for the subject device is substantially equivalent to the compatible implant system predicates listed above.

CONCLUSION

The subject device and the predicate devices have the same intended use, have similar technological characteristics, and are made of similar materials. The subject device and predicate devices encompass the same range of physical dimensions, including diameter and design of the abutments. The subject and predicate devices are packaged in similar materials and are to be sterilized using similar methods.

The data included in this submission demonstrate substantial equivalence to the predicate devices listed above.



June 15, 2018

Terrats Medical SL  
% Linda Schulz  
Regulatory Affairs  
PaxMed International, LLC  
12264 El Camino Real, Suite 400  
San Diego, California 92130

Re: K173908

Trade/Device Name: DESS Dental Smart Solutions  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: Class II  
Product Code: NHA  
Dated: May 16, 2018  
Received: May 17, 2018

Dear Linda Schulz:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

  
**Andrew I. Steen -S**

for Tina Kiang, Ph.D.  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory,  
Infection Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)

K173908

Device Name

DESS Dental Smart Solutions

Indications for Use (Describe)

DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.

All digitally designed custom abutments for use with Aurum™ Abutment or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.

Compatible Implant Systems

Implant System Compatibility	Implant Body	Implant Platform
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5
NobelActive®	3.5, 4.3, 5.0	NP, RP
NobelReplace® Conical	3.5, 4.3, 5.0	NP, RP
NobelReplace® Trilobe	3.5, 4.3, 5.0	NP, RP, WP
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRASStaff@fda.hhs.gov

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

**510(k) Summary**  
**Terrats Medical SL**  
**DESS Dental Smart Solutions**

June 15, 2018

**ADMINISTRATIVE INFORMATION**

Manufacturer Name	Terrats Medical SL Carrer Mogoda 75-99 Barberà del Vallès 08210 Barcelona, Spain Telephone +34 93 564 60 06 Fax +34 93 564 73 17
Official Contact	Roger Terrats, COO
Representative/Consultant	Linda Schulz, BSDH, RDH Floyd G. Larson, MS, MBA PaxMed International, LLC 12264 El Camino Real, Suite 400 San Diego, CA 92130 Telephone: +1-858-792-1235 Fax: +1-858-792-1236 Email: LSchulz@paxmed.com FLarson@paxmed.com

**DEVICE NAME AND CLASSIFICATION**

Trade/Proprietary Name	DESS Dental Smart Solutions
Common Name	Dental implant abutment
Classification Name	Endosseous dental implant abutment
Classification Regulations	21 CFR 872.3630, Class II
Product Code	NHA
Classification Panel	Dental Products Panel
Reviewing Branch	Dental Devices Branch

**PREDICATE DEVICE INFORMATION**

Primary Predicate  
K170588, DESS Dental Smart Solutions, Terrats Medical SL

**Reference Device**

K092341	Low Profile Abutment	Biomet 3i, Inc.
K150669	Neoss TiBase and CoCr Abutments	Neoss Ltd.
K120414	OsseoSpeed™ Plus	Astra Tech AB
K160784	CAM Titanium Blanks	Altatec GmbH

Compatible Implant Systems

K063341	3i OSSEOTITE Certain® Dental Implants	Implant Innovations, Inc.
K063286	OSSEOTITE® Dental Implants	Implant Innovations, Inc.
K101732	OsseoSpeed™	Astra Tech AB
K073075	FRIADENT Implant Systems	DENTSPLY International, Inc.
K142260	NobelActive®	Nobel Biocare AB
K073142	NobelReplace Hexagonal Implants	Nobel Biocare AB
K050705	TiUnite® Implants	Nobel Biocare AB
K050406	NOBELSPEEDY™ Implants	Nobel Biocare USA LLC
K022562	Various Brånemark System Implants – Immediate Function Indication	Nobel Biocare AB
K140878	Straumann® Bone Level Tapered Implants	Straumann USA, LLC
K062129	P.004 Implants	Institut Straumann AG
K130222	Straumann® Dental Implant System SLActive and Roxolid Product Families	Straumann USA, LLC
K112160	Tapered Screw-Vent® X Implant	Zimmer Dental, Incorporated

INDICATIONS FOR USE

DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations.

All digitally designed custom abutments for use with Aurum™ Abutment or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.

Compatible Implant Systems

<b>Implant System Compatibility</b>	<b>Implant Body</b>	<b>Implant Platform</b>
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5
NobelActive®	3.5, 4.3, 5.0	NP, RP
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7

DEVICE DESCRIPTION

DESS Dental Smart Solutions subject devices include four abutment design types (Aurum Base, Pre-milled Blank, CoCr Pre-milled Blank, CoCr Abutment) and one screw type (Aurum Base Screw). Abutments are provided in ten abutment connections compatible with eleven implant systems. Implant platform diameters range from 3.3 mm to 6.5 mm. Corresponding implant body diameters range from 3.25 mm to 6.0 mm. All abutments are provided non-sterile. The following table outlines the body and platform diameters by abutment design and corresponding compatible implant line.

Summary of Abutment Designs

DESS Abutment System	CAD / CAM				Cast-To		OEM Implant System		Connection
	AURUM base™ Non-Engaging Titanium Alloy ASTM F136	AURUM base™ Engaging Titanium Alloy ASTM F136	CoCr Pre-milled Blank Engaging CoCr ASTM 1537	Ti Pre-milled Blank Engaging Titanium Alloy ASTM F136	CoCr Abutment Non-Engaging CoCr ASTM 1537	CoCr Abutment Engaging CoCr ASTM 1537			
Internal Hex "Click"	NP RP WP	NP RP WP	NP RP WP		NP RP WP	NP RP WP	3i Certain®	3.4 (NP) 4.1 (RP) 5.0 (WP)	Internal
External Hex USA	NP RP WP	NP RP WP	NP RP WP		NP RP WP	NP RP WP	3i OSSEOTITE®	3.4 (NP) 4.1 (RP) 5.0 (WP)	External
Internal Hex Conic	RP WP	RP WP	RP WP		RP WP	RP WP	OsseoSpeed™	3.5/4.0 (RP) 4.5/5.0 (WP)	Internal
Internal Hex FD	NP RP WP	NP RP WP	NP RP WP	NP RP WP	NP RP WP	NP RP WP	FRIADENT XiVE®	3.4 (NP) 3.8 (RP) 4.5 (WP)	Internal
Active Hex	NP RP	NP RP	NP RP		NP RP	NP RP	NobelActive®	3.5 (NP) 3.9 (RP)	Internal
							NobelReplace® Conical	3.5 (NP) 3.9 (RP)	Internal
Tri-Lobe	NP RP WP	NP RP WP	NP RP WP		NP RP WP	NP RP WP	NobelReplace® Trilobe	3.5 (NP) 4.3 (RP) 5.0 (WP)	Internal
External Hex Universal	NP RP	NP RP	NP RP WP	NP	NP RP WP	NP RP WP	Brånemark System®	3.5 (NP) 4.1 (RP) 5.1 (WP)	External
Conical BL	NC RC	NC RC	NC RC		NC RC	NC RC	Straumann® Bone Level	3.3 (NC) 4.1 (RC)	Internal
Octagon	RN WN	RN WN	RN WN		RN WN	RN WN	Straumann® Tissue Level	4.8 (RN) 6.5 (WN)	Internal
Internal Hex USA	NP RP WP	NP RP WP	NP RP WP		NP RP WP	NP RP WP	Tapered Screw-Vent	3.5 (NP) 4.5 (RP) 5.7 (WP)	Internal

Aurum Abutment is a titanium abutment that can be used for a direct multi-unit restoration or to support a zirconia superstructure plus a single-unit or multi-unit restoration. Design parameters for the zirconia superstructure are a minimum wall thickness of 0.4 mm, a minimum post height for single-unit restorations of 4.0 mm, and a maximum gingival height of 6.0 mm. All zirconia superstructures are for straight abutments only.

Pre-milled Blank abutments are cylindrical abutments designed for custom abutment fabrication by a CAD/CAM process. All patient-specific abutment fabrication is by prescription on the order of the clinician. Pre-milled Blanks are made of cobalt-chromium alloy or titanium alloy. Design parameters for the Pre-Milled Blanks are, a minimum wall thickness of 0.45 mm, a minimum post height for single-unit restorations of 4.0 mm, a maximum gingival height of 6.0 mm, and a maximum total abutment height of 19 mm. All Pre-Milled Blanks are for straight abutments only.

CoCr Abutments are designed as a cast-to abutment for support of a single-unit or a multi-unit restoration. They are made of cobalt-chromium alloy. Design parameters for CoCr Abutments are, a minimum wall thickness of 0.4 mm, a minimum post height for single-unit restorations of 4.0 mm, and a maximum gingival height of 6.0 mm. All CoCr Abutments are for straight abutments only.

**PERFORMANCE DATA**

Non-clinical testing data submitted to demonstrate substantial equivalence included: sterilization validation to an SAL of 10<sup>-6</sup> according to ISO 17665-1 and ISO 17665-2 to ensure sterilization of the final finished device; biological evaluation according to ISO 10993-1 and cytotoxicity testing according to ISO 10993-5 to demonstrate that all devices are non-cytotoxic, and compatibility analysis by reference to K170588.

No clinical data were included in this submission.

**EQUIVALENCE TO MARKETED DEVICE**

The subject device is substantially equivalent in indications and design principles to the predicate devices shown above. Below are summary tables comparing the Indications for Use and the technological characteristics of the subject device and the predicate devices.

Comparison of Indications for Use Statements

	Indications for Use Statement																																				
<b>Subject Device</b>																																					
DESS Dental Smart Solutions Terrats Medical SL	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations. All digitally designed custom abutments for use with Aurum Base or Pre-Milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p style="text-align: center;">Compatible Implant Systems</p> <table border="1"> <thead> <tr> <th>Implant System Compatibility</th> <th>Implant Body</th> <th>Implant Platform</th> </tr> </thead> <tbody> <tr> <td>3i Certain®</td> <td>3.25, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>3i OSSEOTITE®</td> <td>3.25, 3.75, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>OsseoSpeed™</td> <td>3.5, 4.0, 5.0</td> <td>3.5/4.0, 4.5/5.0</td> </tr> <tr> <td>FRIADENT XiVE</td> <td>3.4, 3.8, 4.5</td> <td>3.4, 3.8, 4.5</td> </tr> <tr> <td>NobelActive®</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>NobelReplace Conical</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>Nobel Replace Trilobe</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Brånemark</td> <td>3.5, 3.75/4.0, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Straumann® Bone Level</td> <td>3.3, 4.1, 4.8</td> <td>NC, RC</td> </tr> <tr> <td>Straumann® Tissue Level</td> <td>3.3, 4.1, 4.8</td> <td>RN, WN</td> </tr> <tr> <td>Tapered Screw-Vent®</td> <td>3.7, 4.1, 4.7, 6.0</td> <td>3.5, 4.5, 5.7</td> </tr> </tbody> </table>	Implant System Compatibility	Implant Body	Implant Platform	3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0	3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0	OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0	FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5	NobelActive®	3.5, 4.3, 5.0	NP, RP	NobelReplace Conical	3.5, 4.3, 5.0	NP, RP	Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP	Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP	Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC	Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN	Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7
Implant System Compatibility	Implant Body	Implant Platform																																			
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0																																			
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0																																			
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0																																			
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5																																			
NobelActive®	3.5, 4.3, 5.0	NP, RP																																			
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP																																			
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP																																			
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP																																			
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC																																			
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN																																			
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7																																			
<b>Primary Predicate Devices</b>																																					
K170588 DESS Dental Smart Solutions Terrats Medical SL	<p>DESS Dental Smart Solutions abutments are intended to be used in conjunction with endosseous dental implants in the maxillary or mandibular arch to provide support for prosthetic restorations. All digitally designed custom abutments for use with TiBase or Pre-milled Blank are to be sent to a Terrats Medical validated milling center for manufacture.</p> <p style="text-align: center;">Compatible Implant Systems</p> <table border="1"> <thead> <tr> <th>Implant System Compatibility</th> <th>Implant Diameter (mm)</th> <th>Platform Diameter (mm)</th> </tr> </thead> <tbody> <tr> <td>3i Certain®</td> <td>3.25, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>3i OSSEOTITE®</td> <td>3.25, 3.75, 4.0, 5.0</td> <td>3.4, 4.1, 5.0</td> </tr> <tr> <td>OsseoSpeed™</td> <td>3.5, 4.0, 5.0</td> <td>3.5/4.0, 4.5/5.0</td> </tr> <tr> <td>FRIADENT XiVE</td> <td>3.4, 3.8, 4.5</td> <td>3.4, 3.8, 4.5</td> </tr> <tr> <td>NobelActive®</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>NobelReplace Conical</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP</td> </tr> <tr> <td>Nobel Replace Trilobe</td> <td>3.5, 4.3, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Brånemark</td> <td>3.5, 3.75/4.0, 5.0</td> <td>NP, RP, WP</td> </tr> <tr> <td>Straumann® Bone Level</td> <td>3.3, 4.1, 4.8</td> <td>NC, RC</td> </tr> <tr> <td>Straumann® Tissue Level</td> <td>3.3, 4.1, 4.8</td> <td>RN, WN</td> </tr> <tr> <td>Tapered Screw-Vent®</td> <td>3.7, 4.1, 4.7, 6.0</td> <td>3.5, 4.5, 5.7</td> </tr> </tbody> </table>	Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)	3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0	3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0	OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0	FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5	NobelActive®	3.5, 4.3, 5.0	NP, RP	NobelReplace Conical	3.5, 4.3, 5.0	NP, RP	Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP	Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP	Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC	Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN	Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7
Implant System Compatibility	Implant Diameter (mm)	Platform Diameter (mm)																																			
3i Certain®	3.25, 4.0, 5.0	3.4, 4.1, 5.0																																			
3i OSSEOTITE®	3.25, 3.75, 4.0, 5.0	3.4, 4.1, 5.0																																			
OsseoSpeed™	3.5, 4.0, 5.0	3.5/4.0, 4.5/5.0																																			
FRIADENT XiVE	3.4, 3.8, 4.5	3.4, 3.8, 4.5																																			
NobelActive®	3.5, 4.3, 5.0	NP, RP																																			
NobelReplace Conical	3.5, 4.3, 5.0	NP, RP																																			
Nobel Replace Trilobe	3.5, 4.3, 5.0	NP, RP, WP																																			
Brånemark	3.5, 3.75/4.0, 5.0	NP, RP, WP																																			
Straumann® Bone Level	3.3, 4.1, 4.8	NC, RC																																			
Straumann® Tissue Level	3.3, 4.1, 4.8	RN, WN																																			
Tapered Screw-Vent®	3.7, 4.1, 4.7, 6.0	3.5, 4.5, 5.7																																			



Reference Devices	
K092341 Low Profile Abutment Biomet 3i, Inc.	BIOMET <i>3i</i> Low Profile Abutments are intended for use as accessories to endosseous dental implants to support a prosthetic device in a partially or completely edentulous patient. A dental abutment is intended for use to support single and multiple tooth prosthesis, in the mandible or maxilla. The prosthesis is screw retained to the abutment.
K150669 Neoss TiBase and CoCr Abutments Neoss Ltd.	<p>Neoss TiBase: Neoss Abutments are designed to be connected to the Neoss Implants and intended for use as an aid in prosthetic rehabilitation. The Neoss TiBase is compatible with the Sirona Dental System inCoris ZI Meso L. All digitally designed copings and/or crowns for use with the Neoss TiBase Abutments are to be designed and milled using the Sirona Dental CAD/CAM System.</p> <p>Neoss CoCr Abutments: Neoss abutments are designed to be connected to the Neoss Implants and intended for use as an aid in prosthetic rehabilitation.</p>
K120414 OsseoSpeed™ Plus Astra Tech AB	<p>Implants: The Astra Tech Dental Implants are intended for both one- and two-stage surgical procedures in the following situations and with the following clinical protocols:</p> <ul style="list-style-type: none"> <li>• replacing single and multiple missing teeth in the mandible and maxilla,</li> <li>• immediate placement in extraction sites and in situations with a partially or completely healed alveolar ridge,</li> <li>• especially indicated for use in soft bone applications where implants with other implant surface treatments may be less effective,</li> <li>• immediate loading in all indications, except in single tooth situations on implants shorter than 8 mm, or in soft bone (type IV) where implant stability may be difficult to obtain and immediate loading may not be appropriate.</li> </ul> <p>The intended use for OsseoSpeed™ Plus 3.0S is limited to replacement of maxillary lateral incisors and mandibular incisors.</p> <p>Abutments: Astra Tech Implant System Plus abutments are intended to be used in conjunction with Astra Tech Implant System Plus in fully edentulous or partially edentulous maxillary and/or mandibular arches to provide support for crowns, bridges or overdentures.</p> <p>Atlantis Abutments: The Atlantis™ Abutment is intended for use with an endosseous implant to support a prosthetic device in a partially or completely edentulous; patient. It is intended for use to support single and multiple tooth prostheses, in the mandible or maxilla. The prosthesis can be cemented, screw retained or friction fit to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.</p> <p>The Atlantis™ Crown Abutment in Zirconia is intended for use with an endosseous; implant to function as a substructure that also serves as the final restoration, in partially or completely edentulous; patients. The prosthesis is screw retained. The abutment screw is intended to secure the crown abutment to the endosseous implant.</p>
K160784 CAM Titanium Blanks Altatec GmbH	<p>CAM Titanium Blanks are intended for the fabrication of abutments and healing caps on CAMLOG® SCREW LINE and CAMLOG® ROOT-LINE implants in the maxilla and mandible. Digitally designed abutments fabricated with CAD/CAM techniques for use with CAM Titanium Blanks are intended to be sent to a CAMLOG validated milling center for manufacture.</p> <p>CAM Titanium Blanks are intended for the fabrication of abutments and healing caps on CONELOG® SCREW-LINE implants in the maxilla and mandible. Digitally designed abutments fabricated with CAD/CAM techniques for use with CAM Titanium Blanks are intended to be sent to a CAMLOG validated milling center for manufacture.</p> <p>CAM Titanium Blanks are intended for the fabrication of abutments and healing caps/gingiva former on iSy® implants in the maxilla and mandible. Digitally designed abutments fabricated with CAD/CAM techniques for use with CAM Titanium Blanks are intended to be sent to a CAMLOG validated milling center for manufacture.</p>

Comparison of Technological Characteristics

Comparison	Subject Device	Primary Predicate Device	Reference Devices			
	DESS Dental Smart Solutions  Terrats Medical SL	K170588  DESS Dental Smart Solutions  Terrats Medical SL	K092341  Low Profile Abutment  Biomet 3i, Inc.	K150669  Neoss TiBase and CoCr Abutment  Neoss Ltd.	K120414  OsseoSpeed Plus  Astra Tech AB	K160784  CAM Titanium Blank  Altatec GmbH
<b>Design</b>						
Abutment Design	CAD/CAM Blank CAD/CAM TiBase Castable Abutment Aurum Abutment	CAD/CAM Blank CAD/CAM TiBase Abutment	Titanium Abutment	CAD/CAM TiBase Castable Abutment	Castable Abutment	CAD/CAM Blank
Prosthesis Attachment	Cement-retained Screw-retained	Cement-retained Screw-retained	Screw-retained	Cement-retained Screw-retained	Cement-retained Screw-retained	Cement-retained
Restoration	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit	Single-unit, Multi-unit
Abutment/Implant Platform Diameter* (mm)	3.4 – 6.5	3.4 – 6.5	3.4 – 5.0	4.1	3.0 – 5.4	3.3 – 6.0
Abutment Angle	0°	0°	0°- 30°	0°- 20°	0°- 30°	0°- 30°
Abutment/Implant Interface	Internal and External	Internal and External	Internal and External	Internal	Internal	Internal
<b>Material</b>						
Abutment	Titanium Alloy CoCr, Zirconia	Titanium Alloy Zirconia	Titanium Alloy	Titanium Alloy, CoCr, Zirconia	Titanium Alloy Zirconia, Gold, PEEK	Titanium Alloy
Screw	Titanium Alloy	Titanium Alloy	Titanium Alloy	Titanium Alloy	Titanium Alloy	Titanium Alloy

\*See Summary of Abutment Designs (p.3) for individual OEM platform diameters

DESS Dental Smart Solutions subject device abutments are substantially equivalent in intended use, design, function, material, size, and Indications for Use to DESS Dental Smart Solutions abutments (K170588). The only change in the Indications for Use statement is the name of the abutment. The Aurum Abutment is substantially equivalent to the TiBase abutment in K170588. Design parameters for minimum wall thickness, maximum gingival height and no angulation correction remain the same as the primary predicate K170588. The minimum post height for single-unit restorations is now 4.0 mm. Titanium and CoCr Pre-milled Blanks are substantially equivalent to the titanium Pre-milled Blanks in K170588. Design parameters for minimum wall thickness, minimum post height for single-unit restorations, maximum gingival height and no angulation correction remain the same as the primary predicate K170588. The CoCr Abutment is substantially equivalent in design and function to OsseoSpeed Plus abutments (K120414) and in material and function to Neoss CoCr abutments (K150669). Design parameters for minimum wall thickness, minimum post height for single-unit restorations, maximum gingival height and no angulation correction remain the same as the primary predicate K170588 and reference device K160784 and are substantially equivalent in intended use, design, function, material, size and final design parameters. The Aurum Abutment has a substantially equivalent titanium post height to the titanium post height of the Low Profile Abutment in K092341. Both are used for single-unit and multi-unit restorations. When used for a single-unit restoration the Aurum Abutment and the Low Profile Abutment are to be used with an additional castable component to create a minimum post height of 4 mm). All are intended for use with endosseous dental implants in the maxilla and mandible to provide prosthetic support. Digital files for subject device DESS Dental Smart Solutions abutments and DESS Dental Smart Solutions abutments cleared in K170588 are to be sent to a validated milling center for manufacture. Implant/abutment interface compatibility for the subject device is substantially equivalent to DESS Dental Smart Solutions abutments (K170588) and to the compatible implant system predicates listed above.

## CONCLUSION

The subject device and the predicate devices have the same intended use, have similar technological characteristics, and are made of similar materials. The subject device and predicate devices encompass the same range of physical dimensions, including diameter and design of the abutments. The subject and predicate devices are packaged in similar materials and are to be sterilized using similar methods.

The data included in this submission demonstrate substantial equivalence to the predicate devices listed above.