

Illumina DRAGEN Bio-IT Platform

Getting Started Guide



This document and its contents are proprietary to Illumina, Inc. and its affiliates ("Illumina"), and are intended solely for the contractual use of its customer in connection with the use of the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose and/or otherwise communicated, disclosed, or reproduced in any way whatsoever without the prior written consent of Illumina. Illumina does not convey any license under its patent, trademark, copyright, or common-law rights nor similar rights of any third parties by this document.

The instructions in this document must be strictly and explicitly followed by qualified and properly trained personnel in order to ensure the proper and safe use of the product(s) described herein. All of the contents of this document must be fully read and understood prior to using such product(s).

FAILURE TO COMPLETELY READ AND EXPLICITLY FOLLOW ALL OF THE INSTRUCTIONS CONTAINED HEREIN MAY RESULT IN DAMAGE TO THE PRODUCT(S), INJURY TO PERSONS, INCLUDING TO USERS OR OTHERS, AND DAMAGE TO OTHER PROPERTY, AND WILL VOID ANY WARRANTY APPLICABLE TO THE PRODUCT(S).

ILLUMINA DOES NOT ASSUME ANY LIABILITY ARISING OUT OF THE IMPROPER USE OF THE PRODUCT(S) DESCRIBED HEREIN (INCLUDING PARTS THEREOF OR SOFTWARE).

© 2019 Illumina, Inc. All rights reserved.

All trademarks are the property of Illumina, Inc. or their respective owners. For specific trademark information, see www.illumina.com/company/legal.html.

Revision History

Document	Date	Description of Change
Document # 1000000076675 v01	April 2019	• Updated pedigree file option. • Added two sections on de novo calling.
Document # 1000000076675 v00	January 2019	Initial release.

Table of Contents

Revision History	iii
Introduction	1
Hardware and Software Installation	1
Run the Self Test	1
Run Your Own Test	2
Generate a Reference	2
Generate an HG19 Reference	2
Load a Reference	3
Process FASTQ Data	3
End-to-End Aligning and Variant Calling Examples	4
Alignment Only Examples	6
RNA Map and Align Only Examples	8
Epigenome Map and Align Examples	9
Variant Calling Only Examples	10
Somatic Examples	11
gVCF and Joint Calling Examples	12
Coverage Metrics Report Example	13
CNV Examples	13
Structural Variant Calling Examples	15
BCL to FASTQ Conversion with Minimal Settings	16
S3 and HTTP Streaming Input Examples	16
Troubleshooting	17
Technical Assistance	18

Introduction

This Getting Started Guide helps you to start processing data as quickly as possible. Make sure that the Illumina® DRAGEN™ Bio-IT Platform server is turned on and that you are logged in.

Hardware and Software Installation

If you are already running the latest version of the DRAGEN software and hardware, you can skip ahead to *Run the Self Test* on page 1.

You can query the current version of software and hardware with the following command:

```
dragen_info -b
```

You can query only the software version with the following command:

```
dragen --version
```

To install a new version of software or hardware, first download the package from the [DRAGEN Bio-IT Platform support pages](#) on the Illumina website to your DRAGEN server. The preferred installation method is to use the self-extracting .run file, as follows:

```
sudo sh dragen-3.2.8.el7.x86_64.run
```

During installation, if you are prompted to switch to a new hardware version, enter 'y'. It is important that the hardware upgrade process is not interrupted. When it is complete, you must halt and power cycle the server. A reboot command does not update the hardware version. You must issue a halt command and power the server off and on.

DRAGEN periodically checks for license renewal by communicating with the license server at lus.edicogenome.com. For servers that are behind a firewall, a proxy can be configured by the network administrator in `/etc/environment`. For example:

```
http_proxy="http://proxy.customer.com:80/"
https_proxy="https://proxy.customer.com:80/"
ftp_proxy="http://proxy.customer.com:80/"
rsync_proxy="http://proxy.customer.com:80/"
no_
proxy="localhost,127.0.0.1,localaddress,.localdomain.com,.customer.com"
```

Run the Self Test

To run the self test, run the command the following command:

```
/opt/edico/self_test/self_test.sh
```

This command performs a thorough test of the hardware and takes about 25 minutes. When complete, it should output the following message:

```
SELF TEST RESULT : PASS
```

If a failure occurs, please contact Illumina Technical Support. You can ignore any tests that indicate NON MANDATORY TEST SKIPPED.

Run Your Own Test

In addition to the self-test, you can test with your own data as follows:

- Generate a reference—See [Generate a Reference](#) on page 2.
- Load a reference—[Load a Reference](#) on page 3.
- Process your data—See [Process FASTQ Data](#) on page 3.

Generate a Reference

If you do not have a reference, you can generate one by using the `dragen --build-hash-table` command and passing in the location of the reference FASTA file. You can specify a set of parameters when building your hash table (see the *DRAGEN BioIT Platform User Guide* (1000000070494)).

For testing purposes, you can run the example shell script or the one of commands shown in the examples in this guide. For these examples, the FASTA file is located in `/staging/human/reference/hg19/hg19.fa`.

Run the example shell script as follows:

```
/opt/edico/examples/build_hash_table.sh
```

Run the `dragen` command as follows:

```
mkdir -p /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
cd /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
dragen --build-hash-table true --ht-reference
/staging/human/reference/hg19/hg19.fa
--output-dir /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_
149
```

The `dragen --build-hash-table` command is multithreaded and defaults to eight threads. This command takes about 15 minutes to run. You can use the `--ht-num-threads` option with a value up to 32 (depending on the number of threads your server supports) to reduce the run time.

The hash table directory name lists key default option values that were used during the hash table build. Illumina recommends following this best practice when you generate your own hash tables and change the directory name accordingly.

Generate an HG19 Reference

If you do not have a FASTA reference, you can get the hg19 FASTA files from UCSC and concatenate them into a single `hg19.fa` file as follows:

```
mkdir /staging/hg19fa
cd /staging/hg19fa
wget hgdownload.cse.ucsc.edu/goldenPath/hg19/bigZips/chromFa.tar.gz
tar -zxvf chromFa.tar.gz
cat chr*.fa > hg19.fa
```

Generate the DRAGEN hashtable reference using the following commands.

```
mkdir /staging/hg19/
/opt/edico/bin/dragen --ht-reference /staging/hg19fa/hg19.fa --output-
directory /staging/hg19/ --build-hash-table true
```

Load a Reference

After the binary reference is loaded into memory on the DRAGEN board, it can be used for processing any number of input data sets. You do not need to reload the reference unless you restart the system, or need to use a different reference hash table.

The reference is loaded automatically the first time you process data with it. You can manually load the reference genome onto the board by using the following shell script or command. The reference directory in this example is /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149.

```
/opt/edico/examples/load_reference.sh
```

OR

```
dragen -l -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_
149
```

The dragen command should take less than one minute, and should display the following message:

```
DRAGEN finished normally
```

If a manual or automatic system reset occurs, then the next time you try to process data, the reference you specify on the command line is automatically reloaded. The reference is also reloaded if you reboot the system.

Process FASTQ Data

After you have loaded your reference, you can process input FASTQ data. Choose the example that best matches your data sets. These commands can take up to approximately 30 minutes to run on a 24 core server with SSD drives on a 30x coverage whole human genome when running end-to-end (FASTQ input to VCF output). The speed scales with input size, so a 60x coverage genome would take twice as long. Exome data takes a fraction of the time. A successful result is indicated by the following message (an application exit code of 0 when run from a script):

```
DRAGEN finished normally
```

This message is followed by a block of metrics such as read count and performance. If there is a problem with the command line options, an error is displayed, followed by help usage. You may need to scroll up to see the error.

The DRAGEN log can be redirected to a file, to keep the record for future reference.

To get help on dragen command line options, run the following command:

```
dragen -h
```

The example commands shown in this document are formatted for visual display and include line feed characters. To avoid copy-paste errors, each example command is contained in an individual shell script in /opt/edico/examples/. These examples have the following requirements:

- ▶ All commands accept either FASTQ or gzipped FASTQ (fastq.gz). DRAGEN automatically determines the file type.
- ▶ All commands include the -f option, which forces the output file to be overwritten if it exists.

- All commands assume that your DRAGEN reference (hash table) directory is /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149, and your FASTA reference file is /staging/human/reference/hg19/hg19.fa. Replace those with the correct references, if needed.
- All command examples assume that the example data package is in /staging/examples (in particular, the .fastq and fastq.gz files are expected to be in /staging/examples/reads).

End-to-End Aligning and Variant Calling Examples

Paired-End BAM Input, VCF Output

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-b /staging/human/unsorted_SRA056922_30x_e10_50M.bam
--enable-map-align true
--enable-map-align-output true
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true
```

Paired-End FASTQ Input, VCF Output (Default)

```
/opt/edico/examples/paired_fastq_in_vcf_out.sh
```

OR

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
```

This command should take approximately 6 minutes to run on a 24-core server.

This example shows the minimum options that must be specified to perform an end-to-end run. By default, duplicate-marking is not performed and no BAM output is produced.

Paired-End Fastq Input, Sorted and Duplicate-Marked, VCF Output

```
/opt/edico/examples/paired_fastq_in_dupmark_vcf_out.sh
```

OR

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
```

```

2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true

```

Paired-End FASTQ Input, Sorted BAM and VCF Output

```
/opt/edico/examples/paired_fastq_in_dupmark_bam_and_vcf_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true
--enable-map-align-output true

```

Paired-End FASTQ Input, Sorted SAM and VCF Output

```
/opt/edico/examples/paired_fastq_in_dupmark_sam_and_vcf_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true
--enable-map-align-output true
--output-format SAM

```

Paired-End FASTQ Input, Sorted CRAM and VCF Output

```
/opt/edico/examples/paired_fastq_in_dupmark_cram_and_vcf_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_

```

```

1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true
--enable-map-align-output true
--output-format CRAM

```

Paired-End FASTQ Input, Sorted BAM and VCF Output, plus Repeat Genotyping VCF

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--enable-duplicate-marking true
--enable-map-align-output true
--repeat-genotype-enable true
--repeat-genotype-specs /opt/edico/repeat-specs/hg19
--repeat-genotype-sex female
--repeat-genotype-ref-fasta
/staging/human/reference/h19/hg19.fa

```

Alignment Only Examples

All the variations for performing alignment shown in these examples can be used in the end-to-end case as well.

Map/Align Single-Ended FASTQ Input, Sorted BAM output (Default)

```
/opt/edico/examples/single_fastq_in_bam_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-1 /staging/examples/reads/SRA056922_30x_rand1_100K.fastq
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_rand1_100K

```

Map/Align Single-ended FASTQ input, Sorted, Duplicate-Marked BAM Output

```
/opt/edico/examples/single_fastq_in_dupmark_bam_out.sh
```

OR

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_rand1_100K.fastq
  --output-directory /staging/examples/
  --output-file-prefix SRA056922_30x_rand1_100K_dup_marked
  --enable-duplicate-marking true
```

Map/Align Paired-End FASTQ Input, Sorted BAM Output (Default)

```
/opt/edico/examples/paired_fastq_in_bam_out.sh
```

OR

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix SRA056922_30x_e10_50M
```

Map/Align Paired-End FASTQ Input, Sorted CRAM Output

```
/opt/edico/examples/paired_fastq_in_cram_out.sh
```

OR

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix SRA056922_30x_e10_50M
  --output-format CRAM
```

Map/Align Paired-End FASTQ Input, Sorted Uncompressed BAM Output

```
/opt/edico/examples/paired_fastq_in_uncompressed_bam_out.sh
```

OR

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix uncompressed_SRA
```

Map/Align Paired-End FASTQ Input, Sorted SAM Output

```
/opt/edico/examples/paired_fastq_in_sam_out.sh
```

OR

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
--output-format SAM
```

Map/Align Paired -End FASTQ Input, UN-Sorted BAM output

```
/opt/edico/examples/paired_fastq_in_unsorted_bam_out.sh
```

OR

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-l /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--output-directory /staging/examples/
--output-file-prefix unsorted_SRA056922_30x_e10_50M
--enable-sort false
```

Map/Align Interleaved Paired-Ended FASTQ Input, BAM Output

```
/opt/edico/examples/interleaved_fastq_in_bam_out.sh
```

OR

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-l /staging/examples/reads/SRA056922_PE_30x_rand1_10K_
interleaved.fastq
--interleaved
--output-directory /staging/examples/
--output-file-prefix SRA056922_PE_30x_rand1_10K_interleaved
```

RNA Map and Align Only Examples

Any of the Map/Align Only examples can be used for RNA. The only difference in the command is to set the `--enable-rna` option to true. DRAGEN automatically picks up the RNA-specific hash tables and uses the RNA spliced aligner in its processing.

RNA Map/Align Paired-Ended FASTQ Input, BAM Output

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix SRA056922_30x_e10_50M
  --enable-rna true
```

Epigenome Map and Align Examples

Prior to performing an epigenome (methylation) map and align run with bisulfite sequencing data, you must first create methylation-specific reference hash tables, as follows:

```
mkdir -p /staging/human/reference/hg19_epigenome

dragen --build-hash-table true --ht-reference
  /staging/human/reference/hg19/hg19.fa --ht-max-seed-freq 64 --
  ht-seed-len 27 --ht-methylated true --output-directory
  /staging/human/reference/hg19_epigenome
```

The above dragen command produces two hash table directories under `/staging/human/reference/hg19_epigenome`: `GA_converted` and `CT_converted`. The `CT_converted` hash table is produced by converting each C base to T in the reference sequences. Similarly, the `GA_converted` hash table is produced from the G->A base-converted reference sequences. The base-converted references have less complexity, and to compensate, the hash table seed length argument (`--ht-seed-len`) is typically increased to 27 for mammalian genomes (default seed length is 21).

Epigenome Map/Align, Directional-protocol, Single-Ended FASTQ Input, BAM Output

The directional (Lister) protocol produces reads from two of the four possible bisulfite sequencing strands. Therefore, when the `--methylation-protocol=directional` option is used, DRAGEN aligns each read or read pair twice with different constraints corresponding to the two possible strands. The following DRAGEN command produces two separate BAM files:

```
mkdir -p /staging/epigenome/directional

dragen --output-directory /staging/epigenome/directional
  --methylation-protocol=directional
  -r /staging/human/reference/hg19_epigenome
  --fastq-file1=/staging/epigenome/reads/sample_1_R1.fastq.gz
  --RGID=rg1
  --RGSM=samp1
  --RGPL=illumina
  --output-file-prefix=sample_1
```

Epigenome Map/Align, Nondirectional-protocol, Paired-Ended FASTQ Input, BAM Output

The nondirectional protocol produces reads from all four possible bisulfite sequencing strands. Therefore, when the `--methylation-protocol=non-directional` argument is used, DRAGEN aligns each read four times and produces four separate BAM files.

```
mkdir -p /staging/epigenome/non-directional
dragen --output-directory /staging/epigenome/non-directional
      --methylation-protocol=non-directional
      -r /staging/human/reference/hg19_epigenome
      --fastq-file1=/staging/epigenome/reads/sample_10_R1.fastq.gz
      --fastq-file2=/staging/epigenome/reads/sample_10_R2.fastq.gz
      --RGID=rg10
      --RGSM=samp10
      --RGPL=illumina
      --output-file-prefix=sample_10
```

Variant Calling Only Examples

The variant calling only examples show how you can pass an existing aligned BAM or CRAM file directly to the DRAGEN Variant Caller. By default, the BAM/CRAM file is passed through the sorting stage prior to variant calling.

To duplicate mark your BAM file before running the DRAGEN Variant Caller, you need to use a separate tool. The DRAGEN Duplicate Marker depends on information provided by the Mapper/Aligner that does not exist in BAM files. To take advantage of the DRAGEN Duplicate Marker, use DRAGEN in end-to-end mode.

The BAM/CRAM files that are used as input to these example commands are not included in the example data set. They are generated by example commands in [Alignment Only Examples on page 6](#).

Unsorted BAM Input, VCF Output (Default)

```
/opt/edico/examples/unsorted_bam_in_vcf_out.sh
```

OR

```
dragen -f
      -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
      -b /staging/human/unsorted_SRA056922_30x_e10_50M.bam
      --enable-variant-caller true
      --vc-sample-name RGSM
      --output-directory /staging/examples/
      --output-file-prefix unsorted_output_SRA056922_30x_e10_50M
      --enable-map-align false
```

Sorted BAM Input, VCF Output

```
/opt/edico/examples/sorted_bam_in_vcf_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-b /staging/human/SRA056922_30x_e10_50M.bam
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix sorted_output_SRA056922_30x_e10_50M
--enable-map-align false
--enable-sort false

```

Sorted CRAM Input, VCF Output

```
/opt/edico/examples/sorted_cram_in_vcf_out.sh
```

OR

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix sorted_output_SRA056922_30x_e10_50M
--enable-sort false
--enable-map-align false
--cram-input /staging/human/SRA056922_30x_e10_50M.cram

```

Somatic Examples

If you are using the Tumor-Normal BAM input option, and your BAM read groups have a shared RGID, DRAGEN cannot determine which read group the reads belong to. Ideally, you should have different RGIDs for each read group, but you can work around the problem by setting the *--prepend-filename-to-rgid* to true.

Paired-End FASTQ Input

```

dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
--tumor-fastq1 /staging/examples/reads/SRA056922_30x_
shuffle16k_e10_50M_1.fastq.gz
--tumor-fastq2 /staging/examples/reads/SRA056922_30x_
shuffle16k_e10_50M_2.fastq.gz
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M

```

Sorted BAM Input

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
--tumor-bam-input /staging/human/SRA056922_30x_e10_50M.bam
--enable-variant-caller true
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix sorted_output_SRA056922_30x_e10_50M
--enable-map-align false
--prepend-filename-to-rgid true
```

gVCF and Joint Calling Examples

Paired-End FASTQ Input, gVCF Output

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
-1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
1.fastq.gz
-2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
2.fastq.gz
--enable-variant-caller true
--vc-emit-ref-confidence GVCF
--vc-sample-name RGSM
--output-directory /staging/examples/
--output-file-prefix SRA056922_30x_e10_50M
```

Joint Calling with gVCF Input

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
--enable-joint-genotyping true
--output-directory /staging/examples/
--output-file-prefix Joint_SRA056922_30x_e10_50M
--variant /staging/examples/SRA056922_30x_e10_50M.gvcf
```

Joint Calling with Three gVCF Files and a Pedigree File Input, VCF Output

```
dragen -f
-r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
--enable-joint-genotyping true
--output-directory /staging/examples/
--output-file-prefix Joint_SRA056922_30x_e10_50M
--variant /staging/examples/mother.gvcf.gz
--variant /staging/examples/father.gvcf.gz
--variant /staging/examples/child.gvcf.gz
--pedigree-file <pedigree file>
```

Table 1 Joint-calling modes, with associated input files and command line options.

VCF to generate	Population joint-called multisample gVCF	Family joint-called multisample gVCF	Population joint-called multisample VCF	Family joint-called multisample VCF
Input file	Multisample combined gVCF file	Multisample combined gVCF file	Multi-sample combined gVCF file or X individual gVCF files	Multi-sample combined gVCF file or X individual gVCF files
Use pedigree file	No	Yes	No	Yes
Command line option	--enable-joint-genotyping true --enable-multi-sample-gvcf=TRUE	--enable-joint-genotyping true --enable-multi-sample-gvcf=TRUE --pedigree-file file.ped	--enable-joint-genotyping true	--enable-joint-genotyping true --pedigree-file file.ped

Coverage Metrics Report Example

By default, DRAGEN reports read coverage over the whole genome, and if available also for the target bed. You can specify up to three additional regions over which coverages are reported. In the following example, DRAGEN generates two coverage reports, **cov_report** and **full_res**, for the additional coverage region 1.

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix SRA056922_30x_e10_50M
  --qc-coverage-region-1 /staging/examples/reads/vc_smoke.callable.bed
  --qc-coverage-report-1 cov_report full_res
```

The **full_res** report corresponds to the Bedtools coverage option, and includes a per base resolution read depth. The **cov_report** includes read depth summaries per region including mean, median, min and max depths.

CNV Examples

These examples show how to use DRAGEN CNV to process already mapped and aligned BAM files using the two modes of normalization supported by the DRAGEN CNV pipeline: Self Normalization and Panel of Normals.

Self Normalization requires that the DRAGEN hashtable be generated with the *enable-cnv=true* option. It is recommended that you always generate a CNV compatible hashtable if you frequently run CNV.

The *enable-map-align option* is set to true by default in the configuration file. Set it to false if you do not need to map and align the input BAM.

Running with Self Normalization

Self normalization is the preferred method for single sample WGS processing. The BAM goes through the entire CNV pipeline with a single command line.

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -b /staging/human/SRA056922_30x_e10_50M.bam
  --output-directory /staging/examples/
  --output-file-prefix dragen_cnv1
  --enable-map-align false
  --enable-cnv true
  --cnv-enable-self-normalization true
```

Running with a Panel of Normals

The Panel of Normals approach requires pregenerating the target.counts file for each sample to be used, and then executing one final command to perform the normalization and copy number variant calling.

To calculate target counts with BAM Input, this example command extracts the signals, including read counts, from the alignments of the BAM file. It generates a *.target.counts file to be used in the normalization step. Target counts should be calculated for each input BAM file, including the case sample under analysis and the normals samples.

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -b /staging/human/SRA056922_30x_e10_50M.bam
  --output-directory /staging/examples/
  --output-file-prefix dragen_cnv1
  --enable-map-align false
  --enable-cnv true
```

The following command performs the normalization and generates the CNV calls. The normals samples should be listed in a text file (**normal.txt** in this example) that provides the path to the *.target.counts files of the normal samples. The case sample *.target.counts file is specified with the **--cnv-input** option. In this example, gcbias correction of the input is disabled.

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  --output-directory /staging/examples/
  --output-file-prefix dragen_cnv2
  --enable-cnv true
  --cnv-input /staging/examples/dragen_cnv1.target.counts
  --cnv-normals-list normal.txt
  --cnv-enable-gcbias-correction false
```

FASTQ Processing

This example runs the DRAGEN CNV caller in Self Normalization mode directly from FASTQ samples. It first maps and aligns the FASTQ and continues directly to CNV calling. This step can be combined with variant calling.

```

dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    1.fastq.gz
  -2 /staging/examples/reads/SRA056922_30x_shuffle16k_e10_50M_
    2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix dragen_cnv
  --enable-map-align true
  --enable-cnv true
  --cnv-enable-self-normalization true

```

Running De Novo CNV Calling

De novo calling requires previously generated normalized signal files (*.tn.tsv) from the single sample analysis. If a pedigree file is supplied, then a de novo state and a de novo quality score will be annotated for the proband sample's records.

```

dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  --cnv-input father.tn.tsv
  --cnv-input mother.tn.tsv
  --cnv-input child.tn.tsv
  --output-directory /staging/examples/
  --output-file-prefix trio_cnv
  --pedigree-file trio.ped
  --enable-cnv true

```

Structural Variant Calling Examples

Structural Variant calling can run in the following modes:

- Standalone—Runs from mapped BAM/CRAM input files. Requires the `--enable-map-align=false` and `--enable-sv=true` options.
- Integrated—Automatically runs on the output of the DRAGEN mapper/aligner. Requires the `--enable-map-align=true`, `--enable-sv=true`, `--enable-map-align-output=true`, and `--output-format=bam` options.

Structural Variant calling can be enabled along with any other caller as well.

Integrated Execution Example

```

dragen -f
  --ref-dir=<HASH_TABLE>
  --sv-reference=<FASTA>
  --enable-map-align=true
  --enable-map-align-output=true
  --output-format=BAM
  --enable-sv=true
  --output-directory=<OUT_DIR>
  --output-file-prefix=<PREFIX>
  -1 <FASTQ1> -2 <FASTQ2>

```

Standalone Joint Diploid Calling Example

```
dragen -f
  --ref-dir=<HASH_TABLE>
  --sv-reference=<FASTA>
  --enable-map-align=false
  --enable-sv=true
  --bam-input=<BAM1>
  --bam-input=<BAM2>
  --bam-input=<BAM3>
  --output-directory=<OUT_DIR>
  --output-file-prefix=<PREFIX>
```

Standalone De Novo Quality Scoring Example

```
dragen -f
  --variant=<TRIO_VCF_FILE>
  --pedigree-file=<PED_FILE>
  --enable-map-align=false
  --sv-scoring=true
  --output-directory=<OUT_DIR>
  --output-file-prefix=<PREFIX>
```

BCL to FASTQ Conversion with Minimal Settings

This example shows how to use DRAGEN to process Illumina BCL format files.

The BCL directory used in the example is not included in the example data package. Please replace `/mnt/san/131022_hsxten008_0123_FC543` with your own BCL directory.

DRAGEN might produce multiple files per sample with the following names: `<SampleName>_001.fastq`, `<SampleName>_002.fastq`, etc. There is no need to concatenate these files before performing Map-Align using DRAGEN. When you specify the first file in the series, DRAGEN read all files in the series as if they were concatenated into one file.

```
dragen --bcl-conversion-only=true
  --bcl-input-directory /mnt/san/131022_hsxten008_0123_FC543
  -- output-directory /staging/examples/
```

S3 and HTTP Streaming Input Examples

DRAGEN can process input files directly from an S3 bucket, or using HTTP presigned URLs, which is known as input streaming. The input files do not need to be downloaded to a local disk prior to being processed. Instead, the files are streamed over the network directly into the DRAGEN processor.

Streaming is supported for compressed FASTQ (*.fastq.gz) files. A future version of DRAGEN will also support streaming from BAM (*.bam) files. Input streaming can be used in all the configurations that use single-end FASTQs, paired-end FASTQs, and FASTQ lists. The following examples show some of the ways to use input streaming.

Streaming FASTQ Input using S3

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 s3://s3-bucket-name/path/to/object_1.fastq.gz
  -2 s3://s3-bucket-name/path/to/object_2.fastq.gz
  --output-directory /staging/examples/
  --output-file-prefix streaming
```

Streaming FASTQ Input using HTTP

```
dragen -f
  -r /staging/human/reference/hg19/hg19.fa.k_21.f_16.m_149
  -1 https://bucket-name.amazonaws.com/path/to/object_
1.fastq.gz?querystring
  -2 https://bucket-name.amazonaws.com/path/to/object_
2.fastq.gz$querystring
  --output-directory /staging/examples/
  --output-file-prefix streaming
```

You need to have permission to access the remote files. If you have access to the file, then DRAGEN is capable of streaming the remote file. The S3 object requires AWS authentication and credentials. The authentication should already be set up on the instance you are running, for example, via IAM policies. An HTTP URL most likely has a query string attached to it, which provides the authentication credentials or necessary tokens to grant permission. The security method may be present in other parts of the URL, for example:

```
https://stagingdl.dnanex.us/security/string/sample_1.fastq.gz
```

Troubleshooting

The DRAGEN software automatically resets the board if any problems are encountered. In the rare case that reset does not occur automatically, you can use the following command to reset:

```
dragen_reset
```

If resetting does not resolve the issue, contact Illumina Technical Support. Run the following command to collect diagnostic and configuration information for Technical Support.

```
sudo sosreport --batch
```

This tool takes several minutes to run and reports the location where it has saved the diagnostic information in `/tmp`. Please include the report when you submit a support ticket for Illumina Technical Support.

For more information, refer to the *DRAGEN Bio-IT Platform User Guide*(1000000070494) on the Illumina Support Site.

Technical Assistance

For technical assistance, contact Illumina Technical Support.

Website: www.illumina.com
Email: techsupport@illumina.com

Illumina Customer Support Telephone Numbers

Region	Toll Free	Regional
North America	+1.800.809.4566	
Australia	+1.800.775.688	
Austria	+43 800006249	+43 19286540
Belgium	+32 80077160	+32 34002973
China	400.066.5835	
Denmark	+45 80820183	+45 89871156
Finland	+358 800918363	+358 974790110
France	+33 805102193	+33 170770446
Germany	+49 8001014940	+49 8938035677
Hong Kong	800960230	
Ireland	+353 1800936608	+353 016950506
Italy	+39 800985513	+39 236003759
Japan	0800.111.5011	
Netherlands	+31 8000222493	+31 207132960
New Zealand	0800.451.650	
Norway	+47 800 16836	+47 21939693
Singapore	+1.800.579.2745	
Spain	+34 911899417	+34 800300143
Sweden	+46 850619671	+46 200883979
Switzerland	+41 565800000	+41 800200442
Taiwan	00806651752	
United Kingdom	+44 8000126019	+44 2073057197
Other countries	+44.1799.534000	

Safety data sheets (SDSs)—Available on the Illumina website at support.illumina.com/sds.html.

Product documentation—Available for download in PDF from the Illumina website. Go to support.illumina.com, select a product, then select **Documentation & Literature**.



Illumina

5200 Illumina Way

San Diego, California 92122 U.S.A.

+1.800.809.ILMN (4566)

+1.858.202.4566 (outside North America)

techsupport@illumina.com

www.illumina.com

For Research Use Only. Not for use in diagnostic procedures.

FOR IVD PERFORMANCE EVALUATION ONLY.

© 2019 Illumina, Inc. All rights reserved.

illumina®