

Auscultation Simulator with Heart and Lung Sounds

S315.200



Gaumard[®]
Simulators for Health Care Education

This simulator is an interactive educational system developed to assist a certified instructor. It is not a substitute for a comprehensive understanding of the subject matter and not intended for clinical decision making.

User Guide 15.8.2
© Gaumard Scientific Company, 2015
All Rights Reserved
www.Gaumard.com

Care and Cautions

Overall Warnings	5
General	5
Operating Conditions	5
Storage Conditions	5
Cleaning	5

Heart and Lung Sounds

Overview	6
Contents	6
Setup	6

Appendix

Warranty	8
Exclusive One-Year Limited Warranty	8

Care and Cautions

Overall Warnings

Remember that damage caused by misuse is not covered by your warranty. It is critical to understand and comply with the following guidelines:

GENERAL

- The simulator is constructed of material that approximates skin texture. Therefore, in handling the model, use the same gentle techniques as you would in working with a patient.
- Ball point pens, ink and markers permanently stain the skin.
- Do not wrap this or any other Gaumard product in newsprint.
- Indelible marks made with ballpoint pens, ink or marker cannot be removed.
- Do not use alcohol, acetone, Betadine® or any other antiseptic which contains iodine in this or any Gaumard® simulator. These products could damage or stain the skin of the simulator.
- Replacement parts are available from Gaumard Scientific or your Distributor.

OPERATING CONDITIONS

Operating the simulator outside these ranges may affect performance:

- Operating temperature: 50°-95° F (10°-35° C).
- Humidity: 5%-95% (non-condensing).

STORAGE CONDITIONS

- Storage temperature: 32°-85° F (0°-29° C).
- Humidity: 40%-60% (non-condensing).
- Do not stack or store heavy materials on top of the carton.

CLEANING

- The simulator should be cleaned with a cloth dampened with diluted liquid dishwashing soap. If medical adhesives remain on the skin, clean with alcohol wipes.
- Apply talcum powder to the face to reduce the silicone tackiness.
- The simulator is “splash-proof” but not water-proof. Do not submerge or allow fluid to enter the interior of the simulator.

Heart and Lung Sounds

Overview

This simulator features a full-size adult upper torso and head with palpable anatomic landmarks.

Heart and Lung Sounds teaching system is a tool used for auscultation training. The system is composed of a Virtual Stethoscope™ and sound ID tags located beneath the skin in a total of 13 locations, 9 on the front and 4 on the back.



Incorporated in the Virtual Stethoscope™ is:

- Power switch
- Sound menu switch



The skin of the simulator has 13 ID Tags located under the skin where each of the heart and lung sounds are normally heard.

The simulator is shipped with removable blue dots attached to the skin to show the location of the sound ID tags and at any time



You may remove the front and back skins of the torso to locate and note the positions of the tags.

Each ID tag pertains to a different area and will signal the Virtual Stethoscope™ to produce a different sound when in use.



Contents

- 1 Adult upper body simulator
- Overlay torso with ID Tag network located beneath the skin
- VS100 Virtual Stethoscope™ Kit
- Set of External Speakers with AC Plug

Setup

1. The Virtual Stethoscope™ comes with two batteries installed.
2. Locate the small stereo jack on the bell and attach the speakers provided.



3. Plug the speakers into a conventional wall outlet and turn the speakers on.

The speakers provide external feedback of the auscultation sounds. If the speakers are disconnected, only the student will hear the heart and lung sounds through the ear pieces on the Virtual Stethoscope™.

4. Follow the Virtual Stethoscope™ User Guide to use the system.

Appendix

Warranty

EXCLUSIVE ONE-YEAR LIMITED WARRANTY

Gaumard warrants that if the accompanying Gaumard product proves to be defective in material or workmanship within one year from the date on which the product is shipped from Gaumard to the customer, Gaumard will, at Gaumard's option, repair or replace the Gaumard product.

This limited warranty covers all defects in material and workmanship in the Gaumard product, except:

Damage resulting from accident, misuse, abuse, neglect, or unintended use of the Gaumard product;

Damage resulting from failure to properly maintain the Gaumard product in accordance with Gaumard product instructions, including failure to properly clean the Gaumard product; and

Damage resulting from a repair or attempted repair of the Gaumard product by anyone other than Gaumard or a Gaumard representative.

This one-year limited warranty is the sole and exclusive warranty provided by Gaumard for the accompanying Gaumard product, and Gaumard hereby explicitly disclaims the implied warranties of merchantability, satisfactory quality, and fitness for a particular purpose. Except for the limited obligations specifically set forth in this one-year limited warranty, Gaumard will not be liable for any direct, indirect, special, incidental, or consequential damages, whether based on contract, tort, or any other legal theory regardless of whether Gaumard has been advised of the possibilities of such damages. Some jurisdictions do not allow disclaimers of implied warranties or the exclusion or limitation of consequential damages, so the above disclaimers and exclusions may not apply and the first purchaser may have other legal rights.

This limited warranty applies only to the first purchaser of the product and is not transferable. Any subsequent purchasers or users of the product acquire the product "as is" and this limited warranty does not apply.

This limited warranty applies only to the products manufactured and produced by Gaumard. This limited warranty does not apply to any products provided along with the Gaumard product that are manufactured by third parties. For example, third-party products such as computers (desktop, laptop, tablet, or handheld) and monitors (standard or touch-screen) are not covered by this limited warranty. Gaumard does not provide any warranty, express or implied, with respect to any third-party products. Defects in third-party products are covered exclusively by the warranty, if any, provided by the third-party.

Any waiver or amendment of this warranty must be in writing and signed by an officer of Gaumard.

In the event of a perceived defect in material or workmanship of the Gaumard product, the first purchaser must:

Contact Gaumard and request authorization to return the Gaumard product. Do NOT return the Gaumard product to Gaumard without prior authorization.

Upon receiving authorization from Gaumard, send the Gaumard product along with copies of (1) the original bill of sale or receipt and (2) this limited warranty document to Gaumard at 14700 SW 136 Street, Miami, FL, 33196-5691 USA.

If the necessary repairs to the Gaumard product are covered by this limited warranty, then the first purchaser will pay only the incidental expenses associated with the repair, including any shipping, handling, and related costs for sending the product to Gaumard and for sending the product back to the first purchaser. However, if the repairs are not covered by this limited warranty, then the first purchaser will be liable for all repair costs in addition to costs of shipping and handling.

Extended Warranty In addition to the standard one year of coverage, the following support plans are available: Two-Year Extension (covers second and third years)

Call for pricing (USA only)



Gaumard®
Simulators for Health Care Education

Contact Us

E-mail Technical Support: support@gaumard.com

Before contacting Tech Support you must:

1. Have the simulator's Serial Number
2. Be next to the simulator if troubleshooting is needed.

E-mail Sales and Customer Service: sales@gaumard.com

Phone: Toll-free in the USA: (800) 882-6655

Worldwide: 01 (305) 971-3790

Fax: (305) 667-6085

Post: Gaumard Scientific

14700 SW 136 Street

Miami, FL 33196-5691

USA

Office hours: Monday-Friday, 8:30am - 4:30pm EST (GMT-5, -4 Summer Time)

General Information

Gaumard®, ZOE®, Michelle®, Mike®, PEDI®, Susie Simon®, Susie®, Simon® Code Blue®, SIMA Models®, SIMA GYN/AID®, Virtual Instruments®, Codemaker®, Code Blue®, NOELLE®, Simulation Made Easy™, HAL®, CPRLink™, Zack™, RITA™, Chloe™, Seatbelt Susie™, Krash Kids™, Premie™, GIGA®, Omni®, SmartSkin™, Victoria®, UNI® are Trademarks of

© Gaumard Scientific Company, 2015

All rights reserved.

Gaumard®, NOELLE®, HAL® are trademarks of Gaumard Scientific Company, Inc.

Patented; Other Patents Pending

Always dispose of this product and its components in compliance with local laws and regulations.