KBC-144 Reference Lab Request Form R05/24 User's Guide



Reference Laboratory Request Form

3121 Beaumont Centre Circle Lexington, KY 40513

Main: (859)276-2534 & Reference: (859)519-3816 & Fax: (859)514-0128

Please follow these instructions carefully

IMPROPERLY LABELED SPECIMENS WILL NOT BE PROCESSED

- Notify Reference Laboratory in advance of specimen arrival.
- Send 24mL anticoagulated (EDTA) OR 10mL anticoagulated (EDTA) and 20mL fresh dotted blood in sterfie, well stoppered, sealed test tubes. For suspected Warm Autoantibodies, 24-30mL EDTA is preferred.

DO NOT Send specimens collected in gel separator tubes.

- 3. Pack specimens carefully to avoid leakage and/or breakage. Specimens should be packed in a leak-proof container and a box.
- 4. Specimens should be packaged and transported according to Department of Transportation Regulations for clinical specimens.
- 5. If specimens will be in transit for more than 4 hours they should be packed as you would red blood cell components.
- If ABO/Rh history is unknown (i.e. due to an ABO discrepancy), an ABO confirmation tube with a different collection date/time will be required prior to issue of crossmatched products.
- 7. THE INFORMATION ON THIS FORM MUST MATCH THE TUBES.
- 8. MISLABELED SPECIMENS WILL NOT BE TESTED -

Specimens will be used for crossmatch if necessary, and MUST be legibly and clearly labeled with:

- □ The full name of the patient
- ☐ Hospital ID number, SSN or DOB.
- ☐ The DATE (Month, day, and year) and TIME the specimen was collected.
- Phlebotomist's initials, name, or identifier.
- For Translusion Reactions: send dearly labeled pre- and post-translusion specimens along with a sample of the unit(s) translused.
 Additional form required: KBC-106 Translusion Reaction Investigation OR KBC-396 TRAIT Investigation Report
- For Crossmatch problems: send patient specimen along with sample of incompatible units of blood, labeled appropriately.
- For homolytic disease of the newborn (HDN): send appropriately labeled specimens from baby, mother, and father (if available).

PRINT OR TYPE ALL INFORMATION - BOTH SIDES of this form MUST BE COMPLETED.

FOR EMERGENCY RELEASE OF UNCROSSMATCHED BLOOD SEND TUBES WITH KBC-264 AND KBC-144

_		OF HOSPITAL,		RY:	1			ORDERING PHYSICIAN:		2			
ATION	PATIENT'S NAME (Last, First):				3			DOB:		4			
PATIENT'S NAME (Last, First): HOSPITAL ID: DATE/TIME SPECIMEN				5	H/H or Pl Count 6		CLINICAL DIAGNOSIS:		7				
DATE/TIME SPECIMEN COLLECTED:					[28]		PHIEBOTOMIST'S INITIALS/IE			. 9			
PATIENT RACE: 10 S			SSN:	11			PATIENT GENDER:		12				
FORM COMPLETED BY/DATE:					13	13		HOSPITAL PHONE NUMBER:		14			
	COMPLETE ADDITIONAL PATIENT HISTORY ON THE BACK OF THIS FORM												
SAME	PLE COMP	LETION PRI	ORITY (-/):		TING REQUESTED (<-):			MPONENT NEEDED	ECIAL NEEDS (~):				
STAT 15 (Ex. Ordus for IMMEDIANE handring platent is certical)					Antbody Identification Pre-DARA work-up Antbody Titration 16			Leukocyte Reduced Packed Red Blood Cells	0	Irradiation CMV Seron	ogative		
								☐ Crossmatched OR	0	☐ HIA Selected-			
Who Approxis the Call-in-Feet?					IXEN IXIN Investigation			☐ Compatbility Screened		IRRADIATED			
15a					Direct Antiglobulin Test			☐ Antigen negative ☐		1 Washed			
D ASAB					ABO/Rh Typing Resolution			Translusable Plasma	٥	Designated	Donor		
☐ ASAP					Molecular Genotype			Cryoprecipitate Pool		Ofter:			
(Ex. Ordurs to transluse soon, but patient is stable)					HDN Investigation			Platelet, Pherests	Н		18		
ROUTINE C					Ofter:			mponent Quantity (~):17a	IJ	ر ت			
(Ex. patient is stable/outpatient, no translusion orders)								02 03 04 0					

KBC-144 R 05/24 COMPLETION SOP: 15-99-144

To be completed by Transfusion Facility (Front)

- 1. Name of the hospital/laboratory
- 2. Name of the ordering Physician
- 3. Patient's last name, first name (middle initial is optional)
- 4. Patient's date of birth (DOB)
- 5. Patient's hospital number/ID
- 6. Patient's hemoglobin/hematocrit (H/H) or Platelet Count
- 7. Patient's clinical diagnosis
- 8. Date and time specimen collected
- 9. Phlebotomist initials/ID
- 10. Patient's race
- 11. Patient's Social Security Number (SSN)
- 12. Patient's gender Male or Female
- 13. Form completed by and date of completion
- 14. Phone number of hospital/laboratory
- 15. Check the Sample Completion Priority for the status of patient sample STAT, ASAP, Routine. Examples are provided
 - a. Document the name of the person approving the call-in-fee
- 16. Check the Testing Requested. If other is marked document the type of work-up on the blank line.
- 17. Indicate the Component Needed by marking the quantity desired, if applicable.
 - a. Check the component quantity needed. If other is marked document the quantity on the blank line.
- 18. Check Special Needs, if applicable. If other is marked document the special need on the blank line.

KBC-144 Reference Lab Request Form R05/24 User's Guide

Patient Name	(Last, First)	19a	$\overline{}$	Par	tient ID:	19b				
			PATIEN	THISTORY	20a Sort	to KBC Before? U YES U NO				
A complete trans										
			and/or patient's			ory				
			no recently transfused o			blank if not applicable.				
	ATIENT/FAMILY		$\overline{\frown}$			OR TRANSFUSIONS?				
	YES		20b		☐ YES ☐ NO	☐ Unknown 20c				
			ANSFUSION HISTO	RY AT CURRENT						
NUMBER OF UNIT		<u> </u>	MOST RECENT	date(s) transfu	JSED: 20e	20f				
	(-	Dd J				ping® ☐ YES ☐ NO(201				
ABO/RH OF UNIT		20g	REACTION(S), I	FANY:	20h					
(If within the last 3 i	inchiency .	,	ENT HAS RECEIVED T	neAtMeast. /						
IISI ANY OTHER I	ACITITES IN W	HICH IHE PAII	EMI LINZ KECTIAED I	KEAIMENI:	20i					
	ANTIBOD	ES IDENTIFIED:	20i							
ANTIBODIES:		Leadles								
	DATE OF L	AST WORKUP:	20k							
DRUGS/MEDICAT	IONS Ploate let	All moderations	the potion) has receive	within the last A av	onths. If nocossary attac	h extra sheets. Please include if patient				
			or similar IgG immune							
See attached list		7	<u></u>							
	21	J								
			F	NANCIES:						
	patient ever bee		22a 🗅			Currently Prognant/Delivering				
	of Homolytic Disc	case of the Newt	× 2m3 □.	Yes 22b [⊒No #ofPre	gnancias: 22c				
Pregnancy Dates:	22d									
Has the patter	Yes 1	h Immuno Globul	in∂ Whon∂	22f						
	J 165 J 1	40	SEROLOG	ICAL RESULTS						
Please give us YOU	R serological findi	ngs. Indicate stre	ingth of reactions seen		ILY (1+, 2+, etc.) or mic	croscopically. 23a				
ABO/Rh (D) TESTII	VG: Istho	re difficulty in AE	3O grouping≎ □	Yes □No Is	s there difficulty in Rh t	yping? ☐ Yes ☐ No [
		Forward Type	1 4 - 6	l el el	Reverse Type	Direct Antiglobolin Tus				
Anit-A Anit	B Anti-A,B	Anti-D	Rh Ctrl Anti-D		A1 Cells B Cells	(DAT) 23c				
\rightarrow	_	+ +	PVIO	7 (7410)		314,349036				
23b		1 1				23d				
	ANTIBOD	Y SCREEN								
(Please attach	copies al your an		n/papal rasultsl	CROSSMATCHES:						
Screen Cel	23e	Lot #				D peG Other:				
Manufacturer:	_ Solid			Method:	□ Gel □ Sofid	23: D LISS D				
Method: 🗆 Gel	□ Solid (_r23g _Jus	Other:		Phase	231				
		37°C A	GEL/Solid	IS	37	AHG GEL/Solid				
	12 3	A	Phase Phase	15	3/	AHG Phase				
SCI		23h			23i					
SCII					[23]					
SCIII			\neg	NUMBER O	F UNITS CROSSMA	TCHED: 23k				
Auto	-		$\neg \neg$	NUMBER O	F UNITS INCOMPA	TIBLE: 231				
COMMENTS OR	QUESTIONS:					KBC USE ONLY				
					Date/Time Received	! 26				
	24					25 26				
Į.	24				Sample Acceptable	PD You DaNo Fedron				
_						27a No ledi27b				
	KI	BC-144 R 05	/24	COMPLET	TION SOP: 15-99-	144				
			7	COMME						

- 19. Document patient information:
 - a. Document the Patient Name (Last, First)
 - b. Document the Patient ID
- 20. For patient history document as follows:
 - a. Check (✓) YES or NO to indicate if the patient has been sent to KBC before.
 - **NOTE:** If the patient has been sent to KBC before, and there have been no changes to the history (transfusions, drugs or pregnancies), skip to number 22.
 - b. Check (✓) YES or NO to indicate if the patient or their family was consulted for history.
 - c. Check (✓) YES, NO or UNKNOWN to indicate if the patient has received prior transfusions.
 - d. Document the number of units transfused
 - e. Document the most recent date(s) of transfusion
 - If the patient has been transfused in the last 3 months, check
 (✓) YES or NO to indicate if a sample can be sent for molecular genotyping.
 - g. Document the ABO/Rh of any units transfused in the last 3 months.
 - h. Document any known transfusion reactions.
 - i. Document any facilities that have treated the patient and may have given transfusions.
 - i. List any known antibodies identified
 - k. Document date of last workup
- 21. List current medications, if applicable.
 - a. If a list is attached, check (✓) see attached.
- 22. Complete the patient's pregnancy information, if applicable.
 - a. Check (✓)YES or NO to indicate if the patient has ever been pregnant. Check (✓) the box if they are currently pregnant or delivering.
 - b. Check (✓) YES or NO to indicate if there is a history of Hemolytic Disease of the Newborn (HDN).
 - c. Document the number of pregnancies.
 - d. Document the date(s) of pregnancy(ies).
 - e. Check (✓) YES or NO to indicate if the patient has received Rh Immune Globulin (RhIG).
 - f. Document the date of RhIG administration, if applicable.

KBC-144 Reference Lab Request Form R05/24 User's Guide

Patient Name	(Last, Fi	rst)	19a	\neg		P	atient ID:			.9b
				\rightarrow	PATIENT	HISTORY		20a Sorti		U YES U NO
A complete tran										
Please consult w									ry	
Patients with unknow Please include infor									dank if not ann	kontido
			CONSULTED	: (OR TRANSFUS	
		YES D N		20b					☐ Unknown	20c
		12.5		ANSFUSIO	N HISTOR	Y AT CURREN			- Olikilowii	
NUMBER OF UNI	TS:				T RECENT D	ATE(S) TRANSE	FUSED:	20e		
		20d	J	(If we	hin the last 3	months, approv	ral for mala	ular geontys	oing₽ □ YES	□ NO 20f
ABO/RH OF UNI			1	REAC	TION(S), IF	ANY:	20	\neg		•
(If within the last 3		20g	')					"		
LIST ANY OTHER I	ACILITIES	IMMHIC	H THE PAT	IENT HAS R	ECEIVED TR	EATMENT:	20i)		
						Ĺ		J		
	LANI	TRODIES I	DENTIFIED	- (2	0i					
	174	ibobics i	DEI (III EE		9					
ANTIBODIES:	DAT	e of l a st	WORKUP	: 20	0k					
										lease include if patient
is receiving Darzale See attached its		umab (DAR	A, Win Kh	o or similar lg	G Immune g	iobulin medicalio	on. Altoch sh	eets If neede	d.	
→ See affected its	' [21								
				_	PDFGN	IANCIES:				
Hasti	e patient e	ver been pr	oanant?	22a	□ Ye		□ No		Currently Proc	nant/Deli yering
Is there a history					→ <u> </u>		□ No		gnancles:	22c
Pregnancy Dates:	22d	+-				$\overline{}$				
Has the patte	nt over reco	atvod Rh Im	muno (Flob	lin?	Whon?	226				
,	☐ Yes	□ No		22e		22f				
		92 B				CAL RESULTS		8		
Please give us YOU										23a
ABO/Rh (D) TESTI	NG:			BO groupin	g≎ □Y	′as □No			ping? 🗆 Yes	□N ₀ □
			orward Type		Anit-D	Rh Ctrl		se Type	- ⊢ '	Direct Antig lobelin Te. (DAT) 23c
Anit-A Anit	-B A	Anti-A,B Anti-D		D Rh Chrl (AHG)		(AHG)	A1 Cells B Cells			□ Poly □NgG □ G3
\rightarrow	_	-				, , , , , , , , , , , , , , , , , , ,			\vdash	
23b										23d
	ANI	IBODY S	CREEN			_				
(Please attach				en/pa <u>nel re</u> :	subs/			CROSSN	AATCHES:	
Screen Cell		23e	Lot f		3f				□ PeG	Other:
Manufacturer:	Manufacturer:					Method:	☐ Gel	Solid		uss 🗆
Method: □ Ge	l □ Sc	olid (□Pa	23g ⊔µs	s 🗆 Orli	er:			Phase 2	231] -	
		T	一	1.6	EL/Solid					GEL/Solid
	IS	37*(c A	HG	Phase	IS		37	AHG	Phase
SCI		(. h							
SCII		1 23	n 	-			- Il	23j		
SCIII		+		-+		NUMBER (OF UNITS	ROSSMAT	CHED:	23k
SP SHITTER		+-	$-\!$	+				NCOMPAT		\rightarrow
Auto		15				INOMBER (J. OHIISI		(BC USE ONI	231
Auto OtuA	OHESTIC						Date/Te	ne Received:		
Auto COMMENTS OR	QUESTIO	M2:								
	QUESTIO	No:							ſ	26
	QUESTION 24	No:					^		25	
		No:					^	Acceptable 2	25 □ Y∞ □ γN	ledi27b
		No.:					^		25	
)	144 R 0				Sample /		25 1 Yes - 1 Na 27a	ledi27b

- 23. Serological Results:
 - a. Check (✓) YES or NO to indicate if there is difficulty in ABO group or Rh type testing.
 - b. Document the serologic results for ABO/Rh
 - c. Check (✓) Poly, IgG or C3 to indicate which reagent was used for the DAT testing, if applicable.
 - d. Document the results of the DAT, if applicable.
 - e. Document the manufacturer of the screen cells being used.
 - f. Document the lot number of the screen cells being used.
 - g. Check (✓) Gel, Solid Phase, PeG, LISS or Other to document the method of testing for antibody screen. If Other, document the method of testing.

<u>NOTE</u>: A copy of the antigram for the cells tested should be included.

- h. Document the serologic results for the antibody screen.
- i. Check (✓) Gel, Solid Phase, PeG, LISS or Other to document the method of testing for crossmatches, if applicable.
- j. Document the serologic results for crossmatched units, if applicable.
- k. Document the number of units crossmatched, if applicable.
- 1. Document the number of units incompatible, if applicable.
- 24. Document any additional Comments or questions.