MEASURING THE BOILING POINT

	STUDENT BOOK	Chapter 1, page 23
1	TOOLBOX	Page 37

Goal

Apply a technique for measuring boiling point.

Materials

- 100-mL beaker
- container of liquid
- hot plate
- ring stand
- thermometer

- thermometer clamp or universal clamp and perforated cork stopper
- beaker tongs
- · ceramic plate

Procedure



- 1. Pour 40 mL of liquid into the beaker.
- 2. Place the beaker on the hot plate.
- **3.** Insert the thermometer into the test tube and clamp it so the bulb is submerged completely and not touching the test tube.
- 4. Heat at medium temperature.
- 5. Record the temperature at which bubbles form.
- 6. Turn off the hot plate.
- 7. Remove the beaker from the hot plate using the beaker tongs and place it on the ceramic plate.
- 8. Clean up and put away materials.

Results

Record your results in the table below. Give the table a title.

Title:

Liquid	Experimental boiling point (°C)	Reference boiling point (°C)

Reflecting on the lab technique

Compare the boiling point obtained during the experiment with that of tables of characteristic properties. Are they similar? If not, explain why.