TECHNICAL INFORMATION

**2020/10/22 (ver. 8)

*2019/11/13 (ver.7)

Approval number etc.: 23000BZI00008000

Category: Program 01 Disease diagnostic program

Highly Controlled Medical Device

Germline mutation analysis program (For determine the risk of acquiring disease) (71058013)

Germline mutation analysis program (Eligibility screening for anticancer drug) (61777003)

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*【Shape, Structure and Principle, etc.】

Description for diagnostic tests using this system (device)

The BRACAnalysis 診断システム diagnostic system provides a Health Care Provider (HCP) with the classification of a patient's clinically significant mutations. Whole blood is collected from the patient in an EDTA blood collection tube. The patient specimen is packaged together with a test request form and held for pickup by a designated distributor in Japan. The specimen is collected and shipped by the distributor to the Myriad laboratory in Salt Lake City, UT, USA. Genomic DNA is extracted, and PCR and Sanger Sequencing is used to detect sequence variants in the *BRCA1* and *BRCA2* genes. Multiplex PCR analysis is used for the detection of large rearrangements in the *BRCA1* and *BRCA2* genes.

The BRACAnalysis 診断システム diagnostic program checks each variant against a database of previously classified variants. Variants not previously classified in the variant database, are classified by applying a defined set of objective criteria from multiple sources. The process and criteria (evidence) used to classify variants have been established in accordance with ACMG guidelines/standards and are consistently applied to all variants identified by Myriad. *BRCA1* and *BRCA2* variants are classified into one of the following five categories:

- Deleterious
- Suspected Deleterious
- Variant of Uncertain Significance
- Favor Polymorphism
- Polymorphism

When the classification for each variant has been loaded in the variant database, a notification is sent to the HCP. Through secure access the HCP reviews a list of variants identified during the processing of the patient's sample. Through the web-based application, the HCP will review the Interim Variant List, and then request generation and presentation of a report on the clinical significance of the variants based on the results of the analysis performed by this system. Variants classified as Deleterious or Suspected Deleterious are clinically significant and are reported as a "Positive" result. Variants classified as a Variant of Uncertain Significance (VUS) do not have sufficient data to determine whether or not they are clinically significant and unless a clinically significant variant is also identified, the test will be reported as "Negative". Variants classified as Polymorphism or Favor Polymorphism are not clinically significant and will not be included on the report. Patients carrying only Polymorphisms or Favor Polymorphisms will be reported as "Negative."

When a portion of the BRACAnalysis 診断システム test is unable to be analyzed, a new sample will be requested. The patient's Test Report will be provided to the HCP indicating the cancellation of the test. A Customer Service representative will contact the HCP to request an additional blood draw and explain the reason for the cancellation and the need for a new sample.

****【Intended Use or Indication】**

This device is used as an aid for detecting germline *BRCA1* or *BRCA2* gene mutations in genomic DNA extracted from whole blood and for determining the eligibility of patients with breast, ovarian, pancreatic, or prostate cancer for olaparib treatment.

This device is also used to identify individuals with a high risk of *BRCA*-Related Breast and/or Ovarian Cancer (HBOC) Syndrome and may be used as an aid for informing medical management decisions.

****【Usage Method, etc.】**

After testing has been completed, a secure email portal will be used for viewing the list of variants identified through testing, as well as the patient's Test Report. An email notification containing a link to a secure email portal will be sent to the Health Care Provider (HCP) indicated on the test request form. The HCP will be required to follow the following steps:

- 1) Click on the link to the secure email portal provided in the email
- 2) Enter the required login information (Registration for access to the secure email portal will be required when the secure email portal is accessed for the first time)
- 3) Within the secure email portal will be a secure email from Myriad containing a list of all variants identified from testing, the Interim Variant List
- 4) Review the list of variants and request the patient's Test Report by selecting, "View Final Test Report"
- 5) View and print the patient's Test Report

The Interim Variant List contains the list of all variants identified through the testing of the patient's sample, except for those variants for which definitive evidence indicates a lack of clinical significance. The patient's Test Report will indicate the "Positive" or "Negative" depending on the clinical significance of the variants identified during analyses performed using this system. The patient's Test Report of "Positive" will include the classification for each variant that is clinically significant (Deleterious or Suspected Deleterious germline *BRCA1* or *BRCA2* mutation). The classification information for these variants is provided on the report. The patient's Test Report of "Negative" will only display variants classified as a Variant of Uncertain Significance (VUS), if identified. Variants classified as Polymorphism or Favor Polymorphism are not clinically significant and will not be included on the report. Patients carrying only Polymorphisms or Favor Polymorphisms will be reported as "Negative." When multiple variants are provided on the report, the overall interpretation of the test results is based on the variant with the highest clinical significance.

Patients evaluated with the BRACAnalysis 診断システム test that are determined to carry a deleterious or suspected deleterious germline *BRCA1* or *BRCA2* mutation can be considered for treatment with olaparib under the supervision of a HCP. In addition, the test may identify patients at a high risk for HBOC Syndrome associated with *BRCA1* and *BRCA2* deleterious or suspected deleterious mutations. It is recommended that these test results be communicated to the patient in a setting that includes appropriate counselling due to potential implications for the patient and their family members. Variants of Uncertain Significance (VUS) are variants whose clinical significance have not yet been determined. Patients determined to carry a VUS (variant whose clinical significance has not yet been determined)

cannot be considered for treatment with olaparib without a Deleterious or Suspected Deleterious germline *BRCA1* or *BRCA2* mutation.

****【Precautions for Use】**

<Important Precautions>

- It should be confirmed that medical sites satisfy the conditions specified by the related societies before requesting a genetic test.
- Refer to the guideline related to the genetic testing for the unaffected patients issued by academic society, when the BRACAnalysis 診断システム test will be conducted with based on patient family history.
- Patients who have undergone a previous allogeneic bone marrow transplant should not be tested with BRACAnalysis 診断システム.
- Patients under consideration for testing who have been diagnosed with a hematologic malignancy, such as leukemia, could generate a positive (deleterious or suspected deleterious) result that does not reflect mutation status in the germline and therefore should not receive testing with BRACAnalysis 診断システム.
- The classification and interpretation of all variants identified reflects the current state of scientific understanding at the time the result report is issued. In some instances, the classification and interpretation of variants may change as scientific information becomes available.
- This system is designed to detect variants and genomic rearrangements (deletions or duplications) involving the promoter and coding exons of *BRCA1* and *BRCA2*. As such, use in the following cases may lead to false negative results.
- When there are variants that may impact patient care outside of the assessed regions.
- When unequal allele amplification is resulted from rare polymorphisms under primer sites.
- When mutations result in some types of errors in RNA transcript processing
- When mutations are insertions that do not result in duplications

Also, in some instances, this system cannot clearly differentiate a duplication from a triplication.

- Whole blood collected in an EDTA blood collection tube is to be held at ambient temperature for pickup which should occur within 5 days of the blood draw.

【Other Precautions】

- When this product is used, the latest Technical Information for olaparib in Japan should be referenced.

【Clinical Significance – Olaparib Treatment in Breast Cancer – D0819C00003 (OLYMPIAD)】

The clinical utility of the BRACAnalysis 診断システム test has been evaluated using breast cancer patient samples collected from the OlympiAD clinical study. AstraZeneca partnered with Myriad to develop BRACAnalysis CDx as a companion diagnostic for olaparib. Myriad submitted BRACAnalysis CDx, for use as a companion diagnostic program, in Japan under the brand name BRACAnalysis 診断システム.

(1) Overview of OlympiAD Study

Clinical results are based on a Phase III, randomised, controlled, open-label, multicentre study D0819C00003 (OlympiAD), in which 302 *gBRCAm* patients who had received no more than 2 prior lines of chemotherapy for metastatic breast cancer and were HER2-negative were randomised 2:1 to receive either olaparib (300 mg bd, tablet formulation) or physician’s choice of chemotherapy (capecitabine, eribulin, or vinorelbine). The primary endpoint was Progression Free Survival (PFS) as determined by blinded independent central review (BICR).

Overall, of the 302 patients randomized onto OlympiAD, 299 were tested, either prospectively or retrospectively, using the BRACAnalysis CDx test. This represented 99% of the full analysis set in OlympiAD. 297 patients were confirmed to carry a deleterious or suspected deleterious *gBRCA* mutation.

(2) Concordance between local *BRCA* results used for enrolment of patients and the BRACAnalysis CDx test

Within the OlympiAD clinical study, patients were required to have documented evidence of a deleterious or suspected deleterious mutation in *BRCA1* or *BRCA2*. Evidence of a qualifying *BRCA* mutation could be from either *BRCA* mutation result generated by local testing or from prospective testing performed by Myriad. As the to-be-marketed BRACAnalysis CDx test, was not available at the time of the start of screening in OlympiAD, prospective testing at Myriad began by using Myriad’s *BRCAAnalysis* test, a *BRCA* test provided as a CLIA (Clinical Laboratory Improvement Amendments) regulated test. The *BRCAAnalysis* test can be considered equivalent to the to-be-marketed BRACAnalysis CDx test as the two tests use the same assay methodology for the primary analysis and the same variant classification procedure.

All appropriate patients that were not prospectively tested at Myriad, with sufficient material available, were subsequently retrospectively tested using BRACAnalysis CDx.

Concordance analyses were performed between the BRACAnalysis CDx results, and local germline *BRCAm* results in OlympiAD. Comparisons were also made between the BRACAnalysis CDx test and the *BRCAAnalysis* test.

Positive percent agreement and Negative percent agreement and overall percent agreement for patient eligibility (*gBRCAm* and non *gBRCAm*) are calculated and presented in the table below.

Percent Agreement between BRACAnalysis CDx and Tests used for Enrolment onto OlympiAD

Comparison	Positive percent agreement (95%CI)	Negative percent agreement (95%CI)	Overall percent agreement (95%CI)
BRACAnalysis CDx vs <i>BRCAAnalysis</i>	228/229 99.6% (97.6, 99.3)	229/229 100% (98.4, 100)	457/458 99.8% (98.8, 99.9)
BRACAnalysis CDx vs local testing	226/227 99.6% (97.6, 100.0)	409/411 99.5% ¹ (98.3, 99.9)	635/638 99.5% (98.6, 99.9)

1: Negative percent agreement with local laboratories could only be calculated using results from one laboratory, as only one local laboratory performed prospective testing.

Overall, the coincidence rate between the BRACAnalysis CDx diagnostic system and other tests in the OlympiAD test was very high. The overall concordance rate of the BRACAnalysis CDx diagnostic system and the laboratory results of each facility exceeded 99%, and this result was in agreement with the overall concordance rate of the results of the BRACAnalysis CDx Diagnostic System and BRACAnalysis®.

Discordant Results

BRACAnalysis and BRACAnalysis CDx

One patient had a discordant result between the BRACAnalysis CDx and the BRACAnalysis test which would have impacted eligibility for OlympiAD. This patient reported *BRCA1*:IVS19+2insT using both tests. At the time of the BRACAnalysis test this variant was classified as a deleterious mutation. However this classification had been revised at the time of the BRACAnalysis CDx to a variant of uncertain significance based on revised ACMG guidelines (Richards 2015).

Local Testing and BRACAnalysis CDx

Two patients were ineligible for enrolment based on local test results that were later determined to be eligible based on the to-be-marketed BRACAnalysis CDx test. One patient had no mutation detected using the local test but had an Alu insertion detected by the BRACAnalysis CDx test, classified as a deleterious mutation. The second patient was reported to carry *BRCA1*: M18R (172T>G) by both the local test and BRACAnalysis CDx, the local laboratory classified this as a variant of uncertain significance whereas Myriad classified the variant as a suspected deleterious mutation.

Additionally, two patients were enrolled based on local test result that were not confirmed as *gBRCAm* by Myriad. One patient reported by a local test to carry *BRCA1*:5385insC, was determined to have no mutation detected by the BRACAnalysis test and therefore is not included in the comparison table above. The second patient was confirmed by both the local test and BRACAnalysis CDx to have the *BRCA1*:IVS9-2A>C variant, classified as suspected deleterious by the local test but classified as a variant of uncertain significance by Myriad.

There were 2 patients who entered the OlympiAD clinical study based on a local *BRCA* test result but no specimen was available for confirmation by the to-be-marketed BRACAnalysis CDx test.

(3) Comparison of efficacy between OlympiAD full analysis set and the population who have been determined to be *gBRCAm* by the BRACAnalysis CDx

The clinical outcome data for the Full Analysis Set and the subset of patients determined to be *gBRCAm* using the BRACAnalysis CDx is shown in the following table.

Clinical Outcome Data of Clinical Study D0819C00003 (OlympiAD)

	Full Analysis Set (FAS)		BRACAnalysis CDx <i>gBRCAm</i> confirmed	
	Olaparib 300 mg bd ¹	Physicians choice chemotherapy ²	Olaparib 300 mg bd ¹	Physicians choice chemotherapy ²
PFS				
Number of events: total	163:205	71:97	160:202	71:95
number of patients (%)	(79.5)	(73.2)	(79.2)	(74.7)
Median PFS (months)	7.0	4.2	7.4	4.2
HR (95% CI)	0.58 (0.43-0.80)		0.57 (0.41-0.78)	
P-value (2-sided)	P=0.0009		P=0.0005	

1: tablet formulation

2: Physician's choice of chemotherapy consisting of either capecitabine, eribulin or vinorelbine

The clinical outcome data for the 297 patients with confirmed *gBRCAm* status by BRACAnalysis CDx was as follows: Median progression-free survival was significantly longer in the olaparib group than in the standard-therapy group (7.4 months vs. 4.2 months; hazard ratio for disease progression or death, 0.57; 95% confidence interval, 0.41 to 0.78; $P < 0.001$). Taken together, the results in the subset of *gBRCAm* patients tested with the BRACAnalysis CDx were very similar to those observed in the 302 patients in the overall OlympiAD study, which supports the effectiveness of the device.

【Clinical Significance – Olaparib Treatment in Ovarian Cancer – D0818C00001 (SOLO1)】

The clinical utility of the BRACAnalysis 診断システム test has been evaluated using ovarian cancer patient samples collected from the SOLO1 clinical trial.

1. Overview of SOLO1 Study

Clinical results are based on a Phase III, randomised, double blind, placebo controlled, multicentre study D0818C00001 (SOLO1), to assess the efficacy of olaparib maintenance monotherapy in newly diagnosed advanced ovarian cancer patients (including patients with primary peritoneal and/or fallopian tube cancer) with *BRCA* mutations (documented mutation in *BRCA1* or *BRCA2*) that were predicted to be deleterious or suspected deleterious (known or predicted to be detrimental/lead to loss of function) who had responded following first-line platinum based chemotherapy.

Overall, of the 391 patients randomized onto SOLO1, 386 were tested, either prospectively or retrospectively, using the Myriad BRACAnalysis or BRACAnalysis CDx tests. 5 additional patients were enrolled based on results from a laboratory in China. Due to sample export restrictions, these samples could not be retrospectively tested by Myriad.

2. Concordance between local *BRCA* results used for enrolment of patients and the Myriad *gBRCA* test results in SOLO1

Within the SOLO1 clinical study, patients were required to have documented evidence of a deleterious or suspected deleterious mutation in *BRCA1* or *BRCA2*. Evidence of a qualifying *BRCA* mutation could be from either a *BRCA* mutation result generated by local testing or from prospective testing performed by Myriad. As the BRACAnalysis CDx test (Brand Name: BRACAnalysis 診断システム in Japan) was not available at the time of the start of screening in SOLO1, prospective testing at Myriad began by using Myriad's BRACAnalysis test, a *BRCA* test provided as a CLIA (Clinical Laboratory Improvement Amendments) regulated test. The BRACAnalysis test can be considered equivalent to the BRACAnalysis 診断システム test as the two tests use the same assay methodology for the primary analysis and the same variant classification procedure.

Concordance of treatment eligibility reported by the local test and the Myriad *gBRCA* tests was calculated for 208 patients where results were available. Of these, 205 were determined as carrying an eligible mutation by the Myriad *gBRCA* tests. Positive percent agreement is calculated as 98.6% (95% CI 95.6-99.5).

The only local laboratory to perform prospective testing (n=3) was in China. As sample export to Myriad's laboratory was not granted, the negative percent agreement and overall percent agreement with local laboratories could not be calculated.

Discordant Results between local testing and the Myriad *gBRCA* tests

In total, there were 3 patients where the local test results and the Myriad *gBRCA* test results were discordant with respect to treatment eligibility:

For 2 patients the local *BRCA* test was reported to have been generated from tumour testing. The mutations identified in these two samples were not reported by the Myriad *BRCA* test, a germline mutation test. These mutations were subsequently confirmed using a tumour based *BRCA* mutation

test and are therefore considered to be somatic deleterious mutations, thus meeting the *BRCAm* eligibility criteria for SOLO1.

The third patient was classified by a local laboratory as having a deleterious mutation. When tested by the Myriad *gBRCA* test and classified using Myriad's variant classification criteria, the case was classified as VUS. The variant was confirmed to be classified as deleterious according to the classification criteria for the local site, and therefore, the patient was randomised in line with the protocol.

3. Comparison of efficacy between SOLO1 full analysis set and the population determined to be *gBRCAm* by Myriad testing

The effectiveness of the BRACAnalysis 診断システム test was based on a subset of 383 ovarian cancer patients confirmed to carry a deleterious or suspected deleterious germline *BRCA1* and/or *BRCA2* mutations following prospective or retrospective testing performed with either the Myriad BRACAnalysis test or the BRACAnalysis CDx test. The clinical outcome data for the Full Analysis Set and the patients determined to be *gBRCAm* by Myriad is shown in the following table.

Clinical Outcome Data of Clinical Study D0818C00001 (SOLO1)

	Full Analysis Set (FAS)		Myriad <i>gBRCAm</i> confirmed	
	Olaparib 300 mg bd ^a	Placebo	Olaparib 300 mg bd ^a	Placebo
PFS				
Number of events: total number of patients (%)	102:260 (39)	96:131 (73)	99:253 (39)	95:130 (73)
Median PFS (months)	Not reached	13.8	Not reached	13.8
HR (95% CI)	0.30 (0.23-0.41)		0.30 (0.22-0.40)	
P-value (2-sided)	<0.0001		<0.0001	

a- tablet formulation

The clinical outcome data for the 383 patients with a confirmed germline *BRCA1/2* mutation was as follows: a 70% reduction in the risk of disease progression or death at any point in time for olaparib vs placebo treated patients (HR 0.30; 95% CI 0.22-0.40; p<0.0001). After a median follow-up of 41 months, median PFS was not reached on the olaparib arm vs 13.8 months for placebo. Taken together, these results are very similar to those observed in the 391 patients in the SOLO1 study, which supports the effectiveness of the BRACAnalysis 診断システム device.

**【Clinical Significance – Olaparib Treatment in Pancreatic Cancer – D081FC00001 (POLO)】

The clinical utility of the BRACAnalysis 診断システム test has been evaluated using pancreatic cancer patient samples collected from the POLO clinical trial.

1. Overview of POLO Study

The olaparib clinical study D081FC00001 (POLO) was a Phase III, randomized, double-blind, placebo-controlled, multicenter trial to assess the efficacy of Lynparza maintenance treatment in patients with metastatic adenocarcinoma of the pancreas who have a deleterious or suspected deleterious germline *BRCA* mutation (*gBRCAm*) and whose disease had not progressed after at least 16 weeks of first-line platinum-based chemotherapy.

Overall, of the 154 patients randomized in the global POLO cohort, 150 were tested with either the Myriad BRACAnalysis[®] test (n=9) or the Myriad BRACAnalysis CDx[®] test (n=141) and all 150 patients

were confirmed to carry a deleterious or suspected deleterious germline *BRCA* mutation. The PMA cohort represented 97.4% of the full analysis set in the global POLO cohort.

2. Concordance between local *BRCA* results and the Myriad *gBRCA* test results in POLO

A total of 51 patients screened as part of the POLO study joined the study as they had a prior known local *gBRCAm* status. Out of those 51 patients, 44 had both a local *BRCA* result and a Myriad *gBRCA* test. Concordance of treatment eligibility reported by the local test and the Myriad *gBRCA* tests was calculated for 44 patients enrolled locally where Myriad results were available. Of these, 44 (100%) were determined as carrying an eligible mutation by the Myriad *gBRCA* tests.

3. Comparison of efficacy between POLO full analysis set and the population determined to be *gBRCAm* by Myriad testing

The effectiveness analysis for BRACAnalysis 診断システム test was based on a subset of 150/154 metastatic pancreatic adenocarcinoma patients who were confirmed with a deleterious or suspected deleterious germline *BRCA1/2* mutation by either the Myriad BRACAnalysis® test or the Myriad BRACAnalysis CDx® test. Four patients were not confirmed to have a germline *BRCA* mutation by the Myriad test as no sample was submitted for testing. The clinical outcome data for the Full Analysis Set and the patients determined to be *gBRCAm* by Myriad is shown in the following table.

Clinical Outcome Data of Clinical Study D081FC00001 (POLO)

	Full Analysis Set (FAS)		Myriad <i>gBRCAm</i> confirmed	
	Olaparib 300 mg bd	Placebo	Olaparib 300 mg bd	Placebo
PFS				
Number of events: total number of patients (%)	60:92 (65)	44:62 (71)	59:89 (66)	44:61 (72)
Median PFS (months)	7.4	3.8	7.4	3.8
HR (95% CI)	0.53 (0.35-0.81)		0.55 (0.36-0.84)	
P-value (2-sided)	p=0.0038		p=0.0060	

bd: twice daily

The PFS data for the 150 patients in the confirmed Myriad *gBRCAm* subset was as follows: HR of 0.55 (95% CI 0.36-0.84; p=0.0060; the median was 7.4 months for olaparib vs 3.8 months for placebo). These results are almost identical to those observed in the 154 patients in the Full Analysis Set (FAS) from the POLO study, which supports the effectiveness of the BRACAnalysis 診断システム device.

**【Clinical Significance – Olaparib Treatment in Prostate Cancer – D081DC00007 (PROfound)】

The clinical utility of the BRACAnalysis 診断システム test has been evaluated using metastatic castration-resistant prostate cancer patient samples collected from the PROfound clinical trial.

1. Overview of PROfound Study

The olaparib clinical study D081DC00007 (PROfound) was a pivotal Phase III, randomized, open-label, controlled, multicenter study, in which patients with metastatic castration-resistant prostate cancer (mCRPC) and a HRR gene mutation who had progressed after prior treatment with an NHA were randomized 2:1 to receive either olaparib (300 mg twice daily, tablet formulation) or investigators choice of NHA (either abiraterone acetate or enzalutamide). Patients with mutations in *BRCA1*, *BRCA2* or *ATM* were randomized in Cohort A (irrespective of co-occurring mutations in one of the 12 other HRR genes), whereas patients with mutations among 12 other genes involved in the HRR pathway (*BARD1*, *BRIP1*, *CDK12*, *CHEK1*, *CHEK2*, *FANCL*, *PALB2*, *PPP2R2A*, *RAD51B*, *RAD51C*, *RAD51D* or *RAD54L*) were randomized in Cohort B.

Of the 387 randomized patients, 288 patients had a blood sample successfully tested, with a conclusive result, for *gBRCA* mutations using the Myriad BRACAnalysis CDx[®] test. In total, 21.5% (62/288) of the patients were confirmed to carry a loss of function (deleterious or suspected deleterious) mutation in either *gBRCA1* or *gBRCA2*.

2. Comparison of efficacy between PROfound full analysis set and the population determined to be *gBRCAm* by Myriad testing

The effectiveness analysis for BRACAnalysis 診断システム test was based on a subset of 62 mCRPC patients who were confirmed with a deleterious or suspected deleterious germline *BRCA1/2* mutation by the Myriad BRACAnalysis CDx[®] test. The clinical outcome data for the Full Analysis Set and the patients determined to be *gBRCAm* by Myriad is shown in the following table.

Clinical Outcome Data of Clinical Study D081DC00007 (PROfound)

	FAS		Myriad <i>gBRCAm</i> subset	
	Olaparib 300 mg bd ^a	Investigators choice of NHA	Olaparib 300 mg bd ^a	Investigators choice of NHA
rPFS by BICR (maturity 72%)				
Number of events: total	180: 256	99: 131	25: 43	17: 19
number of patients (%)	(70.3)	(75.6)	(58.1)	(89.5)
Median rPFS [months]	5.82	3.52	10.12	1.87
HR (95% CI)	0.49 (0.38, 0.63)		0.08 (0.03, 0.18)	
p-value (2-sided)	<0.0001		<0.0001	

a-twice daily

The rPFS data for the 62 patients in the confirmed Myriad *gBRCAm* subset was as follows: HR of 0.08 (95% CI 0.03-0.18; p<0.0001; the median was 10.1 months for olaparib vs 1.9 months for investigators choice of NHA) in Cohort A+B. These results are consistent with the results of the FAS, demonstrating a statistically significant and clinically meaningful improvement in rPFS for olaparib compared with investigators choice of NHA, which supports the effectiveness of the Myriad BRACAnalysis CDx[®] device (Brand Name: BRACAnalysis 診断システム in Japan).

**【Clinical Significance – HBOC Syndrome】

Inherited germline mutations in the *BRCA1* and *BRCA2* genes are the most important cause of Hereditary Breast and Ovarian Cancer (HBOC). Women who inherit a mutation in either *BRCA1* or *BRCA2* have a 43% to 87% lifetime risk for breast cancer, and a 16% to 63% risk for ovarian cancer.¹⁻⁴ Male mutation carriers have significantly increased risks for breast and prostate cancer^{5,6} and carriers of both genders

have an increased risk for pancreatic cancer and possibly melanoma.⁷⁻¹¹ It is estimated that between 1:300 and 1:500 individuals worldwide carry mutations in *BRCA1* or *BRCA2*, and current data indicates that these estimates are similar for individuals of Asian¹² and specifically Japanese ancestry.¹³ A growing body of published evidence has demonstrated that the identification of mutations in at-risk patients has significant clinical utility in guiding medical management that can lead to better patient outcomes through cancer prevention, early detection and treatment of malignancies. As a result, there are numerous professional society guidelines in place that recommend clinical analysis of *BRCA1* and *BRCA2* for patients meeting specified personal and/or family cancer history criteria.

Prevalence of Mutations in *BRCA1* and *BRCA2*

Approximately 1:300 to 1:500 individuals in the general population carry mutations in *BRCA1* or *BRCA2*.¹⁴⁻¹⁶ Approximately 7% of female breast cancer cases¹⁷ and 11-15% of ovarian cancer cases¹⁸⁻²⁰ are caused by mutations in *BRCA1* or *BRCA2*, which are inherited in an autosomal dominant pattern. Although these estimates are derived primarily from European populations, current data indicates that the numbers are similar in individuals of Asian¹² and Japanese ancestry.¹³

The following clinical features increase the likelihood of a *BRCA* mutation^{2, 21-31} and form the basis for the *BRCA1* and *BRCA2* testing criteria outlined in numerous societal guidelines and technology assessments^{9,52,63,64,68,69}:

- Breast cancer diagnosed before age 50
- Triple negative breast cancer (estrogen receptor negative, progesterone receptor negative, and HER-2 negative)
- Epithelial ovarian, primary peritoneal, or fallopian tube cancer diagnosed at any age
- Two primary breast cancers in an individual
- Male breast cancer
- Multiple cases of breast cancer among biological relatives
- The occurrence of both breast cancer and ovarian/primary peritoneal/fallopian tube cancer among biological relatives
- Pancreatic cancer diagnosed at any age
- Prostate cancer at any age, if intraductal, high grade, or metastatic, or if diagnosed in a man with a family history of other HBOC-related cancers

Mutations in the *BRCA1* and *BRCA2* genes account for approximately 84% of families with a clinical history consistent with hereditary breast and ovarian cancer.² While other genes are known to be linked to an increased breast and/or ovarian cancer risk, *BRCA1* and *BRCA2* are widely considered to be the most clinically significant genes for hereditary risk of these cancers based on their prevalence, associated high cancer risks, and the resulting medical management implications.

Two independent studies have concluded that 1 in every 5 patients with breast cancer has a personal or family history that warrants consideration of *BRCA1* and *BRCA2* testing.³²⁻³³ Among unaffected women seen in a primary care or screening mammography setting, approximately 6% report family histories of breast and/or ovarian cancer that warrant further evaluation and consideration of *BRCA1* and *BRCA2* testing.³⁴⁻³⁵ Finally, the 2005 National Health Interview Survey found that 8.9% of male and female respondents had family histories of breast and ovarian cancer that meet the US Preventive Services Task Force criteria for evaluation and consideration of *BRCA1* and *BRCA2* testing.³⁶

The clinical validity of *BRCA1* and *BRCA2* testing is based on its ability to accurately identify individuals who carry a mutation, and therefore, have significantly increased risks of certain cancers. This is documented in the extensive body of published literature, summarized below. In addition, published data confirms that individuals who test negative for a *BRCA1* or *BRCA2* mutation, known to be

associated with breast and/or ovarian cancer in their family, are NOT at significantly increased risk for these cancers.

***BRCA* Mutations Increase Breast and Ovarian Cancer Risk**

The range of risks for breast and ovarian cancer associated with mutations in the *BRCA1* and *BRCA2* genes has been characterized through numerous studies. Without intervention, the majority of women with inherited mutations in *BRCA1* and *BRCA2* will develop breast and/or ovarian cancer. The lower estimates of risk are derived from analyses of mutations in an unselected general population of individuals, whereas higher estimates of mutation-associated cancer risk are derived from families with a strong history of cancer.

Generally, mutations in *BRCA1* and *BRCA2* are associated with a 43% to 87% risk of breast cancer by age 70¹⁻⁴, however, variability of this risk has been demonstrated among families, most likely due to other modifying factors.^{1,6,37,39,44} Most importantly, hereditary breast cancer occurs, on average, at an earlier age than the nonhereditary (sporadic) form. Women in the general population have only a 2% chance of developing breast cancer before age 50.⁴² However, a woman with a mutation in *BRCA1* or *BRCA2* has a 23% to 51% likelihood of developing breast cancer before reaching 50 years of age.^{2,3,4,40}

The risk of ovarian cancer due to inherited *BRCA1* mutations is 13 to 23% to age 50, and 39% to 63%^{3,4,40} by age 70, compared to the general population risk of <1%. Mutations in *BRCA2* confer a risk of ovarian cancer of 0.4 to 4% to age 50, and 16.5% to 27%²⁻⁴ by age 70.

A *BRCA* Mutation Increases the Risk of a Second Cancer

Women who have developed breast cancer are at greatly increased risk of a second cancer if they carry a mutation in *BRCA1* or *BRCA2* compared to women in the general population.^{1,43,45-49} The risk for a second primary breast cancer within five years of the initial diagnosis in women with *BRCA1* mutations is 20%,⁴⁷ and 12% for women who are positive for *BRCA2* mutations.⁴⁸

One recent study⁵⁰ suggested that the age of onset of the first breast cancer may play a role in the level of second breast cancer risk for a woman who carries a *BRCA* mutation, with a younger age of onset (under 40 years of age) correlating with a higher risk of a second breast cancer occurrence. Women diagnosed with initial breast cancer younger than age 40 had a significantly higher risk of contralateral breast cancer in *BRCA1* families (and a similar trend was seen in *BRCA2* families but not statistically significant). 62.9% of women from *BRCA1* mutation families with their first breast cancer diagnosed before age 40 developed contralateral breast cancer within 25 years compared to 19.6% of women whose first breast cancer was over age 40. Another recent study by Malone et al. found that, compared with noncarriers, *BRCA1* and *BRCA2* mutation carriers had a 4.5-fold and 3.4-fold increased risk of contralateral breast cancer, respectively.⁵¹ The relative risk of developing contralateral breast cancer in *BRCA1* mutation carriers compared to noncarriers increased with younger age at first diagnosis with an 11-fold increased contralateral breast cancer risk among women diagnosed before 35, a 4-fold increased risk for women diagnosed from 35-44, and a 2.6-fold increased risk for women diagnosed 45-54 years. Data regarding *BRCA2* mutation carriers did not show a clear age-related trend.

Data indicate that the 10-year risk for ovarian cancer AFTER a breast cancer in *BRCA1* and *BRCA2* mutation carriers is 12.7% and 6.8%, respectively.⁴⁸ Other studies have shown that the risk for ovarian cancer AFTER a breast cancer diagnosis in *BRCA2* mutation carriers is 16%⁴³ and women with a *BRCA1* or *BRCA2* mutation have a 10-fold increased risk compared to women who do not have a *BRCA1/2* mutation.⁴⁶ Due to the ineffectiveness of screening for ovarian cancer, the National Comprehensive Cancer Network (NCCN) and the Society of Gynecologic Oncology (SGO) have recommended that mutation positive women have a prophylactic bilateral salpingo-oophorectomy after child-bearing.⁵²⁻⁵³

Risks in Men with a *BRCA* Mutation

To date, fewer men than women have undergone genetic testing for *BRCA1* and *BRCA2* mutations. However, men who carry *BRCA1* or *BRCA2* mutations have a significantly increased risk of cancer.⁵⁴ In a study of *BRCA2* families, the cumulative risk of all *BRCA*-associated cancers in men was 32% by age 70, compared with 90% for women.⁴³

For *BRCA2* mutation carriers, the lifetime risk of male breast cancer is approximately 6.8%.⁵ This is compared to a risk of 0.1% (1 in 1000) for male breast cancer in the general population.⁴² The risk of male breast cancer in *BRCA1* families is less characterized, but may be up to 1.2% by age 70.⁵

Mutations in *BRCA1* and *BRCA2* are also associated with an increased risk of prostate cancer. Although precise estimates of the risk vary, mutations in *BRCA2* are consistently found to confer significantly higher risks than those in *BRCA1*. Across various studies the increase in risk has been estimated to be up to 8.6-fold in men with *BRCA2* mutations, and up to 4-fold in men with mutations in *BRCA1*.^{66,70} A recent prospective analysis of men from the United Kingdom and Ireland estimated risk to age 85 as 60% for *BRCA2* and 30% for *BRCA1*.⁶⁷ A large analysis of Japanese prostate cancer patients, compared to unaffected Japanese controls, found a 5.6-fold increased risk in men with *BRCA2* mutations. A 2.3-fold increased risk found in men with *BRCA1* mutations was not statistically significant.⁶⁵

Multiple studies have found that men with aggressive or metastatic prostate cancers are more likely to have *BRCA2* mutations compared to men with less advanced disease.^{72,73} Additional studies have established that *BRCA2* mutations are associated with earlier ages and more advanced stages of diagnosis, as well as higher overall mortality.^{63,66,68,70,71} It is not yet known to what extent these same features are also associated with mutations in *BRCA1*. Based on this evidence, professional societies recommend that *BRCA2* status, and possibly *BRCA1* status, should be factored into prostate cancer screening strategies, as well as treatment decisions for men with a prostate cancer diagnosis.^{62,63,68}

Risks of Other Cancers

Both men and women with mutations in *BRCA1* or *BRCA2* have an increased risk of pancreatic cancer compared to the general population.^{7-9,26,30,31,43,56,57} However, *BRCA2* mutation carriers, with a risk of pancreatic cancer of up to 7% by age 80, appear to be at higher risk than *BRCA1* mutation carriers.^{8,9,76} Across multiple studies of pancreatic cancer patients tested with panels of hereditary cancer genes, *BRCA2* mutations are the most, or among the most, common findings, suggesting that this gene may be the most significant contributor to hereditary risk for the disease.^{74-77,81} In addition to the treatment implications for individuals already diagnosed with pancreatic cancer, *BRCA2* and *BRCA1* mutation carriers with a family history of pancreatic cancer are recommended for consideration of pancreatic cancer screening within clinical trials using modalities such as magnetic resonance cholangiopancreatography (MRCP) and endoscopic ultrasound (EUS).^{9,52,69}

An increased risk of colorectal cancer among *BRCA* mutation carriers, compared to the general population, has been described. Specifically, a two-fold increased risk was observed among *BRCA1* carriers,^{56,58} and a three-fold increased risk was observed among *BRCA2* carriers.¹⁹ However, additional studies have not observed this increased risk, so a true association remains unclear.^{8,43,59} At this time, there are no professional/societal guidelines for colorectal cancer screening in mutation carriers.

An increased risk of melanoma has been described in *BRCA2* mutation carriers with a 2.5-fold increase in risk compared to the general population.^{10,11} An increased risk of melanoma in *BRCA1* mutation carriers has not been observed. NCCN guidelines recommend that mutation carriers be counseled about general melanoma risk management, such as reducing exposure to UV radiation and annual full-body skin exams.

A small number of studies, including two of Japanese patients, have found *BRCA1* and *BRCA2* mutations in individuals with biliary tract cancer.⁷⁸⁻⁸⁰ This data is preliminary and at this time, there are no

professional/societal recommendations related to the management of risk for this cancer in mutation carriers.

Cancer Risks among Individuals Testing Negative for a Known Family Mutation

The power of *BRCA1* and *BRCA2* genetic testing is most fully realized when a mutation can first be identified as the cause of cancer in a family, allowing other family members to determine whether they inherited that mutation. Once a *BRCA* mutation has been identified, relatives who test negative for that mutation do not have an increased risk for breast and ovarian cancer, despite the strong family history of cancer. In a prospective study by Domchek et al., women who were negative for their known familial mutation in *BRCA1/2* had no excess risk of either invasive breast cancer or ovarian cancer; the authors suggest that these women should adhere to population based guidelines for breast cancer screening, unless they have other factors that indicate the need for increased screening, such as a personal history of atypical hyperplasia or a significant history of cancer on the side of the family that is NOT carrying the mutation.⁶⁰ A subsequently published similar study also found no statistically significant overall increase in breast cancer risk among women who tested negative for a known familial mutation.⁶¹ One of the most important benefits of genetic testing for cancer risk, therefore, is that it enables an individual to avoid unnecessary medical interventions if it can be shown that he/she does not carry the mutation in the family.

Referenced Literature

1. Ford D, et al. Risks of cancer in *BRCA1*-mutation carriers. Breast Cancer Linkage Consortium. Lancet. 1994 343:692-5. PMID: 7907678.
2. Ford D, Easton DF, Stratton M, et al. Genetic heterogeneity and penetrance analysis of the *BRCA1* and *BRCA2* genes in breast cancer families. *Am J Hum Genet.* 1998;62:676-89.
3. Chen S, et al. Characterization of *BRCA1* and *BRCA2* mutations in a large United States sample. *J Clin Oncol.* 2006 24:863-71. PMID: 16484695.
4. Mavaddat N, et al. Cancer risks for *BRCA1* and *BRCA2* mutation carriers: results from prospective analysis of EMBRACE. *J Natl Cancer Inst.* 2013 105:812-22. PMID: 23628597.
5. Tai YC, et al. Breast cancer risk among male *BRCA1* and *BRCA2* mutation carriers. *J Natl Cancer Inst.* 2007 99:1811-4. PMID: 18042939.
6. Struwing JP, et al. The risk of cancer associated with specific mutations of *BRCA1* and *BRCA2* among Ashkenazi Jews. *N Engl J Med.* 1997 336:1401-8. PMID: 9145676.
7. Lynch HT, et al. *BRCA1* and pancreatic cancer: pedigree findings and their causal relationships. *Cancer Genet Cytogenet.* 2005 158:119-25. PMID: 15796958.
8. van Asperen CJ, et al. Netherlands Collaborative Group on Hereditary Breast Cancer (HEBON) . Cancer risks in *BRCA2* families: estimates for sites other than breast and ovary. *J Med Genet.* 2005 42:711-9. PMID: 16141007.
9. Goggins M, et al. International Cancer of the Pancreas Screening (CAPS) consortium. Management of patients with increased risk for familial pancreatic cancer: updated recommendations from the International Cancer of the Pancreas Screening (CAPS) Consortium. *Gut.* 2020 Jan;69(1):7-17. doi: 10.1136/gutjnl-2019-319352. Epub 2019 Oct 31. PubMed PMID: 31672839.
10. Gumaste PV, et al. Skin cancer risk in *BRCA1/2* mutation carriers. *Br J Dermatol.* 2015 172:1498-506. Epub 2015 Apr 29. PMID: 25524463.
11. Moran A, et al. Risk of cancer other than breast or ovarian in individuals with *BRCA1* and *BRCA2* mutations. *Fam Cancer.* 2012 11:235-42. PMID: 22187320.
12. Hall MJ, et al. *BRCA1* and *BRCA2* mutations in women of different ethnicities undergoing testing for hereditary breast-ovarian cancer. *Cancer.* 2009 115(10):2222-33. PMID:19241424.
13. Nakamura S et al. Prevalence and differentiation of hereditary breast and ovarian cancers in Japan. *Breast Cancer.* 2015 22(5) 462-8. 24249303.
14. Antoniou AC, Gayther SA, Stratton JF, Ponder BA, Easton DF Risk models for familial ovarian and breast cancer. *Genet Epidemiol.* 2000;18:173-90.

15. Coughlin SS, Khoury MJ, Steinberg KK. BRCA1 and BRCA2 gene mutations and risk of breast cancer: public health perspectives. *Am J Prev Med.* 1999;16(2):91-98.
16. Anglican Breast Cancer Study Group. Prevalence and penetrance of BRCA1 and BRCA2 mutations in a population-based series of breast cancer cases. *Br J Cancer.* 2000;83:1301-8.
17. Claus EB, Schildkraut JM, Thompson WD, Risch NJ, et al. The genetic attributable risk of breast and ovarian cancer. *Cancer.* 1990;77:2318-2324.
18. Zhang S, et al. Frequencies of BRCA1 and BRCA2 mutations among 1,342 unselected patients with invasive ovarian cancer. *Gynecol Oncol.* 2011 May 1;121(2):352-7.
19. Risch HA, McLaughlin JR, Cole DEC, et al. Prevalence and penetrance of germline BRCA1 and BRCA2 mutations in a population series of 649 women with ovarian cancer. *Am J Hum Genet.* 2001;68:700-710.
20. Pal T, et al. BRCA1 and BRCA2 mutations account for a large proportion of ovarian carcinoma cases. *Cancer* 2005;104(12):2807-2816.
21. Frank TS, Deffenbaugh AM, Reid JE et al. Clinical characteristics of individuals with germline mutations in BRCA1 and BRCA2: Analysis of 10,000 individuals. *J Clin Oncol.* 2002 20:1480-1490.
22. Gonzalez-Angolo AM, et al. Incidence and outcome of BRCA mutations in unselected patients with triple-receptor negative breast cancer. *Clin Cancer Res.* 2011 Mar 1;17(5):1082-9.
23. Hartman, A.-R., Kaldate, R. R., Sailer, L. M., et al. (2011), Prevalence of BRCA mutations in an unselected population of triple-negative breast cancer. *Cancer.* doi: 10.1002/cncr.2657
24. Comen E, et al. Relative contributions of BRCA1 and BRCA2 mutations to "triple-negative" breast cancer in Ashkenazi women. *Breast Cancer Res Treat.* 2011 Mar 11. [Epub ahead of print]
25. Young SR, et al. The prevalence of BRCA1 mutations among young women with triple negative breast cancer. *BMC Cancer.* 2009 Mar 19;9:86.
26. Ferrone CR, et al. BRCA germline mutations in Jewish patients with pancreatic adenocarcinoma. *J Clin Oncol.* 2009;27(3):433-8.
27. Couch FJ, et al. The prevalence of BRCA2 mutations in familial pancreatic cancer. *Cancer Epidemiol Biomarkers Prev.* 2007 Feb;16(2):342-6.
28. Hall MJ, et al. Family history of pancreatic cancer in a high-risk cancer clinic: implications for risk assessment. *J Genet Couns.* 2008 Aug;17(4):365-72.
29. Slater EP, et al. Prevalence of BRCA2 and CDKN2a mutations in German familial pancreatic cancer families. *Fam Cancer.* 2010 Sep;9(3):335-43.
30. Hahn SA, et al. BRCA2 germline mutations in familial pancreatic carcinoma. *J Natl Cancer Inst.* 2003 Feb 5;95(3):214-21.
31. Murphy MM, et al. Evaluation of candidate genes MAP2K4, MADH4, ACVR1B, and BRCA2 in familial pancreatic cancer: deleterious BRCA2 mutations in 17%. *Cancer Res.* 2002 Jul 1;62(13):3789-93.
32. Shannon KM, et al. Model-based predictions of BRCA1/2 mutation status in breast carcinoma patients treated at an academic medical center. *Cancer.* 2002 Jan 15;94(2):305-13.
33. Dominguez FJ, et al. Prevalence of hereditary breast/ovarian carcinoma risk in patients with a personal history of breast or ovarian carcinoma in a mammography population. *Cancer.* 2005 Nov 1;104(9):1849-53.
34. Hughes KS, et al. Prevalence of family history of breast and ovarian cancer in a single primary care practice using a self-administered questionnaire. *Breast J.* 2003 Jan-Feb;9(1):19-25.
35. Bellcross CA et al. Evaluation of a breast/ovarian cancer genetics referral screening tool in a mammography population. *Genetics in Medicine.* 2009;11:783-789.
36. Hall IJ, Middlebrooks A and Coughlin SS. Population prevalence of first-degree family history of breast and ovarian cancer in the United States: implications for genetic testing. *The Open Health Services and Policy Journal.* 2008;1:4-37.
37. Antoniou A, et al. Average risks of breast and ovarian cancer associated with BRCA1 or BRCA2 mutations detected in case series unselected for family history: a combined analysis of 22 studies. *Am J Hum Genet.* 2003;72:1117-1130.
38. King MC, et al. Breast and Ovarian Cancer Risks Due to Inherited Mutations in BRCA1 and BRCA2. *Science.* Oct 24 2003:643-646.
39. Chen S, Parmigiani G. Meta-analysis of BRCA1 and BRCA2 penetrance. *J Clin Oncol.* 2007 Apr 10;25(11):1329-33.
40. Easton DF, et al. Breast and ovarian cancer incidence in BRCA1-mutation carriers. *Am J Hum Genet.* 1995;56:265-271.

41. Whittemore AS, Gong G, Itnyre J. Prevalence and contribution of BRCA1 mutations in breast cancer and ovarian cancer: results from three U.S. population-based case-control studies of ovarian cancer. *Am J Hum Genet.* 1997;60:496-504.
42. DevCan: Probability of Developing or Dying of Cancer Software, Version 6.0. Statistical Research and Applications Branch, National Cancer Institute, 2005. <http://srab.cancer.gov/devcan>.
43. The Breast Cancer Linkage Consortium: Cancer Risks in BRCA2 Mutation Carriers. *J Natl Cancer Inst.* 1999;15:1310-1316.
44. Begg CB, et al. Variation of Breast Cancer Risk Among BRCA1/2 Carriers. *JAMA.* 2008 Jan 9;299(2):194-201.
45. Metcalfe K, et al. Contralateral breast cancer in BRCA1 and BRCA2 mutation carriers. *J Clin Oncol.* 2004;22:2328-2335.
46. Frank TS, et al. Sequence analysis of BRCA1 and BRCA2: correlation of mutations with family history and ovarian cancer risk. *J Clin Oncol.* 1998;16:2417-242.
47. Verhoog LC, et al. Survival and tumour characteristics of breast-cancer patients with germline mutations of BRCA1. *Lancet.* 1998;351:316-321.
48. Metcalfe KA, et al. The risk of ovarian cancer after breast cancer in BRCA1 and BRCA2 carriers. *Gynecol Oncol.* 2005 Jan;96(1):222-6.
49. Chen Y, et al. Epidemiology of contralateral breast cancer. *Cancer Epidemiol Biomarkers Prev.* 1999 Oct;8(10):855-61.
50. Graeser MK, et al. Contralateral Breast Cancer Risk in BRCA1 and BRCA2 Mutation Carriers. *J Clin Oncol.* 2009 Dec 10;27(35):5862-4.
51. Malone et al. Population-Based Study of the Risk of Second Primary Contralateral Breast Cancer Associated With Carrying a Mutation in BRCA1 or BRCA2. *J Clin Oncol.* 2010 May 10;28(14):2404
52. NCCN Clinical Practice Guidelines in Oncology v.1.2020 Genetic/Familial High-Risk Assessment: Breast and Ovarian. Accessed at www.nccn.org
53. Committee on Practice Bulletins—Gynecology, Committee on Genetics, Society of Gynecologic Oncology. Practice Bulletin No 182: Hereditary Breast and Ovarian Cancer Syndrome. *Obstet Gynecol.* 2017;130:e110–e126.
54. Liede A, et al. Cancer risks for male carriers of germline mutations in BRCA1 or BRCA2: A review of the literature. *J Clin Oncol.* 2004;22:735-42.
55. Kirchoff T, et al. BRCA mutations and risk of prostate cancer in Ashkenazi Jews. *Clinical Cancer Research.* 2004;10:2918-2.
56. Thompson D, et al. Cancer Incidence in BRCA1 mutation carriers. *J Natl Cancer Inst.* 2002;94:135.
57. Lal G, et al. Inherited predisposition to pancreatic adenocarcinoma: role of family history and germ-line p16, BRCA1, and BRCA2 mutations. *Cancer Res.* 2000 Jan 15;60(2):409-16.
58. Brose MS, et al. Cancer Risk Estimates for BRCA1 Mutation Carriers Identified in a Risk Evaluation Program. *J Natl Cancer Inst.* 2002;94:1365-72.
59. Lin KM, et al. Colorectal cancer in hereditary breast cancer kindreds. *Diseases of the Colon and Rectum* 1999;42:1041-1045.
60. Domchek et al. Breast cancer risks in individuals testing negative for a known family mutation in BRCA1 or BRCA2. *Breast Cancer Res Treat.* 2010 Jan;119(2):409-14.
61. Korde LA, Mueller CM, Loud JT, et al.: No evidence of excess breast cancer risk among mutation-negative women from BRCA mutation-positive families. *Breast Cancer Res Treat* 125 (1): 169-73, 2011.
62. NCCN Clinical Practice Guidelines in Oncology v.2.2019 Prostate Cancer Early Detection Accessed at www.nccn.org
63. NCCN Clinical Practice Guidelines in Oncology v.4.2019 Prostate Cancer Accessed at www.nccn.org
64. NCCN Clinical Practice Guidelines in Oncology v.1.2020 Pancreatic Cancer Accessed at www.nccn.org
65. Momozawa Y, et al. Germline pathogenic variants in 7,636 Japanese patients with prostate cancer and 12,366 controls. *J Natl Cancer Inst.* 2019 Jun 19. pii: djz124. doi: 10.1093/jnci/djz124. [Epub ahead of print] PubMed PMID: 31214711.
66. Pilarski R. The Role of BRCA Testing in Hereditary Pancreatic and Prostate Cancer Families. *Am Soc Clin Oncol Educ Book.* 2019 Jan;39:79-86. doi: 10.1200/EDBK_238977. Epub 2019 May 17. Review. PubMed PMID: 31099688.

67. Nyberg T, et al. Prostate Cancer Risks for Male BRCA1 and BRCA2 Mutation Carriers: A Prospective Cohort Study. *Eur Urol.* 2020 Jan;77(1):24-35. doi: 10.1016/j.eururo.2019.08.025. Epub 2019 Sep 6. PubMed PMID: 31495749; PubMed Central PMCID: PMC6926480.
68. Giri VN, et al. Role of Genetic Testing for Inherited Prostate Cancer Risk: Philadelphia Prostate Cancer Consensus Conference 2017. *J Clin Oncol.* 2018 Feb 1;36(4):414-424. doi: 10.1200/JCO.2017.74.1173. Epub 2017 Dec 13. Review. PubMed PMID: 29236593; PubMed Central PMCID: PMC6075860. (Giri 2018A)
69. Stoffel EM, et al. Evaluating Susceptibility to Pancreatic Cancer: ASCO Provisional Clinical Opinion. *J Clin Oncol.* 2019 Jan 10;37(2):153-164. doi: 10.1200/JCO.18.01489. Epub 2018 Nov 20. PubMed PMID: 30457921.
70. Page EC, et al. Interim Results from the IMPACT Study: Evidence for Prostate-specific Antigen Screening in BRCA2 Mutation Carriers. *Eur Urol.* 2019 Dec;76(6):831-842. doi: 10.1016/j.eururo.2019.08.019. Epub 2019 Sep 16. PubMed PMID: 31537406; PubMed Central PMCID: PMC6880781.
71. Mitra A, et al. Prostate cancer in male BRCA1 and BRCA2 mutation carriers has a more aggressive phenotype. *Br J Cancer.* 2008 Jan 29;98(2):502-7. doi: 10.1038/sj.bjc.6604132. Epub 2008 Jan 8. PubMed PMID: 18182994; PubMed Central PMCID: PMC2361443.
72. Pritchard CC, et al. Inherited DNA-Repair Gene Mutations in Men with Metastatic Prostate Cancer. *N Engl J Med.* 2016 Aug 4;375(5):443-53. doi: 10.1056/NEJMoa1603144. Epub 2016 Jul 6. PubMed PMID: 27433846; PubMed Central PMCID: PMC4986616.
73. Giri VN, et al. Germline genetic testing for inherited prostate cancer in practice: Implications for genetic testing, precision therapy, and cascade testing. *Prostate.* 2018 Mar;79(4):333-339. doi: 10.1002/pros.23739. Epub 2018 Nov 18. PubMed PMID: 30450585. (Giri 2018B)
74. Salo-Mullen EE, O'Reilly EM, Kelsen DP, Ashraf AM, Lowery MA, Yu KH, Reidy DL, Epstein AS, et al. Identification of germline genetic mutations in patients with pancreatic cancer. *Cancer.* 2015 Dec 15;121(24):4382-8. doi: 10.1002/cncr.29664. Epub 2015 Oct 6. PubMed PMID: 26440929; PubMed Central PMCID: PMC5193099.
75. Chaffee KG, et al. Prevalence of germ-line mutations in cancer genes among pancreatic cancer patients with a positive family history. *Genet Med.* 2018 Jan;20(1):119-127. doi: 10.1038/gim.2017.85. Epub 2017 Jul 20. PubMed PMID: 28726808; PubMed Central PMCID: PMC5760284.
76. Hu C, et al. Association Between Inherited Germline Mutations in Cancer Predisposition Genes and Risk of Pancreatic Cancer. *JAMA.* 2018 Jun 19;319(23):2401-2409. doi: 10.1001/jama.2018.6228. PubMed PMID: 29922827; PubMed Central PMCID: PMC6092184. (Hu 2018A)
77. Hu C, et al. Multigene Hereditary Cancer Panels Reveal High-Risk Pancreatic Cancer Susceptibility Genes. *JCO Precis Oncol.* 2018;2. doi: 10.1200/PO.17.00291. Epub 2018 Jul 25. PubMed PMID: 31497750; PubMed Central PMCID: PMC6731034. (Hu 2018B)
78. Maynard H, et al. Germline alterations in patients with biliary tract cancers: A spectrum of significant and previously underappreciated findings. *Cancer.* 2020 Feb 3. doi: 10.1002/cncr.32740. [Epub ahead of print] PubMed PMID: 32012241.
79. Terashima T, et al. Germline mutations in cancer-predisposition genes in patients with biliary tract cancer. *Oncotarget.* 2019 Oct 15;10(57):5949-5957. doi: 10.18632/oncotarget.27224. eCollection 2019 Oct 15. PubMed PMID: 31666926; PubMed Central PMCID: PMC6800267.
80. Wardell CP, et al. Genomic characterization of biliary tract cancers identifies driver genes and predisposing mutations. *J Hepatol.* 2018 May;68(5):959-969. doi: 10.1016/j.jhep.2018.01.009. Epub 2018 Jan 31. PubMed PMID: 29360550.
81. Takai E, et al. Germline mutations in Japanese familial pancreatic cancer patients. *Oncotarget.* 2016 Nov 8;7(45):74227-74235. doi: 10.18632/oncotarget.12490. PubMed PMID: 27732944; PubMed Central PMCID: PMC5342048.

【Approval Conditions】

1. It is required to undergo the document-based or on-site inspection in compliance with the provisions of Article 23-2-5, Paragraph 9 of the Act on Securing Quality, Efficacy and Safety of Pharmaceuticals, Medical Devices, Regenerative and Cellular Therapy Products, Gene Therapy Products, and Cosmetics (Act No.145 of 1960, hereinafter referred to as the “Act”) pursuant to the provisions of Article 23-2-17, Paragraph 5 of the Act, and to submit an operational update report of the variant classification process to PMDA one year after approval and annually thereafter for the time being.

2. The applicant is required to take necessary measures and appropriate management strategies to ensure that blood samples and information obtained from them are not used for the purposes other than the items specified in the attached application form, and also implement the latest measures to ensure security and privacy to prevent unauthorized access.

3. The quality control of the variant information must be conducted as specified in the Remarks section of the attached application. When it is planned to revise the quality control method of the variant information specified in the Remarks section of the attached application (except for a minor revision defined in the MHLW Ordinance specified in Article 23-2-5, Paragraph 15 of the Act), the revision must be approved by the Minister of Health Labor and Welfare in compliance with the provisions of Article 23-2-5, Paragraph 15 of the Act pursuant to the provisions of Article 23-2-17, Paragraph 5 of the Act. The applicant should note that, for the approval of the above revision, the provisions of Article 23-2-5, Paragraph 17, Article 23-2-6 and Article 23-2-7 of the Act are applied pursuant to the provisions of Article 23-2-17, Paragraph 6 of the Act.

【Name or Designation etc. of Marketing Authorization Holder or Manufacturer】

Foreign Special Approval Holder

Myriad Genetic Laboratories, Inc. (*in Japanese*)
(Myriad Genetic Laboratories, Inc.)
US

DMAH

Micren Healthcare Co., Ltd.
TEL: 03-3513-6641

Foreign Manufacturer

Myriad Genetic Laboratories, Inc. (*in Japanese*)
(Myriad Genetic Laboratories, Inc.)
US

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1. Ford D, et al. Risks of cancer in BRCA1-mutation carriers. Breast Cancer Linkage Consortium. Lancet. 1994 343:692-5. PMID: 7907678.
2. Ford D, Easton DF, Stratton M, et al. Genetic heterogeneity and penetrance analysis of the BRCA1 and BRCA2 genes in breast cancer families. Am J Hum Genet. 1998;62:676-89.
3. Chen S, et al. Characterization of BRCA1 and BRCA2 mutations in a large United States sample. J Clin Oncol. 2006 24:863-71. PMID: 16484695.
4. Mavaddat N, et al. Cancer risks for BRCA1 and BRCA2 mutation carriers: results from prospective analysis of EMBRACE. J Natl Cancer Inst. 2013 105:812-22. PMID: 23628597.
5. Tai YC, et al. Breast cancer risk among male BRCA1 and BRCA2 mutation carriers. J Natl Cancer Inst. 2007 99:1811-4. PMID: 18042939.

6. Struwing JP, et al. The risk of cancer associated with specific mutations of BRCA1 and BRCA2 among Ashkenazi Jews. N Engl J Med. 1997 336:1401-8. PMID: 9145676.
7. Lynch HT, et al. BRCA1 and pancreatic cancer: pedigree findings and their causal relationships. Cancer Genet Cytogenet. 2005 158:119-25. PMID: 15796958.
8. van Asperen CJ, et al. Netherlands Collaborative Group on Hereditary Breast Cancer (HEBON). Cancer risks in BRCA2 families: estimates for sites other than breast and ovary. J Med Genet. 2005 42:711-9. PMID: 16141007.
9. Goggins M, et al. International Cancer of the Pancreas Screening (CAPS) consortium. Management of patients with increased risk for familial pancreatic cancer: updated recommendations from the International Cancer of the Pancreas Screening (CAPS) Consortium. Gut. 2020 Jan;69(1):7-17. doi: 10.1136/gutjnl-2019-319352. Epub 2019 Oct 31. PubMed PMID: 31672839.
10. Gumaste PV, et al. Skin cancer risk in BRCA1/2 mutation carriers. Br J Dermatol. 2015 172:1498-506. Epub 2015 Apr 29. PMID: 25524463.
11. Moran A, et al. Risk of cancer other than breast or ovarian in individuals with BRCA1 and BRCA2 mutations. Fam Cancer. 2012 11:235-42. PMID: 22187320.
12. Hall MJ, et al. BRCA1 and BRCA2 mutations in women of different ethnicities undergoing testing for hereditary breast-ovarian cancer. Cancer. 2009 115(10):2222-33. PMID:19241424.
13. Nakamura S et al. Prevalence and differentiation of hereditary breast and ovarian cancers in Japan. Breast Cancer. 2015 22(5) 462-8. 24249303.
14. Antoniou AC, Gayther SA, Stratton JF, Ponder BA, Easton DF Risk models for familial ovarian and breast cancer. Genet Epidemiol. 2000;18:173-90.
15. Coughlin SS, Khoury MJ, Steinberg KK. BRCA1 and BRCA2 gene mutations and risk of breast cancer: public health perspectives. Am J Prev Med. 1999;16(2):91-98.
16. Anglican Breast Cancer Study Group. Prevalence and penetrance of BRCA1 and BRCA2 mutations in a population-based series of breast cancer cases. Br J Cancer. 2000;83:1301-8.
17. Claus EB, Schildkraut JM, Thompson WD, Risch NJ, et al. The genetic attributable risk of breast and ovarian cancer. Cancer. 1996;77:2318-2324.
18. Zhang S, et al. Frequencies of BRCA1 and BRCA2 mutations among 1,342 unselected patients with invasive ovarian cancer. Gynecol Oncol. 2011 May 1;121(2):352-7.
19. Risch HA, McLaughlin JR, Cole DEC, et al. Prevalence and penetrance of germline BRCA1 and BRCA2 mutations in a population series of 649 women with ovarian cancer. Am J Hum Genet. 2001;68:700-710.
20. Pal T, et al. BRCA1 and BRCA2 mutations account for a large proportion of ovarian carcinoma cases. Cancer. 2005;104(12):2807-2816.
21. Frank TS, Deffenbaugh AM, Reid JE et al. Clinical characteristics of individuals with germline mutations in BRCA1 and BRCA2: Analysis of 10,000 individuals. J Clin Oncol. 2002 20:1480-1490.
22. Gonzalez-Angolo AM, et al. Incidence and outcome of BRCA mutations in unselected patients with triple-receptor negative breast cancer. Clin Cancer Res. 2011 Mar 1;17(5):1082-9.
23. Hartman, A.-R., Kaldate, R. R., Sailer, L. M., et al. (2011), Prevalence of BRCA mutations in an unselected population of triple-negative breast cancer. Cancer. doi: 10.1002/cncr.2657
24. Comen E, et al. Relative contributions of BRCA1 and BRCA2 mutations to "triple-negative" breast cancer in Ashkenazi women. Breast Cancer Res Treat. 2011 Mar 11. [Epub ahead of print]

25. Young SR, et al. The prevalence of BRCA1 mutations among young women with triple negative breast cancer. *BMC Cancer*. 2009 Mar 19;9:86.
26. Ferrone CR, et al. BRCA germline mutations in Jewish patients with pancreatic adenocarcinoma. *J Clin Oncol*. 2009;27(3):433-8.
27. Couch FJ, et al. The prevalence of BRCA2 mutations in familial pancreatic cancer. *Cancer Epidemiol Biomarkers Prev*. 2007 Feb;16(2):342-6.
28. Hall MJ, et al. Family history of pancreatic cancer in a high-risk cancer clinic: implications for risk assessment. *J Genet Couns*. 2008 Aug;17(4):365-72.
29. Slater EP, et al. Prevalence of BRCA2 and CDKN2a mutations in German familial pancreatic cancer families. *Fam Cancer*. 2010 Sep;9(3):335-43.
30. Hahn SA, et al. BRCA2 germline mutations in familial pancreatic carcinoma. *J Natl Cancer Inst*. 2003 Feb 5;95(3):214-21.
31. Murphy MM, et al. Evaluation of candidate genes MAP2K4, MADH4, ACVR1B, and BRCA2 in familial pancreatic cancer: deleterious BRCA2 mutations in 17%. *Cancer Res*. 2002 Jul 1;62(13):3789-93.
32. Shannon KM, et al. Model-based predictions of BRCA1/2 mutation status in breast carcinoma patients treated at an academic medical center. *Cancer*. 2002 Jan 15;94(2):305-13.
33. Dominguez FJ, et al. Prevalence of hereditary breast/ovarian carcinoma risk in patients with a personal history of breast or ovarian carcinoma in a mammography population. *Cancer*. 2005 Nov 1;104(9):1849-53.
34. Hughes KS, et al. Prevalence of family history of breast and ovarian cancer in a single primary care practice using a self-administered questionnaire. *Breast J*. 2003 Jan-Feb;9(1):19-25.
35. Bellcross CA et al. Evaluation of a breast/ovarian cancer genetics referral screening tool in a mammography population. *Genetics in Medicine*. 2009;11:783-789.
36. Hall IJ, Middlebrooks A and Coughlin SS. Population prevalence of first-degree family history of breast and ovarian cancer in the United States: implications for genetic testing. *The Open Health Services and Policy Journal*. 2008;1:4-37.
37. Antoniou A, et al. Average risks of breast and ovarian cancer associated with BRCA1 or BRCA2 mutations detected in case series unselected for family history: a combined analysis of 22 studies. *Am J Hum Genet*. 2003;72:1117-1130.
38. King MC, et al. Breast and Ovarian Cancer Risks Due to Inherited Mutations in BRCA1 and BRCA2. *Science*. Oct 24 2003:643-646.
39. Chen S, Parmigiani G. Meta-analysis of BRCA1 and BRCA2 penetrance. *J Clin Oncol*. 2007 Apr 10;25(11):1329-33.
40. Easton DF, et al. Breast and ovarian cancer incidence in BRCA1-mutation carriers. *Am J Hum Genet*. 1995;56:265-271.
41. Whittemore AS, Gong G, Itnyre J. Prevalence and contribution of BRCA1 mutations in breast cancer and ovarian cancer: results from three U.S. population-based case-control studies of ovarian cancer. *Am J Hum Genet*. 1997;60:496-504.
42. DevCan: Probability of Developing or Dying of Cancer Software, Version 6.0. Statistical Research and Applications Branch, National Cancer Institute, 2005. <http://srab.cancer.gov/devcan>.
43. The Breast Cancer Linkage Consortium: Cancer Risks in BRCA2 Mutation Carriers. *J Natl Cancer Inst*. 1999;15:1310-1316.
44. Begg CB, et al. Variation of Breast Cancer Risk Among BRCA1/2 Carriers. *JAMA*. 2008 Jan 9;299(2):194-201.
45. Metcalfe K, et al. Contralateral breast cancer in BRCA1 and BRCA2 mutation carriers. *J Clin Oncol*. 2004;22:2328-2335.
46. Frank TS, et al. Sequence analysis of BRCA1 and BRCA2: correlation of mutations with family history and ovarian cancer risk. *J Clin Oncol*. 1998;16:2417-242.
47. Verhoog LC, et al. Survival and tumour characteristics of breast-cancer patients with germline mutations of BRCA1. *Lancet*. 1998;351:316-321.
48. Metcalfe KA, et al. The risk of ovarian cancer after breast cancer in BRCA1 and BRCA2 carriers. *Gynecol Oncol*. 2005 Jan;96(1):222-6.
49. Chen Y, et al. Epidemiology of contralateral breast cancer. *Cancer Epidemiol Biomarkers Prev*. 1999 Oct;8(10):855-61.
50. Graeser MK, et al. Contralateral Breast Cancer Risk in BRCA1 and BRCA2 Mutation Carriers. *J Clin Oncol*. 2009 Dec 10;27(35):5862-4.
51. Malone et al. Population-Based Study of the Risk of Second Primary Contralateral Breast Cancer Associated With Carrying a Mutation in BRCA1 or BRCA2. *J Clin Oncol*. 2010 May 10;28(14):2404
52. NCCN Clinical Practice Guidelines in Oncology v.1.2020 Genetic/Familial High-Risk Assessment: Breast and Ovarian. Accessed at www.nccn.org
53. Committee on Practice Bulletins—Gynecology, Committee on Genetics, Society of Gynecologic Oncology. Practice Bulletin No 182: Hereditary Breast and Ovarian Cancer Syndrome. *Obstet Gynecol*. 2017;130:e110–e126
54. Liede A, et al. Cancer risks for male carriers of germline mutations in BRCA1 or BRCA2: A review of the literature. *J Clin Oncol*. 2004;22:735-42.
55. Kirchoff T, et al. BRCA mutations and risk of prostate cancer in Ashkenazi Jews. *Clinical Cancer Research*. 2004;10:2918-2.
56. Thompson D, et al. Cancer Incidence in BRCA1 mutation carriers. *J Natl Cancer Inst*. 2002;94:135.
57. Lal G, et al. Inherited predisposition to pancreatic adenocarcinoma: role of family history and germ-line p16, BRCA1, and BRCA2 mutations. *Cancer Res*. 2000 Jan 15;60(2):409-16.
58. Brose MS, et al. Cancer Risk Estimates for BRCA1 Mutation Carriers Identified in a Risk Evaluation Program. *J Natl Cancer Inst*. 2002;94:1365-72.
59. Lin KM, et al. Colorectal cancer in hereditary breast cancer kindreds. *Diseases of the Colon and Rectum* 1999;42:1041-1045.
60. Domchek et al. Breast cancer risks in individuals testing negative for a known family mutation in BRCA1 or BRCA2. *Breast Cancer Res Treat*. 2010 Jan;119(2):409-14.
61. Korde LA, Mueller CM, Loud JT, et al.: No evidence of excess breast cancer risk among mutation-negative women from BRCA mutation-positive families. *Breast Cancer Res Treat* 125 (1): 169-73, 2011.
62. NCCN Clinical Practice Guidelines in Oncology v.2.2019 Prostate Cancer Early Detection Accessed at www.nccn.org
63. NCCN Clinical Practice Guidelines in Oncology v.4.2019 Prostate Cancer Accessed at www.nccn.org
64. NCCN Clinical Practice Guidelines in Oncology v.1.2020 Pancreatic Cancer Accessed at www.nccn.org
65. Momozawa Y, et al. Germline pathogenic variants in 7,636 Japanese patients with prostate cancer and 12,366 controls. *J Natl Cancer Inst*. 2019 Jun 19. pii: djz124. doi: 10.1093/jnci/djz124. [Epub ahead of print] PubMed PMID: 31214711.
66. Pilarski R. The Role of BRCA Testing in Hereditary Pancreatic and Prostate Cancer Families. *Am Soc Clin Oncol Educ Book*. 2019 Jan;39:79-86. doi: 10.1200/EDBK_238977. Epub 2019 May 17. Review. PubMed PMID: 31099688.
67. Nyberg T, et al. Prostate Cancer Risks for Male BRCA1 and BRCA2 Mutation Carriers: A Prospective Cohort Study. *Eur Urol*. 2020 Jan;77(1):24-35. doi: 10.1016/j.eururo.2019.08.025. Epub 2019 Sep 6. PubMed PMID: 31495749; PubMed Central PMCID: PMC6926480.
68. Giri VN, et al. Role of Genetic Testing for Inherited Prostate Cancer Risk: Philadelphia Prostate Cancer Consensus Conference 2017. *J Clin Oncol*. 2018 Feb 1;36(4):414-424. doi: 10.1200/JCO.2017.74.1173. Epub 2017 Dec 13. Review. PubMed PMID: 29236593; PubMed Central PMCID: PMC6075860. (Giri 2018A)

69. Stoffel EM, et al. Evaluating Susceptibility to Pancreatic Cancer: ASCO Provisional Clinical Opinion. J Clin Oncol. 2019 Jan 10;37(2):153-164. doi: 10.1200/JCO.18.01489. Epub 2018 Nov 20. PubMed PMID: 30457921.
70. Page EC, et al. Interim Results from the IMPACT Study: Evidence for Prostate-specific Antigen Screening in BRCA2 Mutation Carriers. Eur Urol. 2019 Dec;76(6):831-842. doi: 10.1016/j.eururo.2019.08.019. Epub 2019 Sep 16. PubMed PMID: 31537406; PubMed Central PMCID: PMC6880781.
71. Mitra A, et al. Prostate cancer in male BRCA1 and BRCA2 mutation carriers has a more aggressive phenotype. Br J Cancer. 2008 Jan 29;98(2):502-7. doi: 10.1038/sj.bjc.6604132. Epub 2008 Jan 8. PubMed PMID: 18182994; PubMed Central PMCID: PMC2361443.
72. Pritchard CC, et al. Inherited DNA-Repair Gene Mutations in Men with Metastatic Prostate Cancer. N Engl J Med. 2016 Aug 4;375(5):443-53. doi: 10.1056/NEJMoa1603144. Epub 2016 Jul 6. PubMed PMID: 27433846; PubMed Central PMCID: PMC4986616.
73. Giri VN, et al. Germline genetic testing for inherited prostate cancer in practice: Implications for genetic testing, precision therapy, and cascade testing. Prostate. 2018 Mar;79(4):333-339. doi: 10.1002/pros.23739. Epub 2018 Nov 18. PubMed PMID: 30450585. (Giri 2018B)
74. Salo-Mullen EE, O'Reilly EM, Kelsen DP, Ashraf AM, Lowery MA, Yu KH, Reidy DL, Epstein AS, et al. Identification of germline genetic mutations in patients with pancreatic cancer. Cancer. 2015 Dec 15;121(24):4382-8. doi: 10.1002/cncr.29664. Epub 2015 Oct 6. PubMed PMID: 26440929; PubMed Central PMCID: PMC5193099.
75. Chaffee KG, et al. Prevalence of germ-line mutations in cancer genes among pancreatic cancer patients with a positive family history. Genet Med. 2018 Jan;20(1):119-127. doi: 10.1038/gim.2017.85. Epub 2017 Jul 20. PubMed PMID: 28726808; PubMed Central PMCID: PMC5760284.
76. Hu C, et al. Association Between Inherited Germline Mutations in Cancer Predisposition Genes and Risk of Pancreatic Cancer. JAMA. 2018 Jun 19;319(23):2401-2409. doi: 10.1001/jama.2018.6228. PubMed PMID: 29922827; PubMed Central PMCID: PMC6092184. (Hu 2018A)
77. Hu C, et al. Multigene Hereditary Cancer Panels Reveal High-Risk Pancreatic Cancer Susceptibility Genes. JCO Precis Oncol. 2018;2. doi: 10.1200/PO.17.00291. Epub 2018 Jul 25. PubMed PMID: 31497750; PubMed Central PMCID: PMC6731034. (Hu 2018B)
78. Maynard H, et al. Germline alterations in patients with biliary tract cancers: A spectrum of significant and previously underappreciated findings. Cancer. 2020 Feb 3. doi: 10.1002/cncr.32740. [Epub ahead of print] PubMed PMID: 32012241.
79. Terashima T, et al. Germline mutations in cancer-predisposition genes in patients with biliary tract cancer. Oncotarget. 2019 Oct 15;10(57):5949-5957. doi: 10.18632/oncotarget.27224. eCollection 2019 Oct 15. PubMed PMID: 31666926; PubMed Central PMCID: PMC6800267.
80. Wardell CP, et al. Genomic characterization of biliary tract cancers identifies driver genes and predisposing mutations. J Hepatol. 2018 May;68(5):959-969. doi: 10.1016/j.jhep.2018.01.009. Epub 2018 Jan 31. PubMed PMID: 29360550.
81. Takai E, et al. Germline mutations in Japanese familial pancreatic cancer patients. Oncotarget. 2016 Nov 8;7(45):74227-74235. doi: 10.18632/oncotarget.12490. PubMed PMID: 27732944; PubMed Central PMCID: PMC5342048.

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