The Contrology® Reformer with Tower

THIS IS A MEDICAL DEVICE


STATEMENT OF INTENDED USE

The Contrology™ Reformer device is intended by the manufacturer to be used by people for the purposes of physiotherapy, rehabilitation, the prevention, treatment and relief of injuries, supportive, functional therapy for a disability or for general strengthening and improving the mobility of the human body. Should injury occur during usage, report to the manufacturer and competent authority of the Member State in which the user and/or patient is established.

Important: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions in document subject to change. Please visit www.pilates.com/support/assembly-use-maintenance for current version.

Visit pilates.com/patent for complete and current information on Balanced Body product patents.

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

REFORMER PARTS LIST

DESCRIPTION	PART NO.	QTY
Reformer Frame	N/A	1
Carriage Assembly	N/A	1
Shoulder Assembly (pair)	17976	1
Mat Conversion	N/A	1
Frame foot spacer pads	16839	4
36" Wooden Pole	16457	1
1/4" Black Kneeling Pad	10248	2
Sitting Box	16147	1
Footstrap	16102	2
Extension Straps	16152	2
Leather Straps (pair)	16153	1
Leather Handles with straight wood (pair)	16154	1
Contrology Reformer Springs	16228	4
Contrology Springbar	16100	1
Wooden Gear Block Stopper	16343	1
Foam Spacer, Sitting Box	N/A	2
Legs, Mat Conversion	16903	4
Mat caster Cover	17644	1

* 86" Reformer: (Leather Straps P/N: 16360, Carriage Block P/N: 16484, Mat Conversion. P/N: 17412). 89" Reformer: (Leather Straps P/N: 16361, Carriage Block P/N: 16505, Mat Conv. P/N: 17414)

INSTALL THE CARRIAGE

Clean all the carriage wheels and the metal rail surfaces of the frame. Carefully place the carriage into the frame with the headrest facing away from the footbar.

INSTALLING/REMOVING THE SHOULDER RESTS

1. On the bottom of the shoulder rest brackets, there are protruding studs. These fit into the plate and secure the shoulder rest. Lower the shoulder rest so the studs fit into the key holes on the plate on the carriage. Slide the shoulder rest away from the carriage pad until it stops. See Figure 1.


Figure 1

2. Secure the shoulder rest by turning the post clockwise until tight. See Figure 2. Overtightening these will make them difficult to remove. They should be snug but not too tight.


Figure 2

3. To remove the shoulder rests, turn the posts counter-clockwise 1-2 turns, gently slide the shoulder rest towards the carriage. Then gently lift the shoulder rest up and out of the key holes.

Note: The shoulder rests can be installed in two positions.

One position is wider and more comfortable for broad shoulders. Swapping the left and right shoulder rests will change between standard and wide configuration.

INSTALL THE STRAPS

Place the handles on the shoulder rests posts. Insert the end of the strap with the O-ring up and into the caster opening from the bottom. The rough side of leather should be facing the caster wheel. As you push the strap up, spin the caster wheel towards the carriage to assist with feeding the strap in the wheel assembly. See Figure 3. Once the O-ring is through the caster bring it towards the carriage and connect it to the handle clip on that side of the reformer.

Figure 3

Guide the strap to the eye bolt on the underside of the carriage. Connect the Quick Link clip to the eye bolt on the same side of the carriage as shown in Figure 4. Repeat on other side.

Figure 4

ATTACH THE SPRINGS TO THE UNDERSIDE OF THE CARRIAGE

Figure 5

Attach one spring to one slot in the bracket on the underside of the carriage as shown in Figure 5. The open end of the spring hook should be facing down. Repeat for other springs.

INSTALL THE SPRINGBAR

Place the springbar with hook ends facing up into the gear angle as shown in Figure 6. Connect at least (1) spring to an eye bolt on the springbar.


Figure 6

CARRIAGE BLOCKS

The Carriage Block at the head end of the Reformer is installed on the frame to limit the distance the carriage can travel, see Figure 7. There are two lengths: the Medium Carriage Block comes with the 86" and Long Carriage Block comes with the 89" Reformer. Not intended for, but can be used on 80" Reformers, must purchase separately.


Figure 7

The Gear Block, used by all Reformer lengths, at the foot end is installed on the left side rail of the gear positions as shown in Figure 8. The Gear Block has two settings, short and long. Figure 8 shows the short position, when spring bar is in the middle or "B" position. Flip the Gear Block over to get the long position, for use when the springbar is in the last or "C" position.


Figure 8

The Carriage Blocks are most commonly used for clients who have a limited range of motion, or to stop the carriage from traveling too far. The Gear Block is used to change the carriage starting position.

INSTALL FOOTSTRAPS


Figure 9

Two footstraps are included with the Reformer. Clip one end of the strap to one eye bolt on the inside foot end of the frame, the other clip to the other eye bolt on the opposite side of frame end. Repeat for second footstrap, see Figure 9.

HOW TO INSTALL EXTENSION STRAP

Make sure that your extension strap has a double ended clip attached at one end. Place the handle and the leather loop end of leather strap together, side by side. See Figure 10. Thread one end of the extension strap through both the leather handle and the loop of the leather strap and attach the double ended clip to the O-ring of the end of the extension strap. See Figure 11.

Figure 10

Figure 11

INSTALL THE OPTIONAL JUMPBOARD (86"/218cm and 89"/226cm models only)

Correct installation of the jumpboard is very important. Lay the footbar flat down on the Reformer. Install the jumpboard down into the installed bracket on the inside of the foot end of the Reformer. The padded side should face the carriage. Properly installed the jumpboard will rest on the foot plate bracket, see Figure 12.


Figure 12

FIXED STRAP HOOK (OPTIONAL)


Figure 13

This Reformer comes with an option to install a Fixed Hook for the straps. This can be installed as shown in Figure 13. Non-Adjustable Strap Kit (16492). If your machine came with a Fixed Strap Hook pre-installed, it may be removed and a Kit for Adjustable Straps (16481) may be purchased separately.

LEVELING

If desired or needed, there are (4) square pads that can be placed one under each leg; if your machine is not sitting level on your floor, or if the machine is moving during use.


TOWER PARTS LIST

DESCRIPTION	PART NO.	QTY
Tower Assembly	N/A	1
Vertical Tower Tubes	16968	2
Horizontal Tower Tube	17323	1
Push Through Bar	17598	1
Push Through Bar Spring (Bottom)	16912	1
Push Through Bar Spring	17780	1
Roll Back Bar	17110	1
Roll Back Bar Springs	16910	2
Light Arm Springs	17109	2
Light Leg Springs	17094	2
Cotton Loop (pair)	16992	1
Storage Chain w/ Quick Link	17648	1
Safety Chain w/ Quick Link	17649	2
Leather Contoured Handle Assembly (pair)	17115	1
Allen Wrench, 5/32"	13778	1
Allen Wrench, 3/16"	13777	1
1/4-20 x 5/8" Button Screw	15223	4
Eye Bolt 2-1/2" long	13744	2
Lock Washer	13633	2
Hex Nut, 1/4"-20	13635	2
Quick Link Clip	17160	2
Spring Sleeve	19026	2

INSTALL THE TOWER

1. Remove the Tower tubes with Push Through Bar from the box. Do not cut the tie between the Push Through Bar and the tubes or else it can swing when lifting and hit you or someone else. On the outside of the two vertical tubes there will be (2) threaded holes, these will face out away from the apparatus when installed.

2. Insert the tubes into the tower brackets on the head end of the apparatus, making sure the threaded holes are facing you when installing. Lower the tubes carefully all the way down the tower brackets. They will touch the ground. See Figure 14.
3. Use the provided (4) 1/4"-20 x 5/8" long screws, to secure the tubes to the brackets. It may be necessary to slightly lift the tube to align the holes. Start all the screws by hand before fully tightening. Use the 5/32" Allen wrench to tighten all (4) screws to the brackets. After screws are tight, remove the tie and packaging from Push Through Bar. Install the 2 1/2" long eye bolts with a nut and a lock washer up against the tube in the bottom hole of brackets. Make sure that the eye bolt is aligned vertically when tight. See figure 14 for reference.
4. Install the Leg Springs to the tower on the desired eye bolts and install the cotton loops to these springs. See 15 and A for reference.
5. Attach the Roll Back Bar Springs to the Roll Back Bar and install the springs to the outer eye bolts on the top cross tube. See Figure 15 and B for reference.
6. Install the Arm Springs to the desired position on the tower and attach the handles to the springs. See Figure 15 and C for reference.

- Figure 15:** A) Leg Springs (the longest springs)
B) Roll Down Bar Springs (have 2 dog clips attached)
C) Arm Springs (only 1 dog clip attached)
7. Attach the Storage Chain (shortest chain) to the middle eye bolt on the top cross tube using the Quick Link.
 8. Attach the Safety Chains, one to each of the outer eye bolts on the top cross tube using the Quick Links.
 9. Attach the Push Through Bar Spring with the quick link to one of the top tube outer eye bolts. The other Push Through Bar Spring will connect to the Reformer frame eye bolts using the open hooks. See Figure 16 and 17 for reference.


Figure 16: D) Push Through Bar Springs
E) Quick Link
F) Double Harness Snap

HOW TO INSTALL YOUR WIRE STOPS ON YOUR CONTROLOGY SPRINGS

PARTS INCLUDED:

DESCRIPTION	PART NO.	QTY
Wire Stop	17922	1
10-24 Set Screw	17931	1
3/32" Allen Wrench	17930	1

Thread a provided set screw into each wire stop, install the set screw into the wire stop so about 1/4 of the set screw sticks out. (Note: the wire stops may already have the screw installed. If so, make sure it doesn't interfere with proper installation).

Go to the bottom mounted spring attached to your apparatus. Install a wire stop onto the end of the open hook. Slide the middle opening of the wire stop onto the wire so that the wire end is barely coming out the underside of the wire stop.

Make sure the wire stop length of the body is perpendicular to the wire. See Figure 18. Use the provided Allen wrench to tighten the set screw. Ensure the set screw is very tight to prevent them from coming loose.


Figure 17


Figure 18

ASSEMBLY MAT CONVERSION

1. You will need to install legs on your Mat Conversion pad before using.
2. Find the (4) feet with a thread sticking out and thread them into the corner inserts on the bottom side of the mat conversion.
3. Turn them clockwise until tight. Hand tight is okay. See Figure 19.


Figure 19

INSTALLING THE MAT CONVERSION

1. Remove all springs that are attached to the spring bar.
2. Remove the shoulder pads and brackets from the carriage.
3. Roll the carriage to the head end of the frame. When close to the end of the frame, drop in the Mat Caster Cover assembly between the frame and the carriage. This will act as a stopper for the carriage and cover the casters. See Figure 20. The flaps go over shoulder plates, when the shoulder rests are removed.


Figure 20

4. Lay the footbar flat on the inside of the frame as shown in Figure 21. If using the wooden gear stopper, remove this so the footbar can lay lower in the frame.


Figure 21

5. Lift and install the Mat Conversion inside the frame. There are (2) cutouts on the underside of the mat to clear the footbar pivots, ensure these are on the foot end of the apparatus when installing. See Figure 22 for reference


Figure 22

CONVERT BACK TO A REFORMER

Remove the Mat Conversion, flip the footbar back into position and ensure the kickstand is secure, remove the Mat Caster cover, roll the carriage to home position and attach at least one spring to the spring bar.

INSTALLING YOUR SITTING BOX SPACERS:

1. The foam spacer reduces the movement between the sitting box and the shoulder posts when installed in Short Box on the carriage.
2. Find the (2) foam spacers, see Figure 23 for what they look like.


Figure 23

3. Flip over the sitting box and place the foam spacer on the bottom with the adhesive paper facing up. Slide this under the wood rib in the middle. See Figure 24.


Figure 24

4. Once the cut in the foam spacer is centered under the wood rib, rotate the foam spacer up 90 degrees and push against the wood side. See Figure 25.


Figure 25

5. Repeat this process of the other spacer.
6. If desired, the adhesive paper can be removed and the adhesive can secure the foam spacer to the side of the sitting box. The foam spacer fits snug in the sitting box area, so it is not necessary to do this, but can be if desired.

Cleaning and Maintenance

CLEANING

You can extend the life of upholstery by keeping it clean and free of dirt, oil and perspiration. After each use, wipe down all skin contact parts, including carriage, shoulder and head rests.

Equipment upholstery is coated with BeautyGard®, which offers antibacterial protection. Balanced Body Clean™ is an all-purpose cleaner that is 100% natural, organic and biodegradable. It cleans, disinfects and removes dirt, grease and oil and will extend the life of your upholstery. A mild soap and water solution may also be used.

- » Be sure to dry thoroughly.
- » For heavier soiling please call Balanced Body Technical Support.
- » Wheel tracks and wheels should be cleaned regularly with a cloth to ensure that wheels roll smoothly and quietly.
- » Wash footstrap regularly. Place in pillowcase, wash on gentle cycle, air dry.
- » Be sure to inspect springs regularly.

If you have any questions on cleaning or maintenance, please call Balanced Body Technical Support.

LUBRICATION AND ADJUSTMENT

As needed, lubricate side wheels with 3-in-1 oil.

It is not recommended to lubricate the carriage rolling wheels, do so at your own risk.

IN CASE OF FREIGHT DAMAGE

If there is freight damage, make certain you keep all packaging material. Please call Balanced Body immediately.

Safety First: A guide to proper maintenance and safe use of your Pilates equipment.

For over 45 years, Balanced Body has been introducing safety-related innovations to Pilates equipment. Many of our improvements are now industry standards, resulting in Pilates equipment that's safer today than ever before.

Safety depends on proper maintenance and safe use, in addition to the quality of the equipment. This guide was created to help you use and maintain your equipment for optimum safety. Please read it through carefully and keep for future reference. If you have any questions, give us a call. **Failure to follow these instructions may result in serious injury.**

ALL EQUIPMENT

Springs

Spring inspections are critical to maintain your equipment in safe operating condition. All Balanced Body springs should be replaced at least every two years. Certain environments and usages can shorten the expected life of the springs and you may need to replace the springs more frequently. Therefore, it is very important to inspect springs on a regular basis since worn or old springs lose resilience and may break during use. Injury may result if a spring breaks during use.

During use, do not allow springs to recoil in an uncontrolled manner. This will damage the spring and shorten its expected life.

Inspect springs for gaps and kinks (weekly or monthly, depending on frequency of use). Look for gaps and kinks between the coils when the spring is at rest. It is not unusual for the spring to have a very small gap on the tapered end (a gap is sometimes created during the manufacturing process). However, there should be no gaps in the body of the spring. If you see any gaps or kinks in the body of the spring, discontinue use and replace the springs immediately. See **Figure 1**. Additionally, corrosion anywhere on the coils will shorten the life of the spring. Discontinue using the spring immediately if you see any rust or oxidation during inspection.


Figure 1

Snaps

Inspect snaps for wear (monthly). First, verify that the snap hook is working properly. If the snap hook does not retract and return properly, discontinue using the spring immediately and replace the snap. Eyebolts can cause excessive wear on snap hooks. If the hook shows a lot of wear, discontinue using the spring immediately and call Balanced Body to replace spring or snap. See **Figure 2**.

Figure 2: Good snap: no wear on hook.

Bad snap: excessive wear on hook.

EYEBOLTS, NUTS AND BOLTS

Tighten all equipment bolts and screws (monthly). Verify that all eyebolts, nuts and bolts are tight. See the section titled "How to inspect and tighten nuts and bolts."

ROPES AND STRAPS

Rope and strap wear (quarterly). Ropes should be replaced if you can see the core of the rope through the outer lining, or if the ropes are flattened. Straps should be replaced as soon as any fraying is noticed. Be sure to check the sections of rope or straps that attach to the clips and run through the pulleys.

REFORMERS

Check springbar hooks or eyebolts (quarterly). Balanced Body makes two different springbar systems:

- » Revo Springbar. Make sure springbar hooks and handle are tight.

Standard Springbar. Verify that the nuts securing the springbar hooks are tight. See section titled "How to inspect and tighten nuts and bolts."

Spring rotation (quarterly). You can prolong Reformer spring life by rotating springs of the same weight each quarter. Unhook and move to another position on the springbar. Rotating springs helps them wear more evenly.

Risers on the outside. Wood risers must be installed on the outside of the frame. Risers can loosen over time, so always make sure they are tight.

Springs hooked downward under carriage. Make sure springs are hooked in a downward position. See **Figure 3**.

Secure the carriage. When your Reformer is not in use, be sure that at least two springs secure the carriage to the springbar.

Figure 3: Springs hooked downward

Default settings. Many users have a “default setting” for Reformers. At the end of a session, the user connects a prescribed number of springs in neutral tension, sets the footbar at a pre-determined height, and sets the ropes at a specified length. This ensures that the equipment is ready for the next use, and that the carriage is secured by the springs.

Footstrap under tension in box work. When using the box and footstrap, be sure the footstrap is under tension (with snaps pulling from the top of the eyebolt) before beginning the exercise. See **Figure 4**.

Figure 4: Foot strap under tension

REFORMER WHEEL AND TRACK MAINTENANCE

Clean the tracks and wheels (weekly). For smooth carriage travel and to maintain the longevity of the wheels, we recommend that you wipe down the tracks once a week.

Disconnect the springs and clean the entire length of the tracks with a soft cloth and Balanced Body Cleaner, mild soap with water or a mild commercial cleaner such as, Fantastik® or 409®. Do not use abrasive cleansers or pads, as they can damage the anodizing on the rails. To clean the wheels, hold the cloth against the wheels while you move the carriage. If you feel a bump in the ride, dirt has adhered to the surface of the rails or wheels. Clean hair and debris out of the rails. Hair can wrap around the wheel axles and eventually build up and cause wheel failure. Use tweezers to remove hair from the wheels.

Lubrication. Never spray silicone near or inside the wheels – this can wash the lubricant out of the bearings and ruin the bearings. You can purchase dry silicone at most hardware and auto parts stores. Pulleys sometimes require lubrication to stop a squeak. Direct a very quick spray of dry silicone or Teflon spray into the pulley axle. “Dry” silicone does not have an oil base. Oil-based (“wet”) silicone and WD40 should not be used as they attract dirt. Be careful not to over spray. You may want to remove ropes to avoid getting silicone on them.

Do not lubricate the Allegro 2 rails.

Footbar supports (quarterly). For all Balanced Body footbars with footbar support brackets, verify that the pivot screw attaching the footbar support bracket to footbar is tight, but not so tight that it prevents the support from rotating freely. For Legacy Reformers, tighten the pivot bolt to secure footbar support.

Headrest (monthly). Make sure the hinge screws and bolts on your headrest are tight.

Under the Reformer (monthly). Move Reformers and make sure you clean the floor space underneath.

Standing Platform Footbar Bumpers (wood Reformers only). If your standing platform footbar bumpers (the small plastic pieces that protect the standing platform from the footbar) are broken or damaged, please call Balanced Body to replace.

TRAPEZE TABLE (CADILLAC) & TOWERS

Cotter pins removed. These pins are located in the vertical tubes that align the canopy to the frame and should be removed as soon as installation is complete. Unremoved cotter pins can tear clothing and lacerate the skin. Use pliers to remove the pins.


Figure 5: Cotter pin before and after removal from Trap Table.

Save the pins in case you need to disassemble and reassemble the table for transportation purposes. See **Figure 5**.

Push-Through Bar (PTB) with Sliders.

The PTB moves vertically to accommodate different users and exercises. The sliders on the tubes allow for this vertical movement. Make sure these sliders are clean and easy to move. Before beginning any exercise, ensure that the sliders are properly aligned with the PTB holes and locked into position. Apply a downward force to ensure. If you notice wear on the slider knob pins, please call Balanced Body to replace. Using the PTB in this condition could cause injury.

Weekly maintenance for Push-Through Bars with Sliders. Check to ensure the pins in the plunger knobs protrude and lock into the vertical tubes. Pull both knobs out and move the sliders to a different position. Release the knobs before the next hole and continue to slide the system into position. Once over a hole the pin of the knob will automatically drop into the opening. Once each slider is in a new position apply a downward force on the PTB. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced. Call Balanced Body for replacement parts.

Push-Through Bar (PTB) with T-pin setting. For bottom sprung exercises, if your client's head is below the PTB, use the T-pin setting in addition to the safety strap or chain. Spotting your client is highly recommended. This is important for safety.

Push-Through Bar (PTB) control. Make sure you have enough room around the trap table to safely use the PTB without fear of hitting other people. The PTB can be dangerous if not properly

used. Only trained, experienced users should use the PTB. A spotter should always maintain control of the bar with one hand. If the user should lose control of the bar, the spotter can maintain control of it.

Correct safety strap attachment. For bottom-sprung exercises, the safety strap or chain should always secure the bar.

The safety strap or chain should wrap around the PTB and the canopy frame, not the eyebolts. The strap or chain is only as strong as the weakest link, and the frame and bar are a great deal stronger than eyebolts. **Figure 6.**

Spotting your client is highly recommended. This is important for safety.

Figure 6: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

Setting the PTB for bottom-sprung exercises. For bottom-sprung exercises, the safety strap should be attached so that the angle of the push-through bar is no lower than the 4 o'clock position. This limits the range of the bar and prevents it from potentially coming into contact with the user.

Using the 4th side on the PTB along with the safety strap is highly recommended to prevent injury.

Figure 7: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

CHAIRS

Dismount with control. When dismounting the chair, release the pedals slowly, with control. Don't let the pedal snap back.

Spot users. When a user is standing, sitting or lying on top of the chair, there is increased risk of falling. Standing exercises, in particular, can be unstable. Spotting users will make these exercises safer.

Hourglass spring mounts. If your chair has hourglass spring mounts and the mounts do not successfully retain the springs, replace the fiber washers (they are reddish-brown in color).

Figure 8: Fiber washer

Figure 8. If your chair is a Balanced Body Split-step Pedal Chair (Combo Chair), please consider upgrading to the Cactus Springtree).

UPHOLSTERY CLEANING & MAINTENANCE

Cleaning. You can extend the life of upholstery by keeping it clean and free of dirt, oil and perspiration. After each use, wipe down the upholstery with a solution of mild soap and water. Then wipe it down with clean water and dry with a soft towel.

Disinfecting. Equipment upholstery is coated with BeautyGard®, which offers antibacterial protection. If you want additional disinfection, Balanced Body offers Balanced Body Clean™ disinfecting solution. Use of any other solution (especially those containing essential oils) will shorten the life of some equipment and is not recommended.

For over molded foam padding it is recommended to clean with mild soap or BB Clean and wipe clean with a water dampened cloth to remove any leftover cleaner.

Cleaning Loops. Loops should be washed using mild soap such as Woolite. Hang to dry.

HOW TO INSPECT AND TIGHTEN NUTS AND BOLTS.

Use your fingers to check nuts and bolts for tightness. If you can turn the nut or bolt with your fingers, it's too loose and should be tightened.

To tighten, first tighten using your fingers. Rotate nuts and bolts clockwise to tighten. Insert a screwdriver through eyebolts to hold them steady while you tighten the nuts. Then use a small wrench to tighten the nuts further. **Figure 9.**

Figure 9: Use two fingers to tighten bolts

It is recommended to check that the pins on the PTB protrude and lock into the vertical tubes appropriately. To verify their function, first pull both knobs out and start moving the sliders to a different position. Release the knobs before the next hole and continue to slide the system. Once over a hole the pin of the knob will drop into the opening. Once each slider is in a new position apply a down force on the PTB directly downwards. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced; call Balanced Body for replacement parts.

EQUIPMENT INSPECTION AND MAINTENANCE LOG

We suggest that you keep a maintenance log for each piece of equipment. The log should include:

1. A description of the machine including the serial number, the date and place of purchase, and the manufacturer. All of this information should appear on the invoice.
2. Date and description of all required maintenance and inspections performed.

EQUIPMENT ATTACHMENT POINTS

Before using equipment, ensure that the springs are securely attached. See photos below for examples of securely attached springs to various attachment points that you may encounter on Balanced Body equipment. Note that adding spring resistance makes some exercises easier and others more difficult. Performing exercises with improper spring resistance increases the risk of falling. If unsure of proper spring resistance, please consult a certified Pilates instructor.

"Hour Glass" spring attachment point found on the Allegro 2 and some versions of the Allegro Reformer, as well as the upper (pre-load) attachment points on a Revo springbar.

"Hook" attachment point found on Revo Springbars, Standard Springbars and some versions of the Allegro Reformer. Note the Hour Glass posts in the upper section of this springbar as a pre-load position, shown here without the springs attached.

"Chair Cactus" attachment point found on various Chairs. It may also be flatter and made of sheet metal.

Eyebolt and loop attachment point found on Cadillacs, Wall Units and Towers.

"Hook" attachment point found on Contrology Reformers.

"Reformer Cactus" spring attachment point found on the IQ, Metro IQ, and one version of the Allegro Reformer. It may or may not have notches and bends at the bottoms of the cactuses.

MAINTENANCE SCHEDULE

All Equipment	Day	Wk.	Mo.	Qtr.
Inspect springs for gaps & kinks		✓	✓	
Inspect snaps for wear			✓	
Inspect nuts & bolts for tightness			✓	
Reformers				
Clean wheels and tracks		✓		
Inspect springbar hooks/eyebolts				✓
Rotate springs				✓
Inspect ropes/straps				✓
Inspect footbar supports				✓
Inspect springbars				✓

REPLACEMENT PARTS

To order replacement parts, or if you have any questions, please call:

U.S. and Canada: 1-800-PILATES (1-800-745-2837)
 United Kingdom: 0800 014 8207
 Other locations: +1 916-388-2838

Fax: 916-379-9277
 Email: info@pilates.com
www.pilates.com
 5909 88th Street, Sacramento, CA 95828 USA

ASSEMBLY AND MAINTENANCE PODCASTS

View our library of assembly and maintenance videos at <https://video.pilates.com/categories/bb-garage>.

Or scan this QR code with your mobile device:

REFORMER SAFETY USAGE

Safety and the Universal Reformer

The Pilates Universal Reformer, is designed to teach the exercise repertoire of Joseph H. Pilates. All utilization must take the unique functionality and design of the Reformer into account as well as the skill, expertise and education of the instructor. **Exercises outside of the intended use may cause damage to the equipment and may be dangerous to the client.** Furthermore, exercises improperly taught may result in loss of balance, coordination and control of the client on the moving carriage and may lead to injury. The Pilates Universal Reformer is only to be used under the supervision of a Pilates instructor specifically trained in the safe use of the Reformer.

Balance. Coordination. Control.

One of the unique features of the Reformer is the utilization of springs as resistance against the travel and movement of the Reformer carriage. In order for exercises to be performed safely, three important conditions must always be met.

- User must always maintain and be aware of balance when on the Reformer.
- User movements must coordinate with the speed, timing and variability of the spring resistance.
- ⚠ User movements must be controlled and organized. There should NEVER be slack in the ropes or bouncing of the springs.
- User should move within their ability, strength, and flexibility.

Control Front Advanced Exercise. Reformer movements are complex and client must always move with Balance, Control and Coordination.

Finding an Instructor

Working with a qualified Pilates instructor who is certified to train on the Pilates Reformer is an important component to learn how to use the Reformer safely and effectively. There are many Pilates education programs and it is important for users to be discerning. Look for instructors who have taken a certification with no less than 200 hours on the Reformer.

Spring Resistance

For many new to Pilates anticipation of the spring resistance needs to be trained and understood for safe utilization of the Reformer. In addition, understanding relative spring weights, how to change springs safely, how to check the springs for wear and how to move the spring bar are an integral part of spring safety specific to exercise transitions and set up. Ensuring that the springs are set up appropriately for the specific client, are set on the hooks to prevent detachment of the spring, and are understood in relationship to the movement of the carriage will decrease potential incidents and increase success for both client and instructor.

Check the springs regularly for wear and replace as needed. We recommend replacing springs every 2 years - 1 year for group fitness classes.

Check springs regularly for wear and replace as needed. We recommend replacing springs every 2 years - 1 year for group fitness classes.

Inexperienced Clients

Pilates exercises are integrated, complex motor learning exercises. Practice and experience are necessary to progress. Instructors, please be cautious with how quickly you progress your clients, if they can maintain balance, coordination and control then adding additional complexity is possible. Clients, please be patient with the learning process and understand the value of balance, coordination and control.

Exercise Set Up. Reformer Mounting and Dismounting ⚠️

Potential injuries may occur during the transitions between exercises as clients mount and dismount the Reformer. Proper instruction for exercise set up, mounting and dismounting as per the balance, control and coordination needed for a particular exercise, and for transitions between exercises, is primary to ensuring user success and to mitigate accidents or incidents of concern. For instructors, it is important, when in doubt, to refer to instructional manual and for clients to understand the importance of seeking out highly qualified instruction.

The Sitting Box and Foot Strap ⚠️

Mounting of the box for exercises seated, (Short Box Series), exercises supine, (Backstroke and Teaser), exercises prone, (Overhead Press and Swan), all require specific mounting and dismounting instructions. When the footstrap is required, please make sure it is properly secured to your Reformer (see individual Reformer setup instructions). The strap must maintain contact with the clients foot at all times. Please review instructional manual for more specifics regarding the sitting box and footstrap. Some general suggestions are included here.*

Short Box correct sitting position. When the box is placed horizontally on the carriage for Short Box exercises, the ends extend beyond the carriage platform. For safe mounting clients must always sit toward the center of the box.

Short Box incorrect sitting position. When the box is placed horizontally do not sit on the edges of the box, the box will tip.

Short Box correct strap position. Make Sure the footstrap is securely attached to the Reformer and the foot is in contact with the strap at all times.


Long box correct box and body placement. For advanced exercises, such as Long Box Teaser, it is important that the box is correctly placed on the Reformer and the body position is optimal to maintain balance, control and coordination throughout the movement.

* ⚠️ Sitting Boxes are designed for sitting, lying and standing. It is not designed for jumping, stomping or step aerobic type activity.

The Ropes

The Rope adjustments and locking mechanisms differ from Reformer to Reformer. Safe utilization of the ropes require they are installed correctly. For all Reformers except the A2, ensure rope is threaded through eye straps and cam cleats as pictured. As for the A2, the locking mechanism is underneath the carriage. Make sure the ropes are secure prior to movement.

Picture of rope cleats

Standing Exercises ⚠

Standing exercises are part of the Pilates repertoire and are functional and integral to a complete exercise experience. However, standing on a moving carriage is complex, potentially dangerous and must be handled with appropriate caution. When possible, instructors are encouraged to spot and provide assistance. All other utilization must be done with caution, after the clients abilities have been assessed, and when instructor and client have created a level of trust and communication to minimize any potential loss of balance, coordination and control. Use of stability props, such as gondola poles, ballet bars or personal assistance, are encouraged as client is learning and based on client ability.

Standing Side Splits with Gondola Pole or spotter

Working off the Back End ⚠


There are some lovely contemporary exercise opportunities from the back end (riser end) of the Reformer. However, extreme caution is required. Do not let go or lose of contact with the carriage when the carriage pulled to the back end and the springs are extended fully. This will result in high velocity impact on the front end (footbar end) of the Reformer which may result in potential damage to the frame, bumpers, spring mechanism and more. Sitting, kneeling and pulling while on the long box prone are all possible from the back end. Caution is required.

The Footbar ⚠

There are many different footbars and depending on which Reformer you have the footbar may or may not have a locking mechanism. Please, make sure to always check that the footbar is in a secure position and if there is a locking mechanism that it is locked and secure.

Set and Check

Whenever making any adjustments to the Reformer (ropes, footbar, risers) always set and check before using.

The Allegro 1 Footbar has several locking positions. Please consult your Reformer's manual as footbars differ between Reformer models and options.


High Kneeling Exercises ⚠️

Chest Expansion - The exercise Chest Expansion is a staple in the Pilates repertoire. It is considered an intermediate to advanced exercise because it requires upright balance on the moving carriage. Clients are facing the riser and the timing as a client moves and pulls the springs may create instability. When there is too little resistance, user may pull with too much force, the carriage will free roll and they may lose balance. Where there is too much resistance a client may find themselves losing balance and being pushed forward. **Please spot and use caution until clients are comfortable with the exercise.** ⚠️ **Standing facing risers with hands in straps is not recommended as there is a potential risk of falling and injury.**

Chest Expansion*

Kneeling Arms facing the footbar - This exercise can quickly become unsafe for clients. ⚠️ **We do not recommend this exercise kneeling or standing as there is a potential risk of falling and injury.** Also use caution when clients are in low kneeling as the position may be uncomfortable leading to clients adjusting and potentially losing balance.

⚠️ **Kneeling Arms Facing Footbar.** We do not recommend as balance, control and coordination may be a challenge to maintain – there is a potential risk of falling and injury.

Ballistic Movements ⚠️

Ballistic movements result in slack in the straps and a momentary floating of the carriage, an almost feeling of momentary weightlessness. When the straps become taught the resistance is reengaged. If not anticipated it may create stress in joints, a loss of control or balance for the user. **Under no circumstance should the carriage ever bang against the bumpers or against the back end.** Controlling movements is key to the success of Pilates conditioning and training.

Group Reformer

When teaching group, or taking a group Reformer class, it is important to ascertain the instructional quality of the teachers, the value that is placed on solid Pilates foundations and repertoire and the organizational philosophy on teaching Pilates. Maintaining Balance, Coordination and Control are key to successful and safe group Reformer class experiences.

Listening, the Client / Instructor relationship

In order to maximize the positive conditioning of the Pilates repertoire and Pilates inspired exercises, it is important for client and instructor to establish a good working relationship based on mutual trust and communication. If client indicates concern, caution or general discomfort, it is important to listen, assess and adapt as needed. Creating an environment for the client to share their experience on the Reformer will be key to making sure they are moving within safe ranges of motion, are moving safely and the exercise choices are effective and appropriate for their body.

Reformer Maintenance

Keeping your Reformer in good condition is key to keeping it working optimally for you and for your clients. The Safety Guide document included with your equipment has recommended schedules for inspection and replacement of wear parts. Balanced Body also offers an optional maintenance log that you may find helpful. Please follow recommended maintenance schedules and practices to optimize your Reformer's performance.

SAFETY & MAINTENANCE

SAFETY

To reduce the risk of injury or damage, read all instructions and the following important precautions before using the Contrology® Reformer with Tower.

- » It is the responsibility of the owner to ensure that all users of the Contrology® Reformer with Tower are adequately informed of all precautions.
- » Use the Contrology® Reformer with Tower only on a level surface.
- » Keep hands and feet away from all moving parts. When the Contrology® Reformer with Tower is not in use, leave at least two springs connected to the frame.
- » Keep children under the age of 12 and pets away from the Contrology® Reformer with Tower at all times.
- » If you feel pain, dizziness, or shortness of breath, stop exercising immediately.
- » Do not stand on the footbar.
- » Before beginning any exercise program, consult your physician.

MAINTENANCE

- » Depending on frequency of use, the springs should be safe to use for one to two years (3000 hours). Replace any spring that is kinked, bent or shows separations while at rest.
- » Remove hair and debris from the axles of the wheels. Debris can create pressure on the wheels, causing them to wear. If you hear noise from the bearings, replace them.
- » Repair any tears in the upholstery with an upholstery repair kit available at auto parts stores, or have the repair done locally.
- » Contact Balanced Body for other upholstery and re-upholstery options. Replacement upholstery is available from Balanced Body.
- » Check and tighten all screws and bolts. Replace any missing screws, retention pins, or other parts.
- » Check the footbar bolts on a monthly basis to make sure they are still secure.

CLEANING

Wipe the carriage pad, headrest, footbar and shoulder rests with a soft cloth and a mild, non-abrasive mixture of soap and water after each use. Keep the carriage track and wheels free from dust and dirt. Clean the frame with a mild, non-abrasive mixture of soap and water. Keep the leather straps and springs clear of dust. Cotton loops can be machine-washed. Hang to dry.

NOTE: DO NOT LUBRICATE OR GREASE ANY COMPONENTS ON THE FOOTBAR, CARRIAGE, WHEELS OR RAILS. By adding lubrication or grease it may not lock the components securely and may become unsafe. If you have any issues or questions please contact Balanced Body Technical Support.

Questions? Call Balanced Body Technical Support at 1-800-PILATES or +1-916-388-2838 or info@pilates.com.

Classical Pilates Studio Best Practices and Safety Guidelines

General Studio Safety

1. Keep pathways clear of equipment and accessories.
2. Inspect all apparatus and maintain equipment on regular weekly/monthly/quarterly basis. See manufacturers guidelines.
3. Reset, check and clean all equipment after use.

Contrology Reformer and Tower (Combo Unit)

1. Reformer Combo Unit Reset

- a. In Reformer Mode
 - i. The shoulder blocks are securely attached to the carriage.
 - ii. The headrest is up in the middle position and secure.
 - iii. The spring bar is in first gear.
 - iv. All 4 springs are attached.
 - v. The extension straps are draped over the back of the frame.
 - vi. The handles are stored on the shoulder block hooks.
 - vii. The box is placed behind the Reformer, or to the studio's specifications.
 - viii. Pole, pads, caster cover, and gear blocks are stored out of the walkways and to the studio's specifications.
 - ix. Cadillac springs are stored on the tower, or out of the walkways and to the studio's specifications.
 - x. PTB is hooked up to the center storage chain and stored to the outside of the tower.
- b. In Cadillac / Mat Mode
 - i. Shoulder blocks, extension straps, leather straps, and handles are stored neatly under the Reformer or out of the walkways.
 - ii. The caster cover is in place.
 - iii. Headrest is down.
 - iv. Carriage is at the back of the frame - all springs are unhooked.
 - v. The foot bar is lowered into the frame and the foot bar support is pointing towards the back of the Reformer.
 - vi. The insert mat is placed in the frame.
 - vii. The rollback bar is attached to the top hooks on the tower uprights.
 - viii. The leg springs are attached to the second-highest hooks.
 - ix. The arm springs are attached to the lowest hooks, or on the standing arm spring hooks on the outside of the tower.
 - x. The push-through bar (PTB) is stored to the outside of the tower using the center chain and clip.
2. Recommended monthly safety check
 - a. Springs
 - i. Inspect for wear, kinks, unevenness, separation, and rust. Immediately replace any parts that show signs of weakness. Refer to Figure A at the end of this document.
 - ii. Reformer - Rotate spring quarterly - move the inner 2 springs to the outside for even wear.
 - b. Clips
 - i. Inspect all clips for signs of wear - check the snaps, openings, and clips. Make sure the closure works correctly. Replace immediately if there's an issue. NOTE: Refer to Figure B at the end of this document.
 - c. Leather Straps
 - i. Inspect leather straps for cracks and signs of wear. Immediately replace any parts that show signs of weakness.
 - ii. Check carriage attachments and QuickLinks for wear and make sure the attachments are tight. Refer to Figure C at the end of this document.
 - iii. Tighten eyelets on Reformer frame if necessary.
 - d. Extension Straps, Foot Straps, and Fuzzies
 - i. Inspect stitching on strap ends.
 - ii. Check the stitching on all the straps and loops for wear and lose stitching and replace immediately.
 - iii. Check the fuzzies and the underlying loops for loose stitching and worn material. Replace immediately if signs of wear.
 - e. Handles
 - i. Check snaps and hooks for wear. Make sure all screws are tight.
 - f. Footbar
 - i. Check the side screws and tighten if necessary.
 - g. Carriage
 - i. Inspect headrest screws and attachments. Replace with any sign of wear.
 - h. Frame
 - i. Check the eyebolts on the Reformer Frame (foot strap attachments, lower PTB attachments, and all tower attachments) NOTE: If you can twist an eyebolt by hand, it is too loose.
 - i. Inspects the PTB pivot points.
 - j. Check the PTB safety chains and eyelets for signs of wear. Tighten or replace as necessary.

3. During a session

a. Reformer

- i. When the footbar is up, make sure the support is fully against the inside of the frame.
- ii. Straps and handles should be inside the frame or on the eyelets (aka eyebolts or hooks) at all times.
- iii. When working in 2nd or 3rd gear, the gear block must be used to ensure proper tension with the springs.
- iv. Make sure the headrest is fully seated and secured in whatever position you use.

b. Cadillac / Mat

- i. If you need to take the springs off during the session, store them out of the walkways. Make sure client is safely out of the way when changing springs.
- ii. Make sure the clips/snaps are seated properly. Make sure client is safely out of the way when checking clips or snaps.
- iii. Double-check the clips on the PTB and that the safety chains are connected properly.
- iv. It is critical that you keep control of the PTB whenever working with students. Always use extreme caution.
- v. Set and check the QuickLinks and wire stops on the springs. Make sure that they are fastened and seated correctly. Refer to Figure C and D at the end of this document.
- vi. Verify that all Cadillac springs with open hooks have wire stops attached.
- vii. Be mindful of the mat handles when walking around the apparatus.

For additional info on the "set and check" procedure for optimal safety, [check out this video](#).


Figure A

NOTE: This is not a Contrology spring but kink will look similar

Figure B

QuickLink not fastened. (first pic) QuickLink fastened. (second pic)

Figure C

Bad snap: excessive wear on hook. (not a Contrology snap but wear will look similar).

Figure D

Wire Stops

Contrology Reformer with Tower Exercises

THE ROLL UP, SKILL LEVEL: FOUNDATION

Reps: 3-5

Apparatus Set Up: Mat; Footbar down, Shoulder blocks removed.

Springs released. Insert mat and caster cover pad

Additional Equipment Recommended: Maple Pole

Precautions: Back issues, forward flexion issues

Starting Position

- » Lying on the mat
- » Feet flexed under the foot strap
- » Pole in hands - Arms overhead shoulder-width apart
- » The back is long, the center is engaged and the rib cage is flat

Focus:

- » Strengthen the center
- » Articulate the vertebrae
- » Stretch the back of the body
- » Warm-up the body

Movement Sequence

- » Inhale slowly and bring arms towards the ceiling
- » Bring the head up and continue to round through the spine
- » Exhale as you stretch forward, keeping the support of the round back as you go
- » Inhale and roll back
- » Exhale as the arms come to the ceiling and then overhead
- » Return to the starting position

Variation:

Beginners may do the roll back with knees bent and a smaller range of motion

Pole is recommended, but optional

1

2

3

THE SINGLE LEG CIRCLE, SKILL LEVEL: FOUNDATION

Reps: 5 circles in each direction

Apparatus Set-Up: Mat Footbar down, Shoulder blocks removed. Springs released. Insert mat and caster cover pad

Precautions: Hip issues

Starting Position

- » Lie on the mat
- » Arms are either - overhead, under the handles, or by the side palms down
- » One foot is under the strap - foot flexed
- » The other leg is long and lifted at 90 degrees
- » The back is long, the center is engaged and the rib cage is flat

Focus

- » Strengthen and stabilize the center
- » Massage the hip
- » Introduce dynamic movement of the lower body
- » Stretch the leg

Movement Sequence

- » With the hips anchored, cross the leg across the body
- » Sweep down towards the ankle and back up making a circle with the leg
- » Pause and repeat
- » After 5 repetitions reverse the circle, starting with the leg in the open position

Variation:

Beginners may soften the support leg

If the student is tight in their lower body they may start the circle at a lower angle and limit the range of motion.

1

2

3

4

THE ELEPHANT, SKILL LEVEL: FOUNDATION

Reps: 8-10

Apparatus Set-Up: Reformer

Footbar: Up Shoulder Block: In place. Springs: 2

Insert mat and caster cover pad: Removed

Additional Equipment Recommended: None

Precautions: Hip issues, Inversion Issues, Low Blood Pressure

Starting Position

- » Place the hands on the footbar thumbs and forefingers on the same side of the bar. (Caution - Do not pull up on the bar).
- » Carefully step onto the reformer and place the heels in front of the shoulder blocks. Allow them to settle into the carriage recesses
- » Bring the back into an even C-Curve or round back position
- » Stretch the feet long and lift the toes

Focus

- » Strengthen and stabilize the center
- » Work the hip
- » Stretch the back of the body

Movement Sequence

- » Stabilize the upper body into the footbar
- » Carefully press the carriage out about three inches from the hip without changing the spine position
- » Pull the center in and up, lengthen through the heels and return the carriage home

Variation:

- » More petite students may have to bring the heels forward on the carriage to obtain the correct position
- » Advanced variation - bring one leg to the back

1

2

3

STANDING CHEST EXPANSION,
SKILL LEVEL: FOUNDATION

Reps: 2-3 sets

Apparatus Set-Up: Cadillac

Footbar: N/A **Shoulder Block:** N/A **Springs:** N/A

Insert mat and caster cover pad: N/A

Additional Set-Up: Roll Back Bar attached to outside hooks

Precautions: None

Skills needed to perform exercise: Up Stretch Combo, Elephant with One Leg, Swan—Mat, Pull-Up—Wunda Chair, holding starting position on 1 spring

Starting Position

- » Stand facing the tower a little more than an arms-length away
- » Place the hands on the rollback bar, thumbs on the opposite side of the fingers

Focus

- » Opens the chest
- » Engages the back
- » Stretches the neck
- » Assisted with the breath

Movement Sequence


- » Inhale and bring the bar to the thighs and lift the upper rib cage
- » Hold the breath - turn the head to the right, turn the head to the left
- » Exhale as you bring the head center and return the bar to the first position
- » On the next repetition, turn the head to the left then to the right

Variation:


None.


1

2

3

4

THE TEASER PREP WITH THE PUSH THROUGH BAR,
SKILL LEVEL: INTERMEDIATE

Reps 2 - 3 Sets

Apparatus Set-Up: Cadillac

Foot bar: Down Shoulder Block: Removed Springs: Disengaged
Insert mat and caster cover pad: Inserted

Additional Set-Up: Push Through Bar attached to the upper spring

Precautions: Shoulder issues, neck issues

Skills needed to perform the basic exercise: Assisted Teaser on the Mat, The Teaser on the Mat

Starting Position

- » Lay on the Cadillac bed with the head under the push-through bar
- » Place the hand on the push-through bar, shoulder-width apart
- » Bring the bar down and overhead to test position. It should just barely clear your nose

Focus

- » Opens the chest
- » Stretches the shoulder
- » Works the center

Movement Sequence

- » Bring the push-through bar towards the chest by bringing the arms wide towards the table
- » Rotate the arms so the palms face the ceiling and the chest stretches
- » Reach the arms overhead
- » Bend the elbows out wide to bring the bar back to the first position
- » Repeat 3 times
- » On the third repetition pull the center in and roll the torso up into the half teaser position
- » Roll down
- » Bring the bar to the chest and then the arms overhead
- » Return the bar to the home position

Variation

- » As you bring the torso up, bring the legs up into a full teaser
- » In the full teaser position - shave the head
- » Lower and lift the legs
- » Shave the head and lower and lift the legs
- » Alternate shaving and lowering and lifting

1

2

3

4

THE PUSH THROUGH, SKILL LEVEL: FOUNDATION

Reps: 3

Apparatus Set-Up: Cadillac

Footbar: Down

Shoulder Block: Removed

Springs: Disengaged

Insert mat and caster cover pad: Inserted

Additional Set-Up: Push Through Bar attached to the upper spring

Precautions: Shoulder issues, neck issues

Skills needed to perform the basic exercise: The roll-up, the roll back

1

Starting Position

- » Sit facing the push-through bar with a tall back, legs long, and feet flexed against the pole
- » Place the hands on the bar, shoulder width apart

Focus

- » Stretches the shoulders and chest
- » Opens the back and the hamstrings
- » Works the center
- » Develops coordination

Movement Sequence

- » Keep the arms straight and roll back towards the mat using your center
- » Press the bar down with straight arms, lift the rib cage up and over, and stretch forward toward the poles
- » Continue to reach into the bar, pull the center in, and roll back towards the mat
- » Lift the arms up and lengthen the back to return to the first position.

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
EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

Balanced Body, Inc.
5909 88th St,
Sacramento, CA 95828 USA

Contact Us
1-800-PILATES | +1-916-388-2838
pilates.com | info@pilates.com
17639D | 2025-10

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