The Balanced Body® Allegro Nextgen Reformer™ and Tower User Guide

THIS IS A MEDICAL DEVICE

STATEMENT OF INTENDED USE

The Allegro Nextgen Reformer™ device is intended by the manufacturer to be used by people for the purposes of physiotherapy, rehabilitation, the prevention, treatment and relief of injuries, supportive, functional therapy for a disability or for general strengthening and improving the mobility of the human body.

Important: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Instructions in document subject to change.
Please visit www.pilates.com/support/assembly-use-maintenance for current version.

Visit pilates.com/patent for complete and current information on Balanced Body product patents.

IMPORTANT PRECAUTIONS: PLEASE READ WARNING

To reduce the risk of serious injury, read the following important precautions before using the Allegro Nextgen.

- » Read all instructions in this manual before using the Allegro Nextgen. Review the Setup and Safety video included with the Allegro Nextgen before using the equipment. Use the Allegro Nextgen Reformer only as described in these instructions and the video.
- » It is the responsibility of the owner to ensure that all users of the Allegro Nextgen are adequately informed of all precautions.
- » Use the Allegro Nextgen only on a level surface. Keep hands and feet away from all moving parts. When the Allegro Nextgen is not in use, leave at least two springs connected to the carriage. Keep children under the age of 12 and pets away from the Allegro Nextgen at all times.
- » If you feel pain, dizziness, or shortness of breath, stop exercising immediately.
- » Before beginning any exercise program, consult your physician.

EMERGO EUROPE
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

For printed instructions contact Balanced Body.

Balanced Body, Inc.
5909 88th St,
Sacramento, CA 95828 USA

Contact Us
1-800-PILATES | +1-916-388-2838
pilates.com | info@pilates.com

19093B | 2025-10

How to Assemble the Allegro Nextgen Reformer™

PARTS LIST

Part number	Description	Qty
19171	Complete footbar	1
12680	Shoulder bolts and washers	1 kit
12907	Shoulder rest with brackets	1 pair
11765	Locking shoulder post	2
18903	Twist lock plunger knobs	2
10536	Soft Touch Ropes	1 pair
10046	Double loops (Soft Touch)	1 pair
19080	Allegro Nextgen riser with pulley	2
13629	5/16" Allen wrench	1
12981	Allegro Nextgen wrench	1
14285	Red spring	3
14286	Blue spring	1
14295	Yellow spring	1
15968	Foam Pad for Standing Platform	1
13778	5/32" Allen wrench	1
13891	1/4"-20 Button head screw	2
13633	Split washer	2
11890	Transport bracket preinstalled with wheel	2

INSTALL THE SHOULDER RESTS

1. Note the two protruding studs on the bottom of the shoulder rests. Lower the shoulder rest studs into the keyholes and slide away from the carriage pad. See Figure A.

Note: The shoulder rests can be installed in two positions. One position is wider, and more comfortable for broad shoulders. Swap the left and right shoulder rests to change between standard and wide configurations.


Figure A

2. Now screw the locking shoulder posts onto the exposed bolt on each shoulder rest. See Figure B.


Figure B

Using the Locking Shoulder Posts:

Turn the posts clockwise to tighten them and lock the shoulder rests into position. Over tightening the posts will make them difficult to remove, they should be snug but not too tight. To remove the shoulder rests, turn the post counter clockwise 1-2 turns and gently slide the shoulder rests toward the carriage.

Figure C shows the storage position for the shoulder rests. Note the key hole slots at the head end of the Allegro Nextgen frame.


Figure C

INSTALL THE BLACK PLUNGER KNOBS

3. Look through the large threaded hole in the trunnion and line it up with one of the holes in the frame. Screw the plunger knob into the large hole until it is completely threaded in. See Figure D.
4. Repeat on other side. Tighten both plunger knobs securely with the included open end wrench.


Figure D

To move the trunnions, pull the ball of the plunger knobs away from the frame and turn 90° to lock out. To re-engage the plunger turn the knob until it snaps into one of the holes in the side of frame.

Note: Make sure both trunnion plates are locked in the same horizontal position before moving on.

INSTALL THE INFINITY FOOTBAR

5. Pick up the footbar so the seam of the footbar cover is away from the carriage pad.
6. While squeezing the footbar plunger levers, guide the footbar over the outsides of the trunnion plates. Align each footbar pin with the top hole in the trunnion plate. While continuing to support the footbar, release the footbar plunger levers to engage the trunnion plate. Be sure the pins on both sides are engaged. See Figure E.


Figure E

7. Pivot the footbar slightly until the large hole at the bottom of the footbar aligns with the remaining hole in the trunnion plate.
8. Slide the shoulder bolt with one washer through the large footbar hole. Then screw the shoulder bolt into the trunnion plate as tightly as possible with the large Allen 13629 wrench. Repeat on the other side. See Figure E.

INSTALL THE RISERS

9. Remove locking pins from the frame. Insert the risers through the holes on the top of the frame. Turn the risers until the attached pulleys point straight down the frame toward the footbar. See Figure F.


Figure F

10. Reinsert locking pins until the metal ring touches the frame. You may need to jiggle or slightly lift or rotate the riser to get the locking pins in all the way.

ATTACH THE ROPES

11. To attach the ropes, first engage at least one spring from the carriage to the frame to keep the carriage in the home position. Unroll the ropes and separate them.
12. The Allegro Nextgen now comes with Soft Touch Ropes that do not have the dog-clips to attach the loops. The clip will no longer drag on or bang into the frame! Follow the images below to install your loops. See Figure G.


Figure G

13. Rest the loops over the shoulder rests. Thread the other end of each rope through a riser pulley and back into the cam cleats on the carriage to adjust the length. Be sure to go through the black plastic eyestraps on both sides of the cam cleats. Always push the rope firmly down into the cleats to ensure a good grip.

INSTALL THE TRANSPORT WHEEL BRACKET

14. Remove the 1/4"-20 button head screw and the split washer attached to the last hole in the frame towards the footend of the Reformer, as shown in Figure H, using the 5/32" Allen key. Slide the preassembled transport bracket into the channel in the rail at the foot end of the reformer. Line up the hole in the transport bracket with the hole in the frame rail. Use the 1/4"-20 button head screw and the split washer to attach the transport bracket on to the frame using the 5/32" Allen key. See Figure H for reference. Note: The wheel goes on the inside of the bracket (towards the frame) when it is mounted.


Figure H

15. Repeat this for the other side of the reformer.

STORING THE ALLEGRO NEXTGEN REFORMER

Stacking

Attach the springs to keep the carriage from moving. Move the footbar to the "down" position. Remove the shoulder rests and store in the slots at the head of the frame. Put the head rest down. Remove the risers and store in the holes on the sides of the frame. Place the ropes and loops inside the Allegro Nextgen frame not on the carriage. Alternate the orientation of the Allegro Nextgens as you stack them. Place the head end feet onto the standing platform of the Allegro Nextgen below, and so on. Stack a maximum of 5 reformers with no legs (3 with legs). Do not stack Reformer with Towers.

Standing

It is absolutely critical that the Allegro Nextgen footbar be set and locked in the correct position for standing storage. Move and lock the footbar into the position closest to the foot end of the Reformer frame. Raise the footbar to its 2nd-highest position. Lift the head end of the Allegro Nextgen frame to waist height. Then raise the frame to vertical. Use your foot to steady the wheeled end of the frame as you lift to vertical. The vertical Allegro Nextgen will come to rest on the footbar and the transport wheels. See Figure I.

Figure I

CLEANING & MAINTENANCE

Cleaning

Wipe the carriage pad, headrest, footbar and shoulder rests with a soft cloth and a mild, non-abrasive cleaner after each use. Keep the carriage track and wheels clean from dust and dirt. Wipe the entire track with a towel regularly. Clean the footbar trunnion track weekly. Clean the frame with a mild, non-abrasive cleaner. Keep the ropes and springs clear of dust. Cotton loops can be machine-washed. Hang to dry.

Lubrication

The foot bar trunnion C-channel can be lubricated with dry silicone spray.

MAINTENANCE

» See the included Safety First Guide.

QUESTIONS?

Call Balanced Body technical support at 1-800-745-2838 (US and Canada), or +1-916-388-2838.

How to Install the Legs Kit

PARTS INCLUDED:

Item No.	Part Number	Description	Qty
1	18833	Allegro Nextgen Leg, A	2
2	18834	Allegro Nextgen Leg, B	2
3	18863	Leg Nut Plate	2
4	18989	Wheel Bracket, Allegro Nextgen Ext Leg	2
5	13569	Allen Key, 3/16" Ball Point	1
6	13613	Screw, Button Head, 5/16-18 x 3/4"	12
7	15221	Screw, 1/4-20 x 3/4"	4
8	13613	Nut, Nylock, 1/4-20	4
9	13778	Allen Key 5/32" Short Arm	1
10	18864	Polyurethane Foot Pad, A	2
10	18865	Polyurethane Foot Pad, B	2
11	18866	Double Sided Tape	4
12	13521	Alcohol Wipe	4
13	19076	1/4" Speed Washer	4
14	13724	Washer, Lock, 5/16"	12
15	19017	Transport wheel bracket	2

TOOLS NEEDED (NOT INCLUDED):

- » Two 1/2" or adjustable wrenches
- » One 7/16" wrench
- » 3/8" wrench or socket

PREPARE THE NEW PADS

1. Remove the legs and other materials from the packaging and lay the legs on the floor with the bottom of the leg facing up, so the "L" shape face with three holes are on the floor.
2. Remove the (4) Polyurethane foot pads, and (1) Alcohol wipe from the packaging. Open one wipe and thoroughly wipe down the top side of each foot pad. This is the side opposite the ribbing. If the wipe is dirty or becomes too dry use another wipe for the remaining pads. Set the pads aside and let these dry.
3. Open an Alcohol wipe to clean the bottom of each metal foot on the reformer. If needed use a secondary wipe.
4. If you have a Tower Of Power installed, remove the ropes, loosen the 4 tower knobs, and remove the Tower.
5. Using the included shorter allen wrench, remove the 4 screws holding the head end cover plates and riser pins in place. Set them aside to be re-installed later. See Figure A.



Figure A

1. Attach all springs to keep the carriage stationary. Remove ropes, shoulder rests, and risers. Place the footbar (if installed) in its lowest position and check that all pins are locked so that the footbar is secure. Then, with the help of a friend, lift the Allegro Nextgen from both ends and carefully turn it upside down and set it on the floor. You may want to place a pad or other floor protection beneath the Reformer.

2. Remove the 4 rubber feet from the bottom of the Reformer frame by grasping them firmly with your hand and turning counter-clockwise. You can discard the feet unless you may want to convert back to a flat Allegro Nextgen configuration.

3. Before proceeding ensure the tops of the new pads from step 2 are completely dry and no more alcohol liquid remains on the part. If the pads are not dry, the pad and tape will not bond. Peel one side of tape off the sticker and install it on a pad. Align the 45 degree corner side of the tape with the angled corner of the foot pad.


Figure B

The center hole will go over the center circle, and the flat side will go against the raised alignment pins. See Figure B.

4. Peel off the top side of the applied tape and install this pad assembly on the bottom of a cleaned leg. See Figure C.


Figure C

5. Use one of each: 1/4" screw, speed washer, flat washer, lock nut and the provided Allen wrench, plus a socket or wrench; Align the angled side of the foot with the angled side of the leg, making sure the alignment pins of the foot go into the base of the leg. Again refer to Figure C. Tighten the screw until the end just barely sticks out the top of the nut. **DO NOT OVER TIGHTEN.**

6. Repeat steps 3-5 for all other legs.

7. Insert the leg nut plates shown in Figure D into the head end of the Reformer frame until the big hole lines up with the riser hole in the frame. When it does, drop your riser through the frame and the block to hold it in place. See Figure D and E.


Figure D


Figure E


Figure F

8. With the risers still in place to locate the leg nut plates, place two of the legs onto the frame rail at the head end of the Reformer, lining up the holes in the leg with the three holes in the frame. The left and right legs are not the same so keep looking for one that matches the hole pattern in the frame. Now start threading all six of the 5/16-18 screws with lock washers into the frame to secure the legs in place. Start with the screw that goes into the block that is held in place by the riser and the frame, and once they are all started, tighten them. See Figure F. Now you can remove the risers and set them aside.

NOTE: the long end of the included 3/16" allen wrench is ball shaped to allow for more maneuverability when working in small areas. Once the screws are finger tight, use the short leg of the allen wrench to tighten the screws.

1. Install the legs onto the foot end of the frame in the same way with the 6 screws and lock washers. No leg nut plates are required at the foot end.
2. Using a 7/16" wrench and the included allen wrench, install the transport wheel brackets with the 1/4-20 screws and nuts as shown in Figure G. Make sure the brackets are oriented as shown. Typically these are installed at the foot end of the Reformer as shown in Figure G.


Figure G (Footbar not shown)

NOTE: The transport wheels can also be installed at the head end which allows you to move the Reformer with Tower more easily. Just make sure the carriage is rolled to the head end of the Reformer before lifting the foot end of the machine in this configuration.

3. Using two 1/2" or adjustable wrenches, remove the transport wheels from their current location and install them onto the brackets you just installed. See Figure H.


Figure H

NOTE: It is critical that the washers and wheel get installed in the same order and orientation as they were on the original brackets. Do one wheel at a time and if the wheel does not spin freely after you are done, double check against the other wheel to make sure everything is in the right order.

4. Using Allen Wrench remove bolts holding footstrap brackets in frame and slide out old brackets. Install new brackets and reinstall bolts. See Figure I


Figure I

5. Re-install the head end rail cover plates on the same side you removed them from.

NOTE: If you will be installing a Tower System next, leave the head end rail cover plates off.

6. Your legs are installed and you may now carefully lift and turn your Allegro Nextgen Reformer back over onto the legs. Re-install your shoulder rests, risers, and ropes.

How to Install and Use the Platform Pad

Your Allegro Nextgen Reformer comes with a platform pad. Place the pad directly on top of the platform to use. This pad will make the standing platform height the same as your carriage pad. See Figure A.


Figure A

Assembling the Tower

Prior to assembly, check the parts and components received against the following list. If you are missing any parts, contact the Balanced Body customer service department at 1-800-PILATES (1-800-745-2837) in the U.S. and Canada, or +1-916-388-2838 (international).


PREPARE YOUR ALLEGRO NEXTGEN REFORMER

If you have purchased an Allegro Nextgen 14" leg set, we strongly recommend installing the legs before installing the Tower.

Remove all ropes and risers from the Allegro Nextgen.

PARTS LIST

Description	Part Number	Qty
Tower	18990	1
Push-through bar (PTB)	14594	1
T-pin	14021	2
Tower Bracket, Right	11907	1
Tower Bracket, Left	11906	1
Adjustable pulley and knob assembly	19173	2
5/32" Allen key, long	10551	1
Blue trap spring	14275	2
Yellow trap spring	14272	2
Long spring, purple	14306	2
Long spring, yellow	14279	2
Roll down bar	12228	1
Single cotton loops	10022	1 pair
Safety strap with carabiner	10499	1
Tower lower spring mounts	12915	1

1. Using the provided 5/32" Allen wrench, remove the plastic rail caps and locking pins at the head end of the Allegro Nextgen frame. They will need to go back on the same side, so keep track of which one goes where and keep the screws with them. See Figure A.


Figure A

2. Installing the Tower Mount Brackets:

- a. With the provided 5/32" Allen wrench, loosen the 6 screws on each bracket 1-2 turns. Do not remove them, you just want some free play between the plates. See Figure B.
- b. Make sure the brackets are positioned as in Figure B (use the large grey knobs and screws as reference) and slide them into the C-channel of the frame.


Figure B

3. Re-install the plastic rail end caps on the same side you removed them from. The lanyard and locking pins are no longer needed so they can be set aside with the risers. See Figure C.


Figure C

4. Make sure both brackets are slid all the way up against the plastic rail caps. Using the provided 5/32" allen wrench, tighten the 6 screws on both Tower brackets. See Figures C and D.


Figure D

5. Install the Tower frame onto the brackets.

- a. Position the Tower frame so that the large threaded holes at the bottom face away from the carriage.
- b. Lift the Tower frame over the brackets and lower it around them so it rests on the ledge at the bottom of the brackets. See Figure E.


Figure E

- c. Slide the frame toward the bracket until it is flush.
- d. With one arm or an assistant holding the Tower frame upright and in position, start threading the 4 large knobs into the Tower frame. You may have to push on the frame to get it into position with the knob, and if one is giving you trouble, just move onto the others and come back to it. See Figure F.

IMPORTANT NOTE: Get all four knobs started 3-4 turns into the Tower frame before tightening any of them. It is also critical that the Tower is parallel and flush against the bracket with the knobs. You may have to push/pull/twist the leg of the Tower to make it line up.

- e. Once all four knobs are started, tighten them one by one.


Figure F

6. Install the adjustable pulleys by assembling the components as shown in Figure G with the pulley on the inside of the frame. Tighten the knob when the eyebolt is horizontal as shown in Figure G. Repeat on the other side.


Figure G

7. Install the lower spring mounts in the head end of the frame with the long screw and the 5/32" allen wrench. Tighten the screw with the eyebolt positioned vertically as shown in Figure H.


Figure H

8. Rethread your ropes through the new pulleys.

REPOSITION THE T-PIN

For shipping purposes, the Tower is packaged with the T-pins installed on the inside of the Tower frame. See Figure I.

9. The T-pin should be installed on the outside of the Tower frame for use. To remove the T-pin, push the button on the handle of the pin while holding the push through bar (PTB) securely. Pull the pin free of the PTB and Tower frame. Reinsert the pin from the outside of the Tower frame. See Figure I.


Figure I

INSTALLING MAT CONVERSION

10. Disconnect all carriage springs
11. Move the carriage to the head end of the frame and clear the ropes out of the way.
12. Place the mat in the open space between the standing platform and the carriage.
13. Remove the shoulder rests from the carriage (optional).

USING THE PRE-INSTALLED PUSH-THROUGH BAR (PTB) SAFETY STRAP

Note: Always use the safety strap for bottom-sprung exercises.

Use only with qualified instructor supervision. The PTB safety strap is used only for bottom-sprung exercises. Figure A at right shows the safety strap correctly positioned to secure the push through bar. The strap includes numerous loops or "contact points" so that you can adjust the "stopping point" of the PTB if it is released during exercise.

Figure A

Note: Use the Velcro portion of the safety strap at the top of the Tower frame to secure the PTB upright when not in use.

MOVING AND STANDING AN ALLEGRO NEXTGEN WITH TOWER

Note: Always exercise caution when moving, lifting, or storing an Allegro Nextgen or an Allegro Nextgen with Tower. Use proper lifting technique to avoid injury.

TO MOVE:

To move a Reformer with Tower a short distance you may be able to lift the head end and roll it on the transport wheels as if you did not have a Tower. Be aware that you will be lifting the Tower and Reformer so it will be heavier than the Reformer by itself. The Tower will also be rising quite high so you will have to be careful if you do not have high ceilings.

If moving the Reformer and Tower is not possible, remove the Tower, then move the Reformer.

1. Pull the ropes back through the pulleys on the riser.
2. Loosen the four gray Tower knobs and remove the Tower from the Reformer frame. If this is too heavy, the springs and push thru bar can be removed first. If you do not remove the push thru bar or roll down bar, be aware that they can swing freely.
3. Now with the carriage held in the home position by at least 2 springs, lift the head end of the Reformer to roll.

Note: If you have the extended legs installed and have the transport wheels at the head end, you can leave the Tower installed and lift the foot end of the machine to roll it around. Make sure the carriage is at the head end of the frame before lifting. You will take up more room as you move around in this configuration but for some circumstances this is acceptable.

TO STAND:

4. Remove the Tower frame as explained in the "To Move" section. **The Tower must be removed before standing.**
5. It is absolutely critical that the Allegro Nextgen footbar be set and locked in the correct position for standing storage. Move and lock the footbar into the position closest to the foot-end of the Reformer frame. Raise the footbar to its 2nd highest position. Lift the head-end of the Allegro Nextgen frame to waist height. Then raise the frame to vertical. Use your foot to steady the wheeled end of the frame as you lift to vertical. The vertical Allegro Nextgen will come to rest on the footbar and the transport wheels.

Note: It is absolutely critical that the Allegro Nextgen footbar be set and locked in the correct position for standing storage.

Installing Your Wall Security Strap

The Allegro Nextgen Wall Security Strap is designed to mount on walls with standard wood framing. If you have a different type of wall construction (brick, concrete, etc.), it may require a different fastener. Check with your local hardware store or contractor.

Parts List (included)

- » 1 adjustable belt strap
- » 1 3-1/2" eye-bolt

Tools Needed (not included)


- Drill with 1/8" bit
- » Phillips or flathead screwdriver

INSTALL THE EYE-BOLT

1. Find a stud in the wall where you would like to stand your Allegro Nextgen. Make sure you have ample space on both sides of the stud as well as above for the Allegro Nextgen. Using a 1/8" drill bit, pre-drill a hole into the stud. The hole should be positioned 5' to 5 1/2' above the floor. Secure the eye-bolt into the pre-drilled hole with the screwdriver until the shaft of the eyebolt is completely within the wall.

USING THE WALL SECURITY STRAP

2. Follow your Allegro Nextgen or Allegro Nextgen Tower instructions for proper standing procedure.
3. Move the Allegro Nextgen near the wall with the footbar pointed away from the wall (perpendicular to the wall).
4. Wrap the strap ends around one rail and connect the ends like a belt, using the adjustable clip.
5. Adjust the strap to remove any slack.

NOTE: Do not try to wrap the strap ends around the entire frame - it will not fit.

IMPORTANT: Balanced Body is not responsible for any damage or injury caused by improper wall mount installation, or the use of incorrect or improper wall mounting equipment.

If you have any questions, please call Balanced Body at 1-800-PILATES (1-800-745-2837) in the U.S. and Canada, or +1-916-388-2838 (international).

Safety First: A guide to proper maintenance and safe use of your Pilates equipment.

For over 45 years, Balanced Body has been introducing safety-related innovations to Pilates equipment. Many of our improvements are now industry standards, resulting in Pilates equipment that's safer today than ever before.

Safety depends on proper maintenance and safe use, in addition to the quality of the equipment. This guide was created to help you use and maintain your equipment for optimum safety. Please read it through carefully and keep for future reference. If you have any questions, give us a call. **Failure to follow these instructions may result in serious injury.**

ALL EQUIPMENT

Springs

Spring inspections are critical to maintain your equipment in safe operating condition. All Balanced Body springs should be replaced at least every two years. Certain environments and usages can shorten the expected life of the springs and you may need to replace the springs more frequently. Therefore, it is very important to inspect springs on a regular basis since worn or old springs lose resilience and may break during use. Injury may result if a spring breaks during use.

During use, do not allow springs to recoil in an uncontrolled manner. This will damage the spring and shorten its expected life.

Inspect springs for gaps and kinks (weekly or monthly, depending on frequency of use). Look for gaps and kinks between the coils when the spring is at rest. It is not unusual for the spring to have a very small gap on the tapered end (a gap is sometimes created during the manufacturing process). However, there should be no gaps in the body of the spring. If you see any gaps or kinks in the body of the spring, discontinue use and replace the springs immediately. See **Figure 1**. Additionally, corrosion anywhere on the coils will shorten the life of the spring. Discontinue using the spring immediately if you see any rust or oxidation during inspection.


Figure 1

Snaps

Inspect snaps for wear (monthly). First, verify that the snap hook is working properly. If the snap hook does not retract and return properly, discontinue using the spring immediately and replace the snap. Eyebolts can cause excessive wear on snap hooks. If the hook shows a lot of wear, discontinue using the spring immediately and call Balanced Body to replace spring or snap. See **Figure 2**.

Figure 2: Good snap: no wear on hook.

Bad snap: excessive wear on hook.

EYEBOLTS, NUTS AND BOLTS

Tighten all equipment bolts and screws (monthly). Verify that all eyebolts, nuts and bolts are tight. See the section titled "How to inspect and tighten nuts and bolts."

ROPES AND STRAPS

Rope and strap wear (quarterly). Ropes should be replaced if you can see the core of the rope through the outer lining, or if the ropes are flattened. Straps should be replaced as soon as any fraying is noticed. Be sure to check the sections of rope or straps that attach to the clips and run through the pulleys.

REFORMERS

Check springbar hooks or eyebolts (quarterly). Balanced Body makes two different springbar systems:

- » Revo Springbar. Make sure springbar hooks and handle are tight.

Standard Springbar. Verify that the nuts securing the springbar hooks are tight. See section titled "How to inspect and tighten nuts and bolts."

Spring rotation (quarterly). You can prolong Reformer spring life by rotating springs of the same weight each quarter. Unhook and move to another position on the springbar. Rotating springs helps them wear more evenly.

Risers on the outside. Wood risers must be installed on the outside of the frame. Risers can loosen over time, so always make sure they are tight.

Springs hooked downward under carriage. Make sure springs are hooked in a downward position. See **Figure 3**.

Secure the carriage. When your Reformer is not in use, be sure that at least two springs secure the carriage to the springbar.

Figure 3: Springs hooked downward

Default settings. Many users have a "default setting" for Reformers. At the end of a session, the user connects a prescribed number of springs in neutral tension, sets the footbar at a pre-determined height, and sets the ropes at a specified length. This ensures that the equipment is ready for the next use, and that the carriage is secured by the springs.

Footstrap under tension in box work. When using the box and footstrap, be sure the footstrap is under tension (with snaps pulling from the top of the eyebolt) before beginning the exercise. See **Figure 4**.

Figure 4: Foot strap under tension

REFORMER WHEEL AND TRACK MAINTENANCE

Clean the tracks and wheels (weekly). For smooth carriage travel and to maintain the longevity of the wheels, we recommend that you wipe down the tracks once a week.

Disconnect the springs and clean the entire length of the tracks with a soft cloth and Balanced Body Cleaner, mild soap with water or a mild commercial cleaner such as, Fantastik® or 409®. Do not use abrasive cleansers or pads, as they can damage the anodizing on the rails. To clean the wheels, hold the cloth against the wheels while you move the carriage. If you feel a bump in the ride, dirt has adhered to the surface of the rails or wheels. Clean hair and debris out of the rails. Hair can wrap around the wheel axles and eventually build up and cause wheel failure. Use tweezers to remove hair from the wheels.

Lubrication. Never spray silicone near or inside the wheels – this can wash the lubricant out of the bearings and ruin the bearings. You can purchase dry silicone at most hardware and auto parts stores. Pulleys sometimes require lubrication to stop a squeak. Direct a very quick spray of dry silicone or Teflon spray into the pulley axle. "Dry" silicone does not have an oil base. Oil-based ("wet") silicone and WD40 should not be used as they attract dirt. Be careful not to over spray. You may want to remove ropes to avoid getting silicone on them.

Do not lubricate the Allegro 2 rails.

Footbar supports (quarterly). For all Balanced Body footbars with footbar support brackets, verify that the pivot screw attaching the footbar support bracket to footbar is tight, but not so tight that it prevents the support from rotating freely. For Legacy Reformers, tighten the pivot bolt to secure footbar support.

Headrest (monthly). Make sure the hinge screws and bolts on your headrest are tight.

Under the Reformer (monthly). Move Reformers and make sure you clean the floor space underneath.

Standing Platform Footbar Bumpers (wood Reformers only). If your standing platform footbar bumpers (the small plastic pieces that protect the standing platform from the footbar) are broken or damaged, please call Balanced Body to replace.

TRAPEZE TABLE (CADILLAC) & TOWERS

Cotter pins removed. These pins are located in the vertical tubes that align the canopy to the frame and should be removed as soon as installation is complete. Unremoved cotter pins can tear clothing and lacerate the skin. Use pliers to remove the pins.

Figure 5: Cotter pin before and after removal from Trap Table.

Save the pins in case you need to disassemble and reassemble the table for transportation purposes. See **Figure 5**.

Push-Through Bar (PTB) with Sliders.

The PTB moves vertically to accommodate different users and exercises. The sliders on the tubes allow for this vertical movement. Make sure these sliders are clean and easy to move. Before beginning any exercise, ensure that the sliders are properly aligned with the PTB holes and locked into position. Apply a downward force to ensure. If you notice wear on the slider knob pins, please call Balanced Body to replace. Using the PTB in this condition could cause injury.

Weekly maintenance for Push-Through Bars with Sliders. Check to ensure the pins in the plunger knobs protrude and lock into the vertical tubes. Pull both knobs out and move the sliders to a different position. Release the knobs before the next hole and continue to slide the system into position. Once over a hole the pin of the knob will automatically drop into the opening. Once each slider is in a new position apply a downward force on the PTB. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced. Call Balanced Body for replacement parts.

Push-Through Bar (PTB) with T-pin setting. For bottom sprung exercises, if your client's head is below the PTB, use the T-pin setting in addition to the safety strap or chain. Spotting your client is highly recommended. This is important for safety.

Push-Through Bar (PTB) control. Make sure you have enough room around the trap table to safely use the PTB without fear of hitting other people. The PTB can be dangerous if not properly

used. Only trained, experienced users should use the PTB. A spotter should always maintain control of the bar with one hand. If the user should lose control of the bar, the spotter can maintain control of it.

Correct safety strap attachment. For bottom-sprung exercises, the safety strap or chain should always secure the bar.

The safety strap or chain should wrap around the PTB and the canopy frame, not the eyebolts. The strap or chain is only as strong as the weakest link, and the frame and bar are a great deal stronger than eyebolts. **Figure 6.**

Figure 6: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

Spotting your client is highly recommended. This is important for safety.

Setting the PTB for bottom-sprung exercises. For bottom-sprung exercises, the safety strap should be attached so that the angle of the push-through bar is no lower than the 4 o'clock position. This limits the range of the bar and prevents it from potentially coming into contact with the user.

Using the 4th side on the PTB along with the safety strap is highly recommended to prevent injury.

Figure 7: Safety strap holding the push-through bar at 4 o'clock. The strap is secured to the PTB and canopy frame, not the eyebolts.

CHAIRS

Dismount with control. When dismounting the chair, release the pedals slowly, with control. Don't let the pedal snap back.

Spot users. When a user is standing, sitting or lying on top of the chair, there is increased risk of falling. Standing exercises, in particular, can be unstable. Spotting users will make these exercises safer.

Hourglass spring mounts. If your chair has hourglass spring mounts and the mounts do not successfully retain the springs, replace the fiber washers (they are reddish-brown in color).

Figure 8: Fiber washer

Figure 8. If your chair is a Balanced Body Split-step Pedal Chair (Combo Chair), please consider upgrading to the Cactus Springtree).

UPHOLSTERY CLEANING & MAINTENANCE

Cleaning. You can extend the life of upholstery by keeping it clean and free of dirt, oil and perspiration. After each use, wipe down the upholstery with a solution of mild soap and water. Then wipe it down with clean water and dry with a soft towel.

Disinfecting. Equipment upholstery is coated with BeautyGard®, which offers antibacterial protection. If you want additional disinfection, Balanced Body offers Balanced Body Clean™ disinfecting solution. Use of any other solution (especially those containing essential oils) will shorten the life of some equipment and is not recommended.

For over molded foam padding it is recommended to clean with mild soap or BB Clean and wipe clean with a water dampened cloth to remove any leftover cleaner.

Cleaning Loops. Loops should be washed using mild soap such as Woolite. Hang to dry.

HOW TO INSPECT AND TIGHTEN NUTS AND BOLTS.

Use your fingers to check nuts and bolts for tightness. If you can turn the nut or bolt with your fingers, it's too loose and should be tightened. To tighten, first

tighten using your fingers. Rotate nuts and bolts clockwise to tighten. Insert a screwdriver through eyebolts to hold them steady while you tighten the nuts. Then use a small wrench to tighten the nuts further. **Figure 9.**

Figure 9: Use two fingers to tighten bolts

It is recommended to check that the pins on the PTB protrude and lock into the vertical tubes appropriately. To verify their function, first pull both knobs out and start moving the sliders to a different position. Release the knobs before the next hole and continue to slide the system. Once over a hole the pin of the knob will drop into the opening. Once each slider is in a new position apply a down force on the PTB directly downwards. The pins of the sliders should not come out of the holes. If the pins do not stay in the tube holes, the pull knobs need to be replaced; call Balanced Body for replacement parts.

EQUIPMENT INSPECTION AND MAINTENANCE LOG

We suggest that you keep a maintenance log for each piece of equipment. The log should include:


1. A description of the machine including the serial number, the date and place of purchase, and the manufacturer. All of this information should appear on the invoice.
2. Date and description of all required maintenance and inspections performed.
3. Date and description of each repair, including name and contact information for person or company performing the repair.

EQUIPMENT ATTACHMENT POINTS

Before using equipment, ensure that the springs are securely attached. See photos below for examples of securely attached springs to various attachment points that you may encounter on Balanced Body equipment. Note that adding spring resistance makes some exercises easier and others more difficult. Performing exercises with improper spring resistance increases the risk of falling. If unsure of proper spring resistance, please consult a certified Pilates instructor.

"Hour Glass" spring attachment point found on the Allegro 2 and some versions of the Allegro Reformer, as well as the upper (pre-load) attachment points on a Revo springbar.

"Hook" attachment point found on Revo Springbars, Standard Springbars and some versions of the Allegro Reformer. Note the Hour Glass posts in the upper section of this springbar as a pre-load position, shown here without the springs attached.

"Chair Cactus" attachment point found on various Chairs. It may also be flatter and made of sheet metal.

[Click here for Balanced Body Patent Data.](#)

Eyebolt and loop attachment point found on Cadillacs, Wall Units and Towers.

"Reformer Cactus" spring attachment point found on the IQ, Metro IQ, and one version of the Allegro Reformer. It may or may not have notches and bends at the bottoms of the cactuses.

MAINTENANCE SCHEDULE

All Equipment	Day	Wk.	Mo.	Qtr.
Inspect springs for gaps & kinks		✓	✓	
Inspect snaps for wear			✓	
Inspect nuts & bolts for tightness			✓	
Reformers				
Clean wheels and tracks		✓		
Inspect springbar hooks/eyebolts				✓
Rotate springs				✓
Inspect ropes/straps				✓
Inspect footbar supports				✓
Inspect springbars				✓

REPLACEMENT PARTS

To order replacement parts, or if you have any questions, please call:

U.S. and Canada: 1-800-PILATES (1-800-745-2837)

United Kingdom: 0800 014 8207

Other locations: +1 916-388-2838

Fax: 916-379-9277

Email: info@pilates.com

www.pilates.com

5909 88th Street, Sacramento, CA 95828 USA

ASSEMBLY AND MAINTENANCE PODCASTS

View our library of assembly and maintenance videos at <https://video.pilates.com/categories/bb-garage>.

Or scan this QR code with your mobile device:

REFORMER SAFETY USAGE

Safety and the Universal Reformer

The Pilates Universal Reformer, is designed to teach the exercise repertoire of Joseph H. Pilates. All utilization must take the unique functionality and design of the Reformer into account as well as the skill, expertise and education of the instructor. **Exercises outside of the intended use may cause damage to the equipment and may be dangerous to the client.** Furthermore, exercises improperly taught may result in loss of balance, coordination and control of the client on the moving carriage and may lead to injury. The Pilates Universal Reformer is only to be used under the supervision of a Pilates instructor specifically trained in the safe use of the Reformer.

Balance. Coordination. Control.

One of the unique features of the Reformer is the utilization of springs as resistance against the travel and movement of the Reformer carriage. In order for exercises to be performed safely, three important conditions must always be met.

- User must always maintain and be aware of balance when on the Reformer.
- User movements must coordinate with the speed, timing and variability of the spring resistance.
- ⚠ User movements must be controlled and organized. There should NEVER be slack in the ropes or bouncing of the springs.
- User should move within their ability, strength, and flexibility.

Control Front Advanced Exercise. Reformer movements are complex and client must always move with Balance, Control and Coordination.

Finding an Instructor

Working with a qualified Pilates instructor who is certified to train on the Pilates Reformer is an important component to learn how to use the Reformer safely and effectively. There are many Pilates education programs and it is important for users to be discerning. Look for instructors who have taken a certification with no less than 200 hours on the Reformer.

Spring Resistance

For many new to Pilates anticipation of the spring resistance needs to be trained and understood for safe utilization of the Reformer. In addition, understanding relative spring weights, how to change springs safely, how to check the springs for wear and how to move the spring bar are an integral part of spring safety specific to exercise transitions and set up. Ensuring that the springs are set up appropriately for the specific client, are set on the hooks to prevent detachment of the spring, and are understood in relationship to the movement of the carriage will decrease potential incidents and increase success for both client and instructor.

Check the springs regularly for wear and replace as needed. We recommend replacing springs every 2 years - 1 year for group fitness classes.

Check springs regularly for wear and replace as needed. We recommend replacing springs every 2 years - 1 year for group fitness classes.

Inexperienced Clients

Pilates exercises are integrated, complex motor learning exercises. Practice and experience are necessary to progress. Instructors, please be cautious with how quickly you progress your clients, if they can maintain balance, coordination and control then adding additional complexity is possible. Clients, please be patient with the learning process and understand the value of balance, coordination and control.

Exercise Set Up. Reformer Mounting and Dismounting ⚠️

Potential injuries may occur during the transitions between exercises as clients mount and dismount the Reformer. Proper instruction for exercise set up, mounting and dismounting as per the balance, control and coordination needed for a particular exercise, and for transitions between exercises, is primary to ensuring user success and to mitigate accidents or incidents of concern. For instructors, it is important, when in doubt, to refer to instructional manual and for clients to understand the importance of seeking out highly qualified instruction.

The Sitting Box and Foot Strap ⚠️

Mounting of the box for exercises seated, (Short Box Series), exercises supine, (Backstroke and Teaser), exercises prone, (Overhead Press and Swan), all require specific mounting and dismounting instructions. When the footstrap is required, please make sure it is properly secured to your Reformer (see individual Reformer setup instructions). The strap must maintain contact with the clients foot at all times. Please review instructional manual for more specifics regarding the sitting box and footstrap. Some general suggestions are included here.*

Short Box correct sitting position. When the box is placed horizontally on the carriage for Short Box exercises, the ends extend beyond the carriage platform. For safe mounting clients must always sit toward the center of the box.

Short Box incorrect sitting position. When the box is placed horizontally do not sit on the edges of the box, the box will tip.

Short Box correct strap position. Make Sure the footstrap is securely attached to the Reformer and the foot is in contact with the strap at all times.


Long box correct box and body placement. For advanced exercises, such as Long Box Teaser, it is important that the box is correctly placed on the Reformer and the body position is optimal to maintain balance, control and coordination throughout the movement.

* ⚠️ Sitting Boxes are designed for sitting, lying and standing. It is not designed for jumping, stomping or step aerobic type activity.

The Ropes

The Rope adjustments and locking mechanisms differ from Reformer to Reformer. Safe utilization of the ropes require they are installed correctly. For all Reformers except the A2, ensure rope is threaded through eye straps and cam cleats as pictured. As for the A2, the locking mechanism is underneath the carriage. Make sure the ropes are secure prior to movement.

Picture of rope cleats

Standing Exercises ⚠

Standing exercises are part of the Pilates repertoire and are functional and integral to a complete exercise experience. However, standing on a moving carriage is complex, potentially dangerous and must be handled with appropriate caution. When possible, instructors are encouraged to spot and provide assistance. All other utilization must be done with caution, after the clients abilities have been assessed, and when instructor and client have created a level of trust and communication to minimize any potential loss of balance, coordination and control. Use of stability props, such as gondola poles, ballet bars or personal assistance, are encouraged as client is learning and based on client ability.

Standing Side Splits with Gondola Pole or spotter

Working off the Back End ⚠


There are some lovely contemporary exercise opportunities from the back end (riser end) of the Reformer. However, extreme caution is required. Do not let go or lose of contact with the carriage when the carriage pulled to the back end and the springs are extended fully. This will result in high velocity impact on the front end (footbar end) of the Reformer which may result in potential damage to the frame, bumpers, spring mechanism and more. Sitting, kneeling and pulling while on the long box prone are all possible from the back end. Caution is required.

The Footbar ⚠

There are many different footbars and depending on which Reformer you have the footbar may or may not have a locking mechanism. Please, make sure to always check that the footbar is in a secure position and if there is a locking mechanism that it is locked and secure.

Set and Check

Whenever making any adjustments to the Reformer (ropes, footbar, risers) always set and check before using.

The Allegro 1 Footbar has several locking positions. Please consult your Reformer's manual as footbars differ between Reformer models and options.

High Kneeling Exercises ⚠️

Chest Expansion - The exercise Chest Expansion is a staple in the Pilates repertoire. It is considered an intermediate to advanced exercise because it requires upright balance on the moving carriage. Clients are facing the riser and the timing as a client moves and pulls the springs may create instability. When there is too little resistance, user may pull with too much force, the carriage will free roll and they may lose balance. Where there is too much resistance a client may find themselves losing balance and being pushed forward. **Please spot and use caution until clients are comfortable with the exercise.** ⚠️ **Standing facing risers with hands in straps is not recommended as there is a potential risk of falling and injury.**

Chest Expansion*

Kneeling Arms facing the footbar - This exercise can quickly become unsafe for clients. ⚠️ **We do not recommend this exercise kneeling or standing as there is a potential risk of falling and injury.** Also use caution when clients are in low kneeling as the position may be uncomfortable leading to clients adjusting and potentially losing balance.

⚠️ **Kneeling Arms Facing Footbar.** We do not recommend as balance, control and coordination may be a challenge to maintain – there is a potential risk of falling and injury.

Ballistic Movements ⚠️

Ballistic movements result in slack in the straps and a momentary floating of the carriage, an almost feeling of momentary weightlessness. When the straps become taught the resistance is reengaged. If not anticipated it may create stress in joints, a loss of control or balance for the user. **Under no circumstance should the carriage ever bang against the bumpers or against the back end.** Controlling movements is key to the success of Pilates conditioning and training.

Group Reformer

When teaching group, or taking a group Reformer class, it is important to ascertain the instructional quality of the teachers, the value that is placed on solid Pilates foundations and repertoire and the organizational philosophy on teaching Pilates. Maintaining Balance, Coordination and Control are key to successful and safe group Reformer class experiences.

Listening, the Client / Instructor relationship

In order to maximize the positive conditioning of the Pilates repertoire and Pilates inspired exercises, it is important for client and instructor to establish a good working relationship based on mutual trust and communication. If client indicates concern, caution or general discomfort, it is important to listen, assess and adapt as needed. Creating an environment for the client to share their experience on the Reformer will be key to making sure they are moving within safe ranges of motion, are moving safely and the exercise choices are effective and appropriate for their body.

Reformer Maintenance

Keeping your Reformer in good condition is key to keeping it working optimally for you and for your clients. The Safety Guide document included with your equipment has recommended schedules for inspection and replacement of wear parts. Balanced Body also offers an optional maintenance log that you may find helpful. Please follow recommended maintenance schedules and practices to optimize your Reformer's performance.

Sample Reformer Exercises

FOOTWORK, ALL LEVELS

10 reps • Springs: 2 – 4 springs • Bar: Middle or High • Head rest: Up

FOCUS

- » Breathing – exhale out/inhale in, or inhale out/exhale in
- » Neutral spine
- » Pelvic stability
- » Hip, leg and ankle alignment
- » Hip, leg and ankle strengthening
- » **Circulation**
- » **Isolation** – release unnecessary tension in the upper body and hips

Precautions

Sensitive to ankle, knee, hip flexion, spinal compression

MOVEMENT

Starting position

Supine on carriage, feet on foot bar, legs hip width apart

Heels

Heels on foot bar, push back and return

Toes

Ball of foot on bar, heels slightly raised, push back and return

Prehensile

Ball of foot wrapped around the bar, push back and return

Pilates V

Ball of foot on bar, turned out, low releve, heels together, push back and return

Flex/Releve

Ball of foot on bar, parallel, legs straight, plantar flex ankle, dorsiflex ankle, plantar flex ankle, bend knees, push back to starting position

2nd position

Heels at ends of bar, slight turn out, push back and return

Running in place

Ball of foot on bar, dorsiflex one heel, bend the other knee, alternate legs 20-50 times

Heels

Toes

Prehensile

Pilates V

HUNDRED, ALL LEVELS

Sets: 10 • Springs: 1 – 3 • Bar: None Loops: Regular • Headrest: Up

FOCUS

- » Percussive breathing - in for 5, out for 5
- » Stable pelvis – imprinted or neutral
- » Hollow abdominals
- » Abdominal strength
- » Shoulders down
- » Neck long

Precautions

Back injuries, neck injuries, hip flexor injuries, osteoporosis

Prerequisites

Hundred on the mat

MOVEMENT

Starting Position

Lying supine on carriage, knees at 90 degrees, hands in loops, arms to ceiling

Level 1

Knees bent at 90 degrees, reach arms to sides as the head and upper body lift off the carriage, pulse arms with breath

Level 2

Legs straight up to ceiling, reach arms to sides as the head and upper body lift off the carriage, pulse arms with breath

Level 3

Reach arms to sides as the head and upper body lift off the carriage, straighten legs to ceiling then lower legs keeping low back on mat, pulse arms with breath

KNEELING ABDOMINALS, LEVELS 1-2

10 reps • Springs: 0 – 2 • Bar: None • Loops: None

FOCUS

- » Breath – Exhale as knees pull in/ inhale to return to starting position
- » Pelvic stability
- » Abdominal strength
- » Pelvic stability
- » Scapular stability

Precaution

Shoulder, elbow and wrist problems, inability to kneel

Prerequisites

Mat all fours hollowing (pregnant cat abdominals)

MOVEMENT

Starting position

Kneel on all fours facing straps with hands on the edges of the frame

Standard Exercise (flat back)

Face straps with knees against shoulder rests, keeping spine neutral and parallel to the ground, pull carriage toward the head by flexing at the hips

Variation – Oblique 1

Place both hands on one frame allowing torso to side bend and rotate while keeping hips as square as possible, pull both legs in increasing the movement of the torso

Variation – Single leg 1

Lift one leg off carriage and extend it in a straight line behind you, keeping the shoulders as square as possible pull one leg in maintaining a neutral back position

FEET IN STRAPS, LEVEL 1

6 reps • Set up: 2 springs • Bar: AnyLoops: Regular

FOCUS

- » Breath - exhale out/inhale in or inhale out/exhale in
- » Spine to mat or neutral spine
- » Hollow abdominals
- » Abdominal strength
- » Pelvic Stability
- » Hamstring, adductor and gluteal strength
- » Hamstring and adductor flexibility
- » Leg and hip alignment
- » Hip range of motion

Precautions

Hip flexor injury, limited hamstring flexibility, back injuries, weak abdominals

Prerequisites

Adequate hamstring flexibility, ability to stabilize the back, Hundred

MOVEMENT

Starting Position

Supine on carriage, loops around arches

Leg lowers

Hips stable, begin with inner thighs together, lower and raise legs

Variations: Legs parallel, turned out or turned in Magic circle or ball between the legs

Scissors

Hips stable, begin with inner thighs together, open legs to sides and return

Variations: Legs stay over hips, legs move down toward the bar while carriage moves (V's)

Circles

Hips stable, begin with inner thighs together, moving legs down and around in circles or D's, reverse directions Variations: Legs parallel, turned out or turned in, or knees in straps

SHORT BOX - ABDOMINALS, ALL LEVELS

6 reps • Springs: 4 - 5 • Box: ShortStrap: Foot strap • Pole in hands

FOCUS

- » Breath – inhale down/exhale back or exhale down/inhale back
- » Abdominal strengthening
- » Abdominal hollowing
- » Back strengthening
- » Back flexibility

Precautions

Back, neck and shoulder injuries, osteoporosis

Prerequisites

Reformer Rollbacks, Cadillac Rollbacks

MOVEMENT

Starting Position

Sitting on short box, facing footbar, feet under foot strap, knees slightly bent

Standard Exercise

Roll down with a long curve and roll back up

Oblique variation

Roll down with a long curve, rotate the torso to each side, roll back up

ARM WORK, LEVEL 1 - 2

4 - 10 reps • Springs: 1-2 • Box: long, short or none • Loops: very short, short or regular

FOCUS

- » Breath – inhale pull/exhale release
- » Biceps, triceps, pectoralis and deltoid strengthening
- » Scapular stabilization
- » Torso stabilization
- » Sitting posture

Precautions

Wrist, arm or shoulder problems, back problems with limited sitting ability

MOVEMENT

Sitting variations for all exercises

Sitting on carriage, cross-legged, legs straight or kneeling

Exercises Facing the Straps

Biceps

Holding very short loops in hands, bend elbows to pull straps to shoulders.

Triceps/Posterior Deltoid

Loops in hands, arms straight, pull straps back level with hips and pulse arms back.

Exercises Facing Footbar

Serve a tray

Sitting facing footbar, regular loops in hands, elbows bent, reach forward, straighten arms, open arms to the side palm up and return.

Hug a tree

Sitting facing footbar, regular loops in hands, arms out to sides, soft elbows, bring fingertips toward each other.

Biceps

Serve a tray

Hug a tree

Triceps

Sample Tower Exercises

DEFINITION OF EXERCISE SET-UP TERMS

- » **Level:** The level of expertise needed to undertake exercise.
- » **Reps:** How many times the exercise is performed.
- » **Springs:** Number and location of springs on the Allegro Nextgen Tower.
- » **Loops:** Which loops should be used during exercise.*
- » **Focus:** What should be emphasized during exercise.
- » **Precautions:** Physical conditions that may limit or exclude a participant. Exercises may need to be modified for people with these conditions.
- » **Prerequisites:** Specific exercises that must be mastered before undertaking a new exercise.
- » **Starting Position:** Where to begin the exercise on the Allegro Nextgen Tower. *If applicable

ROLLBACKS, LEVEL 1

Springs: 2 long yellow or 2 short yellow springs from high position **Reps:** 6-10
Rollback Bar or handles

Focus

- » Breathing – exhale roll down, inhale at the bottom, exhale to roll up
- » Balance between abdominals and lumbar extensors
- » Soft neck and shoulders
- » Maintain C-curve
- » Soft hip flexors

Precautions

Shoulder and neck problems, some low back problems, osteoporosis

Starting position


Sit facing Tower, holding on to bar or handles, knees soft, feet on metal bars.

Standard Exercise

Hold bar with arms straight, roll down, curving back and staying lifted, roll back up maintaining slight flexion in spine

Oblique Variations

- » Wooden bar (Water skiing)
- » Sit diagonally on table, place left foot against pole, cross right foot over ankle, place left hand on bar and reach right arm open while rotating torso to the right.
- » Roll down and up maintaining rotation, then switch sides.

Water skiing

Roll Backs

SEATED PUSH THROUGH, LEVEL 1

6 reps

Springs: 1 short red spring from above on Push-through Bar

FOCUS

- » Breath – exhale stretch forward, inhale stretch up or reverse
- » Spinal flexion and extension
- » Hamstring flexibility
- » Scapula stability/mobility
- » Abdominal hollowing and lift
- » Coordination of breath with spinal mobility

Precautions

Some back injuries, shoulder injuries

Prerequisites

Mat – Spine Stretch

Starting Position

Sit on table facing Tower, with feet against upright bars, and both hands on the Push-through Bar (knees can be bent or straight depending on flexibility).

Standard Exercise

Push bar down, curve spine forward beginning with top of head, reach the bar forward and stretch.

Return by hollowing out abdominals and stacking one vertebra on top of another until sitting up on the sit bones.

Press bar up and lean forward from hips with a flat back before beginning again

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
» Cacciavite Phillips

STATEMENT OF INTENDED USE

The Allegro Nextgen Reformer™ device is intended by the manufacturer to be used by people for the purposes of physiotherapy, rehabilitation, the prevention, treatment and relief of injuries, supportive, functional therapy for a disability or for general strengthening and improving the mobility of the human body. Should injury occur during usage please report to the manufacturer and competent authority of the Member State in which the user and/or patient is established.

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