

Assembly of The Balanced Body® Clinical Reformer®

THIS IS A MEDICAL DEVICE

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

 Balanced Body, Inc.
5909 88th St,
Sacramento, CA 95828 USA

Contact Us
1-800-PILATES | +1-916-388-2838
pilates.com | info@pilates.com

17285 | 2021-06

STATEMENT OF INTENDED USE

The Clinical Reformer® device is intended by the manufacturer to be used by people for the purposes of physiotherapy, rehabilitation, the prevention, treatment and relief of injuries, supportive, functional therapy for a disability or for general strengthening and improving the mobility of the human body. Should injury occur during usage please report to the manufacturer and competent authority of the Member State in which the user and/or patient is established.

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
How to Assemble your Clinical Reformer

PARTS LIST (INCLUDED):

DESCRIPTION	Part Number	Qty.
Ankle Velcro (pair)	101-001	1
Padded Double Loops (pair)	101-009	1
Neoprene Handles (pair)	101-016	1
Thigh Velcro (pair)	101-020	1
Carriage Extension Stopper	200-211	1
Adjustable Ropes w/snaps (pair)	210-026	1
Padded Footstrap	210-031	1
Risers w/pulley	950-184	2
Infinity Foot Plate	707-187	1
Sitting Box	721-002	1
Shoulder Rest Kit, TwistLock	950-251	1
TwistLock Post	616-400	2

DESCRIPTION	Part Number	Qty.
Spring, Reformer, Red	SPR9070	3
Spring, Reformer, Blue	SPR9071	1
Spring, Reformer, Yellow	SPR9241	1
Spring, Reformer, Green	SPR9293	1
Infinity Footbar	950-189	1
Infinity Trunion, Right (installed)	950-112	1
Infinity Trunion, Left (installed)	950-113	1
Plunger lock out knob (installed)	GEN8000	2
Wood Screw	GEN9073	2
Footbar Stopper	REF0702	2
Standing Platform Extender	17873	1
Standing Platform Extender Pad	17878	1

TOOLS NEEDED (NOT INCLUDED):

Phillips Screwdriver

IN CASE OF FREIGHT DAMAGE If there is freight damage, make certain you keep all packaging material. You must call Balanced Body within 3 business days.

ATTACH THE FOOTBAR

This can be done by one person, but it is easier and safer to have a person hold each side.


Figure A

Figure B


Figure C

1. Pull the plunger knob straight out so that the recessed cross-pin comes out of the slots. Then turn the knobs so that they are locked with the retracted pin.
2. Slide the trunion into the aluminum C-channels that run down both sides of the Reformer frame. (Figure A and B)
3. Twist the plunger knob until it locks into the desired position.
4. Attach the black stops to the end of the channels with the provided wood screws. (Figure C)

INSTALL THE CARRIAGE

Turn the carriage over and clean all the wheels of shipping dust, etc. Clean the full length of the rail surfaces inside the frame. Carefully place the carriage into the frame with the headrest facing away from the footbar.

INSTALL THE SHOULDER RESTS

Note the two protruding studs on the bottom of the shoulder rests. Lower the shoulder rest studs into the key holes and slide shoulder rests away from the carriage pad. See Figure D. Note: The shoulder rests can be installed in two positions. One position is wider and more comfortable for broad shoulders. Swap the left and right shoulder rests to change between standard and wide configurations.


Figure D

INSTALL THE LOCKING SHOULDER POSTS

Screw the locking shoulder posts onto the exposed bolt on each shoulder rest. Turn the posts clockwise to tighten them and lock the shoulder rests into position. See Figure E. Over tightening the posts will make them difficult to remove. They should be snug but not too tight. To remove the shoulder rests, turn the post counter-clockwise 1-2 turns and gently slide the shoulder rests towards the carriage, then lift them up and out.


Figure E

ATTACH THE SPRINGS TO THE UNDERSIDE OF THE CARRIAGE

1. Locate the bracket with five eyebolts attached beneath the carriage. Hook the springs onto the eyebolts. The hook should be facing down.
2. The bodies of the springs should rest on the flat, padded brace under the carriage (Figure F).

Figure F

INSTALL THE SPRINGBAR

3. If your Reformer is configured with a no-roll springbar, install it into the "gears" under the footbar with its hooks facing up.

INSTALL THE RISERS

4. Remove the detent pin from the metal riser receiver and insert the riser with the pulley facing the carriage. Re-insert the detent pin to set the desired riser height.

ATTACH THE ROPES

5. To attach the ropes, first engage at least one spring from the carriage to the springbar to stabilize the carriage. Unroll the ropes and separate them.
6. Clip a cotton loop onto the end of each rope and hang the loops on the shoulder rests. (Figure G.) Thread the other end of each rope through a pulley and back into the cam cleats on the carriage to adjust the length. Be sure to go through the chrome eyestraps on both sides of the cam cleats. Always push the rope firmly down into the cleats to ensure a good drip.

Figure G

USING THE EXTENSION STOPPER

To limit the carriage travel, insert the extension stopper into the desired hole in the top of the frame rail.

INSTALLING THE PLATFORM EXTENDER

To install the Studio Reformer Standing Platform Extender on a Reformer with an Infinity Footbar, first loosen the two black plastic knobs beneath the standing platform at the end of the frame. Adjust the footbar so that it is laying horizontally with the Reformer frame (Figure H)

Figure H

Gently lower the tongue of the Standing Platform Extender into the silver bracket on the inside of the frame. Gently tighten the black knobs on the outside of the frame. (Figure I).

Figure I

EASY TO USE PLUNGER LOCKOUT KNOBS

Balanced Body's Allegro® Reformer, and Clinical Reformer® and Studio Reformer® with Infinity footbars come with plunger lockout knobs that make adjusting and moving the Reformer footbar along the frame much easier.

To use, pull the plunger knob straight out so that the recessed cross-pin comes out of the slots as shown in Figure H. Then turn the knob so that the cross-pin rests on the bolt head as shown in Figure I. Repeat on other side.

With the knobs in the "locked out" position, reposition the footbar along the frame of your Reformer as desired. Then turn the knobs so that the cross-pin is once again recessed, firmly locking the footbar into its new position.


Figure J


Figure K

INSTALLING THE JUMPBOARD ON A CLINICAL REFORMER WITH AN INFINITY FOOTBAR®

To install the jumpboard (regular or Jumpus Maximus) on a Reformer with an Infinity Footbar, first loosen the two black plastic knobs beneath the standing platform at the foot end of the frame.

Adjust the footbar to the 4th position from the end of the Reformer. Then adjust the footbar to the second highest vertical position (Figure L). Next, place the tongue of the jumpboard into the silver bracket on the inside of the frame. Finally, re-tighten the black knobs on the outside of the frame.

Figure L

If you have any questions, please call Balanced Body at 1-800-PILATES (1-800-745-2837) or +1-916-388-2838.

ADDITIONAL SAFETY AND USAGE

For additional usage and safety information or to view complete equipment manual translations (**USER GUIDE #15451**), please visit <https://www.pilates.com/store/faq-safety/>. For printed instructions, contact your seller.

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WARNING: This machine is potentially dangerous, and Balanced Body Inc. assumes no liability towards use or misuse. If you are not fully aware of the safety issues and use of this machine, do not use it. Use only with qualified instructor supervision. Please refer to the Additional Usage and Safety section for more information.

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[Click here](#) for **Balanced Body Patent Data**.

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QUESTIONS?

Please call Balanced Body Technical Support at
1-800-PILATES (U.S.) or +1-916-388-2838 (International).