
Ecoline
VC-MS/CA8-6
VC-MS/CA4-12

Ecoline VC-280
Ecoline VC-380

Ecoline VC-360

Ecoline
VC-Easy-Load

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**Tubing pumps
with variable speed drive**

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Please note

We recommend you read this operating manual carefully. When operating a pump, certain hazards cannot be excluded.

ISMATEC® does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC® does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- ▶ The circuit between mains supply and pump has to be connected to earth ground.
- ▶ The pump must not be operated outside the designed operating and environmental conditions.
- ▶ The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes.
- ▶ The pump must be switched OFF when cassettes or tubing are inserted or changed.

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Caution

Rotating pumphead creates a pinch and crush hazard.

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Safety precautions

- ▶ The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when routing tubing in explosion-proof chambers.
- ▶ Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- ▶ Tubing pumps consist of rotating parts (rotor of the pump-head). Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed is completely shut.
- ▶ Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.

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Safety precautions

- ▶ Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the rotating roller head.
- ▶ Only new fuses, according to the specifications stated on Page 11 in this manual, must be used.
- ▶ The fuse-holder must not be short-circuited.
- ▶ Do not open or remove the housing while the pump is operating.
- ▶ Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- ▶ For service and repairs carried out by the customer or by third-party companies ISMATEC® denies any responsibility.

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Warranty

1 year from date of delivery

We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs for instruments under warranty are charged to the customer.

Our warranty becomes invalid in case of:


- ▶ improper operation by the user, or if the pump is diverted from its proper use,
- ▶ unauthorized modification or misuse by the user or by a third-party,
- ▶ improper site preparation and maintenance,
- ▶ operation outside environmental and electrical specifications for the product,
- ▶ use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- ▶ damages caused by contamination or leaks due to torn or burst tubing.

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Product

Package contents

- ▶ Pump (model as ordered)
- ▶ Power Cord
- ▶ Operating Manual

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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Mounting screws only for Easy-Load® (supplied
with every drive)

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a) Are supplied with socket head cap screws and an
allen wrench

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


Rear panel

1. Mains socket
2. Mains switch ON/OFF (Start/Stop)
3. Opening for screw-driver
4. Fuse-holder containing 2 fuses (change always both fuses)
5. Voltage selector (see Page 12)
6. Analog input: Speed control (0–5 or 0–10 V / 0–20 V or 4–20mA)
Digital input (TTL level): Start/Stop, rotation direction

⚠ Ventilation slits (a)

For sufficient ventilation keep the slits at the back and the base slits clear.

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Socket/Power cord

Use exclusively the originally supplied power cord.

The socket must be connected to earth ground (protective conductor contact).

Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage. If necessary, the voltage setting must be changed and the 2 fuses must be replaced.

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
Voltage setting

115V / 230V and changing the fuses

➔ Switch the pump OFF, pull out the mains plug.

1. Pull out the fuse-holder by opening it with a small screw-driver (size 0).
2. Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.

3. Insert 2 new fuses
230 V_{AC}: 2 x T1A/250V
115 V_{AC}: 2 x T2A/250V

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

4. Shut the fuse-holder. The voltage rating is visible in the window.

⚠ Before starting-up

Check that the voltage setting visible in the window of the fuse-holder complies with the local mains voltage.

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
Operating panel

1. Operating LED
on/off / Overload (blinking)
2. Rotation direction
left / stand-by / right
3. Speed setting
2-digit potentiometer, 1–99%, 1% steps

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Starting the pump

- ▶ Check that the voltage indicated on the identification label of the pump complies with the local mains voltage. Change, if necessary, as indicated on Page 12.
- ▶ Insert the type of tubing required for the specific application. (For the tubing handling see Pages 18–24).
- ▶ Connect the pump tubing to the system.
- ▶ Set the rotation direction.
- ▶ Connect the pump to the mains and switch it on.
 - ➔ The roller head will begin to rotate immediately unless the potentiometer is set to «00» or the rotation direction switch is in the stand-by position.

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Caution

Please be aware that siphoning may occur at the moment the cassettes are released or the tube-bed is opened!

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Easy-Load Schlauchbett öffnen

Pumping against pressure

The Ecoline pumps and the tubing are designed for applications at differential pressure conditions not exceeding 1.5 bar.

When the pump is idle

When the pump is idle, it is recommended that pressure be released from the tubing. Releasing the cassette on the right side is sufficient. This protects the tubing from unnecessary strain and prolongs its service life.

VC-MS/CA8-6 and VC-MS/CA4-12

Cassettes can easily be released by pushing up the fixing-tongue

VC-360 Open the tube-bed at the notch lever without removing the tubing

VC-280/1 and VC-380/1 Remove the tubing from the tube-bed

Easy-Load Open the tube bed

Pompage contre pression

The Ecoline pumps and the tubing are destined for applications at differential pressure conditions not exceeding 1.5 bar.

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Typical causes are:

excess backpressure, highly viscous medium, plugging or mechanical blockage, wrong tubing wall thickness, or low mains voltage.

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Overload condition

If the pump can't reach its set speed, the 'POWER' LED starts flashing. The pump continues to deliver but not at the desired flowrate. It can be operated for some time under a small overload.

However, a reduced lifetime of motor and gearbox must be taken into account. As soon as the load decreases and the pump reaches its set speed, the LED stops flashing.

Should the load increase further instead, a motor current limiter stops the device to prevent damage. The pump does not rotate but its 'POWER' LED keeps flashing. To reset, turn the device off. Remove the cause for the overload before operating the pump again.

Overload situations occur mostly in the combined presence of large tubing inner diameters and high speed.

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Please note

Please refer to our detailed tubing documentation.

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(Best.Nr. SC0179) einzusprühen.

Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, it is recommended that new tubing be run-in with water or the medium to be pumped for at least 1 to 3 hours before the application is started.

Tubing life

The service-life of the tubing depends on the application and the tubing material used.

Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog.

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray

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Inserting the tubing

Models VC-280/1, VC-380/1

1. ➡ Switch the pump OFF.
Loosen the screws of the protection
plate and remove the plate
2. Move the upper tubing clamp to
the middle position.
3. Insert the tubing and secure it with
the upper tubing clamp.
4. Move the lower tubing clamp to
the middle position.

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

Inserting the tubing

Models VC-280/1, VC-380/1

5. Turn the rotor counter-clockwise and simultaneously insert the tubing along the tube-bed until the tubing is squeezed by two rollers (tubing must not be twisted).
6. Secure the tubing with the lower tubing clamp.
7. Replace the protection cover and fasten it with the two screws.
8. Switch the pump ON.
➔ Run the pump for a few seconds; depending on the rotation direction, open the upper or lower tubing clamp and restretch the tubing carefully.

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Interchangeable rotor

for VC-280/1 and VC-380/1

2 rollers for higher flow rates
3 rollers for lower pulsation

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


Changing the rotor

1. ➔ Switch the pump OFF.
2. Loosen the screw on the rotor using an Allen key (size 2 mm).
3. Remove the rotor by pulling it towards you.
➔ Be careful that the wedge and spacer remain on the shaft.
4. Insert the new rotor
➔ The screw on the rotor must be aligned with the wedge on the shaft.
5. Fasten the rotor screw tight using the Allen key.

⚠ For manipulating a pump-head also observe the safety precautions listed in the operating manual.

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Inserting the tubing

Model VC-360

- ➔ Switch the pump OFF.
- Press the notch lever upwards.
- Open the tube-bed.
- Push the upper tubing clamp upwards; insert the tubing.
- Loop the tubing round the rotor; press the lower tubing clamp downwards; insert the tubing.
- Close the tube-bed and press it to the right so that the notch lever locks into place.
- ➔ Switch the pump ON.
- Run the pump for a few seconds; depending on the rotation direction open the upper or lower tubing clamp; restretch the tubing carefully.

Important

When the pump is idle, open the tube-bed in order to protect the tubing from unnecessary load.

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
Inserting the tubing

Model Easy-Load®


- ➔ Switch the pump OFF.
- Open the tube bed.
➔ Rotate the lever to the left.
- Load the tubing. Center the tubing between the retainers.
- Close the tube bed.
➔ Rotate the lever to the right.
- To prevent the tubing from moving, fix it with the tubing retainers on the left and right side of the pump-head.
➔ Press the retainer slightly towards the head and simultaneously downwards. The retainer locks in the position according to the tube size.
➔ For optimum tubing life, keep tubing straight where it enters and exits the pump-head.
- ➔ Switch the pump ON.

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Inserting the tubing

Model VC-MS/CA8-6 & 4-12

1. ➔ Switch the pump OFF.
2. Remove the cassette by slightly pressing the fixing-tongue and lifting it simultaneously.
3. Insert the tubing with one collar into the cassette.
4. Let the tubing hang down (prevent it from being twisted).
5. Insert the tubing with the second collar at the other end of the cassette.
6. Reinsert the cassette into the roller-head.
7. ➔ Switch the pump ON.

Important

When the pump is idle, release all cassettes at the fixing-tongue.

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Flow rates

In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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
Digitale Eingänge (TTL-Pegel)

Digital inputs (TTL-level)

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| | |
|--------|--------------|
| Pin 2, | shield |
| Pin 3, | stop/run |
| Pin 4, | direction |
| Pin 6, | speed intern |

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Analog interface

➔ When using the analog interface, the switch for the rotation direction (see Page 13, paragraph 2) must be in the neutral middle position.

Pin 1, GND (ground)

Reference potential for all other inputs.

Pin 2, shield (GND)

This pin couples the shielding of the cable the ground.

Pin 3, stop/run

The pump starts when connected to Pin 1 (GND).

Pin 4, direction

In open position the pump operates clockwise; when connected to Pin 1 (GND) it pumps counter-clockwise.

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Pin 8

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Analog interface

➔ Switch for the rotation direction must be in the neutral middle position.

Pin 5, speed IN

For external speed control (0–5V, 0–10V, 0–20mA und 4–20mA).

Input impedance and input range can be selected via a DIP-switch inside the pump (see Page 29).

Pin 6, speed intern

When connected to the ground, the drive speed can be controlled by the potentiometer integrated in the control panel.

Pin 7, +5V

This pin provides the internal supply voltage of 5 Volts for the supply of a potentiometer (R_I = 330Ω).

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
Inutilisé

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Settings of switch S1

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| Pins | | Impedanz | DIP-Switch 1 | DIP-Switch 2 | DIP-Switch 3 | DIP-Switch 4 | DIP-Switch 5 | DIP-Switch 6 |
|-------------------|--------|----------|--------------|--------------|--------------|--------------|--------------|--------------|
| Pin 5 speed IN | 0-5V | 100kΩ | OFF* | OFF* | OFF* | OFF* | OFF* | OFF* |
| | 0-10V | 20kΩ | OFF | ON | OFF | OFF | OFF | OFF |
| | 0-20mA | 240Ω | ON | OFF | OFF | OFF | OFF | OFF |
| | 4-20mA | 240Ω | ON | OFF | ON | OFF | OFF | OFF |

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➔ Check that the pump is disconnected from the mains.

For reaching the DIP-switch, you need to unscrew the 4 screws on both sides of the pump. Lift the housing carefully (mind the cables) and remove it.

Operation via footswitch

With DIP switch 4 the user can switch between two possibilities:
 Dip-Switch 4 off = On/Off
 Dip-Switch 4 on = On, as long as pressure is maintained on the footswitch.

⚠ The instrument should only be opened by a qualified technician!

Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Footswitch

Order No. IS3572

This footswitch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc. Both hands are free for handling the bottles and tubing.

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
MS/CA Click'n'go

MS/CA Anpresshebel
MS/CA pressure lever
MS/CA levier de pression

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Accessories

Spare-cassettes in POM-C

MS/CA Click'n'go
Order No. IS3510A
MS/CA pressure lever*
Order No. IS0649A

Spare-cassettes in PVDF

MS/CA pressure lever*
Order No. IS3629A

*** The cassettes with pressure lever are available on request. This type of cassette provides better results at elevated differential pressure conditions (>1 bar).**

 **When using new tubing** for the first time, it may occur that, depending on the tubing used (hardness and diameter), the pump cannot be primed and, hence, does not deliver the liquid. If that is the case, it is recommended that the tubing be wet and the pump run with the tubing inserted for about 15 to 30 minutes.

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Disposal

Please retain packaging materials until the product warranty ends. Afterwards please discard packaging materials in an environmentally friendly manner according to local regulations. Once the useful life of the product has ended, please ensure proper disposal according to local laws. Plastic and electronic components should be disposed of at a recycling facility. Please refer to local regulations regarding proper disposal.

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Maintenance

Provided the Ecoline pump is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs, please send the defective Ecoline pump to your ISMATEC® representative.

Please give information on defect and date of purchase.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- ▶ spare-parts
- ▶ parts lists
- ▶ wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Differential pressure

¹⁾ Possible differential pressure with appropriate tubing material, tubing with small I.D.s may enable higher pressures.

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Ecoline VC-280/1, VC-380/1

DC Motor

| | |
|---------------------------|-----------------------------|
| microprocessor controlled | |
| Speed range | 3.5–350 rpm |
| Channels | 1 |
| Pump rollers | 2 or 3 |
| Differential pressure | max. 21.7 psi ¹⁾ |
| Remote control | (see Pages 27–29) |

Mains connection/Fuse rating

| | |
|-------------------|--------------|
| 230 V (50 Hz) | 2 x T1A/250V |
| 115 V (60 Hz) | 2 x T2A/250V |
| Power consumption | 100 Watts |

Operating conditions

| | |
|---|---------------|
| Temp. | +41° to 104°F |
| Rel. humidity | max. 80% |
| not condensing, at normal laboratory conditions | |

| | |
|-------------------|--|
| Dimensions/Weight | |
| DxWxH (in.) | 10 ¹ / ₈ " x 6 ³ / ₄ " x 5 ³ / ₈ " |
| Weight (lb) | 12 / 12 |

CE-compatibility proved according to:
EN 61326-1 / EN 61010-1

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Ecoline VC-360, VC-Easy-Load

DC Motor

| | |
|---------------------------|-----------------------------|
| microprocessor controlled | |
| Speed range | 3.5–350 rpm |
| Channels | 1 |
| Pump rollers | 3 |
| Differential pressure | max. 21.7 psi ¹⁾ |
| Remote control | (see Pages 27–29) |

Mains connection/Fuse rating

| | |
|-------------------|--------------|
| 115 V (60 Hz) | 2 x T2A/250V |
| 230 V (50 Hz) | 2 x T1A/250V |
| Power consumption | 100 Watts |

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CE-compatibility proved according to:

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Please Note

We also recommend you observe our **Warranty Terms as well as our Terms and Conditions of Sale.**

In case of any queries, please contact your local Ismatec® representative.

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Ecoline VC-MS/CA8-6 & 4-12

DC Motor

| | |
|---------------------------|-------------------|
| microprocessor controlled | |
| Speed range | 3.5–350 rpm |
| Channels | 8 4 |
| Pump rollers | 6 12 |
| Differential pressure | max. 21.7 psi |
| Remote control | (see Pages 27–29) |

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CE-compatibility proved according to:

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EN 61326-1 / EN 61010-1

ISMATEC®

Learn more
about the
Ismatec® product
line by visiting:
www.ismatec.com

For ordering and technical support, please contact:

North America sales@ismatec.com | 1-800-323-4340 | 1-847-549-7600

Europe sales.europe@ismatec.com | +49 (0) 9377 9203-0

 **Cole-Parmer®**