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## Enterolert\* 250 Test Kit

**Introduction**

Enterolert\* 250 is designed to detect *Enterococcus* species in bottled water samples. It is based on IDEXX’s proprietary Defined Substrate Technology\* (DST\*). Enterolert 250 provides confirmed presence/absence results in 24 hours. Enterolert 250 uses Ortho-Nitrophenyl-β-D-Glucoside as a nutrient indicator and incorporates a specifically designed blue background colour in its formulation. When the substrate is metabolized by enterococci, the sample turns from blue to green to indicate a positive result. Any change from the original colour to green is considered a positive result. No ultraviolet light source is required. Enterolert 250 detects enterococci at 1 cfu in 250 mL bottled water samples in 24 hours.

**Storage**

Store at 2–25°C away from direct light and humidity.

**Presence/Absence (P/A) Procedure for Noncarbonated Water Samples**

- Add contents of one snap pack to a 250 mL sample in a sterile, transparent vessel.
- Cap vessel and shake until contents are dissolved.
- Incubate at 41 ±1°C for 24 hours; results are valid for up to 28 hours of incubation.
- Read results according to Result Interpretation table, below.

**Presence/Absence (P/A) Procedure for Carbonated Water Samples**

- De-gas sample by vigorously shaking sample in its original container until no gas can be seen or heard (loosen cap to allow gas to escape after shaking).
- Add contents of one snap pack to a 250 mL sample in a sterile, transparent vessel.
- Shake vessel vigorously for 30–60 seconds to release carbonation (loosen cap to allow gas to escape after shaking).
- Secure cap and incubate at 41 ±1°C for 24 hours; results are valid for up to 28 hours of incubation.
- Read results according to Result Interpretation table, below.

**Result Interpretation**

Appearance	Result
Blue	Negative for enterococci
Green	Positive for enterococci

- Enterolert 250 results are definitive at 24–28 hours. In addition, positives for enterococci observed before 24 hours and negatives observed after 28 hours are also valid.

**Procedural Notes**

- Use only sterile, nonbuffered, oxidant-free water for dilutions.
- For comparison, a negative control can be used when interpreting results.
- This insert may not reflect your local regulations. For compliance testing, be sure to follow appropriate regulatory procedures.
- Enterolert 250 is a primary water test. Enterolert 250 performance characteristics do not apply to samples altered by any pre-enrichment or concentration.
- Enterolert 250 has not been validated for use with flavored bottled water.
- The presence of a high mineral content (especially magnesium and/or calcium) can cause Enterolert 250 reagent to become cloudy but does not affect test result outcome.
- Aseptic technique should always be followed when using Enterolert 250. Dispose of in accordance with Good Laboratory Practice.

**Quality Control Procedures**

- One of the following quality control procedures is suggested for each lot of Enterolert 250:
  - IDEXX-QC Enterococci DW<sup>1</sup>: *Enterococcus faecalis*, *Serratia marcescens*
  - For each of the American Type Culture Collection<sup>2</sup> (ATCC) bacterial strains (*Enterococcus faecalis* 29212, *Aerococcus viridans* 11563), streak the culture onto labeled TSA or blood agar plates and incubate at 35 ±2°C for 18–24 hours.
    - For positive controls, prepare proper dilutions to reach a concentration of approximately 10–50 CFU/250 mL.
    - For non-target strains, prepare proper dilutions to a suggested inoculation up to 10<sup>3</sup> CFU/250 mL.
- Follow the P/A Procedure above.
- Results should match the Result Interpretation table above.

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