

SoClean 2 - CPAP cleaning machine



1 Simple one-time setup

The SoClean CPAP cleaning unit connects with all types of PAP equipment including CPAP, Bi-PAP, BPAP, and VPAP.

2 Close the lid and go!

The SoClean sanitizer starts automatically after you place your CPAP mask in the chamber and shut the lid.

3 Kills CPAP germs in minutes

Activated Oxygen cleaning powers through your CPAP reservoir, hose, and mask.

Fits all popular mask types and CPAP machines


Try it Risk Free

Don't use a humidifier? No problem you can still use the SoClean CPAP cleaner.

Try Now Free\Risk-Free for 30-Days Trial

How does the SoClean cleaning unit work?

The SoClean CPAP cleaner and sanitizing machine uses activated oxygen cleaning to eliminate any germs and bacteria that may be lingering in your home CPAP machine. Watch our CPAP cleaning demo video to learn more.



WATCH A DEMO

SoClean the world 1st automated CPAP cleaner machine

SoClean 2 is CPAP Cleaner and sanitizer device that is used to clean and sanitize CPAP mask, reservoir, and hose without the need of separating your CPAP machine apart. So Clean 2 CPAP cleaner and sanitizer doesn't require chemicals or water to work because the tool takes advantage of Ozone known as activated oxygen.

How does So Clean 2 CPAP sanitizer Works?

Soclean CPAP sanitizer has the ability to produce activated oxygen that is pumped all over the system to kill molds, germs and bacteria. Activated oxygen is safe and natural hence does not lead to any danger to the human body. The inside chamber of the cleaner is designed in such a way that a wide range of masks could be used within the CPAP cleaning device. It can certainly fit perhaps even the most full face mask types since the inside chamber is large enough which measuring 7.25 ins height, 5.5 ins. length and 5.5 inches width.

Soclean 2 - CPAP cleaning machine is a powerful automated CPAP sterilizer for all CPAP devices. So Clean 2 - CPAP cleaner and sanitizer is capable of getting rid of 99.9 % of germs, molds and bacteria with any need for chemicals and water. Soclean 2 - CPAP cleaner is registered by FDA and has been tested in a laboratory.

Soclean 2 CPAP cleaning device has a front display, indicator lights all throughout the cycle, right or left hose attachment approach and east filter cartridge kits.

Features of Soclean 2 - CPAP sanitizer

Automatic daily CPAP sanitizing solution - So Clean 2 CPAP sanitizing machine automatically sterilizes your CPAP equipment after just a self-explanatory initial setup. You can also sanitize the equipment by hand by just a pressing a button. Soclean CPAP cleaning device does not require a lot of expertise to utilize and is well designed to make sure that the equipment is definitely safe to put into action.

Quick and Convenient to use - As I mentioned before, you don't have to disassemble the parts apart during the sanitizing process because this machine is developed in such a way that you can sanitize hose, reservoir and mask which still attached to your CPAP equipment.

Clearing up your CPAP equipment only takes couple of minutes while using this machine but it's advised you leave the cover which is on the machine shut for about 2 hours for after the initial cycle for correct and efficient sanitation.

Compatibility with almost all CPAP devices - So Clean 2 CPAP cleaner is suitable for nearly all CPAP devices. For those kinds of devices with a heated humidifier, a unique tube insert is included so that the mask and tubes can be cleaned automatically simultaneously with the chambers.

Very Effective eliminating germs and bacteria - Soclean 2 - CPAP cleaning machine can kill 99.9 % of mold, bacteria and virus. This makes it really effective in getting rid of virus, mold and bacteria from your CPAP mask, reservoir and hose.

SoClean 2 Go Portable CPAP Cleaner and Sanitizer machine

SoClean 2 Go basically utilizes the exact same effective sterilizing procedure as the SoClean sanitizer. SoClean2Go represent the finest choice for frequent travelers. You can take your elegant and compact sterilizing device with you anywhere you go. This is also the very first portable CPAP devices sanitizing machine that is portable.

Soclean 2 - cartridge filter kit

SoClean filter kit consists of one replacement filter with O-ring attached and one replacement check-valve for the SoClean -2. The filter and check valve ought to be changed every 6 months. SoClean 2 cartridge filter kit converts activated oxygen back to regular oxygen when it leaves the SoClean-2 chamber. The check-valve avoids water from your CPAP from damaging your CPAP electronic devices.

The Science Behind The Soclean 2 - CPAP cleaning device

Activated oxygen, also referred to as Ozone (O₃), is a safe disinfecting agent that So Clean 2 - CPAP sanitizer utilizes to eliminate CPAP bacteria, mold or viruses inside your CPAP devices. The generator inside the So Clean CPAP cleaner breaks down the chemical bond of common oxygen (O₂) and allows it to recombine into a brand-new molecule that has three atoms of oxygen instead of two. Activated oxygen has a brief lifecycle. After being generated it automatically breaks down back to oxygen within 2 hours.

How Activated Oxygen Works ?

1. Your So Clean 2 - CPAP cleaner and sanitizer creates and pumps the activated oxygen through the supply tube and into your humidifier tank. This cleans up not only the water, but likewise the inner walls of your reservoir.
2. The activated oxygen pass through your CPAP hose As it takes a trip, the activated oxygen reaches every crevice, removing all threatening germs and bacteria.
3. The activated oxygen reaches your CPAP mask. The activated oxygen passes into and out of your mask during this step, cleaning it simply as it did with your tube and tank.
4. Activated oxygen exits the chamber through an unique filter. SoClean filter converts any excess activated oxygen back into the oxygen that we breathe (O₂).

Conclusion

So Clean 2 - CPAP sanitizer is intelligently developed to make sanitizing your CPAP tool easier. It uses safe activated oxygen as opposed to water and chemical and does not need a lot of skills in operating it. Finally, it is convenient and compatible with almost all CPAP machine. Click here for more info how you can try soclean CPAP cleaner and sanitizing device free of charge for 30-days trial risk-free.