

Philips Nightbalance lunoa Positional Sleep Therapy

Philips NightBalance Lunoa is a sleep position therapy alternative for patients with positional obstructive sleep apnea. NightBalance Lunoa is a palm-sized, ultra-thin device that is worn around the chest in a soft, adjustable belt. Over the course of the night, it provides mild vibrations to encourage a shift from sleeping on your back to sleeping on your side. It monitors your sleep positions and changes the strength of the vibrations as required, without disrupting your sleep. All the data is uploaded straight into an online control panel and can be accessible by both patient and physician.



What is positional obstructive sleep apnea?

When sleeping on you back, positional Obstructive Sleep Apnea (POSA) is a specific diagnosis where breathing obstructions can take place more often. People with POSA require a different treatment than those with OSA-- the goal being to manage sleep positions during the night.

Main Features of Philips NightBalance Lunoa

No mask required; small, lightweight, easy to use and worn in a soft belt around your chest

Clinically validated option to CPAP for positional obstructive sleep apnea patients

Records and monitors your response to treatment

Thirty days cash back guarantee

You can now buy using your HSA or FSA

Try the thirty days refund guarantee

What's consisted of inside the package?

- NightBalance Lunoa
- Chest Strap
- Charging/Docking Station with power supply
- Travel Case
- User Guide

Philips Nightbalance lunoa FAQ

Q: I have a pacemaker/ICD can I utilize the Nightbalance?

A: Philips Nightbalance lunoa can not be utilized by patients who have a pacemaker.

Q: I sleep in an upright position or propped up on more than two pillows.Can I utilize the Nightbalance?

A: a. The Philips Nightbalance lunoa can not be utilized by patients that sleep in the upright position or need more than two pillows during sleep.

Since of a medical condition(back/shoulder surgery or osteoarthritis), b.The Philips Nightbalance is not advised for the treatment of patients who are supposed to sleep inthe supine position.

Q: What is the warranty for the Philips Night balance?

A: The Nightbalance lunoa device and docking station are warrantied for 2 years.The warranty will not cover: Any parts that undergo regular wear and tear. This consists of the chest strap and sensor machine batteries.Any damage or malfunction due to incorrect usage, unauthorized modifications or repair work to the machine, or not storing the machine as described in the Guidelines for Use.

Q: How typically do I have to charge the device?

A: The NightBalance Lunoa should be charged every other day for optimal performance. It is suggested that the device is charged every morning after use by merely putting the machine on the docking station. Once the device is positioned on the docking station, all the information will likewise be moved to the patient portal.

Q: What is the adjustment program?

A: The Adaptation Program of the Nightbalance is created to gradually train your body to react to the Sensor Device's vibrations during sleep. The Adaptation Program is nine nights in overall. You need to use the Lunoa for 9 nights before the device provides full positional treatment.

a. Analysis Phase-- no vibrations are given during the very first two nights of the therapy. The Lunoa gathers standard data on your sleep patterns during nights one and two. This information is utilized to tailor the therapy to your individual sleeping habits later on.

Philips NightBalance Lunoa

b. Buildup phase-- the Lunoa begins to provide vibrations. During nights three through 9 of the therapy, the amount of vibrations when laying on you back build up gradually. This permits you to get used to the treatment.

KEEP IN MIND: It is recommended to utilize the Nightbalance every night, or as agreed with your medical expert for best treatment results.



Resources :

<https://s3.amazonaws.com/snoringcures/Philips-Nightbalance-lunoa.html>

<https://storage.googleapis.com/nightbalance/Philips-Nightbalance-lunoa.html>

<https://snoring.blob.core.windows.net/snoringcures/Philips-Nightbalance-lunoa.html>

<https://s3.amazonaws.com/night.balance/Philips-Nightbalance-lunoa.html>

<https://snoring.blob.core.windows.net/nightbalance/Philips-Nightbalance-lunoa.html>