

User Manual

Note: Please use 5V2A adapter ! !

Quick Start



- ① Brightness selection: 4000-7000-10000 Lux
- ② Power: ON/OFF
- ③ Timer: 30-60-90 min to turn off

Intended Applications:

This therapeutic lamp offers rapid and persistent effect on seasonal affective disorder (SAD), and pregnancy depression, etc. It uses a natural light source to create a mellow aura that cheers you up and eases off your everyday discomfort. If you live in a place where sufficient sunlight is a luxury, try this product. Better yet, it is light and portable, so you can use it at anywhere and anytime.

Usage Advisory:

1. Only people OVER 12 years old are suggested to use this product when necessary.
2. This lamp is equipped with a strong light source, please DO NOT directly look at the light at any time, and the lamp should be placed 14 inches away from your eyesight.

Operational instructions:

1. Use buttons 1-3 to turn on the lamp and select your desired brightness and sunshine time to start the therapy.
2. A total 90 mins daily usage consisted of 30 mins at 10000 lux, then 30 mins each at 4000/7000 lux settings, respectively, is suggested to maximize the therapeutic outcome.

Manufacturer's contact information

Name: Shenzhen Leyu Interconnection Industrial Co.,Ltd

Address: Room 13B, Beiligong innovation centre, No 9 yuexin 5th road, Nanshan district, Shenzhen, PRC-518000

Email: amazonblbs@163.com