

Delie®

A photograph of three women wearing elegant, light-colored dresses, posing in front of a draped curtain. One woman is standing in the background, blowing a kiss. Two women are seated in the foreground, smiling at the camera.

Leading brand
IN COLOMBIA

Your second skin

The first in ecological garments

20 YEARS ANNIVERSARY

WE ARE THE #1

SMART SPECIALTY FABRICS

The combination of our fabrics makes the garments keep their shape for longer.

WE ARE PREFERRED BY SPECIALISTS

Our 20-year history makes us allies of the best surgeons in the country.



Viveltex protec
for your skin care

The Lycra with which our girdles are lined contains microcapsules with Viveltex Protec components that provide multiple benefits to your skin



Spirulina plus seaweed

Contribute to the toning and reduction of cellulite, healing, reduction of measurements, hydration of the skin and anti-aging effect.




Vitamin B and E

Promotes recovery from wounds.



Vitamin K

Helps to recover from bruises, irritations, scars and spider veins.



At Delié we deliver well-being to men and women with excellent quality garments.



We are pioneers in ecological girdles.



DELIÉ IS PRESENT IN

More than
20
countries

In Colombia
we have about
own points of sale

60



WEB EXPERIENCE
WWW.DELIE.COM



Girdles lined 100% in lycra. LYCRA®
For greater comfort and softness.

COCOA COLOR


It neutralizes with your skin tone making the girdle the perfect complement to your daily look.

OUTER FABRIC - POWERNET



The Powernet fabric is responsible for providing compression to your body. At Delié we take care of every detail that is why our fabrics are softer, more delicate and meet the highest quality standards.

We are the experts in comfort molding



4 TYPES OF GLUETES

TRADITIONAL
BUTTOCK HEART
ENHANCEMENT
FIT 360

At Delié we value diversity, so you can choose garments that adapt to your needs and type of silhouette
¡In Delié there is place for everyone!

ONE BY ONE



The supplies and textiles of our girdles are developed with an exact color percentage and this results in aesthetic, ammonic and imperceptible garments.



EXPERIENCE

We provide you with a secure and omnichannel service so that you can live an experience of complicity at any time and in any place

We're all Delié

20 years
being your second skin.



Certifications that guarantee the quality of our garments.

We are your accomplices at all times

We are grateful to life and new beginnings, which have given us the opportunity to continue dreaming, believing and growing for 20 years.

It has been two decades now, weaving a story with threads of love, perseverance and passion, together with you and all those who have been part of this dream. For 20 years we have lived experiences and teachings that have helped us to grow, learn and strengthen our desire to move forward, creating, innovating and leading, because we still believe in the beauty of authenticity, the importance of comfort and the strength of self-love.

Thank you for being part of our inspiration, for allowing us to be your accomplices and for believing in Delié.

• Your second skin •



CHARACTERISTICS OF OUR PRODUCTS

Compression level

It is the pressure that the girdle exerts on your body.

1.



High compression: Keeps your body firmer, more controlled and with better definition of the silhouette.



Medium compression: Gives you a moderate level of control, stylize your body providing greater comfort.

Adjustment system

It is the adjustment system with which the girdle adapts to your body

2.



Buckle: As measurements are lost, a perfect fit is achieved. Ideal for: post-surgical in surgeries without horizontal incision, postpartum in natural childbirth, daily or aesthetic use.



Zipper: Allows easy posture and is less noticeable with outer clothing. Ideal for: Surgeries with horizontal incision, post-caesarean section, daily or aesthetic use.



Smooth: You find it in medium compression garments, they are light and more invisible.

Intimate Zone


3.



Side: The hole is covered by a cloth that covers the area, it moves to the side.



Central: They are two overlapping fabrics that run in the central part.



Zipper: It has greater width and is ideal for postoperative or for people who wear the girdle all day.

CHARACTERISTICS OF OUR PRODUCTS

Types of buttocks

We are #1 because we understand that bodies are different, that's why we design 4 types of buttocks.

1.



Traditional buttock: It is ideal for women with a moderate butt who want to enhance buttock in a natural way.

2.




Heart shape buttock: Perfect for women with pronounced curves. Provides subtle enhancement and natural scooping and centering the buttocks.

3.



Ultra-enhancing buttock: Suitable for all body types or silhouette. Recommended for women who want to highlight their curves and see their buttocks more pronounced.

4.



Fit 360: Ideal for women with big buttock as it has more capacity in the gluteal area and more waist definition.

Back and arms coverage

1.



3/4 Sleeve: Recommended for postoperative use after an intervention in arms, back and chest.

2.



3/4 Sleeve and free breasts: Use after arm and back surgery. This girdle allows you to wear the bra of your choice.

3.




Sleeveless and bra: Perfect to gain control over the entire back and bust area. It is ideal to use after a liposculpture or postpartum process to shape your figure more easily to breastfeed your baby.

4.



Sleeveless, free bust: This girdle is recommended to use after liposuction on the back, helping to control the adipose tissue in this area, its free bust allows you to use the bra of your choice.

5.



Upper back load: Provides control to your back after a cosmetic procedure or for daily use.

6.



Low back hip: It is ideal for daily use with all kinds of garments.

7.



Strapless: Perfect to use with necklines or sleeveless, it provides molding to your figure.

Leg coverage

1.



Ankle: It molds your legs in a uniform way, recommended in postoperative. You can use it with leggings, jeans and slim pants.

2.



Knee: This girdle controls and shapes your thighs. Its use is almost imperceptible with clothes. Outside, wear it with leggings, jeans, long skirts and tight pants.

3.



Mid-leg: It is ideal for people looking to shape the hip area, it is perfect to use with pants, dresses and skirts.

4.



Hipster: Recommended to stylize your hips, use it with short garments such as shorts, dresses and skirts.

POST SURGERY

Delii has specialized products to support the healing process at every stage. The VIVELTEX technology along with compression are the perfect mix to reduce tissue inflammation and pain (due to movement) whilst providing support and improving scarring by acting as a 'Second Skin'

First posture medium compression girdles:

- They are recommended for the first weeks of recovery - use 24/7.
- Ideal for sleeping when you are using the high compression girdle in your day to day.

High compression girdles

- After surgery according to your doctors recommendation
- We recommend the use of our cocoa color to hide body fluid stains.
- Always follow the instructions of your doctor, for use of control garments after your surgery.

Girdles with bra included, all in one garment.




- Recommended for breast surgery, its double adjustment system facilitates cleaning of the area and medical controls.
- Recommended for liposuction of the back, arms * and armpit area

References with sleeves

- Double abdominal reinforcement in all our references, to give greater compression of the area.
- All our high compression girdles have double abdominal reinforcement.
- Recommended for upper back, armpit and arm surgery.
- Recommended to conceal sagging arm



MEASUREMENT GUIDE

-  **A: High Inflammation.** This indicator shows the correct size you should use between day 1 and 8 after the procedure
-  **M: Medium Inflammation.** This indicator shows the correct size you should use between day 9 and 30 after the procedure
-  **B: Low Inflammation.** After 30 days post procedure or daily use

[illegible]

10



09084



Central

Color



Size



2XS, XS, S, M, L, XL, 2XL, 3XL
4XL, 5XL

11



09198



Side

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL

12



09354



Zipper

Color Size



2XS, XS, S, M, L, XL, 2XL, 3XL

13



09036



Central


Color Size





2XS, XS, S, M, L, XL, 2XL, 3XL


14




 09022





 Central

Color





Size


2XS, XS, S, M, L, XL, 2XL, 3XL


15




 09103






 Side


Color





Size

2XS, XS, S, M, L, XL, 2XL, 3XL


 09174





 Central

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL




Post-surgical
High compression-knee

16



09347



Side

Color




Size

2XS, XS, S, M, L, XL, 2XL, 3XL



09008



Zipper

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL



4XL, 5XL

Post-surgical
High compression-knee

17






Post-surgical
High compression-knee

18



09360



Central

Color

Size



2XS, XS, S, M, L, XL, 2XL, 3XL



09052



Side

Color

Size



2XS, XS, S, M, L, XL, 2XL, 3XL

Post-surgical
High compression-knee

19



20



21



22



09175



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

09035



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

09047



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL 4XL, 5XL

09021



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL 4XL, 5XL

23



24



09043



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

09020



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

09215



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL 4XL, 5XL

09216



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

25



26



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

27



09111



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL
4XL, 5XL

09048



Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL




Post-surgical
High cheek compression

28



09086



Color


Size

Side



2XS, XS, S, M, L, XL, 2XL, 3XL

09053



Color

Size

Side



2XS, XS, S, M, L, XL, 2XL, 3XL



29



09353



Zipper

Color

Size



2XS, XS, S, M, L, XL, 2XL, 3XL



30



09010



Central



Color

Size

2XS, XS, S, M, L, XL, 2XL, 3XL
4XL, 5XL



31



09361



Zipper



Color

Size

2XS, XS, S, M, L, XL, 2XL, 3XL
4XL, 5XL

32



POST-SURGICAL BRA

Support is synonymous of relief, that is why our bra provides you with all the security so that you can achieve the results after the procedure:

- Elastic materials that adapt to your body
- Support the bust with a slight compression that helps reduce inflammation.
- Full cups for greater coverage and comfort
- Smooth edges to improve the healing of your body.
- Better coverage in the arm pit area.

SUGGESTED SIZE CHART

UNDER BUST MEASUREMENT (in)	27,5-30	30,3-32,7	33-35,4	36-37,7	38-40,5
BRA SIZE	34	36	38	40	42



33



Post-surgical
Bra



09364



Zipper

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL
4XL, 5XL



09027



Color



Size

32, 34, 36, 38, 40



09348



Color



Size

32, 34, 36, 38, 40

34

Post-surgical
Bra



09055



Color



Size

32, 34, 36, 38, 40

35



09349




Color



Size

32, 34, 36, 38, 40



• Internal opening
for prosthetics



09065



Color



Size

32, 34, 36, 38, 40

DAILY USE

ESTHETIC AND POST-SURGICAL

The girdles for daily use are garments that help to stylize and enhance your figure. We offer you different styles according to your needs, by using them daily it helps you:

- Control and define the desired zones.
- Hide sagging and firm the skin.
- Naturally enhances and defines different parts of the body such as: bust, waist and bottoms.
- Improves posture, contributing to the relief of back pain.

TALLA BASE PANTALÓN / GUIDE SIZE BASED ON PANTS																
MUJER / WOMEN					TALLA SEGÚN LA INFLAMACIÓN / SIZE ACCORDING TO INFLAMMATION LEVEL											
	CADERA	CINTURA		HIPS	WAIST	3XS	2XS	XS	S	M	L	XL	2XL	3XL	4XL	5XL
6	87-92	61-65	CENTÍMETROS (cm)	34-36	24-26											
7	93-97	66-70		36-38	26-28											
8	98-102	71-75		38-40	28-30											
10	103-107	76-81		40-42	30-32											
12	108-112	82-86		42-44	32-34											
14	113-117	87-91		44-46	34-36											
16	118-122	92-96		46-48	36-38											
18	123-127	97-101		48-50	38-40											
20	128-132	102-106		50-52	40-42											
22	133-137	107-111	INCHES (in)	52-54	42-44											
24	138-142	112-116		54-56	44-46											

B: Low inflammation.
This indicator gives us the size you should wear after day 30 of surgery for cosmetic procedures or for daily use



Daily use
Medium Compression

37

09075




Central





Color Size
2XS, XS, S, M, L, XL, 2XL, 3XL

38



Daily use esthetic
High compression ultra enhancement



 09231



Color




Size





2XS, XS, S, M, L, XL, 2XL, 3XL

39



Daily use esthetic
High compression ultra enhancement



 09230



Color






Size


2XS, XS, S, M, L, XL, 2XL, 3XL

40




 09112





Side

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL

41



 09066



Side

Color



Size

2XS, XS, S, M, L, XL, 2XL, 3XL

 09046



Central

Color



Size


2XS, XS, S, M, L, XL, 2XL, 3XL

42



43

09106



Side

Color Size



2XS, XS, S, M, L, XL, 2XL, 3XL

09068



Central

Color Size



2XS, XS, S, M, L, XL, 2XL, 3XL





Daily use - esthetic
High compression
Waistbands



45

Daily use-esthetic
High compression waistbands

 09173 Color Size
 2XS, XS, S, M, L, XL, 2XL, 3XL






46



 09006  Color  Size
2XS, XS, S, M, L, XL, 2XL, 3XL

47



 09001    Color  Size
2XS, XS, S, M, L, XL, 2XL, 3XL
Central

Daily use-esthetict
High compression waistbands

Daily esthetic and post-surgical use
Panties, cacheteros and cyclists



48



49



Daily esthetic and post-surgical use
Panties, cacheteros and cyclists

Daily esthetic and post-surgical use
Panties, cacheteros and cyclists



09346



Zipper

Color

Size



2XS, XS, S, M, L, XL, 2XL, 3XL



09350



Zipper

Color

Size



2XS, XS, S, M, L, XL, 2XL, 3XL

50



09197



Central



Color

Size

2XS, XS, S, M, L, XL, 2XL, 3XL

51



09226



Zipper



Color

Size

2XS, XS, S, M, L, XL, 2XL, 3XL




52





53

Daily esthetic and post-surgical use
Panties, cacheteros and cyclists



Daily esthetic and post-surgical use
Panties, cacheteros and cyclists

 09227





Zipper


Color


 

Size



2XS, XS, S, M, L, XL, 2XL, 3XL

 09116




Central

Color

Size

2XS, XS, S, M, L, XL, 2XL, 3XL


ESSENTIAL BRA

The perfect ally for stylize your figure, enhance your bust and shape your outline.




For confident, classic, sensual and romantic women, our Essential Bra is perfect for all types of breasts, their straps adjustable allow you to adjust it to the perfect size, its wide base gives you greater control and support, preventing them from stubborn fat in the armpits and back. The soft padding provides comfort and the side bars give greater compression perfectly complementing the use of the girdle.

54



 15078

Color

Size

32, 34, 36, 38, 40



Eco
control
eco girdles

56



05021



Zipper



Color

Size


XS, S, M, L, XL, 2XL

* Cocoa color coming soon

57



05066



Zipper



Color

Size

XS, S, M, L, XL, 2XL

* Cocoa color coming soon

MATERNITY

During your pregnancy process we have two ideal products to give you greater comfort, always taking care of the well-being of your baby. Our trusses and support band are designed to support the weight of the belly without generating compression in that area and alleviating back discomfort.



Color Size
09044 XS, S, M, L, XL, 2XL, 3XL

Maternal Support Band



Color Size
09051 XS, S, M, L, XL, 2XL



NURSING BRAS



60

Color Size

09232 32, 34, 36, 38, 40

NURSING BRAS



61

Color Size

09057 32, 34, 36, 38, 40



Color Size

09234 32, 34, 36, 38, 40



Color Size

09233 32, 34, 36, 38, 40

POSTPARTUM


We design garments thought of your comfort and well-being to accompany your maternity process.



62

RECOMMENDED FOR POSTPARTUM

All of our high compression girdles are suitable for your postpartum process, keep in mind that if your delivery was by cesarean section we recommend girdles with a side zipper, if it was natural childbirth references with a central clasp are recommended.



63

More options for you,
choose your
favorite color



—Beige Girdles—


NEW COLOR



09066




09360



09043



09197



09355

BLACK GIRDLES



09047



09360



09230



09226



09116



09197



09215



09348



09355



09021



09231



09346



09363



09066



09227



09350



09362




Enjoy your curves
without
measurements



PLUS SIZES

FROM SIZE 4XL





Control garments

Male

ffm
for men

TALLA BASE PANTALÓN / GUIDE SIZE BASED ON PANTS																		
HOMBRE / MEN				TALLA SEGÚN LA INFLAMACIÓN / SIZE ACCORDING TO INFLAMMATION LEVEL														
	PECHO	CINTURA	CADERA		CHEST	HIPS	WAIST	3XS	2XS	XS	S	M	L	XL	2XL	3XL	4XL	5XL
30	CENTÍMETROS (cm)	84-89	74-79	89-94	INCHES (in)	33-35	29-31	35-37										
32		90-94	80-84	95-99		35-37	31-33	37-39										
34		95-99	85-89	100-104		37-39	33-35	39-41										
36		100-104	90-94	105-109		39-41	35-37	41-43										
38		105-109	95-99	110-114		41-43	37-39	43-45										
40		110-114	100-104	115-119		43-45	39-41	45-47										
42		115-119	105-109	120-124		45-47	41-43	47-49										
44		120-124	110-114	125-129		47-49	43-45	49-51										
46		125-129	115-119	130-134		49-51	45-47	51-53										
48		130-134	120-124	135-139		51-53	47-49	53-55										
50	135-139	125-129	140-144	53-55	49-51	55-57												

- A: High Inflammation.** This indicator shows the correct size you should use between day 1 and 8 after the procedure
- M: Medium Inflammation.** This indicator shows the correct size you should use between day 9 and 30 after the procedure
- B: Low Inflammation.** After 30 days post procedure or daily use

We have designed special products for your needs, understanding that you always want to look good, feeling comfortable and safe at all times. In our men's line we have, post-surgical, everyday shirts ideal to use as a complement to your outerwear, boxers with medium and high control with different coverings.

At all times our garments are functional and with internal structures that control and give lumbar support.

Learn more in www.delié.com

72



POST-SURGICAL WITH SLEEVES



09018



Color

Size

XS, S, M, L, XL, 2XL, 3XL

73



POST-SURGICAL WITHOUT SLEEVES



09016



Color

Size

XS, S, M, L, XL, 2XL, 3XL

74



VEST GIRDLE


09017





Color
Size


XS, S, M, L, XL, 2XL, 3XL

75



ONE PIECE GIRDLE


09238




Color
Size



XS, S, M, L, XL, 2XL, 3XL

76



DAILY USE VEST

09237

Color: ☐ ☒

Size: XS, S, M, L, XL, 2XL, 3XL



DAILY USE SHIRT

DAILY USE TANK SHIRT V NECK

09015

Color: ☐ ☐ ☒

Size: XS, S, M, L, XL, 2XL, 3XL

09236

Color: ☐ ☒

Size: XS, S, M, L, XL, 2XL, 3XL

78



79



BOXER

 09079

* Only on request



Color

Size

XS, S, M, L, XL, 2XL, 3XL

 09235



Color

Size

XS, S, M, L, XL, 2XL, 3XL



 15092




Color

Size

XS, S, M, L, XL, 2XL, 3XL



COMPLEMENTARY



09037
Color Size
S, M, L, XL

81

STABILIZER BAND

BRA STRAP PROTECTORS


09193

One size

REMOVABLE STRAP

09083

Color One size



SPLIT HOOD

09029 XS, S, M, L, XL

09031 XS, S, M, L, XL



SLEEVES



09034

Color



Size

XS, S, M, L, XL



BOARDS:

Our line of boards helps you control and accentuate the figure of your body to the maximum, it also prevents breakage in the girdle and contributes to having a better posture, that is why we recommend them for post-operative or daily use.



09109

Color



One size



Size: 30*20 cm
Rigid Foam





FIRST TIME USE CIRCULAR PAD

 09121
Color Size
 S,M,L
Soft foam and acrylic



SHAPING PEAR PAD

 09101
Color One size
 Size: 23*10 cm
Rigid foam

88



PEAR SHAPED PAD



09146

Color One size



Size: 21,5*25 cm
Soft foam and acrylic



89



FIT BACK PAD



09357

Color One size



Size: 21,5*25 cm
Soft foam and acrylic



90



09356

Color One size



91

MARKING FOAMS

09358

We want to accompany
you at all times, that is
why we evolve with you.
Giving you the security
and confidence you need.

Discover the Delié
universe.

Delié®






In your sport routines



Feel
good
feel
Delié



In your intimacy

Fashion
with invisible
shaping

In your day to day



Swimwear with
UV protection

On your vacation

Delié®



SPORTS LATEX BELT

17205 Color Size XS, S, M, L, XL, 2XL, 3XL



LATEX SPORTS BELT

 17206  Color Size
XS, S, M, L, XL, 2XL, 3XL

INDEX

Country code



International

Color code



CF NG BL BG

How to request your reference

Country + Reference + Color + Size

Ejm/Eg: 09047CFXS

We are your accomplices at all time	Pág 2	09185	Pág 26
We are the # 1	Pág 4	09186	Pág 21
Characteristics of our products	Pág 5	09193	Pág 82
Post surgery	Pág 8	09197	Pág 50
Measurement table	Pág 9	09198	Pág 11
		09215	Pág 25
		09216	Pág 25
		09226	Pág 51
		09227	Pág 52
		09230	Pág 39
		09231	Pág 38
		09232	Pág 61
		09233	Pág 60
		09234	Pág 61
		09235	Pág 79
		09236	Pág 77
		09237	Pág 76
		09238	Pág 75
		09346	Pág 48
		09347	Pág 16
		09348	Pág 34
		09349	Pág 34
		09350	Pág 49
		09353	Pág 29
		09354	Pág 12
		09355	Pág 20
		09356	Pág 90
		09357	Pág 89
		09358	Pág 91
		09360	Pág 18
		09361	Pág 31
		09364	Pág 32
		15078	Pág 54
		15092	Pág 79
		17205	Pág 93
		17206	Pág 94
		Post-surgical	Pág 6
		Post-surgical bra	Pág 31
		Daily use esthetic and post-surgical	Pág 34
		Daily use - esthetic waistbands	Pág 37
		Essential bra	Pág 54
		Eco garments	Pág 55
		Maternity	Pág 58
		Postpartum	Pág 62
		Recommended for postpartum	Pág 63
		Black girdles	Pág 66-67
		Beige girdles	Pág 65
		Plus sizes	Pág 69
		Control garments Male	Pág 68
		Complementary	Pág 70
		Delié	Pág 93

ffm

Delié for men



For the man who likes
to look and feel good