



# Combining PET/CT and RTP

David W. Townsend, PhD

*Departments of Medicine and Radiology  
University of Tennessee, Knoxville, TN*



**Advances in PET-SPECT/CT and RTP**

*Melbourne, Australia December 13th 2005*

# Topics

---

- PET/CT in Radiation Therapy Planning
- Monitoring response to therapy
- Imaging biomarkers: beyond FDG

# PET/CT for image-guided radiation therapy


---

- **define the appropriate target volume**
  - functional changes precede, and do not necessarily match, structural changes
- **define the optimal dose to deliver**
  - PET can quantify basic tumor biological properties that can be used for determining target volumes and dose
- **follow the response to radiation therapy**
  - the same properties can be followed during and after therapy to assess treatment response


# Integrating PET/CT in therapy planning





- RT flat table top
- patient immobilization device
- laser positioning system
- 70 cm patient port




PET/CT scan



CT simulation, contouring




Treatment plan




CT simulator

Fusion



# Defining target volumes with molecular imaging

---



PTV: Planning target volume	physicist
CTV: Clinical target volume	clinician
GTV: Gross tumor volume	CT / MRI
BTV: Biological tumor volume	PET or SPEC

PTV: tumor motion and setup errors

CTV: includes microscopic disease

GTV: from high resolution anatomy

BTV: from molecular imaging

# Issue #1



---

- Defining the Biological Tumor Volume

# Defining BTV based on Standardized Uptake Value

$$\text{SUV} = \frac{a_{\text{ROI}}(\mu\text{Ci/ml})}{A(\text{mCi}) / W(\text{kg})}$$


- Choice of threshold (T) affects BTV
- T varies with size and background



- SUV<sub>max</sub> may depend on B but not T
- SUV<sub>mean</sub> and contour (BTV) depend on T
- SUV depends on BMI

# Factors affecting apparent tumor $^{18}\text{F}$ FDG uptake

---





- lesion size and partial volume effects
  - tumor heterogeneity
  - cellular avidity for glucose
  - image reconstruction parameters
  - definition of region-of-interest
  - patient blood glucose levels
  - time after injection of FDG
  - difference between FDG and glucose (LC)
- **These factors will also affect the SUV**

# PET/CT in RTP

*Courtesy Dr Sara Grehl, University Hospital Essen*

58 year old male with oro/naso-pharyngeal cancer crossing the midline (T4) with unilateral lymph node involvement on CT. Left lymph node negative on CT but positive on PET/CT (N2). Modified boost-PTV (70 Gy) includes left lymph node.





# Issue #2

---

- Monitoring respiration (lung cancer)

# PET respiratory gating




# Results

- 9 out of 22 lung lesion patients exhibited lesions with significant motion
- 6 of those with sufficient statistics and regular breathing pattern to evaluate 4-gate studies


Gates	2 (11 lesions)	4 (8 lesions)
Average increase [SUV]	<b>14%</b>	<b>17%</b>
Max increase [SUV]	65%	39%
Duty cycle per gate	15-80% typical: 30-60%	10-50% typical: 20-30%
SUV range	2.1 – 35	2.3 - 35



CT plan



PET/CT plan



Dose-Volume Histogram


46 year-old female with NSCLC: cT2pN2 diagnosed 01/04 post 3 x CTX; PET/CT for RTP.

Compared to pre-CTX, the bronchial tumor is reduced both in size and metabolic activity ( $SUV_{max}$  reduced from 14 to 4.6). Atelectasis was sufficiently differentiated on PET/CT. PET/CT-based conformal treatment plan yields a cumulative mean lung dose of 8 Gy compared to 15 Gy with an assumed CT-based conventional treatment plan. Thus, **PET/CT allowed for dose escalated treatment of this NSCLC.**





BGO-PET: 320 MBq <sup>18</sup>F-FDG, 80 min pi, 3.5 min per bed, 3 beds

Dual-row CT: 110 mAs, 130 kVp, IV and water-based oral contrast, 4 mm slices

# Change of treatment volume



↓ Increase in PTV ↓



PTV from CT

PTV from PET/CT

Decrease in treatment volume

*Erdi et al:* 11 patients: 7/11 PTV ↑ (19% ave)  
4/11 PTV ↓ (18% ave)

*Schmucking et al:*  
27 patients: 2/27 PTV ↑  
25/27 PTV ↓ (12% ave)

# PET/CT in radiation therapy

---

Total of 39 patients for radiation treatment

	head and neck	lung	pelvis
GTV $\uparrow$ > 25%	17%	17%	33%
GTV $\downarrow$ > 25%	33%	67%	19%

- GTV changed in 56% of cases
- PTV change >20% in 46% of cases
- decreased GTV differences between oncologists

# Response to therapy

---

## Assessment measures:

- Kinetic analysis from dynamic scan: CMRGlu
- Patlak-Gjedde analysis:  $K_i$ , net rate of FDG trapping
- Standardized Uptake Value (SUV)
- Visual inspection

## Issues for assessing response:

- requires at least 20% change in SUV between studies
- valid for tumors with sufficient baseline metabolic activity
- wait 4-6 weeks after completion of radiation therapy
- radiation-induced inflammatory response or tumor

# Response to therapy: Lymphoma


*biograph 16*



Cancer Imaging and  
Tracer Developer



PET/CT fused image



**Pre-therapy: 9/04**

203 lbs  
12.3 mCi  
90 min uptake  
120 kVp; 159 mAs  
7 beds; 4 min per bed



**Post-therapy: 1/05**



219 lbs  
11.7 mCi  
98 min uptake  
120 kVp; 193 mAs  
7 beds; 4 min per bed

31 year-old male with lymphoma. Status pre and post chemotherapy

# Response to therapy: head and neck *biograph 16*



Cancer Imaging and  
Tracer Developer




Multiple 2-8 mm  
hyper-metabolic  
lymph nodes



1 month post-  
chemotherapy  
and radiation

# Quantitation of response to therapy: size and SUV



First cycle of  
chemotherapy

Second cycle of  
chemotherapy

Week 0

Week 1

Week 2




Week 3

Week 4

Week 5

Week 6

Week 10




## Summary: response to therapy

Tumor	No. of Patients	Criterion (% SUV)	Median survival (months)	
			Responders	Non responders
Lymphoma	30	Visual	>24	5
Esophagus	37	35	>48	20
	22	30	>38	18
Stomach	35	35	>48	17
Head and Neck	47	50	>120	40
Lung	57	20	9	5

## Combined PET/CT imaging for RTP: outcome

---



The use of PET(/CT) can result in a change of treatment plan in more than 50% of cases.

*Does this change improve outcome?*

*Grosu AL, Weber WA, Franz M, et al. Re-irradiation of recurrent high-grade gliomas using amino acid PET(SPECT)/CT/MRI image fusion to determine gross tumor volume for stereotactic fractionated radiotherapy. Int J Radiat Oncol Biol 2005 (In press)*

$^{11}\text{C}$ -methionine PET with CT/MRI fusion  
Defining GTV with molecular imaging improved mean survival from 5 months to 9 months.

# Issue #3

---


- Automated synthesis of radiotracers

# PET tracers for tumor imaging


---

- blood flow [15O] water
- glucose transport/utilization [18F] FDG
- tumor hypoxia [64Cu] ATSM
- angiogenesis [18F] RGD
- amino acid synthesis [11C] methionine
- cell replication (analog) [18F] FLT
- apoptosis [18F] annexin-V
- somatostatin receptors [68Ga] DOTATOC
- chemotherapy drugs [18F] fluorouracil
- prostate cancer [18F] fluorocholeline


# Hypoxia ( $^{60}\text{Cu}$ -ATSM)-guided RTP



Head phantom and face mask showing fiducial markers for CT and PET registration; with a PET/CT scanner, registration would be automatic.







hypoxic tumor volume ( $h\text{GTV}$ )



GTV and  $h\text{GTV}$

# PET/CT imaging with $^{68}\text{Ga}$ -DOTATOC



Imaging tumors that express somatostatin receptors

*M Hofmann and A. Bockisch*  
*University of Bern and University of Essen*

130 kV, 100 mAs, 5 mm  
6 min/bed, 8 beds, 6i/16s, 2 mm, 256

# Cancer Imaging and Tracer Development Research Program

*University of Tennessee, Knoxville*



*Cancer Imaging and  
Tracer Development*

## Acknowledgements:

Claude Nahmias, PhD

Karl Hubner, MD

Meixiang Yu, PhD

Fang Gao, PhD

Xukui Wang, MD

Bjoern Jakoby, MS

Misty Long, R.T. (R) (N)

Chris Carr, R.T. (R) CNMT

James Avenell, BS LVMT

Linda Paschal, BS

Khanh Nguyen, MD

*Jonathan Carney, PhD*

**National Cancer Institute**



- Sensation 16 CT
- 3D LSO PET
- PICO-3D electronics
- HI-REZ PET detectors