

PUT A LID ON IT

Mucus Traps Leaking During Transport?

If you are having a problem with accidental tube separation during specimen transport, try Busse's 40cc and 80cc Mucus Specimen Traps packaged with an *extra cap*. The Secured Specimen Transport Caps are designed to seal the trap for pneumatic tube transport. After collection, simply remove the standard tubed cap and replace it with the Transport Cap.

They work with any suction device, have a removable adaptor to prevent contamination, and come with collapse-proof, heavy-walled tubing that can be connected to any catheter for suctioning. Both the Transport Cap and the standard cap screw on the container to help maintain a sterile pouring surface. And the free-standing, see-through plastic container is clearly marked with easy-to-read 1cc graduations and patient ID area.



LATEX FREE

Ordering Information

| | | |
|----------------------------|---|--------------------------------|
| Cat. No. 405 | Specimen Trap with Transport Cap, 40cc | 24/ carton |
| Cat. No. 406 | Specimen Trap, 40cc | 24/ carton |
| Cat. No. 407 | Specimen Trap, 80cc | 100/ carton (4 boxes of 25) |
| NEW Cat. No. 408 | Specimen Trap with Transport Cap, 80cc | 50/ carton |

busse
Hospital Disposables

P.O. Box 11067, Hauppauge, NY 11788-0920
(800) 645-6526 • FAX: (631) 435-4721
www.busseinc.com • mail@busseinc.com

©2011, Busse Hospital Disposables, Hauppauge, NY