

NOTICE TO HOKANSON NIVP3 CUSTOMERS

Summary

In January 2020, Microsoft[®] will end extended support for their Windows[®] 7 operating system. After this date, the operating system will no longer receive free technical support, or software and security updates.

Impact on Hokanson NIVP3 Software

The last version of Hokanson NIVP3 Noninvasive Vascular Program software, version 5.40, only runs on the Windows[®] 7, 32-bit operating system. NIVP3 software will not be upgraded to run on newer operating systems.

Systems running Windows 7, including NIVP3 software will continue to work at their current capacity even after Microsoft support ends. However, Hokanson's ability to support NIVP3 software will be limited, due to the end of Microsoft support for the operating system.

Users should exercise caution to prevent the computer running NIVP3 software from malware and security threats according to your facility's computer and network security policies.

December 2019