

THE BD PORTFOLIO OF IV TECHNOLOGIES — Product Specifications

BD Insyte™ Autoguard™ BC Shielded IV Catheter With Blood Control Technology

BD REF NO. NON-WINGED	BD REF NO. WINGED	GAUGE	CATHETER LENGTH (in.)	CATHETER OD (mm)	CATHETER ID (mm)	UNITS BOX/CASE	FLOW RATE (mL/min)
382512	382612	24	0.75	0.66-0.74	0.48-0.56	50/200	20
382523	382623	22	1.00	0.84-0.91	0.61-0.69	50/200	37
382533	382633	20	1.00	1.04-1.12	0.76-0.84	50/200	63
382534	382634	20	1.16	1.04-1.12	0.76-0.84	50/200	61
382537	382637	20	1.88	1.04-1.12	0.76-0.84	50/200	54
382544	382644	18	1.16	1.27-1.35	0.94-1.02	50/200	95
382547	382647	18	1.88	1.27-1.35	0.94-1.02	50/200	87
382554	382654	16	1.16	1.70-1.78	1.32-1.40	50/200	193
382557	382657	16	1.77	1.70-1.78	1.32-1.40	50/200	185

The 18- to 22-gauge catheters are capable of withstanding high-pressure injection procedures.

BD Nexiva™ Closed IV Catheter System—Single Port

BD REF NO.	GAUGE	CATHETER LENGTH (in.)	CATHETER OD (mm)*	CATHETER ID (mm)*	EXTENSION TUBE ID (mm)	UNITS BOX/ CASE	FLOW RATE H ₂ O (mL/min)
383510	24	0.56	0.71	0.53	1.22	20/80	19
383511	24	0.75	0.71	0.53	1.22	20/80	18
383512	22	1.00	0.90	0.67	1.22	20/80	33
383516	20	1.00	1.10	0.83	1.65	20/80	61
383517	20	1.25	1.10	0.83	1.65	20/80	58
383518	20	1.75	1.10	0.83	1.65	20/80	51
383519	18	1.25	1.31	0.98	1.65	20/80	84
383520	18	1.75	1.31	0.98	1.65	20/80	79

For power injection use, refer to instructions for use. *Average dimension.

BD Nexiva™ Closed IV Catheter System—Dual Port

BD REF NO.	GAUGE	CATHETER LENGTH (in.)	CATHETER OD (mm)*	CATHETER ID (mm)*	EXTENSION TUBE ID (mm)	UNITS BOX/ CASE	FLOW RATE H ₂ O (mL/min)
383530	24	0.56	0.71	0.53	1.22	20/80	19
383531	24	0.75	0.71	0.53	1.22	20/80	18
383532	22	1.00	0.90	0.67	1.22	20/80	33
383536	20	1.00	1.10	0.83	1.65	20/80	61
383537	20	1.25	1.10	0.83	1.65	20/80	58
383538	20	1.75	1.10	0.83	1.65	20/80	51
383539	18	1.25	1.31	0.98	1.65	20/80	84
383540	18	1.75	1.31	0.98	1.65	20/80	79

For power injection use, refer to instructions for use. *Average dimension.

BD Nexiva™ Diffusics™ Closed IV Catheter System

BD REF NO.	GAUGE	CATHETER LENGTH (in.)	EXTENSION TUBE ID (mm)	UNITS BOX/CASE	GRAVITY FLOW RATE (mL/min)	MAX CT FLOW RATE (mL/sec)*
383590	24	0.75	1.22	20/80	20	3.0
383591	22	1.00	1.65	20/80	44	6.5
383592	20	1.00	1.65	20/80	67	10.0
383593	20	1.25	1.65	20/80	64	10.0
383594	18	1.25	1.65	20/80	90	15.0

*See instructions for use for all power injection guidelines, cautions and warnings.

THE BD PORTFOLIO OF IV TECHNOLOGIES — Product Specifications

BD PosiFlush™ SureScrub™ Pre-Filled Saline Syringe

BD REF NO.	DESCRIPTION	CURRENT NDC NO.	NEW NHRIC NO.	UNITS BOX/CASE
306559	10 mL		8290-306559	30/480

BD PosiFlush™ Normal Saline Syringes, All With 10 mL Diameter

306544	3 mL	08290-0910-03	8290-306544	30/480
306545	5 mL	08290-0920-05	8290-306545	30/480
306546	10 mL	08290-0940-10	8290-306546	30/480
306547	10 mL, Generally Pump Compatible		8290-306547	30/480
306549	3 mL with Blunt Plastic Cannula	08290-0911-03	8290-306549	30/120
306550	5 mL with Blunt Plastic Cannula	08290-0921-05	8290-306550	30/120
306551	10 mL Standard Plunger Rod with Blunt Plastic Cannula	08290-0941-10	8290-306551	30/120

BD PosiFlush™ SF Saline Syringe

306553	10 mL Sterile Field Syringe	08290-0950-10	8290-306553	30/240
--------	-----------------------------	---------------	-------------	--------

BD PosiFlush™ Heparin Lock Flush Syringes, 10 USP Units/mL

306413	3 mL		8290-306413	30/480
306414	5 mL		8290-306414	30/480
306509	6 mL Fill in 10 mL Syringe	08290-0360-06	8290-306509	30/120
306510	5 mL Fill in 10 mL Syringe	08290-0360-05	8290-306510	30/120
306511	5 mL Fill in 5 mL Syringe	08290-0350-05	8290-306511	30/120
306512	3 mL Fill in 3 mL Syringe	08290-0340-03	8290-306512	30/120
306521	3 mL Fill in 10 mL Syringe	08290-0360-03	8290-306521	30/120
306525	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	08290-0361-05	8290-306525	30/120
306527	2 mL Fill in 3 mL Syringe	08290-0340-02	8290-306527	30/120

BD PosiFlush™ Heparin Lock Flush Syringes, 100 USP Units/mL

306423	3 mL		8290-306423	30/480
306424	5 mL		8290-306424	30/480
306513	3 mL Fill in 3 mL Syringe	08290-0390-05	8290-306513	30/120
306514	5 mL Fill in 10 mL Syringe	08290-0390-03	8290-306514	30/120
306515	5 mL Fill in 5 mL Syringe	08290-0380-05	8290-306515	30/120
306516	3 mL Fill in 5 mL Syringe	08290-0380-03	8290-306516	30/120
306517	3 mL Fill in 3 mL Syringe	08290-0370-03	8290-306517	30/120
306528	2 mL Fill in 3 mL Syringe	08290-0370-02	8290-306528	30/120
306531	5 mL Fill in 10 mL Syringe with Blunt Plastic Cannula	08290-0391-05	8290-306531	30/120

BD™ SASH Kits: 2 BD Pre-Filled Normal Saline Syringes, 1 Pre-Filled Heparin Lock Flush Syringe

306537	5 mL Fill in 10 mL Syringe	08290-0420-05	8290-306537	10/30
--------	----------------------------	---------------	-------------	-------

Products are not made with latex or preservatives.

For ordering information, consult
your local BD sales representative
or call (888) 237-2762.



BD Medical
9450 South State Street
Sandy, Utah 84070
(888) 237-2762
www.bd.com/ivcaths