

Keep safety injections short & simple



BD SafetyGlide™ Insulin Syringe 6mm Needle

The only safety insulin syringe meeting diabetes needle length recommendations¹

A shorter 6mm needle length may help reduce the chance of patient intramuscular injections – a potential risk for hypoglycemia.

One-handed safety mechanism to help prevent accidental needle sticks.

Improved with a larger pad on the pushrod.

½ unit dosing scale on 0.3mL barrel size

Improve confidence in small dose delivery.

Innovation with no price impact:

Pricing is the same as current
BD SafetyGlide™ insulin syringe.



Offered in 3 barrel sizes

1mL 31G x 6mm (15/64")
NDC/HRI# 08290 (328446)

0.5mL 31G x 6mm (15/64")
NDC/HRI# 08290 (328447)

0.3mL with ½ unit scale 31G x 6mm (15/64")
NDC/HRI# 08290 (328449)

1. Needles < 8 mm are preferable to decrease the chance of intramuscular injection. "Teaching Injection Technique to People with Diabetes." American Association of Diabetes Educators. Diabeteseducator.org, 1 Jan. 2013. Web. 1 Jan. 2015. See also Frid A, Hirsch L, et al. Diabetes Metab, 2010; 6 (suppl): S2-S18.



Helping all people
live healthy lives