

UNIVERSAL DRIVERS

- Long and short-handled drivers work with all major implant systems
- Compatible with 95% of all implant systems sold in the United States
- Long-handled drivers provide outstanding access between teeth and are specifically designed for maxillary cases
- Short drivers are designed for anterior mandible cases
- Drivers are color-coded for easy identification
- Fully autoclavable
- Sold individually and in a set of 7
- Short-handled drivers are 1.75 inches long
- Long-handled drivers are 4.25 inches long



K539

Long-Handled Drivers
Set of 7

Individual Drivers,
Long 4.25"



K540

Short-Handled Drivers
Set of 7

Individual Drivers,
Short 1.75"

TIPS

SHORT

LONG

	TIPS	SHORT	LONG
—	Straight Screwdriver Straight Systems	ID535	ID530
★	SCS Screwdriver Straumann	ID536	ID533
★	UniGrip Driver Replace Select, Implant Direct	ID534	ID527
⬠	Hex 1.25mm/.050 in Lifecore, Biohorizons, Astra, Zimmer, Cortex, Implant Direct, MIS, Sybron, Camlog, Implant Logistics	ID543	ID538
⬠	Hex 1.3mm/.050 in Lifecore, Biohorizons, Astra, Zimmer, Cortex, Implant Direct, MIS, Sybron, Camlog, Implant Logistics	ID524	ID531
⬠	Hex 1.2mm/.048 in 3i System, Megagen	ID537	ID532
⬠	Hex 0.9mm/.035 in Friadent, Ankylos, 3i Cover Screw	ID529	ID528
■	Square 1.3mm/.050 in Gold Screw/Nobel Biocare	ID525	ID526

