





THE SHARPENATOR

Take Scaling to New Heights

- Efficiently sharpens instruments without the need of two stones
- Coarse on one side and medium on the other
- The coarse side quickly re-shapes your blade while the medium side delivers a precision edge
- Made from a synthetic material, the stone is autoclavable, non-loading, requiring no oil
- Perfect for chair-side use
- Takes up no counter space



Disposable Chair Side Sharpening Guide



- Face of the blade parallel to floor
- Stone always at the 110° position
- Toe pointing towards you

K-0300

The Sharpinator Kit

Includes one of each:

Medium (white)/ Coarse (pink) Sharpinator Stone, Ergonomic Test Stick, Sharpinator, Chair-side Disposable Sharpening Guide*

