

Infiltration Pump B89055



VITRUVIAN™ Infiltration Pump B89055

The Black & Black Vitruvian™ Infiltration Pump is a peristaltic pump using three rollers creating pumping action to deliver tumescent fluid to the operative site. Powered by a stepper motor, this machine provides consistent pumping action at the desired rate of your choice. Its dual voltage power supply can be used with 100V/240V electrical service.

Our Infiltration Pump is designed with a rate-adjusting knob, prime button, and reverse switch for ease of use and optimal user settings. Large, bright LEDs clearly display both volume and rate. The rate-adjusting knob can be adjusted 0 to 525 mL/min.

For best performance use Black & Black Softouch Pump Tubing (B89213).



- Peristaltic Pump with three rollers
- Momentary/maintained pneumatic foot-switch functions for ON/OFF action or continuous action, respectively.
- Dual voltage motors 100V/240V
- Bright, easy-to-read display screen
- Volume [mL] and Rate [mL/min] displayed
- Reverse switch for flow direction preference
- Prime button for fast filling of tubing line

B89055	VITRUVIAN™ INFILTRATION PUMP
B89213	13' SOFTOUCH PUMP TUBING SET-SINGLE SPIKE, BOX OF 10



5175 South Royal Atlanta Dr.
Tucker, GA 30084
770.414.4880 phone
1.877.252.2517 toll-free
770.414.4879 fax
info@blackandblacksurgical.com
www.blackandblacksurgical.com