

Bledsoe Achilles Boot



**Bledsoe
Brace
Systems™**



Ultra-breathable
foam wrap for
patient comfort

Strong aluminum alloy
provides strength, yet
allows formability to
change widths

- **Inter-locking wedges:**
Four cast inter-locking 10° urethane wedges provide fast change of foot angle.
- **No casting needed:**
Eliminates casting series required to stretch Achilles tendon during rehab.
- **Reformable:** Aluminum boot shell and upright system that, unlike plastic boots, can be reformed to fit virtually any leg or ankle shape.
- **Unique rocker design:**
Allows patients to closely simulate a near normal gait without abnormal loads on any joints.
- **Support and control:**
Contains three layers of pads that contour to the ankle shape and support the ankle for total control.
- **Availability:**
Are also available in Lo-Top and PTB (patella tendon bearing) versions.
- **Optional Air Ankle / Heel Pad:** Allows for ability to adjust ankle pressure for increased comfort and provides pneumatic ankle support.

Save time by
eliminating
series casting
during rehab



Bledsoe Achilles Boot

Restrictions on use

This device is offered for sale by or on the order of a physician or other qualified medical professional. This device is not intended for reuse on a second patient or for extended use beyond 180 days. This device is generally not intended for patients in excess of 250 lbs. (114 kg.) This device is not intended to prevent injuries or reduce or eliminate the risk of injury.

Warning

Do not leave this device in the trunk of a car on a hot day or in any place where the temperatures may exceed 150°F (65°C) for any extended period of time. If any additional pain or symptoms occur while using this device, discontinue use and seek medical attention.



**Bledsoe
Brace
Systems™**

by Medical Technology, Inc.

2601 Pinewood Drive
Grand Prairie, TX 75051 USA

1-888-BLEDSOE
(1-888-253-3763)

1-800-527-3666

972-647-0884

Local Fax:
972-660-5495

International Fax:
972-606-0649

www.bledsoebrace.com

Indications

Indicated for immediate application following injury to or reconstruction of the achilles tendon.

Contraindications

Contraindicated for unstable fractures or for fractures of the proximal tibia or fibula.

Bledsoe Achilles Boot vs. other methods

Features	Bledsoe Achilles Boot	Other Methods
Skin Comfort	Open-cell, breathable foam laminate increases comfort and reduces possibility of skin problems.	N/A
Rocker Bottom	Computer-designed rocker tread allows a near normal walking gait & low heel prevents pelvic tilt.	N/A
Security & Control	Straps engage and lock to brace uprights holding the leg steady.	N/A
Shell & Uprights	Strong aluminum alloy provides strength yet allows formability to change widths.	N/A
Control	Boot has 4, 10 degree removable wedges allowing the foot to be plantar flexed while maintaining good ground contact.	N/A

Ordering Information Achilles Boot

Note: Leg measure taken from knee joint line to medial malleolus.

Shoe Size			Leg Measure	Boot	
Men	Women	Euro	Length	Size	Number
2.5-5.5	3.5-6	33-37	11"-13" (27.5 - 32.5 cm)	S	AL?36013
5.5-9	6-9.5	37-41	13"-15" (32.5 - 37.5 cm)	M	AL?36015
9-12.5	9.5-13	41-46	15"-17" (37.5 - 42.5 cm)	L	AL?36017
12.5-17	13 +	46-49 +	17"-19" (42.5 - 47.5 cm)	XL	AL?36019

Note: When ordering replace the ? with one of the following: O= ankle / heel pad 2 = air ankle / heel pad

Accessories: Rain, Cold Weather, or Night Cover -Replacement cuffs, plates, strap kits, sock, and tibial plate.
Call for more information.

Bledsoe limited warranty

This device is warranted to be free from defects in material and workmanship for a period of 90 days. These warranties apply to devices that have not been modified or subjected to misuse, abuse or neglect.



For more detailed information on this brace including fitting and adjustment instructions, visit our Web site: www.bledsoebrace.com

For product information or questions pertaining to sales or service, please contact the sales representative in your area or Bledsoe Brace Systems.

