

SoClean 2 CPAP cleaning machine



Simple one-time setup

The SoClean CPAP cleaning unit connects with all types of PAP equipment including CPAP, Bi-PAP, BPAP, and VPAP.

1

2 Close the lid and go!

The SoClean sanitizer starts automatically after you place your CPAP mask in the chamber and shut the lid.

3

Kills CPAP germs in minutes


Activated Oxygen cleaning powers through your CPAP reservoir, hose, and mask.

Try it Risk Free

Fits all popular mask types and CPAP machines

Don't use a humidifier? No problem you can still use the SoClean CPAP cleaner.

Try Now Free\Risk-Free for 30-Days Trial



How does the SoClean cleaning unit work?

The SoClean CPAP cleaner and sanitizing machine uses activated oxygen cleaning to eliminate any germs and bacteria that may be lingering in your home CPAP machine. Watch our CPAP cleaning demo video to learn more.

WATCH A DEMO

SoClean automated CPAP cleaner and sanitizing device

SoClean 2 is automated CPAP Cleaner machine which is used to clean and sterilize CPAP mask, reservoir, and hose without the need of detaching your CPAP device apart. Soclean CPAP cleaner and sanitizer doesn't require water or chemicals to work due to the fact that the tool takes advantage of Ozone known as activated oxygen.

How does So Clean 2 CPAP cleaner and sanitizer Works?

Soclean 2 - CPAP sanitizing machine is able to produce activated oxygen that is supplied all the way through the equipment to eliminate germs, bacteria and molds. Activated oxygen is safe and natural hence doesn't lead to pretty much any danger to the body. The inside chamber of the cleaner is designed in such a way that different type of masks could be used within the CPAP cleaning device. It can fit maybe even the most full face mask kinds given that the inside chamber is big enough that measuring 5.5 inches length, 5.5 inches width and 7.25 ins height.

Soclean 2 - CPAP cleaning machine is a powerful automated CPAP sterilizer for all CPAP devices. So Clean 2 - CPAP sanitizer is able to get rid of 99.9 % of germs, bacteria and molds without using water and chemical. So Clean 2 - CPAP cleaner and sanitizer is registered by FDA and has been checked in a lab.

Soclean 2 CPAP cleaner comes with a front display screen, indicator lights throughout the cycle, right or left hose insertion approach and east filter cartridge system.

Features of Soclean 2 - CPAP sanitizer

Easy daily CPAP sanitizing process - Soclean 2 CPAP cleaner and sanitizer sanitizes your CPAP devices immediately after just a self-explanatory initial setup. You can also clean the tools by hand by just a pushing a button. So Clean 2 CPAP sanitizing device does not need a lot of expertise to use and is well designed to make certain that the equipment is always safe to put into action.

Quick and Convenient to use - As I mentioned above, you don't have to disassemble the parts apart when the self-cleaning process due to the fact that this tool is developed in such a way that you can sanitize reservoir, mask and hose which still attached to your CPAP machine.

Clearing up your CPAP equipment simply spends few minutes while putting into action this tool but it's recommended you leave the lid which is on the device closed for about two hrs for after the initial process for proper and efficient sanitization.

Compatibility with all CPAP devices - So Clean CPAP cleaning device is suitable for almost all CPAP devices. For those kinds of equipment with a heated humidifier, a special tube insert is included so that the CPAP mask and tubing can be cleaned automatically simultaneously with the chambers.

Effective eliminating bacteria and germs - So Clean CPAP sanitizing device able to wipe out 99.9 per cent of virus, mold and bacteria. This makes it really highly effective in eliminating mold, bacteria and virus from your CPAP hose, reservoir and mask.

SoClean 2 Go Portable CPAP Cleaner and Sanitizer device

SoClean 2 Go essentially uses the very same powerful sanitizing procedure as the SoClean sanitizer. SoClean2Go represent the best choice for frequent travelers. You can take your elegant and compact sanitizing machine with you anywhere you go. This is likewise the first portable CPAP devices sterilizing machine that is portable.

Soclean 2 cartridge filter

SoClean 2 - cartridge filter kit consists of one replacement filter with O-ring connected and one replacement check valve for the SoClean . The filter and check-valve should be replaced every six months. So clean filter kit converts activated oxygen back to regular oxygen as it leaves the SoClean chamber. The check valve avoids water from your CPAP from damaging your CPAP electronic devices.

The Science Behind The Soclean CPAP cleaning device

Activated oxygen, also referred to as Ozone (O₃), is a safe disinfecting agent that So Clean 2 CPAP sanitizer uses to get rid of CPAP bacteria, mold or infections inside your CPAP devices. The generator inside the Soclean CPAP sanitizing device breaks down the chemical bond of typical oxygen (O₂) and allows it to recombine into a brand-new molecule that has 3 atoms of oxygen instead of two. Activated oxygen has a short life cycle. After being generated it automatically breaks down back to oxygen within 2 hours.

How Activated-Oxygen Works ?

1. Your So Clean CPAP cleaner and sanitizer creates and pumps the activated oxygen through the supply tube and into CPAP humidifier tank. This cleans not only the water, but likewise the inner walls of your tank.
2. The activated oxygen pass through your CPAP tube As it takes a trip, the activated oxygen reaches every crevice, eliminating all threatening bacteria and germs.
3. The activated oxygen reaches your mask. The activated oxygen passes into and out of your mask throughout this action, cleaning it just as it finished with your hose and tank.
4. Activated oxygen exits the chamber through an unique filter. SoClean filter converts any excess activated oxygen back into the oxygen that we breathe (O₂).

Summary

Soclean CPAP cleaner and sanitizer is wisely designed to make cleaning your CPAP machine much easier. It uses safe and natural activated oxygen as an alternative to water and chemical and does not need a lot of skills in operating it. Finally, it is hassle-free and compatible with most of CPAP machine. Check out here for more details how you can try soclean CPAP cleaner device free for 30-days trial without risking your money.