

Sleep8 CPAP cleaner FAQ

Q. Is Sleep8 covered by insurance?

A. Right now, [Sleep8](#) is not covered by insurance. Nevertheless, you can use your HSA/FSA accounts under code 9274 for payment.

Q. Does my mask have to be wiped with a towel or CPAP wipe before make use of?

A. We strongly encourage you to clean your mask with a CPAP wipe, or damp cloth prior to putting it in the Sleep8 Disinfecting Bag. Specific skin oils as well as cosmetics including lotions and creams can leave deposit on your mask which can respond to ozone. It is suggested to ensure your CPAP mask is free of any type of oils prior to making use of the Sleep8.

Q. Can I charge my Sleep8 and also use it at the very same time?

A. Yes. Your Sleep8's unique charging system enables simultaneous charging and cleaning.

Q. Does Sleep8 work with the unique DreamStation Go 12mm tubing?

A. Sleep8 is presently enhancing our bag design to additionally widely support DreamStation Go 12mm tubes. As a workaround, you can utilize your DreamStation Go mask tubing adapter to safely attach to your Sleep8 bag. You can safely leave your mask inside the bag during the cleaning cycle..

Q. How much time does Sleep8 require to charge?

A. The device takes about one and one-half (1.5) hours to fully charge.

Q. Is the Sleep8 Sterilizing Bag recyclable?

A. Yes, our team believes in the atmosphere. Our bags are reusable and pleasant.

Q. The length of time Sleep8 requires to charge?

A. The device takes approximately one and one-half (1.5) hrs to fully bill.

Q. Does the Sleep8 machine use chemicals to clean up?

A. No, your Sleep8 CPAP cleaner converts room air into activated oxygen (O3) by adding one little molecule of oxygen, an atom. There are no added chemicals, solvents or liquids made use of by your tool to disinfect.

Resources :

<https://s3.amazonaws.com/cpap.cleaner/Sleep8-CPAP-cleaner.pdf>