

Lumin CPAP mask and accessories cleaner and sanitizer

Disinfect your CPAP equipment
in just 5 minutes!

AS EASY AS 1-2-3

1
Open the Lumin sliding tray

2
Place your CPAP mask and humidifier chamber on the Lumin tray

3
Close & push the button to start

UV kills 99% of bacteria. No harmful ozone! Only \$249.

Lumin™

BUY NOW

Lumin CPAP mask and accessories cleaning tool by 3B Medical is the most convenient, fastest, and best way to disinfect a CPAP mask and accessories by using a high-powered UV light, which works as an efficient germicide and get rid of 99% any bacteria, molds, and viruses that comes in contact with your CPAP equipment.

3B Lumin CPAP sanitizer machine will sanitizes your CPAP machine for you within 5 minutes. This way you can rest assured you can utilize them right away as they come out of the cleaner, without having to air dry them. You simply require to pull out the main tray of the Lumin cleaner, put your CPAP machine inside and then press start. Within 5 minutes the cleaner will indicate that your items are sanitized and prepared to use.

How Does 3B Lumin CPAP and Mask Cleaning Device works?

Lumin CPAP cleaner and sanitizer tool runs using a specialized high powered germicidal UV bulb. The drawer is made from polished aluminum which reflects the Ultraviolet light to enhance coverage. Ultraviolet disinfection is effective on kill 99.99% germs, molds and bacteria and is well suited for keeping your CPAP mask and humidifier clean from bacteria that can make you sick.

Daily usage of the 3B Lumin CPAP and mask cleaner and sanitizer tool will help keep you healthy. Because the 3B Lumin CPAP mask cleaner and sanitizer device runs utilizing a type of high energy light, it is extremely efficient in decontaminating your CPAP accessories. To keep your mask, hose pipe and water chamber sanitize and decrease the threat of infections, we suggest utilizing 3B Lumin CPAP mask and accessories cleaner device twice weekly on your mask, water chamber and tube.

Benefits of Lumin CPAP and Mask Cleaner Device

Five Minute Sanitization - The 5 minute ozone FREE sanitizing cycle without any need for water or damaging chemicals.

No Harmful Ozone - The use of UV-C is eco-friendly, leaves no residue or hazardous gases or chemicals. UV-C systems are currently in usage to sanitize ambulances, emergency service and other high touch places.

Effective -UV light will sanitize approximately 99% of damaging bacteria and mold that can cause infection and illness. UV light is likewise the best disinfection option on the market, there is no harmful OZONE.

Budget-friendly -Lumin CPAP and mask cleaning machine budget friendly at less than other brand names. A healthy CPAP experience should be budget friendly for everybody.

No Maintenance is needed -The Lumin CPAP Mask sanitizer Tool needs no special care of upkeep and there are no parts that need to be replaced on a set up basis. The high powered UVC light bulb was designed specifically for the Lumin CPAP mask and accessories cleaner tool and creates a high-powered disinfection wavelength of light and outlast the life of your gadget.

Multipurpose -The Lumin CPAP Mask Cleaner Tool is also a general-purpose disinfecting device and can decontaminate any item placed in the drawer. You can keep it on your bathroom counter and sanitize any important items such as your hearing aid, tooth brush, and small toys .

Conclusion

Lumin CPAP mask and accessories cleaner device is a smart investment for anyone who requires to preserve CPAP equipment clean and sterilized. {Lumin CPAP and mask cleaner device is capable of kill 99.99% germs, molds and bacteria that may live inside your CPAP mask and accessories..

Resources :

<https://s3.amazonaws.com/eliminate-snoring/cpap-cleaner/3B-Lumin-CPAP-cleaner.html>