

# ICIS View: Quick Start Guide

The screenshot displays the ICIS View software interface. At the top, there is a search bar and user information: "PODNIK, PETE | 27 Male 01-Jan-1933". Below this is a study selection dropdown menu showing "CT LUMBAL WZ" and "09-Oct-2001 11:56 (63065)". The main area shows a 3D axial CT scan of the lumbar spine. On the left, a vertical toolbar contains icons for pan, zoom, and other navigation functions. At the bottom, there are playback controls and a page indicator "105 of 249".

Callouts and their corresponding functions:

- Image Toolbar**: Points to the vertical toolbar on the left side of the interface.
- Select a study**: Points to the study selection dropdown menu.
- Change study options**: Points to the menu icon (three horizontal lines) next to the study selection dropdown.
- Cine viewer tools**: Points to the playback controls at the bottom of the interface.
- Reset layout**: Points to the circular arrow icon in the top right toolbar.
- Navigate series**: Points to the left and right arrow icons in the top right toolbar.
- Display series tray**: Points to the icon with a list of three items in the top right toolbar.
- Compare Studies**: Points to the icon with two overlapping document pages in the top right toolbar.