








DIRECT SUPPLY® SLING COMPATIBILITY GUIDE








Direct Supply Sling Compatibility Guide – Floor Lifts

Direct Supply Slings, Multi-Brand Compatible		Direct Supply®	Invacare®	Drive®	Joerns®	
						
Product #	Description	Atlas Floor Lift	Reliant 450 Floor Lift	Drive Patient Lift	Drive Levantar Floor Lift	Hoyer® HML400 Hydraulic Lift
		E7793 GV270 GV269	RPL450-2 RPL450-1	13240	FLP500	HML400
#GV258	Full Body Sling, Standard, Extra Small	•	•	•	•	•
#GV257	Full Body Sling, Standard, Small	•	•	•	•	•
#GV263	Full-Body Sling, Standard, Medium	•	•	•	•	•
#GV262	Full-Body Sling, Standard, Large	•	•	•	•	•
#GV261	Full-Body Sling, Standard, Extra Large	•	•	•	•	•
#GV255	Full Body Sling, Mesh, Extra Small	•	•	•	•	•
#GV256	Full Body Sling, Mesh, Small	•	•	•	•	•
#GV260	Full-Body Sling, Mesh, Medium	•	•	•	•	•
#GV259	Full-Body Sling, Mesh, Large	•	•	•	•	•
#GV254	Full Body Sling with Commode Opening, Mesh, Medium	•	•	•	•	•
#GV253	Full Body Sling with Commode Opening, Mesh, Large	•	•	•	•	•
#GV252	Full Body Sling with Commode Opening, Mesh, Extra Large	•	•	•	•	•
#GV251	Divided-Leg Sling, Small	•	•	•	•	•
#GV250	Divided-Leg Sling, Medium	•	•	•	•	•
#GV249	Divided-Leg Sling, Large	•	•	•	•	•
#GV248	Divided Leg Sling, Extra Large	•	•	•	•	•
#GV275	Toileting Sling, Small	•	•	•	•	•
#GV274	Toileting Sling, Large	•	•	•	•	•
#GV273	Toileting Sling, Extra Large	•	•	•	•	•

Buying compatible slings means no need to spend time searching for the right sling for each

Direct Supply Sling Compatibility Guide – Sit-To-Stand Lifts

Direct Supply Slings, Multi-Brand Compatible		Direct Supply®	Invacare®	Medline®	Tollos®	
						
		Atlas Sit-to-Stand Lift	Reliant 350 Sit-to-Stand Lift	Electric Stand Assist Lift	Powered Base Stand Assist Lift	Steady-Aid 3500/7500 Sit-to-Stand Lifts
Product #	Description	E7794 GV272 GV271	RPS350-1 RPS350-2	MDS400SA	MDS500SA	TSA-350 TSA-750
#GV265	Padded Standing Sling, Extra Small	•	•	•	•	•
#GV268	Padded Standing Sling, Small	•	•	•	•	•
#GV267	Padded Standing Sling, Medium	•	•	•	•	•
#GV266	Padded Standing Sling, Large	•	•	•	•	•
#GV264	Gluteal Strap - Multi-Brand	•	•	•	•	•
#GV278	Transfer Sling, Small	•	•	•	•	•
#GV277	Transfer Sling, Medium	•	•	•	•	•
#GV276	Transfer Sling, Large	•	•	•	•	•

Buying compatible slings means no need to spend time searching for the right sling for each lift in your community; manage your slings inventory with multi-brand compatible slings.

Don't see your lift in this guide?
Contact your Account Manager for assistance

The following spreader bars were used to complete compatibility testing:

Invacare Reliant 450 (RPL450-1, RPL450-2): **1143629**
 Drive Patient Lift (13240): **13058**
 Hoyer Hydraulic (HML400) and Hoyer HPL402: **222-600**
 Hoyer HPL500 (HPL500-S2, HPL500-EL-S2): **31009196**
 GF Lumex LF1090: **DPL450-SBAR**
 GF Lumex LF1050: **DPL450-SBAR**
 Medline: MDS400EL: **MDS450CRD**
 Medline: MDS88200D: **MDS88200DCR**



APTURA / EQUIPMENT & FURNISHINGS / DSSI / TELS

Direct Supply multi-brand slings have been tested with a variety of lifts, as detailed in this guide. When purchasing this sling for use with a non-Direct Supply branded lift, customers must be aware of the following:

1. The use of this sling with a non-Direct Supply branded lift may void the warranty of the lift. Customers are advised to check the warranty period on their lifts before making a sling purchase.
2. State surveyors may require documentation showing that the Direct Supply branded sling is appropriate for use with a non-Direct Supply branded lift. We recommend that you keep a copy of the sling compatibility guide with your other product manuals and materials.