

Phone: 888-558-5227 651-644-8424

888-558-7329 Fax:

Email: getinfo@lktlabs.com

lktlabs.com Web:

Product Injormation	

Product ID N342687

CAS No.

**Chemical Name** 

Synonym

**Formula** Formula Wt. **Melting Point** Purity Solubility

**Pricing and Availability** 

Bulk quanitites available upon request

Product ID Size **List Price** N342687 20 Tests \$377.00

Store Temp -20°C Ship Temp Blue Ice

**Description** Near-infrared (NIR)-FLIVO® Tracers are used in conjunc-tion with NIR-FLIVO® Free Dye Control Assays. The NIR-FLIVO® Free Dye Control Assay uses the NIR-FLIVO® Free Dye Control reagent (\*Dylight® 690 Free Dye). When injected, both the Free Dye Control and the Tracer reagents will gen-erate a fluorescent signal, but only the Tracer should bind to active caspases and remain inside an apoptotic cell, while signal associated with the Free Dye Control represents background noise.

> FLIVO® is a registered trademark of Immunochemistry Technologies, LLC. \*Dylight® is a registered trademark of Thermo Fisher Scientific, Inc. and its subsidiaries.

References