

TEST REQUISITION FORM

Investigator Information

Site #: _____ Physician Name: _____ NPI #: _____
 Site Name: _____

Institution

Street Address: _____
 City: _____ State: _____ ZIP Code: _____ Country: _____
 Phone: _____ Fax: _____ Physician Email: _____
 Requisition Completed by (Print Name): _____
 Physician Signature: _____ Date (DD/MM/YYYY): / /

Patient/Subject Information

Patient/Subject ID: _____ Medical Record #: _____
 Sex: F M Birthdate (DD/MM/YYYY): / /

Disease State/Diagnosis

Disease State/Diagnosis: _____
 Did this patient give consent to use this sample for future research? Yes No A copy of the consent form may be requested by Ignyta.

Trailblaze™ Tests Requested: (Please check test below)

Pharos™ NGS (NTRK1/2/3, ROS1, ALK)

Specimen information (Attach a Copy of the Pathology Report)

Specimen Located at Third-Party Facility
Please request tissue specimen from the facility indicated below.

Location: _____ City: _____ State: _____ Phone: _____

Collection Date: (DD/MM/YYYY): / /

Collection Time (24hr): : Formalin Fixation Time (24hr): :

Specimen Site: _____ Specimen ID: _____ Primary Metastatic

Paraffin Embedded Tissue

Fresh-Frozen Tissue

Slides

Other: _____

Date Slides Were Cut (DD/MM/YYYY): / /

4µm Qty: _____ Block/Slide ID: _____

10µm Qty: _____ Block/Slide ID: _____

For IGNYTA Lab Use Only:

Notes:

Specimen Requirements and Handling

The quality of laboratory results is highly dependent upon proper specimen collection and handling. Listed below are specimen requirements and handling procedures for tests processed by Ignyta.

Specimen Type ¹⁻³	Requirements for NGS	Storage/Transport ⁵
Formalin Fixed Paraffin Embedded Tissue (Blocks)	Preferred specimen.	Room Temperature. Use cold pack for transport.
Formalin Fixed Paraffin Embedded Unstained Slides⁴	10µm thick tissue sections on positively charged slides. 4 slides per assay and 1 H&E stained slide.	Room Temperature. Use cold pack for transport.

H&E slide required: If a stained H&E cannot be provided, please send one unstained slide cut at 4µm thickness for H&E staining.

1. Two (2) patient identifiers are required on all specimens. Please include date of collection.
2. Specimens should be fixed in 10% Neutral Buffered Formalin (NBF) for a minimum of 6 hours. Do not use zinc or B5 fixatives.
3. Paraffin embedded tissue that has been decalcified is not appropriate for IHC or NGS.
4. Unstained slides should be received within 6 weeks of cutting.
5. Cold pack should not be placed in direct contact with the specimen.

Transportation:

Place specimens in an Ignyta Pathology Specimen Shipping kit. Slides should be enclosed in a plastic slide mailer. Enclose a completed test requisition form and a copy of the pathology report.

Shipping kits may be obtained from Client Services

QUESTIONS? CALL IGNYTA CLIENT SERVICES AT:

1-844-IGNYTA1

4 4 6 - 9 8 2 1
OR INTERNATIONAL **+1-858-255-5950**