

A SET FOR EVERY BODY.

A set for every lifestyle.

Find out more at medtronicdiabetes.com/infusionsets

WHICH INFUSION SET IS RIGHT FOR YOU?


STEEL

Reactions to plastic cannula /
History of bent cannula / Small child
/ Pregnant (up to 2nd trimester)

90 Degree Insertion

MiniMed® Sure-T® infusion set

Simple steel needle with kink-free insertion




PLASTIC CANNULA

Less dexterity / Less active /
Alternative site usage

90 Degree Insertion

MiniMed® Quick-set® infusion set

Comfort combined with ease-of-use




Less active / Alternative
site usage

90 Degree Insertion

MiniMed® Mio® infusion set

Convenience and Customizability




Less active/ Alternative
site usage

90 Degree Insertion

MiniMed® Pro-set® infusion set

With BD FlowSmart™ technology




Active / Lean / Muscular /
Pregnant (from 2nd trimester)

30-45 Degree Insertion

MiniMed® Silhouette® infusion set

Flexibly angled for convenience



Recommendation based on Insulin Pump Therapy: Best practices in choosing and using infusion devices. American Association of Diabetes Educators, 2011. Data on file. - BMI stands for Body Mass Index which is calculated based on an individual's weight and height.

Important Safety Information: Infusion sets are indicated for the subcutaneous infusion of medication, including insulin, from infusion pumps. They are neither intended nor indicated for intravenous (IV) infusion or the infusion of blood or blood products. Inaccurate medication, delivery, infection and/or site irritation may result from improper insertion and maintenance of the infusion insertion and maintenance of the infusion site. Replace the infusion set every 48-72 hours, or per your healthcare professional's instructions.

For a listing of indications, contraindications, precautions, warnings, and potential adverse events, please refer to the Instructions for Use.

Legal manufacturer: Unomedical a/s, Aaholmvej 1-3, Osted, 4320 Lejre, Denmark

Legal manufacturer: Becton Dickinson, 1 Becton Drive Franklin Lakes NJ 07417 Medtronic Diabetes MiniMed and Pro-set are registered trademarks of Medtronic MiniMed

MiniMed, Pro-set, Mio, Quick-set, Silhouette and Sure-T are registered trademarks of Medtronic MiniMed, Inc. BD and FlowSmart trademarks are property of Becton, Dickinson and company.

940M14407-011 20160218 ©2016 Medtronic MiniMed, Inc. All rights reserved.

Medtronic Diabetes
18000 Devonshire Street
Northridge, CA 91325
www.medtronicdiabetes.com

Medtronic