




Show or Hide Colors


The sequence colors shown in Map and Sequence views can be shown or hidden.

Open the Map and Sequence Options Dialog



To open the Options dialog, click View → Map & Sequence Options... .


Show the Sequence Colors






To show the sequence colors, click the "Show colors" button, then click OK.
This button can also be used to toggle between showing and hiding colors.

View the Sequence Colors



Sequence colors will now be shown in Map and Sequence views.

Hide the Sequence Colors





To hide the sequence colors, click **View → Map & Sequence Options...** , click the "Show colors" button, then click OK.

Show or Hide Sequence Colors Using the Side Toolbar



Alternatively, click the "Show colors" button in the side toolbar to toggle between showing and hiding sequence colors.