



Logger Pro 3 Installation Guide


Updated: 8/27/2020

The instructions below are for gaining access to Logger Pro 3 for classroom computers. For students using the software at home, the instructor or IT needs to follow the process at the following link: <https://www.vernier.com/my-account/distribute-lp3/>

Distance Learning is also available using Vernier Graphical Analysis 4, which can also be used by students. Details for installation are on the software page at pltw.org/software.

Procedure

1. Go to the following link: <http://www.vernier.com/pltw-lp3>



2. Fill out the information and select Continue to Account Setup.

Instructions

1. Fill out the customer information form below.
2. Create or sign in to your account.
3. Download Logger Pro 3.

What would you like to learn more about?

How would you like to learn more about our solutions for hands-on science?*

Email, go green!
 Mail
 Both email and mail
 None

First name*

Last name*

School Email*

Job Title*

School Name*

School Address*

School City*

School State*

School Zip*

Institution Type*

[Continue to Account Setup](#)

[Chat with us!](#)

3. If you have an account Login. If you do not, got to the Registration process.

We're here to help educators adjust to remote learning. We offer free webinars, experiments, and a number of additional resources that help students learn at home.

[EXPLORE REMOTE SOLUTIONS](#)

Downloads

Vernier [IDEAS](#) [PRODUCTS](#) [EDUCATOR TRAINING](#)

Login

Username or email address *

Password *

Remember me

[LOG IN](#)

[Lost your password?](#)

Register

Username *


Email address *

A password will be sent to your email address.


Your personal data will be used to support your experience throughout this website, to manage access to your account, and for other purposes described in our [privacy policy](#).

[REGISTER](#)


4. After logging in or going through Registration, select Downloads.




5. Select the download you need. This guide will continue through the Windows installation.




6. Select the installer from the download location.




7. Select Ok




8. Wait for the installer to begin.




9. Ensure that USB devices have been removed and select OK.




10. Accept the license agreement and select OK.




11. Select Next.




12. Select Next.



13. Select Install.



14. Select Finish.



For students using the software at home, the instructor or IT needs to follow the process at the following link: <https://www.vernier.com/my-account/distribute-lp3/>