

MODEL: R5K DIGITAL PILLOW SPEAKER, TV ONLY
R5K DIGITAL PILLOW SPEAKER, 2 LIGHTS
R5K DIRECT ACCESS PILLOW SPEAKER, TV ONLY
R5K DIRECT ACCESS PILLOW SPEAKER, 2 LIGHTS



NCPSTDSTV



FEATURES

- The sealed design reduces the potential risk of infectious fluids or solids settling into the unit
- Large, recognizable, easily pressed Nurse Call and TV push buttons
- Intuitive, standardized button positioning. Mute, closed-captioning, sleep, and option are standard on all models
- Clear, high quality sound to increase patient satisfaction
- Streamlined design facilitates quick and easy cleaning
- Can be programmed to support various hospital TV brands within seconds
- Lightweight case designed to resemble the familiarity of a home remote control
- High-impact, fire resistant UL 94V-0 plastic case
- Volume control up-down buttons
- UL 1069 Listed US/Canada

SPECIFICATIONS

1. 8.5" (21.6 cm) long x 3" (7.6 cm) wide x 1" (2.5 cm) deep
2. A 2 1/4" (5.7 cm) diameter liquid resistant Mylar® for both nurse call intercom audio and TV audio. The pillow speaker is designed for superior voice TV audio reproduction.
3. A 9 3/4 foot flexible vinyl-jacketed cable with an 8-pin DIN plug for connection to the associated patient station and molded strain relief at both ends of the cable. when the pillow speaker cord is removed, a call will be placed to the associated Nurse Console(s)
4. An easily distinguished red Nurse Call button for placing patient calls with associated call assurance LED (Responder IV and Responder 5 systems)
5. High Reliability membrane switches for repeated use
6. For TV entertainment, the pillow speakers must be interfaced with listed digital hospital grade TVs only
7. The pillow speaker and associated cord are electrostatic discharge (ESD) shielded for protection of the nurse call system to ensure system reliability
8. Specially designed cases, cords, and buttons allow the entire unit to undergo gas sterilization without discolorization or deterioration
9. A user-friendly design makes buttons easy to find and easy to press
10. A stainless steel non-tear linen clip to position the pillow speaker within easy reach of the patient
11. Low voltage isolation compliance with the ANSI/UL 1069 requirements

MODEL: R5K DIGITAL PILLOW SPEAKER, TV ONLY
R5K DIGITAL PILLOW SPEAKER, 2 LIGHTS
R5K DIRECT ACCESS PILLOW SPEAKER, TV ONLY
R5K DIRECT ACCESS PILLOW SPEAKER, 2 LIGHTS

DESCRIPTION

The Basic series of Pillow Speaker models are designed for use with the Responder 5000 nurse call systems and digital hospital grade televisions. A variety of pillow speakers are available. The pillow speakers combine in one compact unit: up to two (2) controls for room lighting, a speaker for listening to TV and nurse call audio. All models include an

easily identified Nurse Call button embossed in Braille and TV control button.

When the pillow speaker is plugged into the patient station receptacle, nurse call intercom audio can be automatically transferred to the pillow speaker while simultaneously muting the TV audio.



NCPDSTV



NCPKSL2



NCPKPTV



NCPKPL2

MODEL (Note 1)	TV Control (Note 2)	TV Channel	Volume Up/ Down	Numeric Keypad	Light 1 Control (Note 3)	Light 2 Control (Note 3)	Light 3 Control
NCPDSTV	X	X	X				X
NCPDSSL2	X	X	X		X	X	X
NCPKPTV	X	X	X	X			X
NCPKPL2	X	X	X	X	X	X	X

All models include 2 1/4" speaker and Nurse Call button

NOTES

- All digital pillow speakers require the use of the R4KTVR1 Single Digital TV Isolation Module for the Responder 5000 System stations.
- For TV Entertainment, pillow speaker must be interfaced with listed hospital grade TVs.
- Light control features require the use of the R5KPS1LCA patient station or the R4KSLC1/R4KDL2 Lighting Control modules for other Responder 5000 Stations. Refer to *Responder 5000 Component Installation Guide*.
- Option button can control Zenith or RCA TVs with available built-in radio tuner or Philips TVs with available Channel Guide.



Rauland
 A Division of AMETEK, Inc.
www.Rauland.com

Toll Free +1 800 752 7725
 From Outside the U.S. +1 847 590 7100

Architect and Engineer (A&E) Specifications available online at customerconnection.rauland.com Specifications subject to change without notice.