



RE: Cross-Generation Accuracy

Date: 3/9/2015

ActiGraph is committed to providing researchers and clinicians with reliable, accurate, and cost effective solutions to their activity monitoring needs. Further, ActiGraph recognizes the necessity and importance of ensuring subsequent device generations match the accuracy and repeatability of its predecessor generation. To this end, ActiGraph expends significant engineering hours and man power in development and test to ensure this goal is met. It is this rigorous and obsessive pursuit and desire that affords ActiGraph's customers the ability to remain confident that cross-generation accuracy is maintained.

The data provided by the GT3X+ and w-GT3X+, when subjected to the same acceleration environment, will exhibit outputs within 5% of each other.

Jeff Miller

DocuSigned by:

A handwritten signature in black ink that reads "Jeff Miller". The signature is written in a cursive style and is enclosed within a blue rectangular border.

2AF4F25DAB3B4AA...

Senior Vice President of Engineering