

# TEMPO™ REMOVAL OF BRASS ELBOW ON INTAKE SIDE OF PUMP

Follow these instructions to remove the brass elbow on the intake side of the pump in the event the machine's internal solution hose is getting kinked.

## Tools Needed:

- Phillips screwdriver or drill driver
- 5/8", 7/16", and 9/16" wrenches
- Teflon tape

Labor hours allowed: 1

- 1 Remove the three connecting plates (left-front, right-front, and rear) holding the solution tank to the base of the machine. Pull off the recovery and solution tanks as a unit.
- 2 Use a 5/8" wrench to undo the solution hose at the brass elbow connected to the bottom of the solution tank (**Figure A**).
- 3 Use a 7/16" wrench to undo the brass elbow and solution hose from the intake side of the pump (**Figure B**). Use a 9/16" wrench to remove the brass elbow from the solution hose.
- 4 Wrap teflon tape on the threads on the pump's intake and screw the solution hose without the brass elbow back onto the intake side of the pump with the 9/16" wrench (**Figure C**). Screw the other end back into the elbow still on the bottom of the solution tank.
- 5 Put the machine back together.

