

Thormahlen Harps: Swan or Cygnet Lever (Folk) Gut				
String#	Note	Color	Gauge	Material
1	C	Red	.025 (highest)	Nylon
2	B	Clear	.025	Nylon
3	A	Clear	.025	Nylon
4	G	Clear	.025	Nylon
5	F	Blue	.025	Nylon
6	E	Clear	2nd Oct	Lever Gut
7	D	Clear	2nd Oct	Lever Gut
8	C	Red	2nd Oct	Lever Gut
9	B	Clear	2nd Oct	Lever Gut
10	A	Clear	2nd Oct	Lever Gut
11	G	Clear	2nd Oct	Lever Gut
12	F	Blue	2nd Oct	Lever Gut
13	E	Clear	3rd Oct	Lever Gut
14	D	Clear	3rd Oct	Lever Gut
15	C	Red	3rd Oct	Lever Gut
16	B	Clear	3rd Oct	Lever Gut
17	A	Clear	3rd Oct	Lever Gut
18	G	Clear	3rd Oct	Lever Gut
19	F	Blue	3rd Oct	Lever Gut
20	E	Clear	4th Oct	Lever Gut
21	D	Clear	4th Oct	Lever Gut
22	C	Red	4th Oct	Lever Gut
23	B	Clear	4th Oct	Lever Gut
24	A	Clear	4th Oct	Lever Gut
25	G	Clear	4th Oct	Lever Gut
26	F	Blue	4th Oct	Lever Gut
27	E	Clear	5th Oct	Lever Gut
28	D	Clear	5th Oct	Lever Gut
29	C	Red	5th Oct	Lever Gut
30	B	Clear	wound	steel/silk/steel
31	A	Clear	wound	steel/silk/steel
32	G	Clear	wound	steel/silk/steel
33	F	Blue	wound	steel/silk/steel
34	E	Clear	wound	steel/silk/steel
35	D	Clear	wound	steel/silk/steel
36	C	Red	wound	steel/silk/steel

This string list is provided as a courtesy
by the Sylvia Woods Harp Center
www.harpcenter.com
(808) 212-9525

String gauges are subject to change without notice