
1) Follow the model and finger trace M and m, then read the sentence aloud.



Use your pencil to trace the letters in the boxes then start at the dots and write your own.

		•		•			l .		l
•	••	••	••	1170	1000	I		ı	ı
•		• •	• •			I		ı	ı
•			• •			1		I	ı
•									
•.						• • • • • • • • • • • • • • • • • • • •	• • • • • • • • • • • • • • • • • • • •		. •
•								ı	
-						1 ' : '		ı	
						1 ' : '		ı	
						1		ı	
								ı	

4) Copy this sentence on the line below:

M is for milk.

