

THE HERITAGE COLLECTION®



Application

Patients will welcome a touch of elegance while dining. Rich colors, contemporary faceted styling, and a china-like appearance combine to help de-institutionalize your healthcare food service and making patients feel more at home. Not intended for use in microwave or conventional ovens, hot/cold carts.

Material

Attractive, durable, and stackable, the Heritage Collection® offers the best temperature retention in the industry - as proven in hospitals and nursing homes. Like all Dinex® insulated products, Heritage is manufactured with ozone-safe urethane foam insulation. Available in a variety of contemporary color choices. Latex free construction.

Washing Instructions

Weekly, soaking 6-8 hours, using Wash in commercial dishmachine. Do NOT use bleach. Can withstand washing temperatures up to 180°F.

Item Numbers

8 oz. Mug

- ☐ DX400003 - Onyx
- ☐ DX400015 - Teal
- ☐ DX400023* - Gray
- ☐ DX400050 - Midnight Blue
- ☐ DX400056 - Mauve
- ☐ DX400061 - Cranberry

5 oz. Bowl

- ☐ DX420003 - Onyx
- ☐ DX420015 - Teal
- ☐ DX420023* - Gray
- ☐ DX420050 - Midnight Blue
- ☐ DX420056 - Mauve
- ☐ DX420061 - Cranberry

Lids for 8 oz. Mug & 5 oz. Bowl

- ☐ DX40008700 - Clear-View Flat Lid
- ☐ DX40008714 - Translucent Lid

Item Numbers

8 oz. Tumbler

- ☐ DX410015* - Teal
- ☐ DX410023* - Gray
- ☐ DX410056* - Mauve
- ☐ DX410061* - Cranberry

Lids for 8 oz. Tumbler

- ☐ DX40008700 - Clear-View Flat Lid
- ☐ DX40008714 - Translucent Lid

9 oz. Bowl

- ☐ DX430003 - Onyx
- ☐ DX430015 - Teal
- ☐ DX430023* - Gray
- ☐ DX430050 - Midnight Blue
- ☐ DX430056 - Mauve
- ☐ DX430061 - Cranberry

Lids for 9 oz. Bowl

- ☐ DX43008700 - Clear-View Flat Lid
- ☐ DX43008714 - Translucent Lid

12 oz. Bowl

- ☐ DX430050 - Midnight Blue
- ☐ DX430056 - Mauve
- ☐ DX430061 - Cranberry

Lids for 12 oz. Bowl

- ☐ DX11858700 - Clear-View Flat Lid
- ☐ DX11858714 - Translucent Lid

Insulated Dome

- ☐ DX440015 - Teal
- ☐ DX440023* - Gray
- ☐ DX440050 - Midnight Blue
- ☐ DX440056 - Mauve
- ☐ DX440061 - Cranberry


* Indicates discontinued item. Quantities may be limited; call for availability.

See spec sheet #370 for options on Domes, Insul Bases and Underliners.


DINEX®

Building 2, Suite 106
628 Hebron Avenue
Glastonbury, CT 06033
800.523.9752 Fax: 860.652.3135
www.dinex.com


THE HERITAGE COLLECTION®




5 oz. Bowl




8 oz. Mug



8 oz. Tumbler



9 oz. Bowl



Dome

Items and Dimensions

ITEM NUMBER	DIAMETER	OVERALL HEIGHT	PACK	SHIP WEIGHT	SHIP CUBE	CAPACITY	FREIGHT CLASS	MATERIAL
DX4000series (Mug)	3.75" (9.5cm)	3.5" (8.9cm)	48	9.6 lb. (4.4kg)	1.51	8 oz.	77.5	Polypropylene
DX4100series (Tumbler)	3.75" (9.5cm)	3.5" (8.9cm)	48	8.7 lb. (3.9kg)	1.50	8 oz.	77.5	Polypropylene
DX4200series (5 oz. Bowl)	3.75" (9.5cm)	2.5" (6.4cm)	48	7.1 lb. (3.2kg)	1.10	5 oz.	77.5	Polypropylene
DX4300series (9 oz. Bowl)	4.5" (11.4cm)	2.5" (6.4cm)	48	8.8 lb. (4.0kg)	1.51	9 oz.	77.5	Polypropylene
DX4500series (12 oz. Bowl)	3.75" (9.5cm)	3" (7.6cm)	48	11.5 lb. (5.2kg)	1.50	12 oz.	77.5	Polypropylene
DX4400series (Dome)	10" (25.4cm)	2.75" (7.0cm)	12	9.6 lb. (4.4kg)	1.00	-	77.5	Polypropylene
DX40008700 (Lid)	-	-	2,000	6.7 lb. (3.0kg)	2.20	-	150	Polystyrene
DX40008714 (Lid)	-	-	2,000	6.4 lb. (2.9kg)	2.20	-	250	Polystyrene
DX43008700 (Lid)	-	-	1,000	4.8 lb. (2.2kg)	1.50	-	150	Polystyrene
DX43008714 (Lid)	-	-	1,000	4.7 lb. (2.1kg)	1.50	-	150	Polystyrene
DX11858700 (Lid)	-	-	1,000	5.7 lb. (2.6kg)	1.40	-	150	Polystyrene
DX11858714 (Lid)	-	-	1,000	5.7 lb. (2.6kg)	1.40	-	150	Polystyrene

Note: Heritage Insulated Dome serving height (measured on a standard patient tray) with:

- Insul Base 4-1/8"
- Wax Base Underliner 4-7/16"
- Dome only 2-7/8"



Please confirm that you have the most current specification sheet by visiting www.dinex.com.

Dinex® reserves the right to change specifications and product design without notice.

Such revisions do not entitle the buyer to corresponding changes, improvements, additions or replacements for previously purchased equipment.

Printed in USA

DINEX®

Building 2, Suite 106

628 Hebron Avenue

Glastonbury, CT 06033

800.523.9752 Fax: 860.652.3135

www.dinex.com

320-REV 10/09