

## The IQSense mouthpiece for Midmark® spirometers plus the IQGuard Bacterial/Viral filter





Keep your patient/clinician spirometry environment safe from exhaled bacterial or viral pathogens by using the IQSense mouthpiece\* with the IQGuard Filter.

- Highly effective barrier (99.996%) to cross-contamination
- Transparent for 100% visual inspection
- Small, not intimidating for adults as well as pediatric patients
- Recyclable
- Relatively low cost for disease prevention
- Supplied with 3 IQSense mouthpieces per box of 25 IQGuard filters



When performing spirometry testing, the use of a hollow mouthpiece affords no protection from from any pathogens present in the aerosolized forced exhalation of your patient. Your workspace environment is therefore compromised. Using the IQGuard System ensures that any airborne bacteria or viruses will be contained. Practice safe spirometry!

Part Number	IQ GUARD SYSTEM DESCRIPTION	Units	Packaging
29-3202-025	IQGuard System (25 PulmoGuard IQ & 3 IQSense)	25	25
29-3202-100	IQGuard System (100 PulmoGuard IQ & 12 IQSense)	100	4 boxes/25
29-3202-300	IQGuard System (300 PulmoGuard IQ & 30 IQSense)	300	300



Your Dealer is:

Manufactured by:

