

# **Ultra Fragrant Commode Clips**



Take control of your environment with Airlift® Ultra Fragrant Commode Clips. The patented design is packed with fragrance to keep restrooms smelling fresh for 30 days. With a discreet design, Airlift Ultra Fragrant Commode Clips eliminate odor at the source, improving restroom sanitation perception. For best results, use with common scent products like Solid Air Refreshers, Airlift odor eliminators, and Spartan cleaner/disinfectants.



# AIRLIFT ULTRA FRAGRANT COMMODE CLIPS

#### **FEATURES & BENEFITS**

- Freshens the restroom for 30 days
- · Discreet design eliminates odors at the source
- Ensures the best experience for patrons by providing a fresh environment
- Complies with VOC requirements in all 50 states, North America, and Europe

#### **RESPONSIBLE CHEMISTRY**

• 100% recyclable

#### **SPECIFICATION DATA**

Color Blue / Purple

Fragrance NABC® Fresh Scent / Xcelenté® Scent

Stability: Shelf @ 24°C/75°F 1 year minimum

#### **DIRECTIONS FOR USE**

- 1. Clip on the outside of the toilet bowl rim
- 2. Replace every 30 days



## **ORDERING INFORMATION**





#800300 NABC®

#800400 Xcelenté®

#### **ADDITIONAL RESOURCES**

Scan the QR code to access this product's safety information and documentation.



### COMMON SCENTS: FOR BEST RESULTS, CHOOSE SPARTAN PRODUCTS WITH MATCHING FRAGRANCES









**SAFETY: For institutional and industrial use only.** Be sure to read all directions, precautionary and first aid statements on product labels before using this or any other Spartan product. If questions remain, consult your employer or a physician. Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor and at www.spartanchemical.com. Workplace labels are available in English, Spanish, and French.

**GUARANTEE:** Spartan's modern manufacturing and laboratory control ensures uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within 12 months of the date of manufacture.

