



Q Restore OTC

Q Restore OTC hydrates and soothes your skin while promoting natural healing. Perfect for treating minor cuts, lacerations, abrasions, first and second degree burns, and skin irritations, Q Restore OTC has been evaluated in standard tests that show it can reduce the levels of common microorganisms including: Staphylococcus aureus, Pseudomonas aeruginosa, Escherichia coli, MRSA and VRE, as well as fungi such as Candida albicans and Aspergillus niger.*

Why Q Restore OTC?

For years colloidal silver was used as a disinfectant to ward off harmful microorganisms. Through new manufacturing technology with multiple patents, nano-silver solution has become the new standard by which all other silver products are measured. Q Sciences has harnessed this antibacterial strength with our patented 24 PPM Proprietary Silver Solution in Q Restore OTC. With the damage our skin incurs daily, Q Restore OTC ensures natural healing and protects against bacterial infections and illnesses.*

Features:

- Nano-silver technology
- Made with 24 PPM Proprietary Silver Solution
- Ideal to be used daily with both children and adults
- Scientifically tested

Directions: Clean wound. Use clean applicator to liberally cover the wound with Q Restore OTC and secondary dressing.



KEY BENEFITS:

- Designed to promote natural healing*
- Safe for both children and adults*
- Protects against bacterial infections and illnesses*



*These statements have not been evaluated by the food and drug administration. This product is not intended to diagnose, treat, cure or prevent any disease.

Live Life eXponentially!

QSciences.com

©2014 Q Sciences. All Rights Reserved.



Q Restore OTC

Wound Dressing Gel

Product Information

Q Restore OTC has been evaluated in standard tests that show it can reduce the levels of common microorganisms including:

- Staphylococcus aureus
- Pseudomonas aeruginosa
- Escherichia coli
- MRSA and VRE
- Fungi such as Candida albicans and Aspergillus niger

CONTENTS:

Q Restore OTC is a water-based hydro gel containing: 24 PPM
Proprietary Silver Solution (purified water, nano-silver at .01 micron), TEA, Carbomer.

WARNINGS:

Store at room temperature. Keep out of reach of children. Consult a health care professional if you have any concerns or questions about the wound.